BONE MANAGEMENT®

CATALOGUE

www.meisinger.de

Meisinger since 1888

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...mit den MEISINGER
Bone Management® Systemen.

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Augmentation with precision...

...with the MEISINGER Bone Management® Systems. Quality that convinces.

The in cooperation with leading scientists and practitioners developed product line MEISINGER Bone Management® offers perfectly harmonized system solutions for the controlled optimization of the bony implant site – completely independent of the used implant system. Besides prosthetic aspects, the minimally invasive character of the surgical intervention in complicated indications plays an important role in this product line. Clearly structured and standardized procedures offer the possibility of an implant restoration to a greater number of surgeons. Due to controlled and gentle bone preparation, Bone Management® is a key factor for long-term success in oral implantology.

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Bone Spreading

Crest-Control	4
Split-Control Basic	5
Split-Control 12 mm	6
Split-Control 15 mm	7
Split-Control Professional	8-9
Split-Control Plus	10-11
Bone Scraper	12-13
Master-Pin-Control	14-15

Bone Transfer & Osteosynthesis

Transfer-Control Basic	16
Transfer-Control Basic 2	17
Transfer-Control	18-19
Transfer-Ring-Control I	20
Transfer-Ring-Control II	21
Transfer-Control Plus	22-23
Trephine Ejection Kit	24-25
Screw System TX	26
Screw System TX Professional	27
Trephine Basic Kit	28
Trephine Kit	29
Trephine Kit Long	30
Punch Basic Kit	31
Saw Basic Kit	31

Sinuslift

Crestal-Lift-Control Basic	32
Crestal-Lift-Control	33
External-Lift-Control	34
Internal-Lift-Control	35
Lift-Control Plus	36-37
Surgical Kit 1	38
Surgical Kit 2	38
Surgical Kit 3	39
Surgical Kit 4	39

Implant Site Preparation

myplant two Surgical Kit		40-41
Alveoplasty Surgical Kit		42
Drill-Stop-Control		43
Drill-Stop-Control BLT		44
Drill-Stop-Control myplant two		45
Guided-Drill-Stop-Control		46-47
3D-Navigation-Control		48
3D-Navigation-Control Professional		49
Implant Preparation Kit		50
Implant Preparation Kit with Stop		50
Periimplantitis Kit		51

Further Bone Management® Systems

Surgical Kit	52
Parodontology Kit	53

Bone Management® Accessories

Handinstrumente/Manual Instruments	54-55
Ratschen/Ratchets	56
Chirurgieinstrumente/Surgery Instruments	57-70
Edelstahl-Bohrerstände/Stainless Steel Bur Blocks	71-73

Bone Grafting Material

NanoBone®		74-77
-----------	---	-------

Membranes & Sutures

Pro-Fix™		78-81
Vitala®, Cytoplast™		82-90
Anhang/Affix		91-92

Benex[®]-Control

Root Extraction System | developed with Dr. Benno Syfrig

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Benex[®]-Control is outstandingly suitable for the controlled and safe extraction of dental roots and dental root fragments, so the conventional extraction method using forceps and leverage can be avoided. Thanks to the innovative construction of the extractor, the root can be easily removed in a particularly controlled manner without damaging bone or soft tissue. To guarantee an optimal vertical transmission of force, there must be a flat area providing a supporting surface to the extractor. The quadrant support was developed specifically to provide this optimal support while also protecting existing dentition and providing support for the extractor in case of edentulous circumstances.

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
Developed with
Dr. Benno Syfrig
 Luzern, Swiss

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Crest-Control is a Horizontal Bone Splitting System which allows for a predetermined and minimally invasive widening of the horizontally resorbed alveolar ridge in particular in the distal lower jaw. With the well-proven Horizontal Spreaders, the alveolar ridge can be expanded up to a maximum of 5 mm, allowing an implant of every common system to be inserted in the widened segment.

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Split-Control Basic

Horizontal Bone Spreading and Condensing System

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
Split-Control Basic is a system for gentle widening of the alveolar ridge and simultaneous lateral bone condensing in the case of a horizontal bone deficit. With the aid of screw-like condensing and spreading instruments (Spreaders), it is possible to achieve a controlled and gentle dilatation of horizontally resorbed bone (Bone Spreading). Further, cancellous bone is condensed due to the special geometry of the spreaders (Bone Condensing), so that the primary stability of the inserted implants is increased.

Split-Control Basic is an introductory system with reduced instrument set and spreaders with a length of 12 mm.

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
Split-Control 12 mm is a system for gentle widening of the alveolar ridge and simultaneous lateral bone condensing in the case of a horizontal bone deficit. With the aid of screw-like condensing and spreading instruments (Spreaders), it is possible to achieve a controlled and gentle dilatation of horizontally resorbed bone (Bone Spreading). Further, cancellous bone is condensed due to the special geometry of the spreaders (Bone Condensing), so that the primary stability of the inserted implants is increased.

Split-Control 12 mm provides an extensive range of instruments for a controlled dilatation with spreaders with a length of 12 mm.

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Split-Control 15 mm is a system for gentle widening of the alveolar ridge and simultaneous lateral bone condensing in the case of a horizontal bone deficit. With the aid of screw-like condensing and spreading instruments (Spreaders), it is possible to achieve a controlled and gentle dilatation of horizontally resorbed bone (Bone Spreading). Further, cancellous bone is condensed due to the special geometry of the spreaders (Bone Condensing), so that the primary stability of the inserted implants is increased.

Split-Control 15 mm provides an extensive range of instruments for a controlled dilatation with spreaders with a length of 15 mm.

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Split-Control Professional

Horizontal Bone Spreading and Condensing System | developed with Dr. Dr. Martin Bonsmann

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Split-Control Professional is a system for gentle widening of the alveolar ridge and simultaneous lateral bone condensing in the case of a horizontal bone deficit. With the aid of screw-like condensing and spreading instruments (Spreader), it is possible to achieve a controlled and gentle dilatation of horizontally resorbed bone (Bone Spreading). Further, cancellous bone is condensed due to the special geometry of the spreaders (Bone Condensing), so that the primary stability of the inserted implants is increased.

To be able to respond more precisely in specific indications, the Split-Control Professional contains a broad range of spreader sizes.

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Split-Control Plus

Combined Horizontal Bone Splitting, Spreading and Condensing System

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Split-Control Plus is a combination of the proved Split-Control and Crest-Control systems, so it offers an extensive range of instruments for the controlled bone spreading and condensing. It contains both specially designed screw-like condensing and spreading instruments (Spreader) and horizontal spreaders that allow for controlled and standardized spreading of horizontally resorbed bone. Due to the special geometry of the screw-like spreaders, cancellous bone is gently condensed while the Horizontal Spreaders allow for a dilatation up to a width of 5 mm, so the bone is optimally prepared for the insertion of every common implant.

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The use of autologous bone is the gold standard for bone augmentation. With MEISINGER Bone Scrapers, autologous bone chips can be easily and safely collected and integrated into the present bone defect. Due to the immediate mixing with blood, the gained bone material is ready for immediate use and can be placed on the defect. The curved shape of the bone chips creates a particularly large volume of augmentable material from bone and blood. Thanks to the completely detachable blade, the extracted bone can be easily removed and the individual parts can be safely processed after use. Another advantage: The exchangeable blades guarantee an optimal and constant cutting performance.

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
Master-Pin-Control

34x Pins


Pin System for Membrane Fixation | developed with Dr. Istvan Urban

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Master-Pin-Control is designed especially for the fixation of resorbable, non-resorbable and titanium reinforced membranes. Since the pins have an extra mini-thread, they are a hybrid between a screw and a pin. The thread increases the surface of the pins and thus offers much more stability. In addition, due to the thread and a unique head-design, the pins can be removed easily and safely which is especially important after a successful healing period. With their sharp tip and very stable shank, they can be inserted easily and precisely into dense cortical bone.

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

URBAN
REGENERATION INSTITUTE

Urban Regeneration Institute
www.implant.hu

Developed with
Dr. Istvan Urban
Budapest, Hungary

Transfer-Control Basic

Horizontal and Vertical Bone Replacing System

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
Transfer-Control Basic permits precise and standardized transplantation of bone cylinders for horizontal and vertical bone augmentation. The instruments are available in three diameters. They are perfectly adapted to each other so that the outer diameter of the ablative bur and the size of the wheel cutter correlates with the inner diameter of the trephine. So a press fit can be produced between the cylindrical bone and the prepared recipient site. The press fit can then be strengthened with a fixation screw, which is separately available (Screw System TX, Art.-No. BTX00 / Screw System TX Professional, Art.-No. BTXPR). Such precise fitting of bone cylinders lead to accelerated bone revitalization and wound healing. So a transplantable bone site is reached after only 3-4 months.

Transfer-Control Basic is an introductory system with reduced instrument set. It is suitable for the extraction of bone cylinders with diameters of 4.0 mm, 5.0 mm and 6.0 mm.

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Art.-No. CTRB2

Transfer-Control Basic 2 permits precise and standardized transplantation of bone cylinders for horizontal and vertical bone augmentation. The instruments are available in three diameters. They are adapted to each other so that the outer diameter of the ablative bur and the diameter of the wheel cutter correlates with the inner diameter of the trephine. So a press fit can be produced between the cylindrical bone and the prepared recipient site. The press fit can then be strengthened with a fixation screw, which is separately available (Screw System TX, Art.-No. BTX00 / Screw System TX Professional, Art.-No. BTXPR). Such precise fitting of bone cylinders lead to accelerated bone revitalization and wound-healing. So a transplantable bone site is reached after only 3-4 months.

Transfer-Control Basic 2 is an introductory system with reduced instrument set. It is suitable for the extraction of bone cylinders with diameters of 6.0 mm, 7.0 mm and 8.0 mm.

	Diamantierte Säge Diamond coated saw	Trepane Trepines	Radfräser Wheel cutters	Lagerfräser Ablative burs						
Fig.	231DC*	T229L	T229L	T229L	TC084	TC084	TC084	TC21X	TC21X	TC21X
Shank ¹	204	205	205	205	204	204	204	205	205	205
Size ²	100	060	070	080	006	007	008	060	070	080
Length mm	0.3	14.0	14.0	14.0	6.0	7.0	8.0	14.0	14.0	14.0
□	-	6.0	7.0	8.0	-	-	-	-	-	-
□*	10.0	7.0	8.0	9.0	10.0	10.0	10.0	6.0	7.0	8.0

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Transfer-Control permits precise and standardized transplantation of bone cylinders for horizontal and vertical bone augmentation. The instruments are available in five diameters. They are perfectly adapted to each other so that the outer diameter of the ablative bur and the size of the wheel cutter correlates with the inner diameter of the trephine. So a press fit can be produced between the cylindrical bone and the prepared recipient site. The press fit can then be strengthened with a fixation screw, which is separately available (Screw System TX, Art.-No. BTX00 / Screw System TX Professional, Art.-No. BTXPR). Such precise fitting of bone cylinders lead to accelerated bone revitalization and wound healing. So a transplantable bone site is reached after only 3-4 months.

Transfer-Control is suitable for the extraction of bone cylinders with diameters of 4.0 mm, 5.0 mm, 6.0 mm, 7.0 mm and 8.0 mm.

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The most important condition for a safe healing of a bone graft is a congruent and fresh contact surface of the implant area. The Transfer-Ring-Control I System allows for this precondition for the healing of vertical bone grafts in a simple and controlled manner. This system offers a variety of various instrument sizes and provides high flexibility for the extraction of bone rings. Depending on requirements, various sizes of Initial Bur Trephines, Trephines and Ablative Burs are available. With the aid of the additional tungsten carbide bur, sharp bone edges can be smoothed precisely.

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Transfer-Ring-Control II

One-Stage Bone Ring Augmentation System | Technology by Dr. Bernd Giesenhagen

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
Transfer-Ring-Control II is a system for the bone ring technique, developed by Dr. Giesenhagen. It allows for one-stage proceeding of augmentation and implantation in case of resorbed jaw bone, if implant insertion alone is not possible because of inadequate spatial conditions. Due to the perfectly matched instruments, there is created a congruent bone and implant site, which allows for the treatment of three-dimensional defects with the help of prefabricated, allogenic bone rings. Because of the optimal fit between jaw bone, bone ring and inserted implant and the possibility of the one-stage proceeding, high primary stability is accomplished and healing time is significantly reduced.

Art.-No. BBR00

© Dr. Bernd Giesenhagen

	Pilotbohrer Pilot bur	Hartmetall Kugelfräser Tungsten carbide round drill	Diamantierte Säge Diamond coated saw	Trepane mit Führungsstift Trephines with guiding pin	Lagerfräser Ablative burs		
Fig.	E1001	HM141A	231DC*	229FS	229FS	DD207	DD207
Shank ¹	206	206	204	205	205	204	204
Size ²	020	031	100	060	070	060	070
Length mm	18.0	-	0.3	12.0	12.0	3.2	3.2
	-	-	-	5.0	6.0	-	-
	2.0	3.1	10.0	6.0	7.0	6.0	7.0

¹ 204=RA, 205=RA L, 206=RA XL ² Largest working part diameter in 1/10 mm

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Technology by
Dr. Bernd Giesenhagen
Kassel, Germany

Transfer-Control Plus

Combined Bone Replacing System for Horizontal and Vertical Bone Augmentation

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
Transfer-Control Plus permits precise and standardized transplantation of bone cylinders and rings for horizontal and vertical bone augmentation. Different types of bone burs allow for high flexibility for the extraction of bone grafts. The instruments are perfectly adapted to each other often producing a press fit between the cylindrical bone or ring and the prepared recipient site. This press fit can then be strengthened with a fixation screw which is sold separately (Screw System TX, Art.-No. BTX00 / Screw System TX Professional, Art.-No. BTXPR). Such precise fitting of bone cylinders lead to accelerated bone revitalization and wound-healing. So a transplantable bone site is reached after only 3-4 months.

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Trephine Ejection Kit

Trephine System for Extraction of Cylindrical Bone Grafts | developed by Prof. Dr. Fouad Khoury

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The Trephine Ejection Kit combines several unique ideas, which makes this set something very special. Internally cooled, the trephines provide a safe extraction of bone cylinders. The internal cooling protects the bone tissue from damage due to overheating. The removable trephine working parts allow for a particularly easy extraction of the cut out bone cylinders. With the aid of special ejection instruments, which fit to the trephine diameters, the bone can be pushed out easily. Further, the kit contains initial trephines which provide the safe application of the trephines. The Trephine Ejection Kit allows for fast and safe working and provides high flexibility due to various instrument sizes.

Art.-No. BTE00


Bone core are harvested from the implant site minimal invasive with the pre- and main trephine in different diameter.


Bone core augmentation stabilized with microscrews.

© Prof. Dr. Fouad Khoury

Vorkörner Trepane
Initial bur trephines

Schaft
Shank

Trepane
Trephines


Fig.	230KH	230KH	230KH	230KH	229KH*	229KH	229KH	229KH	229KH
Shank ¹	204	204	204	204	204	900	900	900	900
Size ²	021	025	029	033	000	021	025	029	033
Length mm	2.0	2.0	2.0	2.0	20.5	15.5	15.5	15.5	15.5
	2.1	2.5	2.9	3.3	-	2.1	2.5	2.9	3.3
	3.1	3.5	3.9	4.3	-	3.1	3.5	3.9	4.3


¹ 204=RA, 900=Unmounted ² Largest working part diameter in 1/10 mm

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

Bone core can be easily removed from the 2 pieces trephine.

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Developed with
Prof. Dr. Fouad Khoury
Olsberg, Germany

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


The Osteosynthesis System Screw System TX enables the universal use of fixation for autologous bone cylinders, bone blocks and bone chips for the purpose of preimplantological augmentation. The Torx connection guarantees reliable accommodation of the screws in the driving tool and, at the same time, high tensile forces can be transferred. For each 0.9 mm screw available, the kit contains a larger diameter rescue screw in case the transplant does not allow high tensile fixation. The small diameter and the low head size of the screws allow for a comfortable application for patient and user.

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Schraubendreher TX
Screw driver TX

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The Osteosynthesis System Screw System TX Professional enables the universal use of fixation for autologous bone cylinders, bone blocks and bone chips for the purpose of preimplantological augmentation. The Torx connection guarantees reliable accommodation of the screws in the driving tool and, at the same time, high tensile forces can be transferred. For each 0.9 mm screw available, the kit contains a larger diameter rescue screw in case the transplant does not allow high tensile fixation. The small diameter and the low head size of the screws allow for a comfortable application for patient and user.

The Screw System TX Professional contains additional screws in the lengths of 13.0 mm and 16.0 mm.

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
The Trephine Basic Kit provides a selected set of high-quality trephine drills in different sizes. These tools are designed specifically for creating cylindrical bone grafts quickly and accurately and for removing defined bone sections. The trephine drills have a depth of 10 mm and have been developed specifically with mandibular surgery and implantology in mind. Thanks to the well-defined laser depth marking in 2 mm intervals, the drilling depth can accurately be monitored.

Art.-No. 7120



Fig.	229	229	229	229	229	229
Shank ¹	205	205	205	205	205	205
Size ²	020	040	060	080	100	120
Length mm	10.0	10.0	10.0	10.0	10.0	10.0
	2.0	4.0	6.0	8.0	10.0	12.0
	3.0	5.0	7.0	9.0	11.0	13.0

¹ 205=RA L ² Largest working part diameter in 1/10 mm  Innendurchmesser Internal diameter  Außendurchmesser External diameter

Trephine Kit

Trephine System for Extraction of Cylindrical Bone Grafts

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The Trephine Kit provides a broad set of high-quality trephine drills with fine size graduations. These tools are designed specifically for creating cylindrical bone grafts quickly and accurately and for removing defined bone sections. The trephine drills have a depth of 10 mm and have been developed specifically with mandibular surgery and implantology in mind. Thanks to the well-defined laser depth marking in 2 mm intervals, the drilling depth can accurately be monitored.

Art.-No. 7121

Trepane
Trepines

Fig.	229	229	229	229	229	229	229
Shank ¹	205	205	205	205	205	205	205
Size ²	020	025	030	035	040	045	050
Length mm	10.0	10.0	10.0	10.0	10.0	10.0	10.0
	2.0	2.5	3.0	3.5	4.0	4.5	5.0
	3.0	3.5	4.0	4.5	5.0	5.5	6.0

¹ 205=RA L ² Largest working part diameter in 1/10 mm

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
The Trephine Kit provides a broad set of high-quality trephine drills with fine size graduations. These tools are designed specifically for creating cylindrical bone grafts quickly and accurately and for removing defined bone sections. The trephine drills have a depth of 14 mm* and have been developed specifically with mandibular surgery and implantology in mind. Thanks to the well-defined laser depth marking in 2 mm intervals, the drilling depth can accurately be monitored.

Art.-No. 7122

Trepane Trepines									
Fig.	229	229	229L	229L	229L	229L	229L	229L	229L
Shank ¹	205	205	205	205	205	205	205	205	205
Size ²	020	025	030	035	040	045	050		
Length mm	10.0	10.0	14.0	14.0	14.0	14.0	14.0	14.0	14.0
	2.0	2.5	3.0	3.5	4.0	4.5	5.0	5.5	6.0
	3.0	3.5	4.0	4.5	5.0	5.5	6.0	6.5	7.0

¹ 205=RA L ² Largest working part diameter in 1/10 mm

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Punch Basic Kit

Biopsy Punch Set for Extraction of Tissue Grafts

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The Punch Basic Kit contains Biopsy Punches in five different diameters for the performance of tissue removal. They serve for the application in immediate implant placement, for taking mucosa grafts in the scope of soft tissue management or biopsy as well as for uncovering the gingiva after implantation. The Biopsy Punches create clean cut surfaces and care for minimal tissue traumatization.

Art.-No. 7140


Fig.	225	225	225	225	225
Shank ¹	204	204	204	204	204
Size ²	030	040	050	060	070
Length mm	6.0	6.0	6.0	6.0	6.0
□	3.0	4.0	5.0	6.0	7.0
□+	3.7	4.7	5.7	6.7	7.7

¹ 204=RA ² Largest working part diameter in 1/10 mm
 □ Innendurchmesser Internal diameter □+ Außendurchmesser External diameter

Saw Basic Kit

Selection of Various Diamond Coated Saws for Oral Surgery

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The Saw Basic Kit provides a selection of diamond saws in the three most important diameters. They are used for the application in bone spreading or splitting operations. Due to the various sizes, the set offers the clinician the possibility of having a proper selection of these saws during surgery.

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Crestal-Lift-Control Basic

Crestal Approach Sinus Lift System

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
Crestal-Lift-Control Basic is a system for the performance of an easy and safe internal sinus lift. The elevation of the Schneiderian membrane occurs during the transcrestal drilling process. The specially atraumatic design of the Crestal Drill with its four cutting edges and concave head allows for a safe forming of a conical bone flap and is more ideally suited for collecting bone chips. Additionally, the stop sleeve system, which is tailored to the special instruments, prevents the membrane from being injured or punctured.

Crestal-Lift-Control Basic is an introductory system with reduced instrument set.

© Dr. Dr. Florian Bauer


Fig.	CL001	CL004	CL005	CL006	CL007	CL020
Shank ¹	206	206	206	206	206	206
Size ²	020	033	036	038	041	030
Length mm	16.0	16.0	16.0	16.0	16.0	18.5
	-	-	-	-	-	1.5
	2.0	3.3	3.6	3.8	4.1	3.0

¹ 206=RA XL

² Largest working part diameter in 1/10 mm

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Crestal-Lift-Control is a system for the performance of an easy and safe internal sinus lift. The elevation of the Schneiderian membrane occurs during the transcristal drilling process. The specially atraumatic design of the Crestal Drill with its four cutting edges and concave head allows for a safe forming of a conical bone flap and is more ideally suited for collecting bone chips. Additionally, the stop sleeve system, which is tailored to the special instruments, prevents the membrane from being injured or punctured.

Crestal-Lift-Control provides an extensive range of instruments for the performance of a safe internal sinus lift.

Art.-No. BCL00

© Dr. Dr. Florian Bauer

	Werkkörner (2 Stck.) Initial burs (2 pcs.)	Diamantinstrument Diamond instrument	Lagerfräser Ablative bur	Pilotbohrer Pilot bur	Crestal-Bohrer Crestal drills					
Fig.	186RF	859*	TC21X	CL001	CL002	CL003	CL004	CL005	CL006	CL007
Shank ¹	204	204	205	206	206	206	206	206	206	206
Size ²	018	018	040	020	028	031	033	036	038	041
Length mm	12.0	10.0	14.0	16.0	16.0	16.0	16.0	16.0	16.0	16.0
⚠	-	-	-	-	-	-	-	-	-	-
⚠	1.8	1.8	4.0	2.0	2.8	3.1	3.3	3.6	3.8	4.1

¹ 204=RA, 205=RA L, 206=RA XL

² Largest working part diameter in 1/10 mm

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External-Lift-Control is a system for the easy performance of an external sinus lift. After lateral access through the facial wall of the maxillary sinus, the Schneiderian membrane is detached using the various manual instruments and lifted cranially from the floor of the sinus. Then, a suitable augmentation material can be inserted under the membrane. Due to the special CARBOCER®-coating, the hand-held instruments provide a reflection-free working environment.

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
Internal-Lift-Control is a system for the performance of an internal sinus lift. With the aid of the elevators, the bone inside the implant site is displaced in the direction of the maxillary sinus. This procedure is suitable for the augmentation of small defects with a residual bone height of more than 6 mm.

Art.-No. CSL00

	Vorkörner Initial bur	Pilotbohrer Pilot burs		Erweiterer Expansion burs				
Fig.	186RF	B2001	G2001	B2004	E2004	F2004	G2004	H2004
Shank ¹	204	205	205	205	205	205	205	205
Size ²	018	016	034	027	034	039	045	049
Length mm	12.0	15.0	11.0	12.0	12.0	12.0	11.0	11.0
	-	-	-	1.17	2.44	2.96	3.55	4.07
	1.8	1.6	3.4	2.7	3.4	3.9	4.5	4.9

¹ 204=RA, 205=RA L ² Largest working part diameter in 1/10 mm

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Length mm 190.0

Lift-Control Plus

Combined System for Internal and External Sinus Lift

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The combined Lift-Control Plus system provides perfectly matched instruments for the performance of internal as well as external sinus floor elevation. For the internal sinus lift, the bone inside the implant site is displaced in the direction of the maxillary sinus. This procedure is suitable for the augmentation of small defects with a residual bone height of more than 6 mm. After lateral access through the facial wall of the maxillary sinus for the performance of the external sinus lift, the Schneiderian membrane is detached using the various manual instruments and lifted cranially from the floor of the sinus. After this, a suitable augmentation material can be inserted under the membrane.

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Surgical Kit 1

Lateral Approach System for External Sinus Elevation

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The Surgical kit 1 contains instruments with handpiece shank for the lateral approach in an external sinus lift. Carbide cutters in three diameters and diamond instruments in three grains allow for the gradual thinning of the facial bone wall of the maxillary sinus down to the Schneiderian membrane. The large diameter of the instruments protects the membrane from unintentional perforation.

Art.-No. BSKSL

Fig.	Hartmetallinstrumente Tungsten carbide burs				Diamantinstrumente Diamond instruments		
	HM141F	HM141F	HM141F	HM161RX	801*	801G*	801H*
Shank ¹	104	104	104	104	104	104	104
Size ²	018	031	050	018	050	050	050
Length mm	-	-	-	10.5	-	-	-
	-	-	-	-	-	-	-
	1.8	3.1	5.0	1.8	5.0	5.0	5.0

¹ 104=HP ² Largest working part diameter in 1/10 mm

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The Surgical Kit 2 contains instruments with right-angled shank for the lateral approach in an external sinus lift. Carbide cutters in three diameters and diamond instruments in three grains allow for the gradual thinning of the facial bone wall of the maxillary sinus down to the Schneiderian membrane. Especially the large diameter of the instruments protects the membrane from unintentional perforation.

Art.-No. BSK02

Fig.	Hartmetallinstrumente Tungsten carbide burs				Diamantinstrumente Diamond instruments		
	HM141F	HM141F	HM141F	HM161RX	801*	801G*	801H*
Shank ¹	205	205	205	205	204	204	204
Size ²	018	031	050	018	050	050	050
Length mm	-	-	-	10.5	-	-	-
	-	-	-	-	-	-	-
	1.8	3.1	5.0	1.8	5.0	5.0	5.0

¹ 204=RA, 205=RA L ² Largest working part diameter in 1/10 mm

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The Surgical Kit 3 contains conical Tungsten Carbide Burs with friction grip shank of various sizes and toothing which allow for a broad field of application. The particularly easy cutting and efficient cutters have an extra-long shank for an optimal view on the working area and sharp tips for the safe application in bone.

Art.-No. BSK03

Fig.	HM151	HM161	HM162	HM162SX	HM254LE	HM408M
Shank ¹	316	316	316	314	314	316
Size ²	016	018	016	014	012	016
Length mm	10.0	10.5	10.5	8.0	6.0	9.5
	-	-	-	-	-	-
	1.6	1.8	1.6	1.4	1.2	1.6

¹ 314=FG, 316=FG XL ² Largest working part diameter in 1/10 mm

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The Surgical Kit 4 contains a selection of conical Tungsten Carbide Cutters with handpiece shank for a broad field of application in oral surgery. The cutters have a particularly easy cutting and efficient saw toothing and a specifically designed twist angle. Further, the narrow working parts allow for an optimal view on the working area.

Art.-No. BSK04

Fig.	HM162A	HM162SL	HM163A	HM166A	HM254	HM408M
Shank ¹	104	104	104	104	104	104
Size ²	016	014	014	021	012	016
Length mm	9.3	8.0	5.0	10.5	6.0	9.6
	-	-	-	-	-	-
	1.6	1.4	1.4	2.1	1.2	1.6

¹ 104=HP ² Largest working part diameter in 1/10 mm

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myplant two Surgical Kit

Basic Surgical Kit for the implant diameter 4.0 mm

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The myplant two Surgical Kit Ø 4.0 mm is a reduced surgical kit that contains all the necessary instruments for setting the myplant two implants with a diameter of 4.0 mm. The principle of Enhanced Primary Stability (EPS) offers increased primary stability and thus the extended option of immediate restoration of the inserted implants. The included instruments are perfectly matched to each other and thus enable an exact and tissue-preserving preparation of the implant site. The optimal and clear arrangement of the instruments makes the kit easy and intuitive to use. All instruments are available in the myplant two Surgical Kit 0CK01.

CAUTION: To prevent injury of any anatomical structure, the apical extra lengths of the instruments must be considered.

© Dr. Mario Kirste


Fig.	OIB01	OSB01	ATS01	ATS03	MTS01	MTS03
Shank ¹	204	204	204	204	204	204
Size ²	015	020	024	024	029	029
Length mm	14.0	20.0	8.0	20.0	8.0	20.0
Apical extra lengths	-	-	0.4	0.4	0.5	0.5
+ □ +	1.5	2.0	2.4	2.4	2.9	2.9

¹ 204=RA ² Largest working part diameter in 1/10 mm

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
Other available implant systems	Ø 3.5 mm	X	Ø 4.5 mm
myplant two	Ø 3.5 mm	Ø 4.0 mm	Ø 4.5 mm

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© Prof. Dr. med. dent. Georg-Hubertus Nentwig

Developed with
Prof. Dr. G.-H. Nentwig
 Frankfurt, Germany

Alveoplasty Surgical Kit

Surgical Contouring and Modification Kit

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Art.-No. ALV18

The Alveoplasty Surgical Kit contains the most important instruments for surgical contouring and modification of the alveolar ridge when unwanted unevenness of the jawbone occurs after tooth extraction or bone degeneration. In this case, the kit offers the best way to smooth the jawbone. The Alveoplasty Kit contains two initial burs and two Lindemann burs, which allow for the performance of precise, initial osteotomies, even in hard cortical bone. To perform the alveoplasty itself, five different carbide cutters and a round diamond are included. The round working parts optimally protect the soft tissue against injuries. Thus, the alveolar ridge is prepared carefully and precisely for implant insertion or denture restoration.

	Vorkörner (2 Stck.) Initial burs (2 pcs.)	Chirurgiefräser Surgical Cutters	Hartmetallinstrumente Tungsten Carbide Instruments				Diamantinstrument Diamond Instrument		
Fig.	186RF	HM166	165RF	HM71	HM71	HM72GX	HM77GX	HM251GX	801*
Shank ¹	204	205	205	104	104	104	104	104	204
Size ²	018	021	023	040	050	060	060	040	033
Length mm	12.0	11.0	7.0	-	-	12.0	11.5	8.3	-
	1.8	2.1	2.3	4.0	5.0	6.0	6.0	4.0	3.3

¹ 204=RA, 205=RA L, 206=RA XL ² Largest working part diameter in 1/10 mm Außendurchmesser External diameter

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
The Drill-Stop-Control system provides an easy and safe preparation of the implant site. The increasing diameters of the twist drills are perfectly coordinated with each other. The stop sleeves provide the defined depth stop during the pilot drilling and during the expansion of the implant site. For the application without stop sleeves, the instruments are equipped with laser markings.

CAUTION: To prevent injury of any anatomical structure, the apical extra length of the instruments of 0.8 mm must be considered.

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
The Drill-Stop-Control BLT provides an easy and safe preparation of the implant site for Bone Level Tapered implants of the Straumann® or Dental Ratio® implant system. The increasing diameters of the twist drills are perfectly coordinated with each other. The stop sleeves provide the defined depth stop during the pilot drilling and during the expansion of the implant site. For the application without stop sleeves, the instruments are equipped with laser markings.

CAUTION: To prevent injury of any anatomical structure, the apical extra length of the instruments of 0.4 mm must be considered.

Fig.	187RF	TDS0T	TDS0T	TDS0T	TDS0T	TDS0T
Shank ¹	204	204	204	204	204	204
Size ²	018	022	028	033	035	042
Length mm	19.0	19.0	19.0	19.0	19.0	19.0
	1.8	2.2	2.8	3.25	3.5	4.2

¹ 204=RA ² Largest working part diameter in 1/10 mm Außendurchmesser External diameter

Fig.	CL066	CL065	CL064	CL063	CL062	CL061
Length mm	15.0	13.0	11.0	9.0	7.0	5.0
± max. drilling depth mm	4.0	6.0	8.0	10.0	12.0	14.0

Drill-Stop-Control myplant two

Drill System for Safe Implant Site Preparation

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With the specially for myplant two developed Drill-Stop-Control system, the preparation of the implant site is particularly easy and safe, since the diameters of the twist drills are perfectly matched to the myplant two implant system. The stop sleeves are used both during the pilot and the extension drilling for the defined depth stop. The lengths of the stop sleeves are ideal for subcrestal placement of the implant. For application without stop sleeves, the instruments are equipped with laser markings.

CAUTION: To prevent injury of any anatomical structure, the apical extra lengths of the instruments must be considered.

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Guided-Drill-Stop-Control

Guided Drill System for Safe Implant Site Preparation | developed with Dr. Henriette Lerner

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
The Guided-Drill-Stop system is used for the easy and safe preparation of the implant site. The special feature of the system are the specially developed twist drills with guiding tip, with which an intermediate step is added prior to the actual extension drilling. With this step, the upper part of the implant site is expanded to the desired diameter, while the lower part is used to guide the drill. This avoids deviation from the desired axis. In addition, the stop sleeves provide the defined depth stop during all drillings.

CAUTION: To prevent injury of any anatomical structure, the apical extra length of the instruments of 0.8 mm must be considered.

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
HL-DENTCLINIC
www.hl-dentclinic.de

Developed with
Dr. Henriette Lerner

Baden-Baden, Germany


3D-Navigation-Control

Guided Sleeve System | developed with Dr. Hans-Joachim Nickenig

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
3D-Navigation-Control is a guiding sleeve system that is suitable for all commonly used implant systems and planning programs. Due to the included closed standard sleeves and the open guiding sleeves for the molar region, the system permits the use of surgical guide templates in the anterior region and even in tight molar situations. In addition, the system provides special tool requirements of flapless surgery. The 3D-Navigation-Control system contains a separate 3D-Laboratory Kit so the dental technician is able to install the sleeves in the template precisely.

Anwendungsvideo
Application Video


Art.-No. BNA00

© Dr. Hans-Joachim Nickenig

Inhalt | content: 3D-Navigation-Control & 3D-Navigation-Control Professional

Fig.	HN010	HN011	HN012	HN013	HN014	HN015	HN016	HN017	HN018
Shank ¹	205	206	206	206	206	206	206	206	206
Size ²	020	020	025	030	035	040	020	025	030
Length mm	16.0	19.0	19.0	19.0	19.0	19.0	22.0	22.0	22.0
	2.0	2.0	2.5	3.0	3.5	4.0	2.0	2.5	3.0


¹ 205= RA L, 206=RA XL ² Largest working part diameter in 1/10 mm

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3D-Navigation-Control Professional

Guided Sleeve System | developed with Dr. Hans-Joachim Nickenig

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3D-Navigation-Control Professional is a guiding sleeve system that is suitable for all commonly used implant systems and planning programs. Due to the included closed standard sleeves and the open guiding sleeves for the molar region, the system permits the use of surgical guide templates in the anterior region and even in tight molar situations. In addition, the system provides special tool requirements of flapless surgery. The 3D-Navigation-Control Professional system contains a separate 3D-Laboratory Kit so the dental technician is able to install the sleeves in the template precisely.

The Professional system offers an additional, separate 3D-Fixation Kit, which allows for the (optional) fixation of the surgical guide template.

© Dr. Hans-Joachim Nickenig

Inhalt | content: 3D-Navigation-Control & 3D-Navigation-Control Professional

	Standard Führungshülsen, geschlossen (je 4 Stk.) Standard guiding sleeves, closed (each 4 pcs.)					Molare Führungshülsen, offen (je 4 Stk.) Molar guiding sleeves, open (each 4 pcs.)				Schablonenbohrer Template drills	
Fig.	HN001	HN002	HN003	HN004	HN005	HN006	HN007	HN008	HN009	HN021	HN022
Shank ¹	-	-	-	-	-	-	-	-	-	205	205
Size ²	-	-	-	-	-	-	-	-	-	045/055	060/065
Length mm	6.0	6.0	6.0	6.0	6.0	5.5	5.5	5.5	6.0	13.0	13.0
□	2.0	2.5	3.0	3.5	4.0	2.0	3.0	4.0	5.0	-	-
+□+	4.0	4.0	4.0	4.0	4.5	5.0	5.0	5.0	6.0	5.5	6.5
△	-	-	-	-	-	-	-	-	-	4.5	6.0

¹ 205=RA L ² Largest working part diameter in 1/10 mm △ Minimaldurchmesser Minimal diameter □ Innendurchmesser Internal diameter +□+ Außendurchmesser External diameter

Zusätzlicher Inhalt | Additionally content: 3D-Navigation-Control Professional

	Kugelhauptanker (je 10 Stück) Ball-head anchors (each 10 pcs.)		Schraubendreher Screw driver		Spiralbohrer Twist drill
Fig.	HN023*	HN024*	TC0SW	UG00W	HN019
Shank ¹	-	-	204	-	205
Size ²	-	-	-	-	010
Length mm	8.0	11.0	13.0	12.0	12.0
△	1.0	1.0	-	-	-
+□+	1.3	1.3	-	-	1.0

¹ 204=RA, 205=RA L ² Largest working part diameter in 1/10 mm

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Implant Preparation Kit

Initial Bur System for Safe Implant Site Preparation

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The Implant Preparation Kit provides the easy and safe preparation of the implant site. The very sharp initial drills allow for the safe marking of the implant site and prevent all following instruments from slipping. The extra slim pilot bur provides a particularly gentle first drilling for the determination of the implant length and orientation.

Art.-No. IPK01


Fig.	186RF	A1001
Shank ¹	204	206
Size ²	018	010
Length mm	12.0	18.0
□	1.8	1.0

¹ 204=RA, 206=RA XL ² Largest working part diameter in 1/10 mm

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Implant Preparation Kit with Stop

Initial Bur System with Stop for Safe Implant Site Preparation

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The Implant Preparation Kit provides the easy and safe preparation of the implant site. The stop sleeves allow for a defined depth stop during the initial and pilot drilling. For the application without stop sleeves, the instruments are equipped with laser markings.

CAUTION: To prevent injury of any anatomical structure, the apical extra length of the instruments of 0.8 mm must be considered.

Art.-No. IPK02


Fig.	187RF	TDS15
Shank ¹	204	204
Size ²	018	020
Length mm	19.0	19.0
□	1.8	2.0

¹ 204=RA ² Largest working part diameter in 1/10 mm

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The Periimplantitis Kit contains eight perfectly matched tungsten carbide finishing burs in egg and flame shapes for the idea peri-implantitis treatment. Both shapes come in two sizes and each in standard and ultra-fine toothings. The shapes, sizes and toothings of the working parts are perfectly suited for the intraoral treatment of titanium and for different shapes of implant neck and shoulder. All finishers have an extra-long FG shank (total length 32 mm) which allows particularly the treatment of difficult-to-access, deep areas. For the peri-implantitis treatment, the finishers are used clockwise and are guided around the exposed implant counter-clockwise so the implant surface is smoothed.

Art.-No. 2575


© Dr. Georg Boch

Hartmetallfinierer
Tungsten carbide finishing burs

Fig.	HM379	HM379	HM48L	HM48L	HM379U	HM379U	HM48LU	HM48LU
Shank ¹	317	317	317	317	317	317	317	317
Size ²	014	023	014	023	014	023	014	023
Length mm	3.1	4.2	8.0	8.0	3.1	4.2	8.0	8.0
□	1.4	2.3	1.4	2.3	1.4	2.3	1.4	2.3

¹ 317=FG XXL ² Largest working part diameter in 1/10 mm

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
Surgical Kit

Selection of the Most Popular Manual Instruments for Oral Surgery

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The Surgical Kit provides a selection of the most important manual instruments in oral surgery and thus it is the foundation for all surgical procedures.

Art.-No. BCS00

Handinstrumente
Manual instruments

Fig.	HI301	HI302	HI303	HI304	HI305*	HI306	HI307
Name	Mikro Schere gebogen Micro Scissors curved	Mikro Nadelhalter gerade Micro Needle Holder straight	Mikro-chirur. Pinzette gerade Micro Surg. Forceps straight	Mikro-chirur. Pinzette gebogen Micro Surg. Forceps curved	Skalpellgriff Scapel Handle	Periotom Periotome	Parodontalfeile Hirschfeld Periodontal File Hirschfeld
Length mm	174.0	180.0	180.0	180.0	150.0	176.0	176.0

* Für Skalpellklingen der Größe: 10, 11, 12, 13, 14, 15, 16
* For scalpel blades of the sizes: 10, 11, 12, 13, 14, 15, 16

Handinstrumente
Manual instruments

Fig.	HI308	HI309	HI310	HI311	HI312	HI313
Name	Molt Raspatorium fein Molt Raspatorium fine	Freer Raspatorium scharf / stumpf Freer Raspatorium sharp / dulled	Mikro Meißel Micro Chisel	Ansaugkanüle 3.0 mm Suction Cannula 3.0 mm	Spiegelgriff Mirror Handle	Lukas Löffel scharf Lukas Spoon sharp
Length mm	185.5	185.5	180.0	166.0	126.0	173.0

Parodontology Kit

Selection of the Most Popular Manual Instruments for Periodontal Treatment

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The Parodontology Kit provides a selection of the most important manual instruments for periodontal treatment. The kit contains measuring probes for the evaluation of the periodontal chart and a selection of curettes for an optimum cleaning and concrement removal at anterior teeth, premolars and molars.

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

This ratchet was designed to transfer turning moments up to 70 Ncm. May be disassembled for simple and safe cleaning.

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Torque ratchet for working with defined torque. In blocking position, torques for both the insertion and the loosening of screws can be transmitted, e.g. for the insertion of implants.

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MEISINGER offers one of the most extensive programs in the area of the Surgical Round Drills, Surgical Cutters, and Trephines. The optimized special construction of the instruments guarantee a safe application and fast achievement for the desired working result.

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Surgical Cutters of Tungsten Carbide

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The special, easy-cutting and efficient saw toothing and the specifically selected angle of twist allow a wide field of application of the instruments within the field of oral surgery.

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Surgical round drill of tungsten carbide

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


Surgical round drill of stainless Steel

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6 Schneiden, Laserung am Schaft: Stainless | 6 cutting edges, lasered on shank: stainless

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
* im Saw Basic Kit 7150 enthalten * contained in the Saw Basic Kit 7150

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
Stainless Steel Bur Blocks, non-corrosive, sterilizable

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MEISINGER drill stands are ideally suited for the effective cleaning, sterilization and organization of the required instruments.

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NanoBone® contains absolutely no human or animal ingredients so no material-related contamination risks and alleviates ethical concerns for patient consent. High-tech production processes ensure consistent and excellent quality.

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Thanks to its special structure, **NanoBone®** can be used for natural remodelling. Osteoclasts resorb the material as osteoblasts form new bone. Material resorption and bone formation are linked. This means predictable results and a reliable bone basis.

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NanoBone® is osteoconductive and osteoinductive. In other words, it is not simply a scaffold: it actively promotes bone formation.

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Fast and easy application with NanoBone®

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NanoBone® | granulate is available in either fine or coarse granules to suit the indication.

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The granulate can be mixed with the patient's blood or with a sterile saline solution. The manufacturer recommends mixing with blood. The material is extremely hydrophilic and fully absorbs the blood or saline solution quickly.

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When mixed with blood, **NanoBone® | granulate** takes on a paste-like consistency and can be easily applied with a spatula or augmentation spoon.

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NanoBone® | block was developed in partnership with users as an alternative to the autologous bone block. The block comes in a set with two screws for attachment.

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The block can be shaped with rotating or scraping tools to adapt it to the bone. The specific surgery instructions are included with the set.

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Ideal bone healing in combination with perfect handling – NanoBone® QD

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New applicator


- "Ready to use" direct application from the applicator
- Easy, intuitive use
- Precise placement of the **NanoBone® QD** material
- Controlled and exact portioning
- Allows for small approaches to the application site (outer diameter of the applicator just 8 mm)
- Allows for retrograde filling of bone defects
- Difficult to access bone defects treated due to the long applicator with extended reach

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Improved material

- Shapeable and adaptable to the defect
- Strong cohesion of the material
- Sticks to the bone (adhesion)
- Stable also in aqueous environment and profusely bleeding wounds
- Controlled osteoinduction

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
The Pro-Fix™ Precision Fixation System is manufactured to precise tolerances to ensure easy pick-up of screws, stable transfer to the surgical site and quick engagement in cortical bone. All kit components are labeled, organized, stored and sterilized together to simplify things for the surgeon and surgical team.

Tecapro™ Storage Tray, Tecapro™ Storage Tray

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The autoclavable Tecapro™ Storage Tray is cleverly designed to store not only fixation instruments, but also a variety of membrane fixation, bone fixation, and tenting screws needed in dental bone grafting.

Schrauben-Organizer, Screw Organizer Dial

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Up to 100 fixation screws – including bone fixation, membrane fixation, and tenting – may be stored in a built-in screw organizer dial designed and labeled for easy identification, easy storage and simple re-ordering.

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The stainless steel driver handle and cruciform driver blade fit together to make up a drive system that ensures easy pick-up of screws and stable transfer to the surgical site.

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Blades work universally with all Pro-Fix™ screws.

Art-No.	Inhalt/Content	
PFT	Autoklavierbares Tecapro™ Storage Tray mit Schrauben-Organizer Autoclavable Tecapro™ Storage Tray with screw organizer dial	
PFDH	Schraubendreher-Griff aus rostfreiem Stahl Stainless steel driver handle	
PFDB	76 mm Kreuzschraubendreher 76 mm cruciform driver blade	
PFDB56	56 mm Kreuzschraubendreher 56 mm cruciform driver blade	
PFPD	1,2 mm Durchmesser RA-Pilotbohrer 1.2 mm diameter latch type pilot drill	
PFDBCA	RA-Schraubendreher (24 mm lang, 10 mm exponierte distale Länge) Contra angle blade (24 mm long, 10 mm exposed distal length)	

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Pro-Fix™ Membrane Fixation Screws are designed as an attractive alternative to using tacks for membrane stabilization. Easy pick-up, solid stability of the screw during transfer to the surgical site, and easy placement make membrane fixation fast and easy.

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Further self-drilling Membrane Fixation Screws

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Pro-Fix™ Tenting Screws are designed with a self-drilling tip, polished neck and broader head to maintain space under resorbable and non-resorbable membranes in horizontal and vertical bone regeneration procedures.

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Further self-drilling Tenting Screws

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Further self-tapping Bone Fixation Screws

Art-No.	verpackt à packed à	Inhalt Content	
PFB8	1	1,5 mm x 8,0 mm selbstschneidende Osteosyntheseschraube 1.5 mm x 8.0 mm self-tapping Bone Fixation Screw	
PFB8-5	5	1,5 mm x 8,0 mm selbstschneidende Osteosyntheseschrauben 1.5 mm x 8.0 mm self-tapping Bone Fixation Screws	
PFB10	1	1,5 mm x 10,0 mm selbstschneidende Osteosyntheseschraube 1.5 mm x 10.0 mm self-tapping Bone Fixation Screw	
PFB10-5	5	1,5 mm x 10,0 mm selbstschneidende Osteosyntheseschrauben 1.5 mm x 10.0 mm self-tapping Bone Fixation Screws	
PFB12	1	1,5 mm x 12,0 mm selbstschneidende Osteosyntheseschraube 1.5 mm x 12.0 mm self-tapping Bone Fixation Screw	
PFB12-5	5	1,5 mm x 12,0 mm selbstschneidende Osteosyntheseschrauben 1.5 mm x 12.0 mm self-tapping Bone Fixation Screws	
PFB14	1	1,5 mm x 14,0 mm selbstschneidende Osteosyntheseschraube 1.5 mm x 14.0 mm self-tapping Bone Fixation Screw	
PFB14-5	5	1,5 mm x 14,0 mm selbstschneidende Osteosyntheseschrauben 1.5 mm x 14.0 mm self-tapping Bone Fixation Screws	

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
Vitala® is a natural collagen membrane made from porcine pericardium. The special manufacturing process ensures that the microporous, 3-layered architecture is preserved and no cross-linking chemicals are needed. This extraordinary structure leads to a very flexible membrane with particularly high tensile strength.

The Vitala® collagen membrane is characterized by its smooth handling, high flexibility and ideal adaptability. Once hydrated the Vitala® membrane loses its original shape and becomes flexible, so that it optimally adapts to all defect surfaces. As the surface of the membrane is identical on both sides, both sides can be placed against the defect.

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
Cytoplast™ RTM membranes are made from highly purified type 1 collagen derived from bovine achilles tendon. A great advantage of the RTM membranes is their ideal handling. They combine high flexibility with optimum stability. This makes them flexible, adaptable to any shape and easy to place. At the same time they are stable enough to hold their shape and not to collapse into the defect. The multi-layer structure of the membranes allows tissue integration into the outer layer and prevents the direct passage of bacteria and epithelial cells. The unique fiber orientation ensures maximum tensile strength. Thus, the membranes can be safely fixed by sutures or pins without tearing. Due to the long and predictable resorption time of 26 to 38 weeks, there is no risk of particle loss due to early resorption.

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Cytoplast™
Catalogue separately
available

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
The Cytoplast™ Technique: Cosmetic Socket Preservation Without Primary Closure

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
Cytoplast™ TXT-200 and TXT-200 singles are made of high-density PTFE with a pore size of less than 0.3 µm. Thus, the non-absorbable membrane is insensitive to bacteria and allows a maximum of predictability and aesthetics. When using the TXT-200 membranes, no primary wound closure is necessary. The exposed PTFE membrane can be removed with tweezers after at least 21 days without surgical intervention, no additional anesthesia is required. The hexagonal depressions increase the surface, which leads to an increased stability of the membrane and favors cell adhesion without increasing the porosity. For best results, the side with the depressions should be placed in the direction of the soft tissue. The TXT-200 membranes can be cut to size as needed.

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The Cytoplast™ Ti-150 and Ti-250 Titanium Reinforced membranes are made of high density PTFE. In combination with titanium reinforcement, the membrane offers exceptional dimensional stability and optimum flexibility at the same time. Therefore, the membrane is particularly suitable for the treatment of 3- to 4-wall defects in the course of the alveolar ridge augmentation. The wide variety of membrane shapes and sizes covers a variety of different defects and minimizes the need for trimming. If necessary, however, the membranes can be optimally trimmed to the required shape. The titanium-reinforced membranes are available as Cytoplast™ Ti-150 or Ti-250. The Ti-150 membranes are 40 % thinner than the Ti-250 membranes.

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Designed for narrow single-tooth extraction sites, especially where one or more bony walls are missing.

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Designed to fit periodontal defects in the anterior.

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Designed for single-tooth extraction sites, especially where one or more bony walls are missing.

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Designed to fit periodontal defects in the posterior.

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
Designed for large bony defects, including ridge augmentation in the anterior maxilla.

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Designed for grafting posterior extraction sites and limited ridge augmentation. Designed for large buccal defects.

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Designed for bony defects between adjacent teeth, including ridge augmentation.

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



The high quality of the Cytoplast™ needles and the monofilament PTFE suture allows for atraumatic wound closure after dental procedures. The high-quality needles are made of stainless steel and have a very high bending strength with a minimum of wear. Sharp needles are the result, even after repeated piercing. The achieved, high constant cutting performance allows for a particularly tissue-saving work. The suture is made of 100% medical grade bioinert PTFE. Thanks to its smooth surface, the monofilament prevents the adhesion of bacteria to the thread and ensures a resistance-free gliding through the tissue. The PTFE is a soft, non-resorbable material with consistent tensile strength, allowing the user to work very precisely, planning with a defined healing time, and patient-friendly care.

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

General note: 12 Sutures / Box, Suture length: 45 cm

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


Quelle / Source: Osteogenics Biomedical, Inc.

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The following reference values for rotation speeds apply to surgery in general:

Hand piece (HP):

Optimum: 6.000 - 10.000 rpm
Maximum: 40.000 - 50.000 rpm

Angle piece (RA):

Optimum: 6.000 - 10.000 rpm
Maximum: 40.000 - 50.000 rpm

Friction Grip (FG):

Optimum: 80.000 rpm
Maximum: 100.000 - 120.000 rpm

Please also observe adapting the rotation speed in relation to the diameter of the instrument as well as the prevailing indication and enough cooling. As a general rule, the larger the working part of an instrument, the lower the speed should be set.

Please also pay attention to the special preparation specifications for products made of tool steel.

Please also note the recommended and maximum rotary speeds for the individual instruments on the product packaging and in the corresponding instructions for use.

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The symbols give merely suggestions for the possible implementation of the products. The user decides and takes full responsibility about the precise deployment according to existing indications. Please follow general application and safety instructions for MEISINGER products in the medical and dental area and also the advice for processing. Details can be found on the internet under www.meisinger.de or you can request one by mail.

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Cleaning/Disinfection/Sterilization

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Meisinger stands for high quality medical devices since 1888. The quality management system of a company which manufactures medical devices must meet specific special requirements. These extremely high requirements are defined in ISO 13485 and meticulously complied with by our company. A MDSAP certificate according to ISO 13485:2016 confirms compliance with the requirements of international authorities in the USA (FDA), Canada (Health Canada), Australia (TGR), Japan (MHLW) and Brazil (ANVISA). All medical devices which you purchase from us as customer, comply with all applicable requirements of the Medical Device Directive 93/42/EEC. Our company is certified by an independent Notified Body and certification is performed according to the specifications of standards. Current certificates can be found on our homepage www.meisinger.de

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Cytoplast™ Vitala® Pro-Fix™

 HERSTELLER MANUFACTURER

Osteogenics Biomedical, Inc.
4620 71st Street | Building 78-79
Lubbock, TX 79424 | USA
Phone: 1.888.796.1923
Fax: 806.796.0059
E-Mail: sales@osteogenics.com
Internet: www.osteogenics.com

Cytoplast™, Vitala® und Pro-Fix™ sind eingetragene Warenzeichen der Osteogenics Biomedical, Inc..
Cytoplast™, Vitala® and Pro-Fix™ are registered trademarks of Osteogenics Biomedical, Inc..

All rights reserved.

Reproduction, also by extract and reproducing photos, are only permitted with written authorization of Hager & Meisinger GmbH, Neuss.

All offers, orders, and deliveries are subject to Meisinger's "General Sales and Delivery Terms."

We reserve the right to modify our range of products and their design as well as to deviate from the illustrations and data shown.

registered trademarks

Hager & Meisinger GmbH | Hansemannstr. 10 | 41468 Neuss | Germany
Phone: +49 2131 2012-0 | Fax: +49 2131 2012-222 | E-Mail: info@meisinger.de | Internet: www.meisinger.de

Meisinger USA, L.L.C. | 10150 E. Easter Avenue | Centennial, Colorado 80112 | USA
Phone: +1 (303) 268-5400 | Toll free: +1 (866) 634-7464 | Fax: +1 (303) 268-5407 | E-Mail: info@meisingerusa.com | Internet: www.meisingerusa.com

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¹ Teflon® is a registered trademark of DuPont, USA
NanoBone® ist ein eingetragenes Warenzeichen der ARTOSS GmbH. NanoBone® is a registered trademark of ARTOSS GmbH.