

mis[®] | SEVEN[®]

A self-tapping internal hex. implant by *mis*

P.4-5

P.6-7

P.8

P.9

P.10

P.11

P.12

P.13

P.14

P.15

P.16-17

P.18-19

P.20-21

P.22

P.23

P.24-25

MIS Warranty:

MIS exercises great care and effort in maintaining superior quality products. All MIS products are warranted to be free from defects in material and workmanship. However, should a customer find fault with any MIS product while using it according to the instructions, the defective product will be replaced.

Warning: MIS products should be used by licensed dentists only.



Table of Contents.

Introduction

Advantages

Implant Range, Narrow Platform

Ø3.30mm Implants

Implant Range, Standard Platform

Ø3.75mm Implants

Ø4.20mm Implants

Implant Range, Wide Platform

Ø5mm Implants

Ø6mm Implants

Drilling Procedures

Surgical Kit

Success

Insertion Tools

Insertion Options

Packaging

MIS scientists and engineers are continually seeking superior materials, technologies and procedures with the aim of developing quality products designed to make implant dentistry effective, safe and simple. The SEVEN implant system is the result of an extensive research and development process, offering a truly innovative solution with a unique combination of surgical and restorative benefits.




© MIS Corporation. All rights reserved.



6.

Advantages.






Internal hex.

MIS SEVEN implants feature a color-coded internal hex. connection, assuring proper implant abutment seating, anti-rotational engagement, resistance to lateral forces, excellent esthetic results and more.




Micro-rings

Micro-rings on the implant neck improves BIC (Bone-To-Implant-Contact) at the crestal zone, and reduces pressure on the cortical bone to minimize resorption at the implant neck.




Dual thread

The dual thread doubles the implant insertion rate (2.40mm), allowing for a simpler and faster implant placement. The thread design also contributes to high initial stability properties.



Conical shape

Conical, root-shaped geometry and a unique thread design ensure superior primary stability, offering the ultimate choice for a wide range of clinical cases and loading protocols.



Surface

The surface roughness and micro-morphology is a result of sand-blasting and acid-etching. This MIS established surface technology has provided millions of patients and clinicians with excellent osseointegration results and long-lasting clinical success.




Three spiral channels

Three spiral channels at the apical end of the implant support the self-tapping properties. The channels also collect and amass bone chips in the course of insertion, supporting efficient osseointegration and long-term stability.



Domed apex

The dome-shaped apex prevents over-insertion for safer implant placement procedures.

| Length | 6mm | 8mm | 10mm | 11.50mm | 13mm | 16mm |
|----------|-----|-----|--|--|--|--|
| Ø3.30 mm | | | MF7-10330  | MF7-11330  | MF7-13330  | MF7-16330  |

Surgical Tools



MT-HSH10

Short insertion tool, internal hex.



MT-HLI10

Long insertion tool, internal hex.



MT-NRH10

Long ratchet adapter internal hex. narrow platform



MT-SRH20

Short insertion tool internal hex.



MT-LRH20

Long insertion tool internal hex.


Implant Cover Screw and Healing Caps



MC1-00330



MH-N2330
MH-N3330
MH-N4330
MH-N6330
MH-N8330



MH-52330
MH-53330
MH-54330
MH-55330
MH-56330

SEVEN®

NARROW

Ø3.30mm
Narrow Platform

Titanium Alloy Ti 6Al 4V ELI
Sand-Blasted and Acid-Etched


Final Drill

Each implant is supplied with a single-use final drill corresponding to the correct diameter and length, allowing a short and safe drilling procedure.

| Catalog No. | Dimensions | |
|-------------|---------------------------|--|
| MF7-10330 | Ø3.30mm length 10mm |  |
| MF7-11330 | Ø3.30mm length 11.50mm | |
| MF7-13330 | Ø3.30mm length 13mm | |
| MF7-16330 | Ø3.30mm length 16mm |  |

Ø3.30mm Implant Procedure



| | | | | | |
|----------------------|-----------|----------|---------|----------------|-------|
| Drilling Speed (RPM) | 1200-1500 | 900-1200 | 200-400 | 15-25 | |
| Diameter | Ø1.90 | Ø2.40 | Ø2.40 | Ø2.20 Ø3.20 | Ø3.30 |



Do not use the final drill for bone type 3&4

The drilling sequence is illustrated using a 13mm implant.

Procedure recommended by MIS cannot replace the judgment and professional experience of the surgeon.

| Length | 6mm | 8mm | 10mm | 11.50mm | 13mm | 16mm | |
|-----------|-----------|--|--|--|--|--|--|
| Ø 3.75 mm | | MF7-08375  | MF7-10375  | MF7-11375  | MF7-13375  | MF7-16375  | |
| | Ø 4.20 mm | MF7-06420  | MF7-08420  | MF7-10420  | MF7-11420  | MF7-13420  | MF7-16420  |

Surgical Tools



MT-HSI10

Short insertion tool, internal hex.



MT-HLI10

Long insertion tool, internal hex.



MT-SRH20

Short insertion tool internal hex.




MT-LRH20

Long insertion tool internal hex.

Implant Cover Screw and Healing Caps



MC1-00375



MH-03375
MH-04375
MH-05375
MH-06375



MH-A3375
MH-A4375
MH-A5375
MH-A6375



MH-53375
MH-54375
MH-55375
MH-56375

Titanium Alloy Ti 6Al 4V ELI
Sand-Blasted and Acid-Etched


| Catalog No. | Dimensions | |
|-------------|---------------------------|--|
| MF7-08375 | Ø3.75mm length 8mm | |
| MF7-10375 | Ø3.75mm length 10mm | |
| MF7-11375 | Ø3.75mm length 11.50mm | |
| MF7-13375 | Ø3.75mm length 13mm | |
| MF7-16375 | Ø3.75mm length 16mm | |

Final Drill

Each implant is supplied with a single-use final drill corresponding to the correct diameter and length, allowing a short and safe drilling procedure.

Ø3.75mm Implant Procedure

| Drilling Speed (RPM) | 1200-1500 | 900-1200 | 500-700 | 200-400 | 200-500 | 15-25 | | |
|----------------------|-----------|----------|---------|---------|---------|----------------|-------|-------|
| Diameter | Ø1.90 | Ø2.40 | Ø2.40 | Ø2.80 | Ø2.80 | Ø2.80 Ø3.60 | Ø3.75 | Ø3.75 |



Do not use the final drill for bone type 3&4

The drilling sequence is illustrated using a 13mm implant.

Procedure recommended by MIS cannot replace the judgment and professional experience of the surgeon.

| Catalog No. | Dimensions | |
|-------------|---------------------------|--|
| MF7-06420 | Ø4.20mm length 6mm | |
| MF7-08420 | Ø4.20mm length 8mm | |
| MF7-10420 | Ø4.20mm length 10mm | |
| MF7-11420 | Ø4.20mm length 11.50mm | |
| MF7-13420 | Ø4.20mm length 13mm | |
| MF7-16420 | Ø4.20mm length 16mm | |


Titanium Alloy Ti 6Al 4V ELI
Sand-Blasted and Acid-Etched

Final Drill

Each implant is supplied with a single-use final drill corresponding to the correct diameter and length, allowing a short and safe drilling procedure.

Ø4.20mm Implant Procedure

| | | | | | | | | | | |
|----------------------|-----------|----------|-------|---------|---------|-------|----------------|----------------|--|-------|
| Drilling Speed (RPM) | 1200-1500 | 900-1200 | | 500-700 | 400-700 | | 200-400 | 200-500 | | 15-25 |
| Diameter | Ø1.90 | Ø2.40 | Ø2.40 | Ø2.80 | Ø3.20 | Ø3.20 | Ø3.30 Ø4.10 | Ø3.30 Ø4.20 | | Ø4.20 |



Do not use the final drill for bone type 3&4

The drilling sequence is illustrated using a 13mm implant.

Procedure recommended by MIS cannot replace the judgment and professional experience of the surgeon.

| Length | 6mm | 8mm | 10mm | 11.50mm | 13mm | 16mm |
|--------|---|---|---|---|---|---|
| Type | MF7-06500 | MF7-08500 | MF7-10500 | MF7-11500 | MF7-13500 | MF7-16500 |
| Ø5 mm |  |  |  |  |  |  |
| Ø6 mm |  |  |  |  |  | |

Surgical Tools



MT-HSI10

Short insertion tool, internal hex.



MT-HLI10

Long insertion tool, internal hex.



MT-SRH20

Short insertion tool internal hex.



MT-LRH20

Long insertion tool internal hex.

Implant Cover Screw and Healing Caps



MC1-00470



MH-W3500
MH-W4500
MH-W5500



MH-W3630
MH-W4630
MH-W5630

| Catalog No. | Dimensions | |
|-------------|------------------------|--|
| MF7-06500 | Ø5mm length 6mm | |
| MF7-08500 | Ø5mm length 8mm | |
| MF7-10500 | Ø5mm length 10mm | |
| MF7-11500 | Ø5mm length 11.50mm | |
| MF7-13500 | Ø5mm length 13mm | |
| MF7-16500 | Ø5mm length 16mm | |


Titanium Alloy Ti 6Al 4V ELI
Sand-Blasted and Acid-Etched

Final Drill

Each implant is supplied with a single use final drill corresponding to the correct diameter and length, allowing a short and safe drilling procedure.

Ø5mm Implant Procedure



| | | | | | | | | | | |
|----------------------|-----------|----------|-------|---------|---------|---------|----|---------|---------|----|
| Drilling Speed (RPM) | 1200-1500 | 900-1200 | | 500-700 | 400-700 | 400-600 | | 200-400 | 200-500 | |
| Diameter | Ø1.90 | Ø2.40 | Ø2.40 | Ø2.80 | Ø3.20 | Ø4 | Ø4 | Ø4.10 | Ø4.90 | Ø5 |



Do not use the final drill for bone type 3&4

The drilling sequence is illustrated using a 13mm implant.

Procedure recommended by MIS cannot replace the judgment and professional experience of the surgeon.

| Catalog No. | Dimensions | |
|-------------|------------------------|--|
| MF7-06600 | Ø6mm length 6mm |  |
| MF7-08600 | Ø6mm length 8mm | |
| MF7-10600 | Ø6mm length 10mm | |
| MF7-11600 | Ø6mm length 11.50mm |  |
| MF7-13600 | Ø6mm length 13mm | |

Titanium Alloy Ti 6Al 4V ELI
Sand-Blasted and Acid-Etched

Final Drill


Each implant is supplied with a single use final drill corresponding to the correct diameter and length, allowing a short and safe drilling procedure.


Ø6mm Implant Procedure

| | | | | | | | | | | | | | | |
|----------------------|-----------|----------|-------|---------|---------|---------|---------|---------|----|----|---------|---------|----|----|
| Drilling Speed (RPM) | 1200-1500 | 900-1200 | | 500-700 | 400-700 | 400-600 | 300-500 | 300-500 | | | 200-400 | 200-500 | | |
| Diameter | Ø1.90 | Ø2.40 | Ø2.40 | Ø2.80 | Ø3.20 | Ø4 | Ø4.50 | Ø5 | Ø5 | Ø5 | Ø5.10 | Ø5.90 | Ø6 | Ø6 |



Drilling Procedures.





Ø2.40

Ø2.40

Ø2.80

Ø2.80

Ø3.20

Ø3.20

Ø4

Ø4

Ø4.50

Ø5

Ø5



Final drill

Countersink

Body try-in Ø5mm

Twist drill Ø5mm

Twist drill Ø5mm

Body try-in Ø4.50mm

Twist drill Ø4.50mm

Body try-in Ø4mm

Twist drill Ø4mm

Body try-in Ø3.20mm

Twist drill Ø3.20mm

Body try-in Ø2.80mm

Twist drill Ø2.80mm


Body try-in Ø2.40mm

Pilot drill Ø2.40mm

Spade / Marking drill Ø1.90mm

SEVEN Surgical Kit.

MK-EI48 | With external irrigation drills



MARKING DRILLS



1 MT-TDN19
Marking drill Ø1.90mm
external irrigation



2 MT-SMD10
Spade marking drill

PILOT DRILLS



3 MT-P2406
Pilot drill with built in stopper Ø2.40
for 6mm length implants



4 MT-P2408
Pilot drill with built in stopper Ø2.40
for 8mm length implants



5 MT-P2410
Pilot drill with built in stopper Ø2.40
for 10mm length implants



6 MT-P2411
Pilot drill with built in stopper Ø2.40
for 11.5mm length implants



7 MT-P2413
Pilot drill with built in stopper Ø2.40
for 13mm length implants



8 MT-P2416
Pilot drill Ø2.40
for 13mm length implants

TWIST DRILLS



9 MT-TDT28
Twist drill Ø2.80mm external irrigation



10 MT-TDT32
Twist drill Ø3.20mm external irrigation




11 MT-TDT40
Twist drill Ø4mm external irrigation










12 MT-TDT45
Twist drill Ø4.50mm external irrigation









13 MT-TDT50
Twist drill Ø5mm external irrigation



SURGICAL TOOLS

- 28  MT-PP240
Parallel pin, Ø3/2.40mm
- 29  MT-DE001
Drill extender
- 30  MT-RT070
Surgical Torque Ratchet
- 31  MN-PF330
Direct press fit for closed tray, NP
- 32  MD-PF375
Direct press fit for closed tray, SP/WP
- 33  MT-RDL30
Long driver for 0.05 inch hex.
- 34  MT-RDS30
Short driver for 0.05 inch hex.







BODY TRY-INS

- 14  MT-BTT24
Body try-in Ø2.40mm
- 15  MT-BTT28
Body try-in, Ø2.80mm
- 16  MT-BTT32
Body try-in, Ø3.20mm
- 17  MT-BTT40
Body try-in, Ø4mm
- 18  MT-BTT45
Body try-in, Ø4.50mm
- 19  MT-BTT50
Body try-in, Ø5mm

COUNTERSINKS

- 20  MT-GDN33
Countersink, SP
- 21  MT-GDN50
Countersink, WP

INSERTION TOOLS

- 22  MT-HLI10
Long insertion tool, internal hex.
- 23  MT-HSI10
Short insertion tool, internal hex.
- 24  MT-LRH20
Long insertion tool, internal hex.
- 25  MT-SRH20
Short insertion tool, internal hex.
- 26  MT-NRH10
Long ratchet adapter, internal hex., NP
- 27  MT-RMR10
Long direct hand and ratchet key

Surface Quality.

*All MIS implants undergo the same surface treatments; sand-blasting and acid-etching. The research study was done on the SEVEN implant, however the results are valid for all MIS implant surfaces.

Identification Card and Codification of the Chemical and Morphological Characteristics of 62 Dental Implant Surfaces. Part 3: Sand-Blasted/Acid-Etched (SLA Type) and Related Surfaces (Group 2A, main subtractive process).

Background and Objectives

Dental implants are commonly used in dental therapeutics, but dental practitioners only have limited information about the characteristics of the implant materials they take the responsibility to place in their patients. The objective of this work is to describe the chemical and morphological characteristics of 62 implant surfaces available on the market and establish their respective Identification (ID) Card, following the Implant Surface Identification Standard (ISIS). In this third part, surfaces produced through the main subtractive process (sand-blasting/acid-etching, SLA-type and related) were investigated.

Materials and Methods

Eighteen different implant surfaces were characterized: Straumann SLA (ITI Straumann, Basel, Switzerland), Ankylos (Dentsply Friadent, Mannheim, Germany), Xive S (Dentsply Friadent, Mannheim, Germany), Frialit (Dentsply Friadent, Mannheim, Germany), Promote (Camlog, Basel, Switzerland), Dentium Superline (Dentium Co., Seoul, Korea), Osstem SA (Osstem implant Co., Busan, Korea), Genesis (GC Corporation, Tokyo, Japan), Aadvia (GC Corporation, Tokyo, Japan), MIS Seven (MIS Implants Technologies, Bar Lev, Israel), ActivFluor (Blue Sky Bio, Grayslake, IL, USA), Tekka SA2 (Tekka, Brignais, France), Twinkon Ref (Tekka,

Brignais, France), Bredent OCS blueSKY (Bredent Medical, Senden, Germany), Magitech MS2010 (Magitech M2I, Levallois-Perret, France), EVL Plus (SERF, Decines, France), Alpha Bio (Alpha Bio Tec Ltd, Petach Tikva, Israel), Neoporos (Neodent, Curitiba, Brazil). Three samples of each implant were analyzed.

Superficial chemical composition was analyzed using XPS/ESCA (X-Ray Photoelectron Spectroscopy/Electron Spectroscopy for Chemical Analysis) and the 100nm in-depth profile was established using Auger Electron Spectroscopy (AES). The microtopography was quantified using optical profilometry (OP). The general

morphology and the nanotopography were evaluated using a Field Emission-Scanning Electron Microscope (FE-SEM). Finally, the characterization code of each surface was established using the ISIS, and the main characteristics of each surface were summarized in a reader-friendly ID card.

Results

From a chemical standpoint, in the 18 different surfaces of this group, 11 were based on a commercially pure titanium (grade 2 or 4) and 7 on a titanium-aluminium alloy (grade 5 or grade 23 ELI titanium). 4 surfaces presented some chemical impregnation of the titanium core, and 5 surfaces were covered with residual alumina blasting particles. 15

surfaces presented different degrees of inorganic pollutions, and 2 presented a severe organic pollution overcoat. Only 3 surfaces presented no pollution (and also no chemical modification at all): GC Aadva, Genesio, MIS SEVEN®. From a morphological standpoint, all surfaces were microrough, with different microtopographical aspects and values. All surfaces were nanosmooth, and therefore presented no significant and repetitive nanostructures. 14 surfaces were homogeneous and 4 heterogeneous. None of them was fractal.

and accurate ID card. The SLA-type surfaces have specific morphological characteristics (microrough, nanosmooth, with rare and in general accidental chemical modification) and are the most frequent surfaces used in the industry. However they present different designs, and pollutions are often detected (with blasting/etching residues particularly). Users should be aware of these specificities if they decide to use these products.

Discussion and Conclusion

The ISIS systematic approach allowed to gather the main characteristics of these commercially available products in a clear

Identification card of the MIS SEVEN surface, following the implant Surface identification Standard (ISIS) codification




Fig. 1

Identification Card of the MIS SEVEN® surface: MIS Seven (MIS Implants Technologies, Bar Lev, Israel; Figure 1) was a sandblasted/acid-etched surface on a grade 23 ELI (Extra Low Interstitials) titanium core. No pollution or chemical modification was detected, the surface was moderately microrough, nanosmooth, and homogeneous all over the implant.

¹LoB5 unit, Research Center for Biomaterialization Disorders, Chonnam National University, South Korea. ²Department of Stomatology, School of Dental Medicine, University of Geneva, Switzerland. ³Department of Oral Surgery, Faculty of Medicine, University Federico II of Naples, Italy. ⁴Private Practice, Turin, Italy. ⁵Department of Physics, Seoul National University, Seoul, South Korea. ⁶Private Practice, Paris, France. ⁷Private Practice, Ra'anana, Israel. ⁸Department of Periodontology and Implant Dentistry, College of Dentistry, New York University, New York, USA. ⁹Private Practice, Reims, France. ¹⁰Department of Oral and Maxillofacial Surgery, School of Dentistry, Chonnam National University, South Korea. ¹¹Department of Stomatology, Shanghai Sixth People's Hospital, Shanghai Jiao Tong University, China. ¹²Department of Periodontology and Oral Implantology, University of Guarulhos, Sao Paulo, Brazil. ¹³Department of Periodontics and Oral Medicine, School of Dentistry, University of Michigan, Ann Arbor, USA. *Corresponding author: David M. Dohan Ehrenfest.

Insertion Tools.

The SEVEN key system is designed to facilitate quick, reliable implant procedures. Keys are supplied within the advanced SEVEN surgical kit. The keys are suitable for use with SEVEN NARROW connectors.

Short insertion tool, internal hex.



MT-SRH20

Long insertion tool, internal hex.



MT-LRH20

Short insertion tool, internal hex.




MT-HSI10

Long insertion tool, internal hex.



MT-HLI10

MIS offers a line of specially engineered insertion tools suitable for use either manually or with a ratchet, effectively reducing the number of tools required in the armamentarium.



Insertion options:

- ① Insertion tool and hand key adapter
- ② Motor insertion tool
- ③ Ratchet insertion tool

Please note: In order to assure their efficient operation, tools should be fully inserted into the implants. A complete insertion of the tool optimizes the transfer of force during implant placement and enables simple release of the tool from the hex, whenever necessary.

Tool will not hold implant unless it is completely inserted into the hex.



Package Contents.

Each SEVEN implant comes with a cover screw in the implant package.

Following our “Make It Simple” philosophy, MIS is proud to be the first to include a sterile single-use final drill with every SEVEN implant, to ensure a safe and precise surgical procedure.



Packaging.

A double packing system ensures sterilization and safety. Packages are designed for easy handling during surgery and for ease-of-use with surgical gloves.

Prosthetic platform indication

Prosthetic components are marked by specific colors, representing platform diameters.



Implant diameter & platform indication

The outer tube is color-coded, indicating the implant platform. The numeric indication specifies implant diameter and length.



Implant identification

The sticker on top of the box, specifies implant diameter, length and platform size.

Easy pull tab

The pull tab is easily identified and facilitates convenient opening during surgery.





mis[®]

All rights reserved. No part of this publication may be reproduced, transcribed, stored in an electronic retrieval system, translated into any language or computer language, or be transmitted in any form whatsoever, without the prior written consent of the publisher.



MIS Implants Technologies Ltd.
www.mis-implants.com

The MIS Quality System complies with International Quality Standards: ISO 13485:2003 - Quality Management System for Medical Devices, ISO 9001:2008 - Quality Management System and CE Directive for Medical Devices 93/42/EEC. MIS products are cleared for marketing in the USA and CE approved.