



www.V3-implant.com

The MIS Quality System complies with international quality standards: ISO 13485:2003 - Quality Management System for Medical Devices, ISO 9001: 2008 - Quality Management System and CE Directive for Medical Devices 93/42/EEC. MIS products are cleared for marketing in the USA and CE approved.

V3^{By}_{MIS}[®]
Implant System

P. 4-5

P. 6-7

P. 8-9

P. 10-13

P. 14-15

P. 16-17

P. 18

P. 19

P. 20

P. 21

MIS Warranty:

MIS exercises great care and effort in maintaining the superior quality of its products. All MIS products are guaranteed to be free from defects in material and workmanship. However, should a customer find fault with any MIS product after using it according to the directions, the defective product will be replaced.

Warning: Products should be used by licensed dentists only.

Table of contents.

Introduction

Advantages

Narrow Implants (Ø3.30mm)

Standard Implants (Ø3.90, Ø4.30, Ø5mm)

Surgical Kit

Significant Gain of Bone Volume

Keys & Adapters

Insertion Tools



Package Contents

Packaging

The V3 implant system is the outcome of an exceptionally high-level R&D process that has resulted in an implant that is simple, easy-to-use and offers enhanced functionality and performance. The V3 conical connection implant features built-in design characteristics that provide biological benefits for hard and soft tissues and promotes esthetic results.



© MIS Corporation. All rights reserved.



6.

Advantages.

The unique biologic and mechanical features of the V3 implant encourage bone regeneration and greater volume of bone; supporting highly stable surrounding soft tissues and more esthetic restorations.


All V3 implants, superstructures and tools are color-coded for easy identification of platform sizes.

Blue
indicates a Narrow platform

Purple
indicates a Standard platform




V3 conical connection implants




More bone

The compression-free gaps around the coronal area of the V3 provide a reservoir for blood pooling and the formation of blood clots, for faster implant integration and accelerated bone growth.



Stress reduction

The gaps around the sides of the implant neck result in an open, compression-free zone. Crestal bone loss is minimized by reducing stress in the cortical bone.



Implant neck

The triangular shape at the V3 neck provides high immediate crestal stability and maximum bone preservation. Anchorage is achieved at three points without compromising crestal primary stability.



Platform switching

The V3 incorporates platform switching to preserve the crestal bone around an implant for better bone preservation and to provide an excellent gingival seal for high soft tissue volume.




Conical connection

The 12° conical connection creates an ultimate seal and ideal connection between the implant and abutment with built-in platform switching, reducing micro-movements.



Micro-rings

Micro-rings on the neck of the implant assist in reducing bone stress and bone resorption, and increase BIC.



Surface treatment

The surface roughness and micro-morphology is a result of sandblasting and acid etching. This proven MIS surface technology provides excellent osseointegration resulting in long-lasting clinical success.







Flat apex

The flat apex allows good grip into bone, especially in immediate placement procedures.

8.

V3

Screw type implant range
Narrow Platform

Length		10mm	11.5mm	13mm	16mm
Type		V3-10330	V3-11330	V3-13330	V3-16330
Ø3.30 mm					

Insertion Tools



CT-NSM30
V3 con. con. short insertion tool for motor, SP



CT-NLM30
V3 con. con. long insertion tool for motor, SP



CT-NSR30
V3 con. con. short insertion tool for ratchet, NP




CT-NLR30
V3 con. con. short insertion tool for ratchet, NP


Implant cover screw and healing caps



CC1-00277



CN-HS340
CN-HS440
CN-HS540
CN-HS640
CN-HS840



CN-HA248
CN-HA348
CN-HA448
CN-HA548
CN-HA648
CN-HA848

V3


Ø3.30mm
Narrow Platform

Titanium Alloy Ti 6Al 4V ELI
Sand-Blasted and Acid-Etched
















Catalog No.	Dimensions	
V3-10330	Ø3.30mm length 10mm	
V3-11330	Ø3.30mm length 11.5mm	
V3-13330	Ø3.30mm length 13mm	
V3-16330	Ø3.30mm length 16mm	

Ø3.30mm Implant Procedure

Drilling Speed (RPM)	800-1000	600-800	200-400	200-400	Torque Max. 45N-cm
Diameter	Ø1.90	Ø2.40	Ø2.40	Ø3.30	



- The drilling sequence is illustrated using 11.50mm implants.
- Procedure recommended by MIS cannot replace the judgment and professional experience of the surgeon.

Length	8mm	10mm	11.5mm	13mm	16mm
Type	V3-08390	V3-10390	V3-11390	V3-13390	V3-16390
Ø3.90 mm					
Ø4.30 mm					
Ø5 mm					

Insertion Tools



CT-SSM30
V3 con. con. short insertion tool for motor, SP



CT-SLM30
V3 con. con. long insertion tool for motor, SP



CT-SSR30
V3 con. con. short insertion tool for ratchet, SP




CT-SLR30
V3 con. con. long insertion tool for ratchet, SP


Implant cover screw and healing caps



CC1-00315



CS-HS248
CS-HS348
CS-HS448
CS-HS548
CS-HS648
CS-HS848



VS-HS358
VS-HS458
VS-HS558
VS-HS658
VS-HS858

V3


Ø3.90mm
Standard Platform

Catalog No.	Dimensions	
V3-08390	Ø3.90mm length 8mm	
V3-10390	Ø3.90mm length 10mm	
V3-11390	Ø3.90mm length 11.5mm	
V3-13390	Ø3.90mm length 13mm	
V3-16390	Ø3.90mm length 16mm	

Titanium Alloy Ti 6Al 4V ELI
Sand-Blasted and Acid-Etched

Ø3.90mm Implant Procedure

Drilling Speed (RPM)	800-1000	600-800	450-650	200-400	200-400	Torque Max. 60N-cm
Diameter	Ø1.90	Ø2.40	Ø2.40	Ø3	Ø3	Ø3.90




- The drilling sequence is illustrated using 11.50mm implants.
- Procedure recommended by MIS cannot replace the judgment and professional experience of the surgeon.

Titanium Alloy Ti 6Al 4V ELI
Sand-Blasted and Acid-Etched

Catalog No.	Dimensions	
V3-08430	Ø4.30mm length 8mm	
V3-10430	Ø4.30mm length 10mm	
V3-11430	Ø4.30mm length 11.5mm	
V3-13430	Ø4.30mm length 13mm	
V3-16430	Ø4.30mm length 16mm	

Ø4.30mm Implant Procedure

Drilling Speed (RPM)	800-1000	600-800	450-650	350-550	200-400	200-400	Torque Max. 60N-cm
Diameter	Ø1.90	Ø2.40	Ø2.40	Ø3	Ø3.50	Ø3.50	Ø4.30



- The drilling sequence is illustrated using 11.50mm implants.
- Procedure recommended by MIS cannot replace the judgment and professional experience of the surgeon.

V3


Ø5mm
Standard Platform

Catalog No.	Dimensions	
V3-08500	Ø5mm length 8mm	
V3-10500	Ø5mm length 10mm	
V3-11500	Ø5mm length 11.5mm	
V3-13500	Ø5mm length 13mm	
V3-16500	Ø5mm length 16mm	

Titanium Alloy Ti 6Al 4V ELI
Sand-Blasted and Acid-Etched

Ø5mm Implant Procedure

Drilling Speed (RPM)	800-1000	600-800	450-650	350-550	300-500	200-400	200-400	Torque Max. 60N-cm
Diameter	Ø1.90	Ø2.40	Ø2.40	Ø3	Ø3.50	Ø4	Ø4	Ø5



- The drilling sequence is illustrated using 11.50mm implants.
- Procedure recommended by MIS cannot replace the judgment and professional experience of the surgeon.

14.

Surgical Kit.

MK-0051



MT-C2416



MT-C2413



MT-C2411



MT-C2410



MT-C2408



CT-TDN40



CT-TDN35



CT-TDN30



CT-BTT40



CT-BTT35



CT-BTT30



CT-BTT24



MT-SMD10



MT-PDM24



MT-PD440



MT-BT120



MT-RI030



MT-RE172




MT-RE160

1.

MARKING DRILLS



PILOT DRILLS



2.

BODY TRY-IN



SHARED TOOLS

EXTRACTORS

PROBE

RATCHET WRENCH

CT-NLM30

CT-NSM30

CT-SLM30

CT-SSM30

CT-NLR30

CT-NSR30

CT-SLR30

CT-SSR30

3.

V3 TOOLS

C1 TOOLS

PISTON DRILLS



BODY TRY-IN



DRILL EXTENDER



0.05" HEX. DRIVERS



C1 conical connection insertion tools are supplied separately, MK-0054.



MT-RDS30



MT-RDL30

MIS[®] CONICAL CONNECTION SURGICAL KIT

DO NOT EXCEED 134°C(273°F)




MT-DE001

Significant Gain of Bone Volume.

The unique triangular-shape of the V3 implant's coronal portion encourages bone regeneration and furthers the gain of greater volume of bone; supporting highly stable surrounding soft tissues and resulting in more esthetic restorations. This triangular design provides solid anchorage at three points in the crestal area while forming gaps between the remaining sides of the implant neck and the osteotomy, resulting in a compression-free zone, where a stable blood clot can more easily be achieved. The ingenious combination of compression-free gaps with a firm anchorage is critical in order to establish a stable blood clot;


the first step towards a successful osseointegration process: Hemostasis Phase > Proliferative Phase > Remodeling Phase. This initial gain of bone volume is a highly advantageous starting point.




18.

Keys & Adapters.


V3 implant placement tools are specifically designed to enable quick, reliable implant placement procedure.





Areas **A** and **B** indicate the position of the flat surface of the insertion tool to be aligned with the flat surface of the implant **C**.

This allows the dentist to easily identify the flat side of the implant for desired placement.



The insertion tool allows the delivery of a cover screw or a healing cap onto the implant after insertion.

Insertion Tools.

The advanced insertion tool system allows secure implant placement without the use of a mount.

20.

Package Contents.

Each V3 implant comes with a sterile cover screw and single-use final drill, suitable for all drilling protocols.

The sterile inner tube is fitted with a special titanium sleeve that has an anti-rotation grip, to ensure easy engagement between the insertion tool and the implant.



Packaging.

Providing a simple, immediate identification of implant type, length and diameter, the V3 package is well-designed for ease-of-use during surgical procedures.

Implant diameter & platform indication

The outer tube is color-coded indicating the implant platform. The numeric indication specifies the implant diameter and length.



Prosthetic platform indication

Prosthetic components are marked by specific colors, representing platform sizes.

A double packing system ensures sterilization and safety. Packages are designed for convenience during surgery and for use with surgical gloves.



Implant identification markings


Quick identification of implant size and length. Sticker on the box lid, specifies implant diameter, length and platform size

Easy pull tab

The convenient pull tab facilitates quick and easy opening during surgery.

Logical storage

Packages fit perfectly into clinic drawers for space-saving storage and easy identification.



mis[®]

All rights reserved. No part of this publication may be reproduced, transcribed, stored in an electronic retrieval system, translated into any language or computer language, or be transmitted in any form whatsoever, without the prior written consent of the publisher.