

Implant Procedure M4

Drilling Speed (RPM)	1200-1500	900-1200	500-700	200-500	15-25
Diameter	Ø1.90	Ø2.40	Ø2.80	Ø3.00-Ø3.30	Ø3.30

Ø3.30 mm



Optional for long type 162
Countersink for long type 162
MT-GDN33

Drilling Speed (RPM)	1200-1500	900-1200	500-700	400-700	400-600	400-600	200-500	15-25
Diameter	Ø1.90	Ø2.40	Ø2.80	Ø3.20	Ø3.80	Ø4.50	Ø4.50-Ø5.00	Ø5

Ø5 mm



Optional for long type 162
Countersink for long type 162
MT-GDN50

Drilling Speed (RPM)	1200-1500	900-1200	500-700	400-700	200-500	15-25
Diameter	Ø1.90	Ø2.40	Ø2.80	Ø3.20	Ø3.75	Ø3.75

Ø3.75 mm



Optional for long type 162
Countersink for long type 162
MT-GDN33

Drilling Speed (RPM)	1200-1500	900-1200	500-700	400-700	400-600	400-600	300-500	300-500	200-500	15-25
Diameter	Ø1.90	Ø2.40	Ø2.80	Ø3.20	Ø3.80	Ø4.50	Ø5	Ø5.50	Ø6.00	Ø6

Ø6 mm



Optional for long type 162
Countersink for long type 162
MT-GDN50

Drilling Speed (RPM)	1200-1500	900-1200	500-700	400-700	400-600	200-500	15-25
Diameter	Ø1.90	Ø2.40	Ø2.80	Ø3.20	Ø3.80	Ø3.20-Ø4.20	Ø4.20

Ø4.20 mm



Optional for long type 162
Countersink for long type 162
MT-GDN33

Implant Packaging

Implant diameter & platform indication
The outer tube of each implant has a distinctive color-coded label indicating implant diameter, length and platform size.





mis | MAKE IT SIMPLE

mis | M4
INTERNAL HEX. CONNECTION


MIS® M4

MIS M4 implants combine the benefits of cylindrical and conical implant designs, aiming to achieve excellent primary stability in every clinical scenario. The two main features of M4 implants are: Self-tapping, V-shaped thread design with three spiral channels, allowing smooth insertion even in type 1 bone conditions ■ A flat, cutting, tapered apex, enabling instant grip into bone in immediate placement procedures.



Internal hexagon

MIS M4 implants feature an internal hex. connection. This well established connection assures proper abutment seating, anti-rotational engagement, resistance to lateral forces, excellent esthetic results and more.




Combined cylindrical and conical shape with V-shaped threads

The implant body and thread shape is designed for mild bone compression while achieving maximum initial and long-term stability.



Dual thread

A dual thread design enhances the placement procedure while being gentle to the surrounding bone. The overall insertion rate of M4 is 1.6mm per revolution.



Surface

The surface roughness and micro-morphology is achieved by a combination of sand-blasting and acid-etching. MIS proven surface technology has provided millions of patients and clinicians with excellent osseointegration results and long-lasting clinical success.




Three spiral channels

Three spiral channels at the apical end of the implant support self-tapping properties. The channels also collect bone chips in the course of insertion, supporting efficient osseointegration and long-term stability.




A flat cutting apex

A flat cutting apex allows for final adjustments during placement procedures.




Narrow Platform

Ø3.30





Standard Platform

Ø3.75, Ø4.20



Wide Platform

Ø5, Ø6



Surgical Instruments Kit.

The M4 Surgical Kit (MK-0016) contains the complete range of drills and tools required for a full implant placement procedure.

The kit features a convenient tool layout and protective cover with an easy-opening mechanism for quicker access.

