



MISCELLANEOUS ITEMS

DIA-100FS/DIA-100SS

- Optical grade protective shields and glasses for the face and eyes from flying debris or substances.
- Lightweight and very comfortable.
- Shields can be quickly positioned and fits over eyeglasses
- Easy to clean, washable

Dia-100FS

- Full face protection shield
- 1 Frame & 5 shields

1005-5010	DIA-100FS	REFILL PACKAGE	5 shields only
1005-5011	DIA-100FS	REGULAR PACKAGE	CRYSTAL
1005-5012	DIA-100FS	REGULAR PACKAGE	WHITE
1005-5014	DIA-100FS	REGULAR PACKAGE	BLUE



DIA-100FS

DIA-200B

• ITEM #1005-1031

- Multi Purpose Goggles. Used among dentists and surgeons for eyes protection from flying debris. Preventing eye injuries and contamination from fine dust.
- Spherical lens design for superior optics.
- Comfortable lightweight fit.
- Fits over Eyeglasses.
- Designed for coolness and ventilation.
- Anti-fog coating, Anti-static.

Product Features



Superior cushioning and ventilation



Good air circulation



Easily used with correction glasses



Anti-fog

Anti-static



Anti-fog coating

DIA-400D

• ITEM #1005-1051


- The protective glasses for blue curing light (400-510nm) are used to protect your eyes when curing composite resin.
- Excessive exposure to blue lights will damage your vision and your eyes health.
- Lenses color: Amber

Anti-fog

Anti-UV

Anti-static

Anti-scratch



Product Features



Vertical adjustments



Horizontal adjustments



Ultra light weight



SMALL ROUND TYPE A

- ITEM # 803-003 Red
- ITEM # 803-004 Blue



- Has 6 Silicone Color Coded Silicone Cups.
- Ideal for Organizing Gutta Percha Points & Paper Points.
- Removable Insert for Cleaning.
- Autoclavable silicone cups.
- Diameter: 40m/m, Height: 60m/m



SMALL ROUND TYPE B

- ITEM # 802-003 Red
- ITEM # 802-004 Blue



- Ideal for Organizing Hand Files.
- Holds 20 instruments.
- Removable Insert for Cleaning.
- Diameter: 40m/m, Height: 60m/m



SMALL SQUARE TYPE A

- ITEM #807-003 Red
- ITEM #807-004 Blue



- Has 6 Silicone Color Coded Silicone Cups.
- Ideal for Organizing Gutta Percha Points & Paper Points.
- Removable Insert for Cleaning.
- Autoclavable silicone cups.
- Dimensions: 50 x 50 x 45m/m (W x H x D)



SMALL SQUARE TYPE B

- ITEM #806-003 Red
- ITEM #806-004 Blue



- Ideal for Organizing Hand Files.
- Holds 30 Endo Instruments.
- Removable Insert for Cleaning.
- Dimensions: 50 x 50 x 45m/m (W x H x D)



MEDIUM OBLONG TYPE A

- ITEM #801-003 Red
- ITEM #801-004 Blue



- A Versatile & Practical Organizer with a Clean Grip for Hand & Rotary File as well as Cups for Gutta Percha & Paper Points.
- Holds 36 Instruments, 3 Silicone Cups, 1 Clean Grip.
- Removable Insert for Cleaning.
- Dimensions: 130 x 50 x 50m/m (W x H x D)




MEDIUM OBLONG TYPE B

- ITEM #809-003 Red
- ITEM #809-004 Blue



- A Versatile & Practical Organizer for Hand & Rotary File as well as Cups for Gutta Percha & Paper Points.
- Holds 60 Instruments, 4 Silicone Cups.
- Removable Insert for Cleaning.
- Dimensions: 130 x 50 x 50m/m (W x H x D)



MEDIUM ROUND

- ITEM # 821-001 Silver



- Ideal for Organizing Hand Files & Gutta Percha Points & Paper Points.
- Holds 64 Instruments, 4 Color Coded Silicone Cups.
- Removable Insert for Cleaning.
- Diameter: 92m/m, Height: 69m/m



MEDIUM SQUARE

- ITEM #804-003 Red
- ITEM #804-004 Blue



- A Versatile & Practical Organizer with a Clean Grip for Hand & Rotary File as well as Cups for Gutta Percha & Paper Points.
- Features a Endo Meter for measuring Files & Points.
- Holds 5 Burs, 30 Instruments, 6 Silicone Cups, 1 Clean Grip, Endo Meter and 1 location for small instruments.
- Removable Insert for Cleaning.
- Dimensions: 120 x 50 x 80m/m (W x H x D)



DIA-DELUXE

• ITEM #811-003 Red



- DiaDent's Ultimate Box for Organizing your Files & Points. Holds 4 Spreaders/Pluggers, 56(21mm) & 48(25mm) Instruments, 6 Gates Glidens, 10 Barbed Broaches, 28 Cups, Section for Cotton.
- Removable Insert for Cleaning.
- Dimensions:
215 x 50 x 100m/m
(W x H x D)



SMALL OBLONG-6 COMPT.

• ITEM #805-003 Red

• ITEM #805-004 Blue



- Ideal for Holding Burs & Small Instruments.
- Has 6 Compartments.
- Removable Insert for Cleaning.
- Dimensions:
130 x 15 x 50m/m
(W x H x D)



LARGE OBLONG TYPE A

• ITEM #812-004 Blue



- A Simple & Effective Organizer with Sections for your Hand Files.
- Holds 180 Endo Instruments, 8 Color Coded Silicone Cups for Points, perforated bottom case.
- Removable Insert for Cleaning.
- Dimensions:
218 x 53 x 105m/m
(W x H x D)



LARGE OBLONG TYPE B

• ITEM #812-005



- An Organizer for Hand & Rotary Files as well as Points.
- Store your Instruments in One Place.
- Holds 357 Endo Instruments, 21 Silicone Cups for Points.
- Removable Insert for Cleaning.
- Dimensions:
205 x 53 x 105m/m
(W x H x D)



LARGE OBLONG TYPE C

• ITEM #812-006



- A Very Practical Box. Easy to use inserts individually as needed.
- Holds 206 Endo Instruments, 12 Silicone Cups for Points. Interchangeable Inserts.
- Removable Insert for Cleaning.
- Dimensions:
205 x 53 x 105m/m
(W x H x D)




LARGE OBLONG TYPE D

• ITEM #812-007



- Organizer for Hand & Rotary Files.
- Holds 192 Hand Instruments, 80 Rotary Instruments.
- Removable Insert for Cleaning.
- Dimensions: 205 x 53 x 105m/m
(W x H x D)




LARGE COMBI

• ITEM # 818-001



- Versatile & Practical. Allows separation of specific instruments.
- Holds 5 individual Inserts(A+B+C+D+E).
- Insert A : 6 cups + 6 cups
- Insert B : 2 sets of 6 cups & Holds 12 Endo Instruments
- Insert C : 4 cups, Holds 48 Endo Instruments
- Insert D : 2 cups, Holds 48 Endo Instruments
- Insert E : Holds 72 Endo Instruments
- Removable Insert for Cleaning.
- Dimensions: 205 x 53 x 105m/m
(W x H x D)




LARGE OBLONG TYPE E

• ITEM # 812-008



- Store your files in one place. Holds 308 Endo Instruments, 5 Color Coded Silicone Cups, 3 Locations.
- 2 Separate Inserts. Height Adjustment(not visible).
- Removable Insert for Cleaning.
- Dimensions: 205 x 53 x 105m/m (W x H x D)

WET-HEAT
AUTOCLAVE

ALUMINIUM

STAINLESS
STEEL

PLASTIC



SILICONE



GLASS

MEDIUM COMBI

• ITEM #817-001



- A Compact & Practical File & Point Organizer.
- Holds 3 Individual Inserts(C+D+E).
Insert C : 4 cups, Holds 48 Endo Instruments
Insert D : 2 cups, Holds 48 Endo Instruments
Insert E : Holds 72 Endo Instruments
- Removable Insert for Cleaning.
- Dimensions:
122 × 50 × 100m/m (W × H × D)



SMALL COMBI TYPE C

• ITEM #816-030



- A Compact Organizer for Hand & Rotary Files.
- Holds 48 Endo Instruments, 4 Color Coded Silicone Cups for Points.
- Removable Insert for Cleaning.
- Autoclavable silicone cups.
- Dimensions: 105 × 52 × 44m/m (W × H × D)



SMALL COMBI TYPE E

• ITEM #816-050



- A Compact Organizer for Hand & Rotary Files.
- Holds 72 Endo Instruments.
- Removable Insert for Cleaning.
- Dimensions: 105 × 52 × 44m/m (W × H × D)



SMALL COMBI TYPE G

• ITEM #816-070



- A Compact Organizer for Hand & Rotary Files.
- Holds 48 Endo Instruments Medium Sized Perforations.
- Removable Insert for Cleaning.
- Dimensions: 105 × 52 × 44m/m (W × H × D)



ENDO GRIP(MINI)

• ITEM #833-001



- Commonly Used for Cleaning Files.
- Holds sponge or gauze.
- Available in Blue.
- Diameter:
40m/m, Height: 35m/m



ENDO CLEANER

• ITEM #834-001 Open

• ITEM #834-002 Closed



- Commonly Used for Cleaning and sterilizing Files.
- Top unscrews for easy replacement of Gauze.
- Diameter: 40m/m, Height: 30m/m




ENDO METER

• ITEM #835-002 Type A- Silver

• ITEM #835-003 Type B- Red



- Holds & Measures 10 Files.
- Gauges from 1 to 35mm.
- Features a locking spring mechanism inside for easy use.
- Diameter: 35m/m, Height: 48m/m

WET-HEAT
AUTOCLAVE

ALUMINIUM

STAINLESS
STEEL

PLASTIC



SILICONE




GLASS

LARGE OBLONG EMPTY BOX

• ITEM #820-003 Blue



- An all purpose storage container.
- Ideal for holding accessories.
- Dimensions:
210 × 25 × 50m/m (W × H × D)

WET-HEAT
AUTOCLAVE

ALUMINIUM

STAINLESS
STEEL

PLASTIC



SILICONE



GLASS

GERMICIDE TRAY

• ITEM #711-001



- Used for Cold Sterilization.
- Easy & Safe to handle.
- One Large & One Small compartment.
- Easy loading and unloading.
- Dimensions:
263 x 72 x 115 m/m
(W x H x D)



WASTE CONTAINER

• ITEM #814-001 Small

- Diameter: 56m/m, Height: 47m/m
- Color: Silver

• ITEM #815-002 Large

- Diameter: 70m/m, Height: 74m/m
- Color: Gold

Separates into 3 pieces for easy Cleaning.
Equipped with locking system.



CONTAINER(STERILIZING)



• ITEM #824-001 Type A

- Holds 10 Endo Instruments.
- Dimensions : 185 x 35 x 145 m/m
(W x H x D)

• ITEM #824-002 Type B

- Holds 16 Endo Instruments.
- Dimensions : 185 x 35 x 286 m/m
(W x H x D)

- Well Built, Rigid, Sterilizing Container.
- Features Locking System, Fits many sterilizing units.



BUR STERI BATH

• ITEM # 730-001 White

• ITEM # 730-002 Transparent



- Used for Sterilizing Burs & Small Instruments.
- Removable insert for liquid drainage.
- 1 Basin bowl/ 1 Perforated insert/ 1 Cover.
- Diameter : 88 m/m
- Height : 68 m/m



GERMICIDE ROUND BUR BOX

• ITEM #915-001



- Contains perforated insert for drainage.
- Practical & Durable.
- 1 Basin bowl/ 1 Perforated insert/ 1 Cover.
- Diameter: 75m/m,
- Height: 35m/m




MINI INSTRU. MAT

• ITEM #702-106 Blue



- Holds up to 8 instruments
- Dimensions :
130 x 16 x 103m/m(W x H x D)

WET-HEAT
AUTOCLAVE

ALUMINIUM

STAINLESS
STEEL

PLASTIC



SILICONE




GLASS

MAGNETIC BUR BLOCK(OBLONG-14 HOLE)

Type A

- ITEM #717-103 Yellow
 ITEM #717-105 Beige
 ITEM #717-106 Blue
 - Dimensions:
 77 x 13 x 33m/m (W x H x D)



Type A

Type B

- ITEM #717-203 Yellow
 ITEM #717-205 Beige
 ITEM #717-206 Blue
 - Dimensions:
 76 x 12 x 37mm (W x H x D)
 - 14 Holes
 - Fits both FG and RA types.



Type B

MINI BUR STAND

Type A



- Silver ITEM #838-101
 •Red ITEM #838-103
 •Blue ITEM #838-104



- Holds 12 Friction Grip (high speed) Burs
 - Durable Aluminum
 - Dimensions: 45 x 32 x 25m/m (W x H x D)

Type B



- Silver ITEM #838-201
 •Red ITEM #838-203
 •Blue ITEM #838-204



- Holds 12 Right Angle (low speed) Burs
 - Durable Aluminum
 - Dimensions: 45 x 32 x 25m/m (W x H x D)

MAGNETIC BUR BLOCK(SQUARE-8 HOLE)

- ITEM #717-303 Yellow
 ITEM #717-306 Blue



- Fits both FG and RA types.
 - Dimensions : 40 x 13 x 40m/m
 (W x H x D)




STERI BUR GUARD(22 HOLE)

- ITEM #716-103 Yellow
 ITEM #716-106 Blue



- Fits 22 latch & friction grip burs
 - Lid does not allow the burs to fall out of the stand.
 - Diameter: 75m/m, Height: 35m/m



MAGNETIC BUR BLOCK(ROUND 24 HOLE)



8 Holes

- ITEM #717-503 Yellow
 ITEM #717-505 Beige
 ITEM #717-506 Blue
 - Diameter : 43m/m, Height: 10m/m



8 Holes

24 Holes

- ITEM #717-603 Yellow
 ITEM #717-605 Beige
 ITEM #717-606 Blue
 - Diameter : 73m/m, Height: 10m/m
 - Fits both FG & RA types



24 Holes

PROPHY RUBBER CUP RIBBED & WEBBED

- **ITEM #950-201** : Latch Type
- **ITEM #950-202** : Screw Type
- Packed 144 cups



Latch Type (RC - 001)



Screw Type (RC - 002)

PROPHY RUBBER CUP RIBBED

- **ITEM #950-221** : Latch Type
- **ITEM #950-222** : Screw Type
- Packed 144 cups



Latch Type (RC - 005)



Screw Type (RC - 006)

PROPHY RUBBER CUP WEBBED

- **ITEM #950-211** : Latch Type
- **ITEM #950-212** : Screw Type
- Packed 144 cups



Latch Type (RC - 003)



Screw Type (RC - 004)

PROPHY RUBBER CUP SNAP-ON

- **ITEM #950-251** : Ribbed
- **ITEM #950-252** : Webbed
- Packed 144 cups



Ribbed (RC - 010)



Webbed (RC - 011)

PROPHY RUBBER CUP UMBRELLA SHAPE

- **ITEM #950-231** : Latch Type
- **ITEM #950-232** : Screw Type
- Packed 144 cups



Latch Type (RC - 007)



Screw Type (RC - 008)

PROPHY RUBBER CUP SNAP-ON RIBBED & WEBBED

- **ITEM #950-253**
- Packed 144 cups



(RC - 012)

PROPHY RUBBER CUP WEBBED WITH WASHER LATCH TYPE

- **ITEM #950-241**
- Packed 144 cups



(RC - 009)

CANDLE BRUSH

• ITEM #941-104 :

Latch Type, Packed 24 brushes

• ITEM #941-105 :

Screw Type, Packed 24 brushes and 10 mandrels

- For removing stains from occlusal pits and fissures and for preparing tooth surface prior to sealant and orthodontic application.

- Diameter : Latch Type : 4.50mm /
Screw Type : 5.00mm



Candle Brush (Latch Type)



Candle Brush (Screw Type)

POST BRUSH

• ITEM #941-100 :

Mandrels for Screw Type, packed 20 pcs

• ITEM #941-102 :

Screw Type, Packed 30 brushes

- For quick cleaning root canal preparation, under bridge.

- Available to insert sterilizing liquid during turn it left in Contra Angle.

- Available to flow out the stains like cement, liquid during turn it right.

- Recommended rotation: 300-800rpm.

- Single use.

- Diameter : Latch Type : 1.20mm /
Screw Type : 1.03mm



Post Brush (Screw Type)

PROPHY BRUSH

• ITEM #950-101 : Latch Type

• ITEM #950-102 : Screw Type

- Packed 144 brushes

- Diameter

Latch Type : 4.48mm

Screw Type : 4.48mm



Latch Type (PB - 001)



Screw Type (PB - 002)

PROPHY BRUSH POINTED HEAD

• ITEM #950-112 : Screw Type

- Packed 144 brushes

- Diameter

Screw Type : 4.48mm



Screw Type (PB - 004)

PROPHY BRUSH LARGE HEAD LATCH TYPE

• ITEM #950-121

- Packed 144 brushes

- Diameter : 5.98mm



(PB - 005)

PROPHY BRUSH LARGE POINTED HEAD LATCH TYPE

• ITEM #950-131

- Packed 144 brushes

- Diameter : 5.98mm



(PB - 006)

INSTRUMENT TRAY TYPE A (Perforated)

- ITEM #826-001



- Holds 10 Instruments
- Ideal for Instrument Sterilization.
- Perforations Facilitate Sterilization & Drainage.
- Dimensions:
190 x 35 x 150m/m (W x H x D)




INSTRUMENT TRAY TYPE B (Non-Perforated)

- ITEM #825-002



- Holds 16 Instruments
- Ideal for Instrument Storage.
- Dimensions:
280 x 35 x 180m/m (W x H x D)



INSTRUMENT CASSETTE TYPE A

- ITEM # 913-001



- Holds 14 Instruments.
- Has 4 silicone instrument holding units.
- Easy to open & close.
- Locking system.
- Dimensions:
215 x 25 x 105m/m (W x H x D)




INSTRUMENT CASSETTE TYPE C

- ITEM #913-003



- Holds 40 Instruments.
- Has 4 silicone instrument holding units.
- Easy to open & close.
- Locking system.
- Fits Stainless Steel Rack Type D.
- Dimensions:
280 x 25 x 190m/m (W x H x D)



MIXING SLAB



- High quality, competitively priced Mixing Slabs available in 4 sizes.
- Do not Autoclave.



- ITEM #603-005 LT
- Dimensions : 150 x 15 x 80m/m (W x H x D)
- ITEM #603-004 LM
- Dimensions : 150 x 10 x 75m/m (W x H x D)
- ITEM #603-002 L
- Dimensions : 112 x 6 x 76m/m (W x H x D)

CHAIR SIDE TABLE

- ITEM #712-001 White



- Comes with attachments.
- Can be adapted to several types of chairs.
- Hole Diameter : available from 25 to 50m/m.
- Easy access to instruments.
- Dimensions : 390 x 30 x 235m/m (W x H x D)



TISSUE HOLDER

- ITEM #712-101 White



- Can be adapted to several types of chairs.
- Tissue box holder fits on unit chair pole.
- Hole Diameter : available from 25 to 50m/m.
- Dimensions :
245 x 54 x 125m/m(W x H x D)



INSTRUMENT CASSETTE TYPE B

- ITEM #913-002



- Holds 20 Instruments.
- Has 4 silicone instrument holding units.
- Easy to open & close.
- Locking system.
- Fits Stainless Steel Rack Type C.
- Dimensions:
185 x 25 x 145m/m (W x H x D)




X-RAY FILM HOLDER

ITEM #750-001



- Used for Anterior and Posterior X-Ray.
- Angled for patient comfort in Anterior x-rays.




RETRACTOR (TYPE A)

ITEM #740-001 Adult

ITEM #740-002 Child



- Soft Flexible Plastic
- Pkg. of 2



COTTON ROLL DISPENSER

ITEM #770-001 White



- Holds approximately 60 Cotton Rolls.
- Features pull out drawer for easy cotton dispensing.
- Lid flips back for refilling.
- Non-slip bottom.
- Adjustable width for Cotton Roll in 20, 30 and 40m/m.
- Dimensions:
170 x 101 x 52m/m(W x H x D)



ID-RING

- Silicone rubber
- Standard size(inner diameter 3.6mm)
- Pkg. of 120



- ITEM #942-101 White
- ITEM #942-102 Black
- ITEM #942-103 Brown
- ITEM #942-104 Green
- ITEM #942-105 Orange
- ITEM #942-106 Red
- ITEM #942-107 Gray
- ITEM #942-108 Yellow
- ITEM #942-109 Pink
- ITEM #942-110 Light Brown
- ITEM #942-191 Assorted




CHEEK RETRACTOR (TYPE B)

- ITEM #741-001 Adult
- ITEM #741-002 Child



- Soft Flexible Plastic
- Expands the mouth widely
- Pkg. of 2



ENDO TRAINING BLOC

ITEM #943-100



- Clear acrylic blocs.
- Ideal for training bloc for new techniques.
- Used for .04 tapered NiTi Rotary Files(curved type)
- Pkg. of 10