



# NATIONAL PROMOTION POLAND

---

AUTUMN  
2020

Valid from 15<sup>th</sup> of October  
until 31<sup>th</sup> of December 2020

**MEDIF**  
FUTURE MEDICINE  
for dentistry

# Introducing...

# HuFriedyGroup

The Best In Practice

## Who We Are

We are a combination of companies that have come together under one roof to deliver the products, services, and solutions you need to be the best in practice. With our passion for infection prevention solutions, premium instrumentation, and service and support programs, HuFriedyGroup will help you perfect your craft and continue to deliver the highest level of care to each patient.

## What Can You Expect From HuFriedyGroup?

HuFriedyGroup is here to provide the best dental experience to the industry. Our unique combination of products, services, and support enable dental professionals like you to reduce risk, improve efficiencies, drive compliance, and enhance patient and staff safety.

### World-Class Products

Our unwavering passion and unrivaled expertise in the dental industry is evident in the instruments and infection prevention products you trust each and every day.

### Compliance & Protection

We provide solutions that help you to ensure a safe environment for patients and staff.

### Service & Support


Through our community platforms, customer service and educational programs, we provide you with the information and tools you need to perfect your craft.

## Proud members of HuFriedyGroup:



# RESTORATIVE INSTRUMENTS

## CARVERS HU-FRIEDY



### CVWI8

Carver Willand  
#8 hdl #41

**Sugg. Retail**

zł 183

**Promo**

zł 119

### CVIPC

Carver  
hdl #41 Interproximal


**Sugg. Retail**

zł 183

**Promo**

zł 119

## EXCAVATORS HU-FRIEDY



### EXC17WH

Excavator #17WH  
hdl #41

**Sugg. Retail**

zł 186

**Promo**

zł 121

### EXC18WH

Excavator #18WH  
hdl #41 heavy

**Sugg. Retail**

zł 186

**Promo**

zł 121

## PLASTIC FILLING HU-FRIEDY



### CIGFT16

Plastic Filling #1  
Instrument Goldstein  
hdl #6 Flexi-Thin

**Sugg. Retail**

zł 215

**Promo**

zł 129

### CIGFT26

Plastic Filling #2  
Instrument Goldstein  
hdl #6 Flexi-Thin

**Sugg. Retail**

zł 215

**Promo**

zł 129



### GCPS6


Gingival Cord Packer  
#S6 hdl #41, smooth

**Sugg. Retail**

zł 211

**Promo**

zł 127



### BB216

Burnisher #21  
hdl #6

**Sugg. Retail**

zł 190

**Promo**

zł 114

### BB186

Burnisher #18  
hdl #6

**Sugg. Retail**

zł 190

**Promo**


zł 114

# RESTORATIVE INSTRUMENTS

## AKRO-FLEX™ COMPOSITE INSTRUMENT HU-FRIEDY


### CREATE BEAUTIFUL RESTORATIONS

Akro-Flex™ is an extremely flexible composite instrument which allows clinicians to effortlessly manipulate composite materials during aesthetic restorations. **The inclusion of a ductile material allows Akro-Flex™ to perform as a solid brush.**



#### HYPER-THIN PROFILE

Ultra thin working ends reach narrow interproximal spaces with ease. The flexible, versatile working ends allow for better visibility as compared to traditional composite instruments.



#### UNPARALLELED FLEXIBILITY

By incorporating Nickel Titanium, a material known for its super elasticity, Akro-Flex™ acts as a solid brush. The resilient working ends are excellent when creating fine anatomical detail with delicate, artistic strokes.



#### ERGONOMIC HANDLE

The smooth, lightweight handle offers increased control due to the large diameter. It creates an ergonomically friendly option that provides maximum comfort and helps reduce hand fatigue. Reducing hand fatigue can increase the longevity of a clinician's career.

### AKRO-FLEX

Akro-Flex™  
Composite Instrument

Sugg. Retail


zł 326

Promo

zł 235



WATCH OUR  
**TUTORIAL VIDEO**  
FOR AKRO-FLEX™  
ON YOUTUBE



Patent pending.

# RESTORATIVE INSTRUMENTS BY WERESTORE.IT<sup>†</sup>



## KOL CORNER

Instruments made in collaboration with **Dr. Gaetano Paolone** and **Dr. Salvatore Scolavino**, Opinion leaders Hu-Friedy.

## CALIBRA™ COMPOSITE INSTRUMENT HU-FRIEDY

### FOR CALIBRATED RESTORATIONS

Calibra™ is an **XTS (AlTiN)<sup>††</sup>** coated composite instrument which helps to provide the proper anatomical curvature to the dentin material layer of a restoration and leaves the ideal space for the final layer of enamel material. The **Calibra™ “IN”** point is used to define the dentinal skeleton, while the **Calibra™ “OUT”** point helps to reduce the chances of **excess composite** material being placed during the final layer of the restoration. Calibra™ acts as a **guide**, helping to prevent over-contouring at the preparation margin while assisting in reducing excess composite thickness build up as well as **modeling** the final layer of composite material.

### TNCALIBRA

XTS PFI WeRestore.it  
Calibra™

NEW

Sugg. Retail

Promo

zł 317

zł 222

### HOW CALIBRA™ IN WORKS



### HOW CALIBRA™ OUT WORKS



Calibra™ consists of two terminals: “IN” which defines the dentinal skeleton. “OUT” which allows to verify the thickness and shape of the last layer without excess. This instrument is in fact both the first caliper ever produced to check the thickness (a practical alternative to the sagittal silicone index) and the first modeling tool ever produced for the last layer.

<sup>†</sup> WeRestore.it is not owned by or affiliated with Hu-Friedy. <sup>††</sup> Data on file, available upon request.

CALIBRA™  
OUT

CALIBRA™  
IN



**CALIBRA™** BY  
WERESTORE.IT  
WATCH OUR  
PRACTICAL  
VIDEO

# DIAGNOSTIC INSTRUMENTS HU-FRIEDY



## EXS3A6

Explorer #3A  
hdl #6

Sugg. Retail zł 99

**Promo zł 74**



## EXS6A6

Explorer #6A  
hdl #6

Sugg. Retail zł 99

**Promo zł 74**



## EXS236

Explorer Shepherds  
Hook #23 hdl #6

Sugg. Retail zł 99

**Promo zł 74**



## EXTU176

Explorer Wilkins/Tufts  
#TU17 hdl #6

Sugg. Retail zł 99

**Promo zł 74**



## EXD11/12A6

Explorer #11/12 hdl #6  
After Five

Sugg. Retail zł 134

**Promo zł 100**

1-2-3....15



## PCPUNC156

Probe UNC #UNC15  
hdl #6

Sugg. Retail zł 173

**Promo zł 130**

3.5....11.5



## PCP11.5B6

Probe #11.5B  
hdl #6

Sugg. Retail zł 173

**Promo zł 130**



## EXTU17/236

Explorer Wilkins/Tufts  
#TU17/23 hdl #6

Sugg. Retail zł 134

**Promo zł 100**

3-6-9-12



## PCP126

Probe #12  
hdl #6

Sugg. Retail zł 173

**Promo zł 130**

1-2-3-5-7-8-9-10



## PQOW6

Probe Michigan/Williams  
hdl #6

Sugg. Retail zł 155

**Promo zł 116**



## P15/11.5B6

Probe #15/11.5B  
hdl #6

Sugg. Retail zł 281

**Promo zł 211**

3-6-9-12



## PQ2N6

Birfurcation Probe  
Nabers #2N hdl #6

Sugg. Retail zł 281

**Promo zł 211**



# EVEREDGE™ 2.0 INSTRUMENTS

## EVEREDGE™ 2.0 SCALERS AND CURETTES

# NO SHARPENING... DO YOU REALLY BELIEVE IT?



### PROPHY MAINTENANCE KIT DIN HU-FRIEDY


IMS™ Infinity Cassette with 5 EverEdge 2.0 Instruments:

**IMEDPROPHYM**

| DESCRIPTION                                  | CODE                |
|--|---------------------|
| Scaler Neville #1 anterior Sickle            | <b>SCNEVI19E2</b>   |
| Scaler #204SD Sickle                         | <b>S204SD9E2</b>    |
| Curette Columbia #4R/4L                      | <b>SC4R/4L9E2</b>   |
| Curette Gracey Labriola-Cortellini #13/14 DE | <b>SG13/149LCE2</b> |
| Curette Gracey Labriola-Cortellini #11/12    | <b>SG11/129LCE2</b> |

**Sugg. Retail** **Promo**

**zł 1.811** **zł 1.177**



For more information  
on **EVEREDGE™ 2.0**  
visit our site [hu-friedy.eu](http://hu-friedy.eu)



### STANDARD GRACEY CURETTE HU-FRIEDY

| Description           | Handle 9    | Resin 8 Color |
|-----------------------|-------------|---------------|
| 1/2 Standard Gracey   | SG1/29IE2   | SG1/2C8E2     |
| 7/8 Standard Gracey   | SG7/897E2   | SG7/8C8E2     |
| 11/12 Standard Gracey | SG11/1293E2 | SG11/12C8E2   |
| 13/14 Standard Gracey | SG13/1498E2 | SG13/14C8E2   |

### MINI FIVE GRACEY CURETTES HU-FRIEDY

| Description            | Handle 9     | Resin 8 Color |
|------------------------|--------------|---------------|
| 1/2 Mini Five Gracey   | SAS1/29IE2   | SAS1/2C8E2    |
| 7/8 Mini Five Gracey   | SAS7/897E2   | SAS7/8C8E2    |
| 11/12 Mini Five Gracey | SAS11/1293E2 | SAS11/12C8E2  |
| 13/14 Mini Five Gracey | SAS13/1498E2 | SAS13/14C8E2  |

### UNIVERSAL CURETTES HU-FRIEDY

| Description  | Handle 9   | Resin 8 Color |
|--------------|------------|---------------|
| 1/2 Langer   | SL1/29E2   | SL1/2C8E2     |
| 3/4 Langer   | SL3/49E2   | SL3/4C8E2     |
| 5/6 Langer   | SL5/69E2   | SL5/6C8E2     |
| 17/18 Langer | SL17/189E2 | SL17/18C8E2   |

### SICKLE SCALERS HU-FRIEDY

| Description             | Handle 9   | Resin 8 Color |
|-------------------------|------------|---------------|
| H6/H7 Scaler            | SH6/79E2   | SH6/7C8E2     |
| 204S Scaler             | S204S9E2   | S204SC8E2     |
| 204SD Scaler            | S204SD9E2  | S204SDC8E2    |
| 1 Nevi Anterior         | SCNEVI19E2 | SCNEVI1C8E2   |
| 2 Nevi Posterior        | SCNEVI29E2 | SCNEVI2C8E2   |
| 3 Nevi Scaler Posterior | SCNEVI39E2 | SCNEVI3C8E2   |
| 4 Nevi Scaler Posterior | SCNEVI49E2 | SCNEVI4C8E2   |

Valid for all **EVEREDGE 2.0** range



#### Handle 9

**Sugg. Retail**

**zł 274**

**Promo**

**zł 206**




#### Handle 8 Color

**Sugg. Retail**

**zł 261**

**Promo**

**zł 195**



\*Products tested: Hu-Friedy EverEdge™ 2.0 H7, American Eagle XP H6-7, P.D.T. H6-7, LM/EM/EN H6-7, and G.M. Hartzell H6-7. For more information, contact David Platt, Hu-Friedy Mfg. Co., LLC, 3232 N. Rockwell St., Chicago, IL 60618.

# SCISSORS & NEEDLE HOLDERS

HU-FRIEDY

**DID YOU KNOW?** All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like Hu-Friedy Penetrating Oil (IPS) and Instrument Lubricant Spray (ILS), will prevent rust, corrosion and stiff joints to ensure smooth operation. All hinged instruments should be sterilised in the open position.



## NH5038

PermaSharp™ Needle Holder Crile-Wood #5038, 15 cm

Sugg. Retail zł 717

Promo zł 502



## NHCW

Needle Holder Crile-Wood grooved, 3-0 - 6-0, 15cm

Sugg. Retail zł 215

Promo zł 150



## NMH

Needle Holder Mayo-Hegar grooved, 1-0 - 4-0 Needle, 16cm

Sugg. Retail zł 236

Promo zł 165



## NH5024CORT

Needle Holder Cortellini Micro Castro #5024 round handle, straight

Sugg. Retail zł 2.082

Promo zł 1.457



## NH5024SIM

PermaSharp™ Needle Holder Simion #5024R round handle, straight, 18cm

Sugg. Retail zł 2.082

Promo zł 1.457



## S17SC

SuperCut™ Scissors Iris, #17, straight, 11,5 cm

Sugg. Retail zł 418

Promo zł 292



## S18SC

SuperCut™ Scissors Iris #18 fine angled, 11,5cm

Sugg. Retail zł 418

Promo zł 292



## S14SC

SuperCut™ Scissors LaGrange #14 double angled, 11,5cm

Sugg. Retail zł 534

Promo zł 374



## S16

SuperCut™ Scissors Iris #18 fine angled, 11,5 cm **Black Line**

Sugg. Retail zł 240

Promo zł 168



## SCORTSC

SuperCut™ Scissors Cortellini curved, 12cm

Sugg. Retail zł 520

Promo zł 364



## S5083

Scissors Iris #5083 angled, PermaSharp™, 11,5 cm

Sugg. Retail zł 942

Promo zł 659

### SUPER CUT™

Cutting blade honed/serrated to „razor“ sharpness with opposing blade finished with fine serrations to avoid tissue slippage.

### PERMA SHARP™

Featuring tungsten carbide inserts bonded to surgical-grade stainless steel on both cutting edges for smooth, controlled cutting and maximum durability.



# PERIOSTEALS

## PERIOSTEAL HU-FRIEDY




**P24G**  
Periosteal #24  
long handle

**Sugg. Retail**

zł 380

**Promo**

zł 285



**PPR3S**  
Periosteal Prichard  
#3 hdl #522 short

**Sugg. Retail**

zł 450

**Promo**

zł 337



**PPBUSER6**  
Periosteal Buser  
hdl #6

**Sugg. Retail**

zł 436

**Promo**

zł 327




**P146**  
Periosteal  
Goldman-Fox  
#14 hdl #6

**Sugg. Retail**

zł 390

**Promo**

zł 293




**CL84**  
Surgical Curette  
Lucas #84, hdl #522  
2,2mm

**Sugg. Retail**

zł 404

**Promo**

zł 283




**CL85**  
Surgical Curette  
Lucas #85, hdl #522  
2,5mm

**Sugg. Retail**

zł 404

**Promo**

zł 283




**CL86**  
Surgical Curette  
Lucas #86, hdl #522  
2,8mm

**Sugg. Retail**

zł 404

**Promo**

zł 283



**CL87**  
Surgical Curette  
Lucas #87, hdl #522  
3,5mm

**Sugg. Retail**

zł 404

**Promo**

zł 283



**CM2/4**  
Surgical Curette Molt  
#2/4, hdl #522

**Sugg. Retail**

zł 408

**Promo**

zł 285

## SURGICAL CURETTES HU-FRIEDY

# ENDODONTIC INSTRUMENTS

## EXCEPTIONAL RESTORATION STARTS WITH EXCELLENT INSULATION



**RDPA**  
Rubber Dam Punch  
Ainsworth  
17 cm

Sugg. Retail zł 1.082

**Promo zł 865**



**RDF4E**  
Rubber Dam Clamp  
Forceps #4  
European

Sugg. Retail zł 643


**Promo zł 514**



**RDCSET7**  
Rubber Dam Clamps  
Starter Kit with 7 clamps

Sugg. Retail zł 524

**Promo zł 419**



**RDCOB**  
IMS™ Rubber Dam Clamp  
Organisation Board  
for 8 clamps

Sugg. Retail zł 372


**Promo zł 298**



**RDAF6**  
Rubber Dam Frame  
Adult, 15cm

Sugg. Retail zł 176

**Promo zł 141**



**RDCF5**  
Rubber Dam Frame  
Child  
13cm

Sugg. Retail zł 176

**Promo zł 141**



**IMS-1273**  
IMS™ Accessory Box  
large

Sugg. Retail zł 271


**Promo zł 217**



**IMS-1271**  
IMS™ Accessory Box  
small

Sugg. Retail zł 200

**Promo zł 160**



**IMEDRDAM**  
IMS™ Infinity DIN Collection Rubber Dam Kit


Sugg. Retail zł 3.686

**Promo zł 2.949**

# INSTRUMENT MANAGEMENT SYSTEM

## IMS™ INFINITY SERIES CASSETTES

IMS Infinity Series Cassettes are made of long-lasting, low-maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle. Available in over 100 different sizes, configurations and color options, cassettes can be specialized to meet your needs.



### 1/4 DIN CASSETTE HU-FRIEDY IME14DIN5\_

IMS Cassette 1/4 DIN for 5 instruments  
Dimensions: 72 x 184 x 34 mm  
(without instruments)

**Sugg. Retail**

zł 615

**Promo**

zł 492



### 1/2 DIN CASSETTE HU-FRIEDY IME12DIN1\_


IMS Cassette 1/2 DIN for 10 instruments  
Dimensions: 152 x 184 x 34 mm  
(without instruments)

**Sugg. Retail**

zł 689

**Promo**

zł 551



### DIN SURGERY CASSETTE HU-FRIEDY IMEDIN18\_S

DIN Cassette for 18 instruments  
with an accessory area + 4 clips  
Dimensions: 290 x 184 x 34 mm  
(without instruments)

**Sugg. Retail**

zł 1.001

**Promo**

zł 801

## IMS SIGNATURE SERIES™ TUBS HU-FRIEDY IMS-140\_

Together, procedure tubs and cassettes create a more efficient and organised system for managing instruments and consumable products.

### Includes:

- Signature Series Tub
- Tub Liner
- Tub Divider (package of 8)
- Single Tub Cup w/Cover
- Double Tub Cup w/Cover
- Divided Slide Tray
- Locking Tub Cover



**Sugg. Retail**

zł 446

**Promo**

zł 357

### CHOOSE THE COLOR.

Just add the color coding number to the last part of the partcode (replace " \_ " with the number)

|        |   |   |            |   |    |
|--------|---|---|------------|---|----|
| Red    | ● | 0 | Lavender   | ● | 6  |
| Grey   | ● | 1 | Black      | ● | 7  |
| White  | ● | 2 | Blue       | ● | 8  |
| Orange | ● | 3 | Green      | ● | 9  |
| Purple | ● | 4 | Ocean Blue | ● | 11 |
| Yellow | ● | 5 |            |   |    |

# INFECTION CONTROL

## DENTAPURE™ CARTRIDGE CROSSTEX\*\*

**REVOLUTIONIZE YOUR WATERLINE MAINTENANCE!  
DAILY WATERLINE MAINTENANCE BECOMES A SIMPLE, ANNUAL ROUTINE\*\***

Benefits for both independent water bottle systems & municipal systems plumbed to city water:

- NO routine shocking
- NO distilled water required
- NO harsh chemicals
- NO silver
- NO special disposal requirements
- NO allergenic iodine proteins<sup>1</sup>
- NO concerns with dental materials and bonding<sup>2</sup>
- NO concerns with dental unit corrosion or etching<sup>3</sup>
- An effective method of treating municipally plumbed dental units

For more information on how the multi-award-winning DentaPure™ Cartridge can help simplify your DUWL maintenance protocols, visit: [DentaPure.Crosstex.com/EU](http://DentaPure.Crosstex.com/EU)

| DESCRIPTION                        | QUANTITY              | CODE      | Sugg. Retail | Promo  |
|------------------------------------|-----------------------|-----------|--------------|--------|
| Independent Water Bottle Cartridge | 1 Each                | DPI365B10 | zł 1.430     | zł 858 |
| Municipal Cartridge                | 1 Each                | DPI365M10 | zł 1.430     | zł 858 |
| Iodine Test Strips                 | 50 Test Strips/Bottle | DPTST     | zł 200       | zł 158 |

### Municipal Cartridge DPI365M10

For waterlines plumbed directly to city water. Installs onto incoming lines with included fittings.



**ONE  
DENTAPURE™  
CARTRIDGE  
= 365 DAYS  
COMPLIANT,<sup>†</sup>  
DENTAL UNIT  
WATER\***

### Independent Water Bottle Cartridge DPI365B10

For independent water bottles. Easily installed by office staff onto existing pickup tube with included fitting.

<sup>†</sup>Compliant = Compliant with EPA potable water standard <500 CFU/mL. <sup>††</sup>Or, 240L of water if usage records are kept. <sup>1</sup>Schabelman E., Witting M. The relationship of radiocontrast, iodine, and seafood allergies: a medical myth exposed, J Emerg Med. 2010 Nov;39(5):701-7. doi: 10.1016/j.jemermed.2009.10.014. Epub 2010 Jan 4. <sup>2</sup>Puttaiah R. Effects of Low Grade Iodine in Dental Unit Waterlines On Shear Bond Strength of a Dentin Bonding Agent, Baylor College of Dentistry. <sup>3</sup>Data on file at Crosstex.

## BACTRYL™ DISINFECTANT CROSSTEX\*\*

**READY-TO-USE SPRAY AND WIPES FOR THE CLEANING AND  
DISINFECTION OF SURFACES AND DENTAL MEDICAL DEVICES**

- Ready-to-use, broad spectrum detergent and disinfectant (One Step)
- Fully tested and compliant to most rigid EU Standards (EN14885)
- Wipes are Compliant to new compulsory EN 16615 in clean and dirty conditions.
- VAH-listed Disinfectant
- Molecule effective against gram positive and gram negative bacteria, mycobacteria, bacteria, viruses (HIV, HBV, HCV) and fungi. Effective on all enveloped viruses.
- Fully compatible with surfaces such as worktops, lamps, trolleys, dental units, trays, sinks, instruments and furnishings of dental surgery with long-lasting residual bacteriostatic action
- Safe for dental professionals and patients
- Available as a spray or wipe

**500046**

4 x 1000 ml Bottles  
+ 2 x Spray  
Triggers/Case

Sugg. Retail zł 326

**Promo zł 196**



**500047**

6 x 110 Wipes/  
Container

Sugg. Retail zł 512

**Promo zł 410**



Cantel Medical Italy Srl is the legal manufacturer



# PROCEDURAL KITS

## BASIC IMPLANT SET OMNIA\*

The Basic Implant Set has been designed for all routine surgical activities where a surgical protocol ensuring an adequate degree of sterility and guaranteeing the safety of operators and patients is needed. Due to its reasonable price, it is especially fit for all surgical procedures with immediate loading or for operations, such as re-opening of implants, flap surgery, operations on minor surgical wounds.



**12.S5621.00**

Basic Implant Set

### DESCRIPTION

|   |   |                    |             |
|---|---|--------------------|-------------|
| 1 | Surgically folded Special Plus model gown (size L) and 2 hand towels wrapped in medical paper | SMS Microfibre     | Light blue  |
| 1 | Gown with cotton cuffs  | NWF                | Light blue  |
| 2 | Surgical caps with fixed laces  | Absorbent NWF      | Light blue  |
| 1 | Round surgical cap for the patient  | Absorbent NWF      | Light blue  |
| 1 | Absorbent/liquid proof drape <b>50x75 cm</b> with self-adhesive <sup>†</sup> side             | NWF + Polyethylene | Light blue  |
| 1 | Absorbent/liquid proof drape 75x90 cm with adhesive <sup>†</sup> <b>U-shaped cut 11x9 cm</b>  | NWF + Polyethylene | Light blue  |
| 1 | Omnisleeve sheath 120x7 cm  | PE                 | Transparent |
| 1 | Disposal bag for contaminated material  | PVC                | Red         |

<sup>†</sup> **Bold numbers** indicate the **adhesive** side

**Sugg. Retail**

**Promo**

zł 126

zł 113

## PERIODONTAL SET OMNIA\*

Sterile set for periodontal surgery or small oral surgery. It contains patient drapes and surface drapes. It includes sheaths for hand pieces, adhesive films for handles, keypads and folded compresses. The set is packed in a double packaging. The practical peel-off label allows easy traceability of the operation.



**12.T0013.00**

Periodontal Set

### DESCRIPTION

|   |   |                         |             |
|---|---|-------------------------|-------------|
| 1 | Adhesive <sup>†</sup> <b>U shaped 11 x 9 cm</b> absorbent/water-proof drape cm 75 x 90 cm | Non woven +Polyethylene | Blue        |
| 1 | Absorbent/water proof drape 75 x 90 cm  | Non woven +Polyethylene | Blue        |
| 2 | Absorbent/water proof drapes 50 x 75 cm   | Non woven +Polyethylene | Blue        |
| 2 | Tubing sleeves 120 x 7 cm with fixing adhesives   | LDPE                    | Transparent |
| 2 | Folded compresses 5 x 5 cm, 4-ply   | Non woven               | White       |

<sup>†</sup> **Bold numbers** indicate the **adhesive** side

**Sugg. Retail**

**Promo**

zł 80

zł 72

## IMPLANTOLOGY SET OMNIA\*

Complete sterile set conceived for all Implantology, Parodontology and Maxillo Facial surgeries. The set is composed of patient drape, surface drapes and tubing sleeves. The set is packed in sterile wrapping.

**The box contains 10 sets.**



**12.T0202.00**

Implantology Set

### DESCRIPTION

|   |  |         |             |
|---|--|---------|-------------|
| 1 | Surgical drape 60 x 60 cm                      | SMS     | Light blue  |
| 1 | Surgical Drape 150 x 150 cm with adhesive hole | NWF/TNT | Light blue  |
| 2 | Surgical Drape 50 x 75 cm                      | PP+PE   | Light blue  |
| 2 | Adesivi medicali 9 x 51 cm                     | SPPB    | Light blue  |
| 2 | Adhesive Velcro Tape 14 x 2,5 cm               | Velcro  | White       |
| 2 | 2 Tube Cover 6 x 120 cm                        | PE      | Transparent |

**Sugg. Retail**

**Promo**

zł 887


zł 798



# SURGERY: THE ESSENTIALS

## SUCTION TIPS OMNIA\*

Omnia disposable suction tips for surgical aspirators are made of medical grade plastic material. Suction tips ensure a targeted and safe aspiration during surgical operations. They are available with differently shaped hand-grips and terminal ends according to the particular needs of the dental surgeon, and they are supplied without connection tubing.



| CODE               | DESCRIPTION   | Sugg. Retail | Promo  |
|--------------------|---|--------------|--------|
| <b>32.F5051.00</b> | Aspirator, 2.50 m, with suction tip with ergonomic hand-grip, and that is also enabled for the autologous bone harvesting trap (each pack contains 1 adapter Adat. 05), PVC, 5 pcs. | zł 229       | zł 183 |
| <b>32.F6565.00</b> | Omniasurg model sterile suction tip (each pack contains 1 adapter Adat. 08) (tip length: 20 cm), PVC, 20 pcs.   | zł 143       | zł 114 |
| <b>32.F2065.00</b> | Surgical aspirator, 2.20 m, with Omnia surg suction tip (each pack contains 1 adapter Adat. 05), PVC, 10 pcs.   | zł 257       | zł 206 |



## VISUAL OMNIA\*

The Visual is an individual protection device available to ensure your safety. Its wide polyethylene fog preventing shield and the polyurethane insert make this device also fit for use by people wearing glasses, 12 pcs.

**10.M0200.00**

Visual

| Sugg. Retail | Promo  |
|--------------|--------|
| zł 172       | zł 154 |



## SHELL COVER OMNIA\*

The SHELL COVER gown cover is perfect for routine practice that create high aerosol dispersion levels. Material: Medical grade LDPE waterproof, with high insulation standard. The additional Thumb-Up system helps keeping the gown fixed while worn.

**30.D1002.00**

Shell Cover - Not Sterile 100 pcs

| Sugg. Retail | Promo  |
|--------------|--------|
| zł 1.072     | zł 965 |

**22.D1002.00**

Shell Cover - Sterile 50 pcs

| Sugg. Retail | Promo    |
|--------------|----------|
| zł 1.373     | zł 1.235 |

# IRRIGATION SETS OMNIA\*

## PHYSIODISPENSER AND PIEZOELECTRIC SURGERY

- Proper irrigation during the creation of the implant site is a fundamental element for successful bone integration of the implant and the overall outcome of the operation.
- Wide range of mechanical irrigation lines which are compatible with most implant motor models. Irrigation systems include a vented spike, a roller or clamp for flow control, fittings, and an integrated peristaltic pump section.
- Convenient terminal cap makes the connection of the irrigation line to the hand piece easy to achieve. The cap allows you to connect tubes to the wide variety of cooling needle models available.
- The mechanical irrigation systems offer the following options:
  - a. Internal irrigation only
  - b. External irrigation only
  - c. Internal and external irrigation



| CODE<br>(10 Pcs, par boîte) | DESCRIPTION   | Sugg. Retail | Promo         |
|-----------------------------|---|--------------|---------------|
| <b>32.F1201.00</b>          | Irrigation Tubing for NSK Surgic XT, Surgic XT Plus, Osseocision, Vario Surg, Kyocera Implantor Neo, Aseptico Aeu 6000 series, Aseptico Aeu 7000 series, Bonart Piezo surgical system   | zł 629       | <b>zł 503</b> |
| <b>32.F0134.00</b>          | NSK Surgic XT, Surgic XT Plus, Osseocision®, VarioSurg, SurgicPro, Kyocera Implantor Neo, Bonart Piezo surgical system, Dentium iCTmotor, NSK 3i Osseocision, Delfine Intra-Lock, 1/2Y  | zł 440       | <b>zł 352</b> |
| <b>32.F0139.00</b>          | Nobel Biocare80/ W&H Implantmed/ Implant Med new design 2005/ Implantmed SI 923/ Astratech Implant SI 923/ Surgical Motor System/Ace/ Straumann/InnovaCorp. 3i/Centerpulse/Lifecore/ DU 900/FRIOS® UNIT S/ FRIOS® UNIT S/I, Osseoset 100/Osseoset 200/Osseocare, Elcomed SA 310/SA2006, Piezomed, W&H SI 1023, 1/2Y | zł 440       | <b>zł 352</b> |
| <b>32.F0170.00</b>          | Bien Air Chiropro 980 / EMS PIEZON® MASTER surgery / Esacrom Surgisonic I / UBS-Ultrasonic bone surgery/ Surgybone - Acrobone - Easybone Silfradent, BioSAFin Easy Surgery, MyTUTOR, Heracle Sweden&Martina, Cefla dental Chair, Osseodoc BienAir, 1/2Y   | zł 440       | <b>zł 352</b> |
| <b>32.F1204.00</b>          | Nobel Biocare80/ W&H Implantmed/ Implant Med new design 2005/ Implantmed SI 923/ Astratech Implant SI 923/ Surgical Motor System/Ace/ Straumann/InnovaCorp. 3i/Centerpulse/Lifecore/ DU 900/FRIOS® UNIT S/ FRIOS® UNIT S/I, Drilltech/Osseoset 100/Osseoset 200/Osseocare, 1/2Y                                     | zł 929       | <b>zł 503</b> |



**FOR MORE  
INFORMATION  
ON OUR  
IRRIGATION  
SET.**



For more information on all available irrigation lines, see:

<https://www.omniaspa.eu/en/products/oral-irrigation-kits.html>



@HuFriedyEU



@hu\_friedy\_europe



@Hu-Friedy Mfg. Co., LLC.EU



@OmniaSPA



@Omnia S.p.A.

\*Products provided by Omnia. \*\* Products provided by Crosstex. All other products provided by Hu-Friedy.

For more information on Hu-Friedy products or to place an order, call at 00800 48 37 43 39 or visit [Hu-Friedy.eu](https://Hu-Friedy.eu)

For more information on Crosstex products, call +39 05 24 52 74 53 or visit [crosstex.com](https://crosstex.com)

For more information on Omnia products or to place an order, call +39 05 24 52 74 53 or visit [omniaspa.eu](https://omniaspa.eu)



Biuro Obsługi Klienta  
+48 22 338 70 50

[www.medif.com](https://www.medif.com)

As some of the products included in this catalogue might not be available in some countries, we invite you to verify the availability of the product with your local authorised dealer or to contact our customer care at the free phone number 00800 48 37 43 39. All prices indicated are intended VAT included. All company and product names are trademarks of Hu-Friedy Mfg. Co. LLC, its affiliates or related companies, unless otherwise noted.

©2020 Hu-Friedy Mfg. Co., LLC All rights reserved. HFL-070/1020