



**Conical Connection**  
Product Catalog

**2021**

The MIS vision, to be the preferred choice of dentists worldwide by simplifying their world, has guided the company on its way. Through creative thinking, passion and persistence, driven to be better while maintaining its core values, MIS Implants Technologies Ltd. has evolved into the innovation hub of expertise that it is today, and is building its future through excellence.

In this catalog, you'll find the comprehensive range of MIS implants and their corresponding prosthetic components and tools, as well as regenerative solutions.

The MIS Quality System complies with international quality standards: ISO 13485:2016 - Quality Management System for Medical Devices and Medical Device Directive 93/42/EEC. Please note, not all products are registered or available in every country/region.

IFUs for MIS products may be found at: <https://ifu.mis-implants.com>. Adobe Acrobat is required to view the IFU file on the website. This software may be freely downloaded from the Adobe website



Additional information may be added to the product descriptions without prior notice.  
© MIS Implants Technologies Ltd. All rights reserved.



# TABLE OF CONTENTS



4 - 24.

## IMPLANT SYSTEMS

|             |    |
|-------------|----|
| V3          | 6  |
| C1          | 14 |
| B+ Implants | 22 |
| LINK        | 24 |



25 - 60.

## TOOLS AND KITS

|                                       |    |
|---------------------------------------|----|
| Marking Drills                        | 28 |
| Pilot Drills, External Irrigation     | 29 |
| Step Drills, External Irrigation      | 30 |
| Countersinks                          | 31 |
| V3 Insertion Tools                    | 32 |
| C1 Insertion Tools                    | 34 |
| Surgical Tools                        | 36 |
| Prosthetic Instruments                | 38 |
| Conical Connection Surgical Kit       | 41 |
| Drill Replacement kit                 | 46 |
| Bone Profiler Kit                     | 47 |
| Drill stoppers kit, standard platform | 48 |
| Drill Stoppers Kit, Depth-Based       | 49 |
| Bone Compression Kit                  | 50 |
| Sinus Lift Osteotome Set              | 52 |
| Convex Osteotome Set                  | 53 |
| Prosthetic Tools Kit                  | 54 |
| Abutment Holding System Kit           | 57 |
| Crown Set                             | 58 |
| Zircon Diamond Burr Set               | 59 |
| Control Units                         | 60 |





62 - 97.

#### PROSTHETIC OPTIONS

|                                |    |
|--------------------------------|----|
| Analogs and Impression Copings | 64 |
| Temporary Cylinders            | 66 |
| CONNECT System                 | 68 |
| Multi-Unit System              | 73 |
| Gold Plastic Cylinders         | 76 |
| EZ-Base System                 | 77 |
| Ti-Bases                       | 78 |
| Scan Post Kit                  | 80 |
| Digital Model Analogs          | 81 |
| Titanium Blanks                | 82 |
| Cementable Abutments           | 83 |
| Cementing Posts                | 85 |
| CPK - Complete Prosthetic Kit  | 86 |
| LOCKit Attachment System       | 89 |
| OT-Equator Kits                | 92 |
| Ball Attachment System         | 94 |
| Bars                           | 96 |



98 - 111.

#### MGUIDE

|                                    |     |
|------------------------------------|-----|
| MGUIDE Custom Surgical Template    | 100 |
| MGUIDE Surgical Set                | 101 |
| MGUIDE Surgical Set, Narrow Sleeve | 106 |
| MGUIDE 16mm Drills Kit             | 110 |



112 - 117.

#### REGENERATIVE SOLUTIONS

|           |     |
|-----------|-----|
| 4MATRIX+  | 114 |
| BONDBONE  | 115 |
| 4BONE BCH | 116 |
| 4BONE RCM | 117 |



118 - 125.


#### ADDITIONAL INFORMATION

|                                   |     |
|-----------------------------------|-----|
| Pilot Drills, External Irrigation | 120 |
| Step Drills, External Irrigation  | 121 |
| Restorative Components            | 122 |
| Alloy Composition                 | 123 |



126 - 133.

#### INDEX



## IMPLANT SYSTEMS

MIS proudly offers an extensive range of implant systems, providing simple solutions for a wide variety of clinical scenarios, from basic to complex. Our implants are the outcome of an exceptionally high-level R&D process resulting in implant systems that are simple, yet offer enhanced functionality and performance features. MIS implants are made of biocompatible, medical grade titanium alloy, with a surface of the highest quality. They are treated by both sand-blasting and acid-etching, to promote extraordinary cleanliness and facilitate long-lasting osseointegration.






Implant Systems

## V3 Implant Conical Connection Ø3.30 NP



Implant Cover Screw  
and Healing Caps:




Cover screw

|                                                                                   | Cat. No. | Dimensions (mm) | Description                                       | Material                        | Item includes              |
|-----------------------------------------------------------------------------------|----------|-----------------|---------------------------------------------------|---------------------------------|----------------------------|
|  | V3-10330 | Ø3.30<br>L=10   | V3 conical connection implant,<br>narrow platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill |
|  | V3-11330 | Ø3.30<br>L=11.5 | V3 conical connection implant,<br>narrow platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill |
|  | V3-13330 | Ø3.30<br>L=13   | V3 conical connection implant,<br>narrow platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill |
|  | V3-16330 | Ø3.30<br>L=16   | V3 conical connection implant,<br>narrow platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill |

L = Length

Ø3.30mm Implant Procedure

| Drilling Speed<br>(RPM) | OR           |             |             |                                          |                                          | Max.<br>Torque:<br>45Ncm |
|-------------------------|--------------|-------------|-------------|------------------------------------------|------------------------------------------|--------------------------|
|                         | 800-<br>1000 | 600-<br>800 | 200-<br>400 | 200-<br>400                              | Ø3.30                                    |                          |
| Diameter                | Ø1.90        | Ø2.40       | Ø2.40       | Final drill for<br>bone types<br>1 and 2 | Final drill for<br>bone types<br>3 and 4 | Ø3.30                    |



The drilling sequence is illustrated using 11.50mm implants.

Procedure recommended by MlS cannot replace the judgment and professional experience of the surgeon.



Implant Systems

V3 Implant  
Conical Connection  
Ø3.90 SP



Implant Cover Screw  
and Healing Caps:



CS-H0339  
CS-H0439  
CS-H0539  
CS-H0639  
CS-H0839



CS-HS248  
CS-HS348  
CS-HS448  
CS-HS548  
CS-HS648  
CS-HS848




VS-HS358  
VS-HS458  
VS-HS558  
VS-HS658  
VS-HS858



|                                                                                   | Cat. No. | Dimensions (mm) | Description                                         | Material                        | Item includes              |
|-----------------------------------------------------------------------------------|----------|-----------------|-----------------------------------------------------|---------------------------------|----------------------------|
|  | V3-08390 | Ø3.90<br>L=8    | V3 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill |
|  | V3-10390 | Ø3.90<br>L=10   | V3 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill |
|  | V3-11390 | Ø3.90<br>L=11.5 | V3 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill |
|  | V3-13390 | Ø3.90<br>L=13   | V3 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill |
|  | V3-16390 | Ø3.90<br>L=16   | V3 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill |

L = Length

### Ø3.90mm Implant Procedure



The drilling sequence is illustrated using 11.50mm implants.

Procedure recommended by MIS cannot replace the judgment and professional experience of the surgeon.





Implant Systems

V3 Implant  
Conical Connection  
Ø4.30 SP



Implant Cover Screw  
and Healing Caps:




Cover screw



|                                                                                   | Cat. No. | Dimensions (mm) | Description                                         | Material                        | Item includes              |
|-----------------------------------------------------------------------------------|----------|-----------------|-----------------------------------------------------|---------------------------------|----------------------------|
|  | V3-08430 | Ø4.30<br>L=8    | V3 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill |
|  | V3-10430 | Ø4.30<br>L=10   | V3 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill |
|  | V3-11430 | Ø4.30<br>L=11.5 | V3 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill |
|  | V3-13430 | Ø4.30<br>L=13   | V3 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill |
|  | V3-16430 | Ø4.30<br>L=16   | V3 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill |

L = Length

### Ø4.30mm Implant Procedure



The drilling sequence is illustrated using 11.50mm implants.

Procedure recommended by MIS cannot replace the judgment and professional experience of the surgeon.




Implant Systems

V3 Implant  
Conical Connection  
Ø5 SP



Implant Cover Screw  
and Healing Caps:



|                                                                                    | Cat. No. | Dimensions (mm) | Description                                         | Material                        | Item includes              |
|------------------------------------------------------------------------------------|----------|-----------------|-----------------------------------------------------|---------------------------------|----------------------------|
|   | V3-08500 | Ø5<br>L=8       | V3 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill |
|   | V3-10500 | Ø5<br>L=10      | V3 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill |
|   | V3-11500 | Ø5<br>L=11.5    | V3 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill |
|   | V3-13500 | Ø5<br>L=13      | V3 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill |
|  | V3-16500 | Ø5<br>L=16      | V3 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill |

L = Length

## Ø5mm Implant Procedure

OR

| Drilling Speed<br>(RPM) | 800-<br>1000 | 600-<br>800 |       | 450-<br>650 | 350-<br>550 | 300-<br>500 |    | 200-<br>400                              | 200-<br>400                              | Max.<br>Torque:<br>60Ncm |
|-------------------------|--------------|-------------|-------|-------------|-------------|-------------|----|------------------------------------------|------------------------------------------|--------------------------|
| Diameter                | Ø1.90        | Ø2.40       | Ø2.40 | Ø3          | Ø3.50       | Ø4          | Ø4 | Final drill for<br>bone types<br>1 and 2 | Final drill for<br>bone types<br>3 and 4 | Ø5                       |





The drilling sequence is illustrated using 11.50mm implants.

Procedure recommended by MIS cannot replace the judgment and professional experience of the surgeon.



Implant Systems

## C1 Implant Conical Connection Ø3.30 NP



Implant Cover Screw  
and Healing Caps:



CN-H0333  
CN-H0433  
CN-H0533  
CN-H0633  
CN-H0833



CN-HS340  
CN-HS440  
CN-HS540  
CN-HS640  
CN-HS840



CN-HA248  
CN-HA348  
CN-HA448  
CN-HA548  
CN-HA648  
CN-HA848



Cover screw

|                                                                                   | Cat. No. | Dimensions (mm) | Description                                       | Material                        | Item includes                                         |
|-----------------------------------------------------------------------------------|----------|-----------------|---------------------------------------------------|---------------------------------|-------------------------------------------------------|
|  | C1-10330 | Ø3.30<br>L=10   | C1 conical connection implant,<br>narrow platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |
|  | C1-11330 | Ø3.30<br>L=11.5 | C1 conical connection implant,<br>narrow platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |
|  | C1-13330 | Ø3.30<br>L=13   | C1 conical connection implant,<br>narrow platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |
|  | C1-16330 | Ø3.30<br>L=16   | C1 conical connection implant,<br>narrow platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |

L = Length

Ø3.30mm Implant Procedure

| Drilling Speed<br>(RPM) | 1200-<br>1500 | 900-<br>1200 |       | 200-<br>400 | 15-25 |
|-------------------------|---------------|--------------|-------|-------------|-------|
| Diameter                | Ø1.90         | Ø2.40        | Ø2.40 | Ø3<br>Ø3.20 | Ø3.30 |

Final drill for  
bone types  
1 and 2



Do not use the final drill for bone type 4.

The drilling sequence is illustrated using a 13mm implant.

Procedures recommended by MIS cannot replace the judgment and professional experience of the surgeon.






Implant Systems

C1 Implant  
Conical Connection  
Ø3.75 SP



Implant Cover Screw  
and Healing Caps:



|                                                                                   | Cat. No. | Dimensions (mm) | Description                                         | Material                        | Item includes                                         |
|-----------------------------------------------------------------------------------|----------|-----------------|-----------------------------------------------------|---------------------------------|-------------------------------------------------------|
|  | C1-08375 | Ø3.75<br>L=8    | C1 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |
|  | C1-10375 | Ø3.75<br>L=10   | C1 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |
|  | C1-11375 | Ø3.75<br>L=11.5 | C1 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |
|  | C1-13375 | Ø3.75<br>L=13   | C1 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |
|  | C1-16375 | Ø3.75<br>L=16   | C1 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |

L = Length

## Ø3.75mm Implant Procedure

| Drilling Speed<br>(RPM) | 1200-<br>1500 | 900-<br>1200 |       | 500-<br>700 |    | 200-<br>400 | 15-25 |
|-------------------------|---------------|--------------|-------|-------------|----|-------------|-------|
| Diameter                | Ø1.90         | Ø2.40        | Ø2.40 | Ø3          | Ø3 | Ø3<br>Ø3.60 | Ø3.75 |

Final drill for  
bone types  
1,2 and 3



Do not use the final drill for bone type 4.


The drilling sequence is illustrated using a 13mm implant.

Procedures recommended by MIS cannot replace the judgment and professional experience of the surgeon.



Implant Systems

C1 Implant  
Conical Connection  
Ø4.20 SP



Implant Cover Screw  
and Healing Caps:






|                                                                                    | Cat. No. | Dimensions (mm) | Description                                         | Material                        | Item includes                                         |
|------------------------------------------------------------------------------------|----------|-----------------|-----------------------------------------------------|---------------------------------|-------------------------------------------------------|
|   | C1-08420 | Ø4.20<br>L=8    | C1 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |
|   | C1-10420 | Ø4.20<br>L=10   | C1 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |
|   | C1-11420 | Ø4.20<br>L=11.5 | C1 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |
|   | C1-13420 | Ø4.20<br>L=13   | C1 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |
|  | C1-16420 | Ø4.20<br>L=16   | C1 conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |

L = Length

## Ø4.20mm Implant Procedure



Do not use the final drill for bone type 4.

The drilling sequence is illustrated using a 13mm implant.

Procedures recommended by MIS cannot replace the judgment and professional experience of the surgeon.



Implant Systems

C1 Implant  
Conical Connection  
Ø5 WP



Implant Cover Screw  
and Healing Caps:



CW-HA263  
CW-HA363  
CW-HA563



CW-HS355  
CW-HS455  
CW-HS555  
CW-HS655



Cover screw


|                                                                                   | Cat. No. | Dimensions (mm) | Description                                  | Material                        | Item includes                                         |
|-----------------------------------------------------------------------------------|----------|-----------------|----------------------------------------------|---------------------------------|-------------------------------------------------------|
|  | C1-08500 | Ø5<br>L=8       | C1 conical connection implant, wide platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |
|  | C1-10500 | Ø5<br>L=10      | C1 conical connection implant, wide platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |
|  | C1-11500 | Ø5<br>L=11.5    | C1 conical connection implant, wide platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |
|  | C1-13500 | Ø5<br>L=13      | C1 conical connection implant, wide platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |
|  | C1-16500 | Ø5<br>L=16      | C1 conical connection implant, wide platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |

L = Length

## Ø5mm Implant Procedure

| Drilling Speed (RPM) | 1200-1500 | 900-1200 |       | 500-700 | 400-700 | 400-600 |    | 200-400        | 15-25 |
|----------------------|-----------|----------|-------|---------|---------|---------|----|----------------|-------|
| Diameter             | Ø1.90     | Ø2.40    | Ø2.40 | Ø3      | Ø3.50   | Ø4      | Ø4 | Ø4.10<br>Ø4.90 | Ø5    |

Final drill for  
bone types  
1,2 and 3



Do not use the final drill for bone type 4.

The drilling sequence is illustrated using a 13mm implant.

Procedures recommended by MIS cannot replace the judgment and professional experience of the surgeon.



Implant Systems

## B+ Implants





|                                                                                     | Cat. No.  | Dimensions (mm) | Description                                            | Material                        | Item includes                                         |
|-------------------------------------------------------------------------------------|-----------|-----------------|--------------------------------------------------------|---------------------------------|-------------------------------------------------------|
|    | V3B-10330 | Ø3.30, L=10     | V3 B+ conical connection implant,<br>narrow platform   | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill                            |
|                                                                                     | V3B-11330 | Ø3.30, L=11.5   |                                                        |                                 |                                                       |
|                                                                                     | V3B-13330 | Ø3.30, L=13     |                                                        |                                 |                                                       |
|                                                                                     | V3B-16330 | Ø3.30, L=16     |                                                        |                                 |                                                       |
|    | V3B-08390 | Ø3.90, L=8      | V3 B+ conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill                            |
|                                                                                     | V3B-10390 | Ø3.90, L=10     |                                                        |                                 |                                                       |
|                                                                                     | V3B-11390 | Ø3.90, L=11.5   |                                                        |                                 |                                                       |
|                                                                                     | V3B-13390 | Ø3.90, L=13     |                                                        |                                 |                                                       |
|                                                                                     | V3B-16390 | Ø3.90, L=16     |                                                        |                                 |                                                       |
|    | V3B-08430 | Ø4.30, L=8      | V3 B+ conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill                            |
|                                                                                     | V3B-10430 | Ø4.30, L=10     |                                                        |                                 |                                                       |
|                                                                                     | V3B-11430 | Ø4.30, L=11.5   |                                                        |                                 |                                                       |
|                                                                                     | V3B-13430 | Ø4.30, L=13     |                                                        |                                 |                                                       |
|                                                                                     | V3B-16430 | Ø4.30, L=16     |                                                        |                                 |                                                       |
|    | V3B-08500 | Ø5, L=8         | V3 B+ conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill                            |
|                                                                                     | V3B-10500 | Ø5, L=10        |                                                        |                                 |                                                       |
|                                                                                     | V3B-11500 | Ø5, L=11.5      |                                                        |                                 |                                                       |
|                                                                                     | V3B-13500 | Ø5, L=13        |                                                        |                                 |                                                       |
|                                                                                     | V3B-16500 | Ø5, L=16        |                                                        |                                 |                                                       |
|   | C1B-10330 | Ø3.30, L=10     | C1 B+ conical connection implant,<br>narrow platform   | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |
|                                                                                     | C1B-11330 | Ø3.30, L=11.5   |                                                        |                                 |                                                       |
|                                                                                     | C1B-13330 | Ø3.30, L=13     |                                                        |                                 |                                                       |
|                                                                                     | C1B-16330 | Ø3.30, L=16     |                                                        |                                 |                                                       |
|  | C1B-08375 | Ø3.75, L=8      | C1 B+ conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |
|                                                                                     | C1B-10375 | Ø3.75, L=10     |                                                        |                                 |                                                       |
|                                                                                     | C1B-11375 | Ø3.75, L=11.5   |                                                        |                                 |                                                       |
|                                                                                     | C1B-13375 | Ø3.75, L=13     |                                                        |                                 |                                                       |
|                                                                                     | C1B-16375 | Ø3.75, L=16     |                                                        |                                 |                                                       |
|  | C1B-08420 | Ø4.20, L=8      | C1 B+ conical connection implant,<br>standard platform | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |
|                                                                                     | C1B-10420 | Ø4.20, L=10     |                                                        |                                 |                                                       |
|                                                                                     | C1B-11420 | Ø4.20, L=11.5   |                                                        |                                 |                                                       |
|                                                                                     | C1B-13420 | Ø4.20, L=13     |                                                        |                                 |                                                       |
|                                                                                     | C1B-16420 | Ø4.20, L=16     |                                                        |                                 |                                                       |
|  | C1B-08500 | Ø5, L=8         | C1 B+ conical connection implant,<br>wide platform     | Titanium Alloy<br>Ti 6Al 4V ELI | Cover screw<br>Final drill<br>Temporary PEEK cylinder |
|                                                                                     | C1B-10500 | Ø5, L=10        |                                                        |                                 |                                                       |
|                                                                                     | C1B-11500 | Ø5, L=11.5      |                                                        |                                 |                                                       |
|                                                                                     | C1B-13500 | Ø5, L=13        |                                                        |                                 |                                                       |
|                                                                                     | C1B-16500 | Ø5, L=16        |                                                        |                                 |                                                       |

L = Length



Implant Systems

## LINK - Orthodontic Anchor Screw



|                                                                                     | Cat. No.  | Dimensions (mm) | Description                               | Material                        | Item includes |
|-------------------------------------------------------------------------------------|-----------|-----------------|-------------------------------------------|---------------------------------|---------------|
|    | MG2-06160 | Ø1.60<br>L=6    | LINK - Orthodontic anchor screw           | Titanium Alloy<br>Ti 6Al 4V ELI |               |
|    | MG2-08160 | Ø1.60<br>L=8    | LINK - Orthodontic anchor screw           | Titanium Alloy<br>Ti 6Al 4V ELI |               |
|    | MG2-10160 | Ø1.60<br>L=10   | LINK - Orthodontic anchor screw           | Titanium Alloy<br>Ti 6Al 4V ELI |               |
|    | MG2-12160 | Ø1.60<br>L=12   | LINK - Orthodontic anchor screw           | Titanium Alloy<br>Ti 6Al 4V ELI |               |
|    | MT-TDO01  | Ø1<br>L=28      | LINK - Orthodontic anchor pilot drill     | Stainless steel                 |               |
|  | MT-SCD20  | L=28            | Short cross driver for orthodontic screws | Stainless steel                 |               |
|  | MT-LCD20  | L=54            | Long cross driver for orthodontic screws  | Stainless steel                 |               |
|  | MT-MSD20  | L=140           | Surgical adapter driver                   | Stainless steel +<br>Aluminum   |               |

L = Length





## TOOLS AND KITS

MIS offers a wide range of surgical instruments and innovative tools to assure smooth and accurate surgical procedures. The range incorporates a wide variety of drills, body try-ins, keys, ratchets and insertion tools.

MIS Surgical Kits bring all the RIGHT tools for each and every implant procedure you may need into convenient, accurate and quality packages, making every job simple, and every surgery as safe and comfortable for the patient as possible. All kits feature an ergonomical kit design, for easy, single-handed opening.

MIS Prosthetic Kits, and all of the complementary products included, help to ensure the final results by way of providing solutions for patient requirements, both functional and esthetic.





Tools and Kits

Marking Drills



|                                                                                     | Cat. No. | Dimensions (mm) | Description         | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|-----------------|---------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | MT-SMD10 | L=30            | Spade marking drill | Stainless steel |   |
|  | MT-PDM24 | Ø2.40<br>L=32   | Position drill mill | Stainless steel |   |
|  | MT-PD440 | Ø4<br>L=32.7    | Position drill      | Stainless steel |   |

L = Length



Tools and Kits

## Pilot Drills, External Irrigation



|                                                                                     | Cat. No. | Dimensions (mm)     | Description                                                       | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|---------------------|-------------------------------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | CT-P2408 | Ø2.40/2.0<br>L=28   | Step pilot drill with built-in stopper for 8mm length implants    | Stainless steel |   |
|  | CT-P2410 | Ø2.40/2.0<br>L=28   | Step pilot drill with built-in stopper for 10mm length implants   | Stainless steel |   |
|  | CT-P2411 | Ø2.40/2.0<br>L=33   | Step pilot drill with built-in stopper for 11.5mm length implants | Stainless steel |   |
|  | CT-P2413 | Ø2.40/2.0<br>L=33   | Step pilot drill with built-in stopper for 13mm length implants   | Stainless steel |   |
|  | CT-P2416 | Ø2.40/2.0<br>L=37.5 | Step pilot drill for 16mm length implants                         | Stainless steel |   |

L = Length



Tools and Kits

Step Drills,  
External Irrigation



|                                                                                     | Cat. No. | Dimensions (mm)   | Description                     | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|-------------------|---------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | CT-TDC30 | Ø3/2.40<br>L=37.5 | Step drill, external irrigation | Stainless steel |   |
|  | CT-TDC35 | Ø3.50/3<br>L=37.5 | Step drill, external irrigation | Stainless steel |   |
|  | CT-TDC40 | Ø4/3.50<br>L=37.5 | Step drill, external irrigation | Stainless steel |   |
|  | CT-TDC45 | Ø4.50/4<br>L=37.5 | Step drill, external irrigation | Stainless steel |   |

L = Length



Tools and Kits

## Countersinks



|                                                                                     | Cat. No. | Dimensions (mm)    | Description                    | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|--------------------|--------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | MT-CSN33 | Ø3.30<br>L=31      | Countersink, narrow platform   | Stainless steel |   |
|  | MT-GDN33 | Ø3.75/4.20<br>L=26 | Countersink, standard platform | Stainless steel |   |
|  | MT-GDN50 | Ø5/6<br>L=26.8     | Countersink, wide platform     | Stainless steel |   |

L = Length



Tools and Kits

## V3 Insertion Tools





|                                                                                     | Cat. No. | Dimensions (mm) | Description                                                          | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|-----------------|----------------------------------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | CT-NLM30 | L=29            | V3 long motor insertion tool, conical connection, narrow platform    | Stainless steel |   |
|  | CT-NSM30 | L=25            | V3 short motor insertion tool, conical connection, narrow platform   | Stainless steel |   |
|  | CT-NLR30 | L=32.4          | V3 long ratchet insertion tool, conical connection, narrow platform  | Stainless steel |   |
|  | CT-NSR30 | L=25            | V3 short ratchet insertion tool, conical connection, narrow platform | Stainless steel |   |
|  | CT-SLM30 | L=29            | V3 long motor insertion tool, conical connection, standard platform  | Stainless steel |   |

L = Length

|                                                                                   | Cat. No. | Dimensions (mm) | Description                                                            | Material        | Sold                                                                                                                                                                    |
|-----------------------------------------------------------------------------------|----------|-----------------|------------------------------------------------------------------------|-----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | CT-SSM30 | L=25            | V3 short motor insertion tool, conical connection, standard platform   | Stainless steel |   |
|  | CT-SLR30 | L=32.4          | V3 long ratchet insertion tool, conical connection, standard platform  | Stainless steel |   |
|  | CT-SSR30 | L=22.3          | V3 short ratchet insertion tool, conical connection, standard platform | Stainless steel |   |

L = Length

## 2 V3 Platforms





Tools and Kits

C1 Insertion Tools



|                                                                                     | Cat. No. | Dimensions (mm) | Description                                                            | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|-----------------|------------------------------------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | CT-NLR10 | L=32            | C1 long ratchet insertion tool, conical connection, narrow platform    | Stainless steel |   |
|  | CT-NSR10 | L=24.5          | C1 short ratchet insertion tool, conical connection, narrow platform   | Stainless steel |   |
|  | CT-NLI10 | L=32            | C1 long motor insertion tool, conical connection, narrow platform      | Stainless steel |   |
|  | CT-NSI10 | L=24.5          | C1 short motor insertion tool, conical connection, narrow platform     | Stainless steel |   |
|  | CT-SLR10 | L=33            | C1 long ratchet insertion tool, conical connection, standard platform  | Stainless steel |   |
|  | CT-SSR10 | L=23            | C1 short ratchet insertion tool, conical connection, standard platform | Stainless steel |   |

L = Length



|                                                                                     | Cat. No. | Dimensions (mm) | Description                                                          | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|-----------------|----------------------------------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | CT-SLI10 | L=33            | C1 long motor insertion tool, conical connection, standard platform  | Stainless Steel |       |
|    | CT-SSI10 | L=25.5          | C1 short motor insertion tool, conical connection, standard platform | Stainless Steel |       |
|    | CT-WLR10 | L=33            | C1 long ratchet insertion tool, conical connection, wide platform    | Stainless steel |       |
|    | CT-WSR10 | L=23            | C1 short ratchet insertion tool, conical connection, wide platform   | Stainless steel |       |
|   | CT-WLI10 | L=33            | C1 long motor insertion tool, conical connection, wide platform      | Stainless steel |       |
|  | CT-WSI10 | L=25.5          | C1 short motor insertion tool, conical connection, wide platform     | Stainless steel |   |

L = Length

## 3 C1 Platforms



Narrow platform



Standard platform





Wide platform




Tools and Kits

Surgical Tools



|                                                                                     | Cat. No. | Dimensions (mm) | Description                           | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|-----------------|---------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | MT-BTI20 | Ø1.40<br>L=100  | Implant site depth probe              | Stainless steel |   |
|  | MT-RDE30 | L=32            | Extra-long driver for 0.05 inch hex.  | Stainless steel |                                                                                        |
|  | MT-RDL30 | L=27            | Long driver for 0.05 inch hex.        | Stainless steel |   |
|  | MT-RDS30 | L=21.3          | Short driver for 0.05 inch hex.       | Stainless steel |   |
|  | MT-RDM30 | L=11.5          | Extra-short driver for 0.05 inch hex. | Stainless steel |                                                                                        |

L = Length

|                                                                                     | Cat. No. | Dimensions (mm) | Description                                         | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|-----------------|-----------------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | MT-DE001 | L=24            | Drill extender                                      | Stainless steel |       |
|    | MT-RE172 | L=28.5          | Int. connection abutment extractor                  | Titanium        |       |
|    | MT-RE160 | L=28.5          | Int. connection abutment extractor, narrow platform | Titanium        |       |
|    | MT-PP240 | Ø2.40-3<br>L=24 | Parallel pin                                        | Stainless steel |                                                                                          |
|   | MT-LM005 | L=30            | Long motor adapter for 0.05 inch hex.               | Stainless steel |     |
|  | MT-SM005 | L=24            | Short motor adapter for 0.05 inch hex.              | Stainless steel |                                                                                        |
|  | MT-SRA10 | L=15.5          | Square connection to ratchet adapter                | Stainless steel |                                                                                        |
|  | MT-TP500 | Ø4.50<br>L=28   | Tissue punch                                        | Stainless steel |                                                                                        |
|  | MT-RI030 | L=84            | Monoblock ratchet                                   |                 |   |
|  | MT-RT070 | L=95            | Surgical torque ratchet                             |                 |   |
|  | MT-RA480 | L=7.5           | Ratchet to hex. adapter                             | Stainless steel |                                                                                        |
|  | MT-MFXGS |                 | MIS MULTIFIX Guide Strip                            | Stainless steel |                                                                                        |

L = Length



Tools and Kits

## Prosthetic Instruments



|                                                                                     | Cat. No. | Dimensions (mm) | Description                            | Material                              | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|-----------------|----------------------------------------|---------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | MT-IT100 | L=16.5          | Impression coping key for closed tray  | Stainless steel                       |                                                                                          |
|   | MT-LM005 | L=30            | Long motor adapter for 0.05 inch hex.  | Stainless steel                       |     |
|  | MT-SM005 | L=24            | Short motor adapter for 0.05 inch hex. | Stainless steel                       |                                                                                        |
|  | MT-MSD20 | L=140           | Surgical adapter driver                | Stainless steel<br>+ Aluminum         |   |
|  | MT-HTS10 | L=142.5         | Technical hex. screwdriver             | Stainless steel+<br>plastic           |   |
|  | MM-LCT10 | L=83            | Locator core tool                      | Steel / Titanium<br>Nitride coated Ti |   |
|  | MT-RB225 | L=22            | Ball attachment anchor key             | Stainless steel                       |   |
|  | MT-NB225 | L=22            | Ball attachment anchor key             | Stainless steel                       |                                                                                        |
|  | MT-RDE30 | L=32            | Extra-long driver for 0.05 inch hex.   | Stainless steel                       |                                                                                        |

L = Length

|                                                                                     | Cat. No. | Dimensions (mm) | Description                                         | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|-----------------|-----------------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | MT-RDL30 | L=27            | Long driver for 0.05 inch hex.                      | Stainless steel |       |
|    | MT-RDS30 | L=21.3          | Short driver for 0.05 inch hex.                     | Stainless steel |       |
|    | MT-RDM30 | L=11.5          | Extra-short driver for 0.05 inch hex.               | Stainless steel |                                                                                          |
|    | MT-RE172 | L=28.5          | Int. connection abutment extractor                  | Titanium        |       |
|   | MT-RE160 | L=28.5          | Int. connection abutment extractor, narrow platform | Titanium        |       |
|  | MT-MURL2 | L=21            | Straight Multi-Unit long ratchet key                | Stainless steel |   |
|  | MT-MURS2 | L=16            | Straight Multi-Unit short ratchet key               | Stainless steel |                                                                                        |
|  | MT-MUML2 | L=31            | Straight Multi-Unit long motor key                  | Stainless steel |   |
|  | MT-MUMS2 | L=24            | Straight Multi-Unit short motor key                 | Stainless steel |                                                                                        |
|  | MT-RA480 |                 | Ratchet to hex. adapter                             | Stainless steel |                                                                                        |
|  | MT-CLM21 | L=27.5          | CONNECT long motor insertion tool                   | Stainless steel |                                                                                        |
|  | MT-CSM21 | L=23.7          | CONNECT short motor insertion tool                  | Stainless steel |                                                                                        |

L = Length

## Prosthetic Instruments

|                                                                                     | Cat. No. | Dimensions (mm) | Description                                | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|-----------------|--------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | MT-CLR21 | L=31.6          | CONNECT long ratchet insertion tool        | Stainless steel |                                                                                          |
|    | MT-CSR21 | L=21.6          | CONNECT short ratchet insertion tool       | Stainless steel |                                                                                          |
|    | MT-RSD19 | L=19            | LOCKiT short torque ratchet driver         | Stainless steel |       |
|    | MT-RSD24 | L=24            | LOCKiT long torque ratchet driver          | Stainless steel |       |
|   | MT-ELM10 | L=30            | Long motor key for angled screw channel    | Stainless steel |                                                                                         |
|  | MT-ESM10 | L=24            | Short motor key for angled screw channel   | Stainless steel |                                                                                        |
|  | MT-ELR10 | L=22.5          | Long ratchet key for angled screw channel  | Stainless steel |   |
|  | MT-ESR10 | L=18.7          | Short ratchet key for angled screw channel | Stainless steel |                                                                                        |
|  | MT-LOC10 | Ø7.8<br>L=86    | LOCKiT core tool                           | Stainless steel |   |
|  | MT-TSD10 | L=50            | Screwdriver for titanium ball cap          | Stainless steel |                                                                                        |
|  | ET-IT002 |                 | OT-Equator handling tool                   |                 |   |
|  | MT-RI040 | L=95            | Torque ratchet for prosthetic screws       |                 |   |



Tools and Kits

## Conical Connection Surgical Kit

MK-T051  
Conical connection surgical kit



|                                                                                     | Cat. No. | Dimensions (mm)  | Description                                                    | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|------------------|----------------------------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | MT-SMD10 | L=30             | Spade marking drill                                            | Stainless steel |   |
|  | MT-PDM24 | Ø2.4<br>L=32     | Position drill mill                                            | Stainless steel |   |
|  | MT-PD440 | Ø4<br>L=32.7     | Position drill                                                 | Stainless steel |   |
|  | CT-P2408 | Ø2.4/2.0<br>L=28 | Step pilot drill with built-in stopper for 8mm length implants | Stainless steel |   |

L = Length

## Conical Connection Surgical Kit

|                                                                                     | Cat. No. | Dimensions (mm)    | Description                                                       | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|--------------------|-------------------------------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | CT-P2410 | Ø2.4/2.0<br>L=28   | Step pilot drill with built-in stopper for 10mm length implants   | Stainless steel |       |
|    | CT-P2411 | Ø2.4/2.0<br>L=33   | Step pilot drill with built-in stopper for 11.5mm length implants | Stainless steel |       |
|    | CT-P2413 | Ø2.4/2.0<br>L=33   | Step pilot drill with built-in stopper for 13mm length implants   | Stainless steel |       |
|    | CT-P2416 | Ø2.4/2.0<br>L=37.5 | Step pilot drill for 16mm length implants                         | Stainless steel |       |
|    | CT-BTC24 | Ø2.4<br>L=28.5     | Body try-in                                                       | Titanium        |       |
|  | CT-BTC30 | Ø3<br>L=28.5       | Body try-in                                                       | Titanium        |   |
|  | CT-BTC35 | Ø3.5<br>L=28.5     | Body try-in                                                       | Titanium        |   |
|  | CT-BTC40 | Ø4<br>L=28.5       | Body try-in                                                       | Titanium        |   |
|  | CT-TDC30 | Ø3/2.4<br>L=37.5   | Step drill, external irrigation                                   | Stainless steel |   |
|  | CT-TDC35 | Ø3.5/3<br>L=37.5   | Step drill, external irrigation                                   | Stainless steel |   |
|  | CT-TDC40 | Ø4/3.5<br>L=37.5   | Step drill, external irrigation                                   | Stainless steel |   |
|  | MT-CSN33 | Ø3.3<br>L=31       | Countersink, narrow platform                                      | Stainless steel |   |

L = Length



|                                                                                     | Cat. No. | Dimensions (mm)   | Description                                                               | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|-------------------|---------------------------------------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | MT-GDN33 | Ø3.75/4.2<br>L=26 | Countersink, standard platform                                            | Stainless steel |       |
|    | MT-GDN50 | Ø5/6<br>L=26.8    | Countersink, wide platform                                                | Stainless steel |       |
|    | CT-NLM30 | L=29              | V3 long motor insertion tool,<br>conical connection, narrow platform      | Stainless steel |       |
|    | CT-NLR30 | L=32.4            | V3 long ratchet insertion tool,<br>conical connection, narrow platform    | Stainless steel |       |
|   | CT-NSM30 | L=25              | V3 short motor insertion tool,<br>conical connection, narrow platform     | Stainless steel |     |
|  | CT-NSR30 | L=25              | V3 short ratchet insertion tool,<br>conical connection, narrow platform   | Stainless steel |   |
|  | CT-SLM30 | L=29              | V3 long motor insertion tool,<br>conical connection, standard platform    | Stainless steel |   |
|  | CT-SLR30 | L=32.4            | V3 long ratchet insertion tool,<br>conical connection, standard platform  | Stainless steel |   |
|  | CT-SSM30 | L=25              | V3 short motor insertion tool,<br>conical connection, standard platform   | Stainless steel |   |
|  | CT-SSR30 | L=22.3            | V3 short ratchet insertion tool,<br>conical connection, standard platform | Stainless steel |   |
|  | CT-NLI10 | L=32              | C1 long motor insertion tool,<br>conical connection, narrow platform      | Stainless steel |   |
|  | CT-NLR10 | L=32              | C1 long ratchet insertion tool,<br>conical connection, narrow platform    | Stainless steel |   |

L = Length

## Conical Connection Surgical Kit

|                                                                                     | Cat. No. | Dimensions (mm) | Description                                                               | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|-----------------|---------------------------------------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | CT-NSI10 | L=24.5          | C1 short motor insertion tool,<br>conical connection, narrow platform     | Stainless steel |       |
|    | CT-NSR10 | L=24.5          | C1 short ratchet insertion tool,<br>conical connection, narrow platform   | Stainless steel |       |
|    | CT-SLI10 | L=33            | C1 long motor insertion tool,<br>conical connection, standard platform    | Stainless Steel |       |
|    | CT-SLR10 | L=33            | C1 long ratchet insertion tool,<br>conical connection, standard platform  | Stainless steel |       |
|    | CT-SSI10 | L=25.5          | C1 short motor insertion tool,<br>conical connection, standard platform   | Stainless Steel |       |
|  | CT-SSR10 | L=23            | C1 short ratchet insertion tool,<br>conical connection, standard platform | Stainless steel |   |
|  | CT-WLI10 | L=33            | C1 long motor insertion tool,<br>conical connection, wide platform        | Stainless steel |   |
|  | CT-WLR10 | L=33            | C1 long ratchet insertion tool,<br>conical connection, wide platform      | Stainless steel |   |
|  | CT-WSI10 | L=25.5          | C1 short motor insertion tool,<br>conical connection, wide platform       | Stainless steel |   |
|  | CT-WSR10 | L=23            | C1 short ratchet insertion tool,<br>conical connection, wide platform     | Stainless steel |   |
|  | MT-RE160 | L=28.5          | Int. connection abutment extractor,<br>narrow platform                    | Titanium        |   |
|  | MT-RE172 | L=28.5          | Int. connection abutment extractor                                        | Titanium        |   |

L = Length


|                                                                                    | Cat. No. | Dimensions (mm) | Description                     | Material        | Sold                                                                                                                                                                      |
|------------------------------------------------------------------------------------|----------|-----------------|---------------------------------|-----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|   | MT-DE001 | L=24            | Drill extender                  | Stainless steel |     |
|   | MT-RDL30 | L=27            | Long driver for 0.05 inch hex.  | Stainless steel |     |
|   | MT-RDS30 | L=21.3          | Short driver for 0.05 inch hex. | Stainless steel |     |
|   | MT-BTI20 | Ø1.40<br>L=100  | Implant site depth probe        | Stainless steel |     |
|  | MT-RT070 | L=95            | Surgical torque ratchet         |                 |   |

L = Length



6

## Marking Options - Body Try-ins





Tools and Kits

## Drill Replacement kit Conical Connection



**MK-CD01**  
Drill replacement kit

|                                                                                     | Cat. No. | Dimensions (mm)     | Description                                                       | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|---------------------|-------------------------------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|   | CT-P2408 | Ø2.40/2.0<br>L=28   | Step pilot drill with built-in stopper for 8mm length implants    | Stainless steel |       |
|  | CT-P2410 | Ø2.40/2.0<br>L=28   | Step pilot drill with built-in stopper for 10mm length implants   | Stainless steel |   |
|  | CT-P2411 | Ø2.40/2.0<br>L=33   | Step pilot drill with built-in stopper for 11.5mm length implants | Stainless steel |   |
|  | CT-P2413 | Ø2.40/2.0<br>L=33   | Step pilot drill with built-in stopper for 13mm length implants   | Stainless steel |   |
|  | CT-P2416 | Ø2.40/2.0<br>L=37.5 | Step pilot drill for 16mm length implants                         | Stainless steel |   |
|  | CT-TDC30 | Ø3/2.40<br>L=37.5   | Step drill, external irrigation                                   | Stainless steel |   |
|  | CT-TDC35 | Ø3.50/3<br>L=37.5   | Step drill, external irrigation                                   | Stainless steel |   |
|  | CT-TDC40 | Ø4/3.50<br>L=37.5   | Step drill, external irrigation                                   | Stainless steel |   |

L = Length



Tools and Kits

## Bone Profiler Kit Conical Connection

MK-0060  
Bone Profiler Kit



|                                                                                     | Cat. No. | Dimensions (mm) | Description                                                    | Material        | Sold                                                                                  |
|-------------------------------------------------------------------------------------|----------|-----------------|----------------------------------------------------------------|-----------------|---------------------------------------------------------------------------------------|
|  | CT-BP320 | Ø3.20<br>L=28.4 | Bone Profiler, conical connection, narrow platform             | Stainless steel |  |
|  | CT-BP355 | Ø3.55<br>L=28.4 | Bone Profiler, conical connection, standard platform           | Stainless steel |  |
|  | CT-BP460 | Ø4.60<br>L=28.4 | Bone Profiler, conical connection, wide platform               | Stainless steel |  |
|  | CT-BGP32 | Ø3.20<br>L=11.7 | Bone Profiler guide pin, conical connection, narrow platform   | Stainless steel |  |
|  | CT-BGP35 | Ø3.55<br>L=11.7 | Bone Profiler guide pin, conical connection, standard platform | Stainless steel |  |
|  | CT-BGP46 | Ø4.60<br>L=11.7 | Bone Profiler guide pin, conical connection, wide platform     | Stainless steel |  |

L = Length



Tools and Kits

Drill Stoppers Kit,  
Standard Platform

Following the use of MIS length-based pilot drills, the C1 Standard Platform Stoppers Kit ensures drilling to the exact desired depth. The kit includes the most commonly used implant lengths: 8, 10, 11.5 and 13mm. Drill Stoppers Kits are suitable for both V3 and C1 conical connection implant systems.



| Cat. No. | Dimensions (mm) | Description | Material | Sold |
|----------|-----------------|-------------|----------|------|
|----------|-----------------|-------------|----------|------|



MK-BC101

Drill stoppers kit, conical connection,  
standard platform

L = Length




Tools and Kits

## Drill Stoppers Kit, Depth-Based

MK-CDS(06-13)   
Drill stoppers kit



|                                                                                     | Cat. No. | Dimensions (mm) | Description                                                       | Material | Sold                                                                                             |
|-------------------------------------------------------------------------------------|----------|-----------------|-------------------------------------------------------------------|----------|--------------------------------------------------------------------------------------------------|
|  | MK-CDS08 | L=8             | Drill stoppers kit for 8mm length implants, conical connection    |          |  INSIDE KIT |
|  | MK-CDS10 | L=10            | Drill stoppers kit for 10mm length implants, conical connection   |          |  INSIDE KIT |
|  | MK-CDS11 | L=11.5          | Drill stoppers kit for 11.5mm length implants, conical connection |          |  INSIDE KIT |
|  | MK-CDS13 | L=13            | Drill stoppers kit for 13mm length implants, conical connection   |          |  INSIDE KIT |

L = Length



Tools and Kits

## Bone Compression Kit

MK-0020  
Bone Compression Kit



|                                                                                     | Cat. No. | Dimensions (mm)    | Description              | Material        | Sold                                                                                  |
|-------------------------------------------------------------------------------------|----------|--------------------|--------------------------|-----------------|---------------------------------------------------------------------------------------|
|  | MO-CO160 | Ø1.80/2.90<br>L=25 | Convex compression screw | Stainless steel |  |
|  | MO-CO200 | Ø2.30/3.30<br>L=25 | Convex compression screw | Stainless steel |  |
|  | MO-CO240 | Ø2.70/3.70<br>L=25 | Convex compression screw | Stainless steel |  |
|  | MO-CO280 | Ø3.10/4.50<br>L=25 | Convex compression screw | Stainless steel |  |
|  | MO-SO300 | Ø3<br>L=25         | Concave sinus screw      | Stainless steel |  |
|  | MO-SO350 | Ø3.50<br>L=25      | Concave sinus screw      | Stainless steel |  |
|  | MO-SO400 | Ø4<br>L=25         | Concave sinus screw      | Stainless steel |  |


L = Length




|                                                                                     | Cat. No. | Dimensions (mm) | Description                                 | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|-----------------|---------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | MT-SMD10 | L=30            | Spade marking drill                         | Stainless steel |       |
|    | MT-HSI10 | L=24.40         | Short motor insertion tool, internal hex.   | Stainless steel |       |
|    | MT-HLI10 | L=28.20         | Long motor insertion tool, internal hex.    | Stainless steel |       |
|    | MT-SRH20 | L=22.30         | Short ratchet insertion tool, internal hex. | Stainless steel |       |
|   | MT-LRH20 | L=32.30         | Long ratchet insertion tool, internal hex.  | Stainless steel |       |
|  | MT-RI030 | L=84            | Monoblock ratchet                           |                 |   |

L = Length

## 2 Surgical Procedures



Bone Compression



Sinus Lifting



Tools and Kits

Sinus Lift  
Osteotome Set

The MIS Sinus Lift Osteotome Set is used for closed sinus augmentation procedures. The set includes four straight and four bayonet concave tools. All tools are laser marked for depths of: 6, 11 and 15mm, and are equipped with stoppers for convenience and safety.



| Cat. No. | Dimensions (mm) | Description | Material | Sold |
|----------|-----------------|-------------|----------|------|
|----------|-----------------|-------------|----------|------|



|         |  |                          |  |  |
|---------|--|--------------------------|--|--|
| MK-0023 |  | Sinus lift osteotome set |  |  |
|---------|--|--------------------------|--|--|



Tools and Kits

## Convex Osteotome Set

The MIS Convex Osteotome Set is used for closed sinus augmentation procedures. The set includes four straight and four bayonet convex tools. All tools are laser marked for depths of 6, 11 and 15mm, and are equipped with stoppers for convenience and safety.



| Cat. No. | Dimensions (mm) | Description | Material | Sold |
|----------|-----------------|-------------|----------|------|
|----------|-----------------|-------------|----------|------|



MK-0024

Convex osteotome set


L = Length



Tools and Kits

Prosthetic Tools Kit

MK-0039  
Prosthetic Tools Kit



|  | Cat. No. | Dimensions (mm) | Description                            | Material        | Sold |
|--|----------|-----------------|----------------------------------------|-----------------|------|
|  | MT-RDS30 | L=21.3          | Short driver for 0.05 inch hex.        | Stainless steel |      |
|  | MT-RDL30 | L=27            | Long driver for 0.05 inch hex.         | Stainless steel |      |
|  | MT-SM005 | L=24            | Short motor adapter for 0.05 inch hex. | Stainless steel |      |
|  | MT-LM005 | L=30            | Long motor adapter for 0.05 inch hex.  | Stainless steel |      |
|  | MT-MURL2 | L=21            | Straight Multi-Unit long ratchet key   | Stainless steel |      |


L = Length

|                                                                                     | Cat. No. | Dimensions (mm) | Description                                                                     | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|-----------------|---------------------------------------------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | MT-MUML2 | L=31            | Straight Multi-Unit long motor key                                              | Stainless steel |       |
|    | MT-RE172 | L=28.5          | Int. connection abutment extractor                                              | Titanium        |       |
|    | MT-RE160 | L=28.5          | Int. connection abutment extractor, narrow platform                             | Titanium        |       |
|    | MT-RSD19 | L=19            | LOCKiT short torque ratchet driver                                              | Stainless steel |       |
|   | MT-RSD24 | L=24            | LOCKiT long torque ratchet driver                                               | Stainless steel |       |
|  | MT-RB225 | L=22            | Ball attachment anchor key                                                      | Stainless steel |   |
|  | MT-UN002 | L=24            | Implant analog for universal hand driver, internal hex., narrow platform        | Stainless steel |   |
|  | MT-US002 | L=24            | Implant analog for universal hand driver, internal hex., standard platform      | Stainless steel |   |
|  | MT-UW002 | L=24            | Implant analog for universal hand driver, internal hex., wide platform          | Stainless steel |   |
|  | CT-UN001 | L=24            | Implant analog for universal hand driver, conical connection, narrow platform   | Stainless steel |                                                                                        |
|  | CT-US001 | L=24            | Implant analog for universal hand driver, conical connection, standard platform | Stainless steel |                                                                                        |
|  | CT-UW001 | L=24            | Implant analog for universal hand driver, conical connection, wide platform     | Stainless steel |                                                                                        |

L = Length

Prosthetic Tools Kit




|                                                                                    | Cat. No. | Dimensions (mm) | Description                               | Material                   | Sold                                                                                                                                                                      |
|------------------------------------------------------------------------------------|----------|-----------------|-------------------------------------------|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|   | ET-IT002 |                 | OT-Equator handling tool                  |                            |     |
|   | MT-LOC10 | Ø7.8<br>L=86    | LOCKIT core tool                          | Stainless steel            |     |
|   | MT-MSD20 | L=140           | Surgical adapter driver                   | Stainless steel + Aluminum |     |
|   | MT-ELR10 | L=22.5          | Long ratchet key for angled screw channel | Stainless steel            |     |
|  | MT-RI040 | L=95            | Torque ratchet for prosthetic screws      |                            |   |

L = Length



1 Prosthetic Solution

2 Implant Connections







Tools and Kits

## Abutment Holding System Kit

The abutment holding system allows convenient grip of prosthetic components for adjustments and preparation of both ceramic and titanium abutments.



| Cat. No. | Dimensions (mm) | Description | Material | Sold |
|----------|-----------------|-------------|----------|------|
|----------|-----------------|-------------|----------|------|



MK-C101

Abutment holding system kit, conical connection



Tools and Kits

Crown Set

The MIS Crown Set is a delivery system comprised of an auto-mix, easy-to-use dispensing syringe and disposable dispensing tips. The resin-based cement combines high retention properties along with the ability to remove a restoration if needed. The cement sets within 120 -150 seconds and allows easy removal of excess material. Note: Abutments should be sand-blasted prior to cementation to achieve maximum retention.



| Cat. No. | Dimensions (mm) | Description | Material | Sold |
|----------|-----------------|-------------|----------|------|
|----------|-----------------|-------------|----------|------|



MK-0022

Crown set





Tools and Kits

## Zircon Diamond Burr Set

MK-DB01   
Zircon Diamond Burr Set



|                                                                                     | Cat. No. | Dimensions (mm) | Description          | Material | Sold                                                                                  |
|-------------------------------------------------------------------------------------|----------|-----------------|----------------------|----------|---------------------------------------------------------------------------------------|
|  | MT-A4022 |                 | Diamond burr, 2.20mm |          |  |
|  | MT-B4018 |                 | Diamond burr, 1.80mm |          |  |
|  | MT-F2121 |                 | Diamond burr, 2.10mm |          |  |
|  | MT-K2023 |                 | Diamond burr, 2.30mm |          |  |
|  | MT-C3125 |                 | Diamond burr, 2.50mm |          |  |


L = Length




Tools and Kits

Control Units



|                                                                                     | Cat. No.  | Dimensions (mm) | Description                                       | Material | Sold |
|-------------------------------------------------------------------------------------|-----------|-----------------|---------------------------------------------------|----------|------|
|   | MWH-SE236 |                 | MIS basic physio and contra angle set (230 V)     |          |      |
|  | MWH-SE237 |                 | MIS basic physio and contra angle set (115 V)     |          |      |
|  | MWH-SE225 |                 | MIS basic physio and contra angle LED set (230 V) |          |      |
|  | MWH-SE226 |                 | MIS basic physio and contra angle LED set (115 V) |          |      |
|  | MWH-2STUB |                 | Irrigation system 1/2 way                         |          |      |
|  | MWH-MMOTR |                 | Micro motor and cord kit                          |          |      |
|  | MWH-M7559 |                 | Surgical contra-angle 20:1                        |          |      |
|  | MWH-CAL02 |                 | Surgical contra-angle LED 20:1                    |          |      |





## PROSTHETIC OPTIONS

MIS offers a comprehensive range of prosthetic systems, which are suitable for any indication and procedure type, including analog and digital workflows. Combined with a wide variety of supporting tools and digital solutions, our products are designed to provide the biological advantages that are critical for successful treatment.

MIS superstructures are designed to meet all prosthetic needs and fit all accepted working methods. Abutments are available in a wide range of heights, gingival collar heights, angles and margin types, as well as a line of concave emergence profile abutments which are designed to facilitate better soft tissue management and growth. MIS offers solutions for both screw retained and cement retained treatments, for temporary and permanent restorations.





Impression Coping

## Analogs and Impression Copings




|  | Cat. No.                                     | Implant              | Dimensions (mm)         | Description                                          | Material                                                 | Platform |
|--|----------------------------------------------|----------------------|-------------------------|------------------------------------------------------|----------------------------------------------------------|----------|
|  | CN-IO330<br>VN-IO401<br>CS-IO375<br>CW-IO500 | NP<br>NP<br>SP<br>WP | C1<br>V3<br>V3/C1<br>C1 | Ø3.3, L=11<br>Ø3.3, L=11<br>Ø3.75, L=11<br>Ø5, L=11  | Impression coping for open tray, conical connection      | Titanium |
|  | CN-IO405<br>VN-IO405<br>CS-IO485<br>CW-IO555 | NP<br>NP<br>SP<br>WP | C1<br>V3<br>V3/C1<br>C1 | Ø3.3, L=15<br>Ø3.3, L=15<br>Ø3.75, L=15<br>Ø5, L=15  | Impression coping for open tray, conical connection      | Titanium |
|  | CN-IO331<br>VN-IO331<br>CS-IO391<br>CW-IO501 | NP<br>NP<br>SP<br>WP | C1<br>V3<br>V3/C1<br>C1 | Ø2.8, L=11<br>Ø2.8, L=11<br>Ø3.2, L=11<br>Ø4.1, L=11 | Thin impression coping for open tray, conical connection | Titanium |
|  | CN-IO335<br>VN-IO335<br>CS-IO395<br>CW-IO505 | NP<br>NP<br>SP<br>WP | C1<br>V3<br>V3/C1<br>C1 | Ø2.8, L=15<br>Ø2.8, L=15<br>Ø3.2, L=15<br>Ø4.1, L=15 | Thin impression coping for open tray, conical connection | Titanium |
|  | CN-IC401<br>VN-IC401<br>CS-IC481<br>CW-IC551 | NP<br>NP<br>SP<br>WP | C1<br>V3<br>V3/C1<br>C1 | Ø3.9, L=11<br>Ø3.9, L=11<br>Ø4.5, L=11<br>Ø4.5, L=11 | Impression coping for closed tray, conical connection    | Titanium |

H = Gingival height

C = Crown height

L = Length

|                                                                                     | Cat. No. | Implant |       | Dimensions (mm) | Description                                                  | Material | Platform                                                                                                                                                                                                                                                          |
|-------------------------------------------------------------------------------------|----------|---------|-------|-----------------|--------------------------------------------------------------|----------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | CN-IC405 | NP      | C1    | Ø3.9, L=15      | Impression coping for closed tray, conical connection        | Titanium |          |
|                                                                                     | VN-IC405 | NP      | V3    | Ø3.9, L=15      |                                                              |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CS-IC485 | SP      | V3/C1 | Ø4.5, L=15      |                                                              |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CW-IC555 | WP      | C1    | Ø4.5, L=15      |                                                              |          |                                                                                                                                                                                                                                                                   |
|    | CN-IT300 | NP      | C1    | Ø4, H=2, C=6    | Direct impression coping for closed tray, conical connection | Titanium |          |
|                                                                                     | VN-IC400 | NP      | V3    | Ø4, H=2, C=6    |                                                              |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CS-IT300 | SP      | V3/C1 | Ø4.8, H=2, C=6  |                                                              |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CW-IT300 | WP      | C1    | Ø5.5, H=2, C=6  |                                                              |          |                                                                                                                                                                                                                                                                   |
|    | CN-PF330 | NP      | C1    | Ø2.8, L=16      | Direct press fit for closed tray, conical connection         | Titanium |          |
|                                                                                     | VN-PF340 | NP      | V3    | Ø2.8, L=15.6    |                                                              |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CS-PF375 | SP      | V3/C1 | Ø3.2, L=16.3    |                                                              |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CW-PF500 | WP      | C1    | Ø5.5, L=16.3    |                                                              |          |                                                                                                                                                                                                                                                                   |
|    | CN-RSM10 | NP      | C1    | Ø3.3, L=12.7    | Implant analog, conical connection                           | Titanium |          |
|                                                                                     | VN-RSM10 | NP      | V3    | Ø3.3, L=12.7    |                                                              |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CS-RSM10 | SP      | V3/C1 | Ø3.75, L=12     |                                                              |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CW-RSM10 | WP      | C1    | Ø5, L=12        |                                                              |          |                                                                                                                                                                                                                                                                   |
|   | MN-S0160 | NP      | V3/C1 | L=7.6           | Direct prosthetic screw, 30 Ncm                              | Titanium |       |
|                                                                                     | MD-S0220 | SP      | V3/C1 |                 |                                                              |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | MD-S0220 | WP      | V3/C1 |                 |                                                              |          |                                                                                                                                                                                                                                                                   |
|  | MN-G0160 | NP      | V3/C1 | L=23            | Direct impression coping guide pin for open tray             | Titanium |    |
|                                                                                     | MD-G0213 | SP      | V3/C1 |                 |                                                              |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | MD-G0213 | WP      | V3/C1 |                 |                                                              |          |                                                                                                                                                                                                                                                                   |


H = Gingival height

C = Crown height

L = Length



## 3 Platforms





Temporary Restoration

## Temporary Cylinders





|  | Cat. No. | Implant | Dimensions (mm) | Description     | Material | Platform |
|--|----------|---------|-----------------|-----------------|----------|----------|
|  | CN-TC013 | NP      | V3/C1           | Ø4, H=1, C=10   | Titanium | NP SP WP |
|  | CS-TC013 | SP      | V3/C1           | Ø4.8, H=1, C=10 |          |          |
|  | VS-TC013 | SP      | V3/C1           | Ø5.8, H=1, C=10 |          |          |
|  | CW-TC013 | WP      | C1              | Ø5.5, H=1, C=10 |          |          |
|  | CN-TCI13 | NP      | C1              | Ø4, H=1, C=10   | Titanium | NP SP WP |
|  | VN-TCI13 | NP      | V3              | Ø4, H=1, C=10   |          |          |
|  | CS-TCI13 | SP      | V3/C1           | Ø4.8, H=1, C=10 |          |          |
|  | VS-TCI13 | SP      | V3/C1           | Ø5.8, H=1, C=10 |          |          |
|  | CW-TCI13 | WP      | C1              | Ø5.5, H=1, C=10 |          |          |
|  | CN-TC023 | NP      | V3/C1           | Ø4, H=2, C=10   | Titanium | NP SP WP |
|  | CS-TC023 | SP      | V3/C1           | Ø4.8, H=2, C=10 |          |          |
|  | VS-TC023 | SP      | V3/C1           | Ø5.8, H=2, C=10 |          |          |
|  | CW-TC023 | WP      | C1              | Ø5.5, H=2, C=10 |          |          |

H = Gingival height

C = Crown height

L = Length



|                                                                                     | Cat. No. | Implant |       | Dimensions (mm)   | Description                                                  | Material | Platform                                                                                                                                                                                                                                                          |
|-------------------------------------------------------------------------------------|----------|---------|-------|-------------------|--------------------------------------------------------------|----------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | CN-TCI23 | NP      | C1    | Ø4, H=2, C=10     | Direct temporary cylinder, anti-rotation, conical connection | Titanium |          |
|                                                                                     | VN-TCI23 | NP      | V3    | Ø4, H=2, C=10     |                                                              |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CS-TCI23 | SP      | V3/C1 | Ø4.8, H=2, C=10   |                                                              |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | VS-TCI23 | SP      | V3/C1 | Ø5.8, H=2, C=10   |                                                              |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CW-TCI23 | WP      | C1    | Ø5.5, H=2, C=10   |                                                              |          |                                                                                                                                                                                                                                                                   |
|    | CN-TC033 | NP      | V3/C1 | Ø4, H=3, C=10     | Direct temporary cylinder, free rotation, conical connection | Titanium |          |
|                                                                                     | CS-TC033 | SP      | V3/C1 | Ø4.8, H=3, C=10   |                                                              |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | VS-TC033 | SP      | V3/C1 | Ø5.8, H=3, C=10   |                                                              |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CW-TC033 | WP      | C1    | Ø5.5, H=3, C=10   |                                                              |          |                                                                                                                                                                                                                                                                   |
|    | CN-TCI33 | NP      | C1    | Ø4, H=3, C=10     | Direct temporary cylinder, anti-rotation, conical connection | Titanium |          |
|                                                                                     | VN-TCI33 | NP      | V3    | Ø4, H=3, C=10     |                                                              |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CS-TCI33 | SP      | V3/C1 | Ø4.8, H=3, C=10   |                                                              |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | VS-TCI33 | SP      | V3/C1 | Ø5.8, H=3, C=10   |                                                              |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CW-TCI33 | WP      | C1    | Ø5.5, H=3, C=10   |                                                              |          |                                                                                                                                                                                                                                                                   |
|    | CN-TPA10 | NP      | C1    | Ø4, H=1.5, C=10   | Temporary plastic abutment, conical connection               | PEEK     |          |
|                                                                                     | VN-TPA10 | NP      | V3    | Ø4, H=1.5, C=10   |                                                              |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CS-TPA10 | SP      | V3/C1 | Ø4.8, H=1.5, C=10 |                                                              |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CW-TPA10 | WP      | C1    | Ø5.5, H=1.5, C=10 |                                                              |          |                                                                                                                                                                                                                                                                   |
|   | CN-TPH50 | NP      | C1    | Ø3.5              | Direct temporary plastic cylinder, conical connection        | PEEK     |                                                                                           |
|                                                                                     | VN-TPH50 | NP      | V3    |                   |                                                              |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CS-TPH50 | SP      | V3/C1 |                   |                                                              |          |                                                                                                                                                                                                                                                                   |
|  | MN-S0160 | NP      | V3/C1 | L=7.6             | Direct prosthetic screw, 30Ncm                               | Titanium |    |
|                                                                                     | MD-S0220 | SP      |       |                   |                                                              |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | MD-S0220 | WP      |       |                   |                                                              |          |                                                                                                                                                                                                                                                                   |

H = Gingival height

C = Crown height

L = Length

## 2 Abutment Options



Titanium



Polyether ether ketone



Screw Retained Restoration

## CONNECT System

Please note, CONNECT products might not be registered or available in every country/region.










|                                                                                     | Cat. No.                         | Platform       | Implant              | Dimensions (mm) | Description                          | Material | Platform |
|-------------------------------------------------------------------------------------|----------------------------------|----------------|----------------------|-----------------|--------------------------------------|----------|----------|
|   | CS-C1540<br>CW-C1540             | SP<br>WP       | V3/C1<br>C1          | Ø4, H=1.5       | CONNECT abutment, conical connection | Titanium | SP WP    |
|  | CN-C2040<br>CS-C2040<br>CW-C2040 | NP<br>SP<br>WP | V3/C1<br>V3/C1<br>C1 | Ø4, H=2         | CONNECT abutment, conical connection | Titanium | NP SP WP |
|  | CN-C3040<br>CS-C3040<br>CW-C3040 | NP<br>SP<br>WP | V3/C1<br>V3/C1<br>C1 | Ø4, H=3         | CONNECT abutment, conical connection | Titanium | NP SP WP |
|  | CS-C4040<br>CW-C4040             | SP<br>WP       | V3/C1<br>C1          | Ø4, H=4         | CONNECT abutment, conical connection | Titanium | SP WP    |
|  | CS-C1557<br>CW-C1557             | SP<br>WP       | V3/C1<br>C1          | Ø5.7, H=1.5     | CONNECT abutment, conical connection | Titanium | SP WP    |
|  | CS-C2057<br>CW-C2057             | SP<br>WP       | V3/C1<br>C1          | Ø5.7, H=2       | CONNECT abutment, conical connection | Titanium | SP WP    |
|  | CS-C3057<br>CW-C3057             | SP<br>WP       | V3/C1<br>C1          | Ø5.7, H=3       | CONNECT abutment, conical connection | Titanium | SP WP    |
|  | CS-C4057<br>CW-C4057             | SP<br>WP       | V3/C1<br>C1          | Ø5.7, H=4       | CONNECT abutment, conical connection | Titanium | SP WP    |

H = Gingival height

C = Crown height

L = Length


|                                                                                     | Cat. No.                         | Platform | Implant | Dimensions (mm)                         | Description                                               | Material | Platform |
|-------------------------------------------------------------------------------------|----------------------------------|----------|---------|-----------------------------------------|-----------------------------------------------------------|----------|----------|
|    | MM-H0540<br>MM-H1540<br>MM-H3040 |          |         | Ø4, H=0.5<br>Ø4, H=1.5<br>Ø4, H=3       | CONNECT healing cap                                       | Titanium |          |
|    | MM-H0557<br>MM-H1557<br>MM-H3557 |          |         | Ø5.7, H=0.5<br>Ø5.7, H=1.5<br>Ø5.7, H=3 | CONNECT healing cap                                       | Titanium |          |
|    | MM-IO040                         |          |         | Ø4<br>C=7                               | CONNECT short impression coping, open tray, free rotation | Titanium |          |
|    | MM-IOI40                         |          |         | Ø4<br>C=7                               | CONNECT short impression coping, open tray, anti-rotation | Titanium |          |
|    | MM-IO011                         |          |         | Ø4<br>C=11                              | CONNECT long impression coping, open tray, free rotation  | Titanium |          |
|   | MM-IOI11                         |          |         | Ø4<br>C=11                              | CONNECT long impression coping, open tray, anti-rotation  | Titanium |          |
|  | MM-IC040                         |          |         | Ø4<br>C=8.3                             | CONNECT Impression coping, closed tray, free rotation     | Titanium |          |
|  | MM-ICI40                         |          |         | Ø4<br>C=8.3                             | CONNECT Impression coping, closed tray, anti-rotation     | Titanium |          |
|  | MM-IO057                         |          |         | Ø5.7<br>C=7                             | CONNECT Impression coping, open tray, free rotation       | Titanium |          |
|  | MM-IOI57                         |          |         | Ø5.7<br>C=7                             | CONNECT Impression coping, open tray, anti-rotation       | Titanium |          |
|  | MM-IC057                         |          |         | Ø5.7<br>C=8.3                           | CONNECT Impression coping, closed tray, free rotation     | Titanium |          |
|  | MM-ICI57                         |          |         | Ø5.7<br>C=8.8                           | CONNECT Impression coping, closed tray, anti-rotation     | Titanium |          |
|  | MM-RSM41                         |          |         | Ø4<br>L=11                              | CONNECT analog                                            | Titanium |          |

H = Gingival height

C = Crown height

L = Length

## CONNECT System

|                                                                                     | Cat. No. | Platform | Implant | Dimensions (mm) | Description                                       | Material | Platform |
|-------------------------------------------------------------------------------------|----------|----------|---------|-----------------|---------------------------------------------------|----------|----------|
|    | MM-RSM57 |          |         | Ø5.7<br>L=11    | CONNECT analog                                    | Titanium |          |
|    | MM-TC041 |          |         | Ø4<br>C=10      | CONNECT temporary cylinder,<br>free rotation      | Titanium |          |
|    | MM-TCI41 |          |         | Ø4<br>C=10      | CONNECT temporary cylinder,<br>anti-rotation      | Titanium |          |
|    | MM-TC057 |          |         | Ø5.7<br>C=10    | CONNECT temporary cylinder,<br>free rotation      | Titanium |          |
|   | MM-TCI57 |          |         | Ø5.7<br>C=10    | CONNECT temporary cylinder,<br>anti-rotation      | Titanium |          |
|  | MM-CE046 |          |         | Ø4<br>C=6       | CONNECT final esthetic abutment,<br>free rotation | Titanium |          |
|  | MM-CEI46 |          |         | Ø4<br>C=6       | CONNECT final esthetic abutment,<br>anti-rotation | Titanium |          |
|  | MM-CF048 |          |         | Ø4<br>C=8       | CONNECT final esthetic abutment,<br>free rotation | Titanium |          |
|  | MM-CFI48 |          |         | Ø4<br>C=8       | CONNECT final esthetic abutment,<br>anti-rotation | Titanium |          |
|  | MM-CE056 |          |         | Ø5.7<br>C=6     | CONNECT final esthetic abutment,<br>free rotation | Titanium |          |
|  | MM-CEI56 |          |         | Ø5.7<br>C=6     | CONNECT final esthetic abutment,<br>anti-rotation | Titanium |          |
|  | MM-CF054 |          |         | Ø5.7<br>C=4     | CONNECT final esthetic abutment,<br>free rotation | Titanium |          |

H = Gingival height

C = Crown height


L = Length

|                                                                                     | Cat. No. | Platform | Implant | Dimensions (mm) | Description                                               | Material        | Platform |
|-------------------------------------------------------------------------------------|----------|----------|---------|-----------------|-----------------------------------------------------------|-----------------|----------|
|    | MM-CFI54 |          |         | Ø5.7<br>C=4     | CONNECT final esthetic abutment, anti-rotation            | Titanium        |          |
|    | MM-AN204 |          |         | Ø4<br>C=8       | CONNECT cementable angulated abutment, 20°, anti-rotation | Titanium        |          |
|    | MM-CCR04 |          |         | Ø4<br>C=14      | CONNECT CoCr plastic cylinder, free rotation              | Titanium        |          |
|    | MM-CCR14 |          |         | Ø4<br>C=14      | CONNECT CoCr plastic cylinder, anti-rotation              | Titanium        |          |
|    | MM-CCR05 |          |         | Ø5.7<br>C=15.4  | CONNECT CoCr plastic cylinder, free rotation              | Titanium        |          |
|  | MM-CCR15 |          |         | Ø5.7<br>C=15.4  | CONNECT CoCr plastic cylinder, anti-rotation              | Titanium        |          |
|  | MM-SP104 |          |         | Ø4<br>L=10      | CONNECT scan post, anti-rotation                          | Titanium        |          |
|  | MM-SP105 |          |         | Ø5.7<br>L=10    | CONNECT scan post, anti-rotation                          | Titanium        |          |
|  | MM-MAN40 |          |         | Ø4<br>L=10      | CONNECT digital model analog                              | Titanium        |          |
|  | MM-MAN57 |          |         | Ø5.7<br>L=12    | CONNECT digital model analog                              | Titanium        |          |
|  | MT-CLM21 |          |         | L=27.5          | CONNECT long motor insertion tool                         | Stainless steel |          |
|  | MT-CSM21 |          |         | L=23.7          | CONNECT short motor insertion tool                        | Stainless steel |          |

H = Gingival height

C = Crown height

L = Length



|                                                                                    | Cat. No.  | Platform | Implant | Dimensions (mm) | Description                               | Material        | Platform |
|------------------------------------------------------------------------------------|-----------|----------|---------|-----------------|-------------------------------------------|-----------------|----------|
|   | MT-CLR21  |          |         | L=31.6          | CONNECT long ratchet insertion tool       | Stainless steel |          |
|   | MT-CSR21  |          |         | L=21.6          | CONNECT short ratchet insertion tool      | Stainless steel |          |
|   | *MT-IT100 |          |         | L=16.5          | Impression coping key for closed tray     | Stainless steel |          |
|   | MM-S0160  |          |         | L=5.8           | CONNECT prosthetic screw                  | Stainless steel |          |
|  | MM-SA160  |          |         | L=6.2           | Prosthetic screw for angled screw channel | Stainless steel |          |

H = Gingival height

C = Crown height

L = Length

\*MT-IT100 - For use with impression coping for closed tray




2

Diameters

1


Internal Connection





Screw Retained Restoration

## Multi-Unit System



|  | Cat. No.                         | Platform       | Implant              | Dimensions (mm) | Description                                    | Material | Platform |
|--|----------------------------------|----------------|----------------------|-----------------|------------------------------------------------|----------|----------|
|  | CM-N1480<br>CM-S1480<br>CM-W1480 | NP<br>SP<br>WP | V3/C1<br>V3/C1<br>C1 | Ø4.8, H=1       | Multi-Unit abutment, conical connection, 30Ncm | Titanium | NP SP WP |
|  | CM-N2480<br>CM-S2480             | NP<br>SP       | V3/C1                | Ø4.8, H=2       | Multi-Unit abutment, conical connection, 30Ncm | Titanium | NP SP    |
|  | CM-N3480<br>CM-S3480<br>CM-W3480 | NP<br>SP<br>WP | V3/C1<br>V3/C1<br>C1 | Ø4.8, H=3       | Multi-Unit abutment, conical connection, 30Ncm | Titanium | NP SP WP |
|  | CM-S4480                         | SP             | V3/C1                | Ø4.8, H=4       | Multi-Unit abutment, conical connection, 30Ncm | Titanium | SP       |
|  | CM-N5480<br>CM-S5480             | NP<br>SP       | V3/C1                | Ø4.8, H=5       | Multi-Unit abutment, conical connection, 30Ncm | Titanium | NP SP    |

H = Gingival height

C = Crown height

L = Length

## Multi-Unit System



|                                                                                     | Cat. No.             | Platform | Implant     | Dimensions (mm) | Description                                                         | Material        | Platform |
|-------------------------------------------------------------------------------------|----------------------|----------|-------------|-----------------|---------------------------------------------------------------------|-----------------|----------|
|    | CM-SO171<br>CM-WO171 | SP<br>WP | V3/C1<br>C1 | Ø4.8, H=1       | Solid angulated Multi-Unit abutment, 17°, conical connection, 20Ncm | Titanium        | SP WP    |
|    | CM-SO172<br>CM-WO172 | SP<br>WP | V3/C1<br>C1 | Ø4.8, H=2       | Solid Multi-Unit angulated abutment, 17°, conical connection, 20Ncm | Titanium        | SP WP    |
|    | CM-SO301<br>CM-WO301 | SP<br>WP | V3/C1<br>C1 | Ø4.8, H=1       | Solid Multi-Unit angulated abutment, 30°, conical connection, 20Ncm | Titanium        | SP WP    |
|    | CM-SO302<br>CM-WO302 | SP<br>WP | V3/C1<br>C1 | Ø4.8, H=2       | Solid Multi-Unit angulated abutment, 30°, conical connection, 20Ncm | Titanium        | SP WP    |
|    | MU-BGC48             |          |             | L=3.3           | Gold bar coping abutment for Multi-Unit                             | Gold-e          |          |
|  | MU-RSM48             |          |             | L=12            | Analog for Multi-Unit                                               | Titanium        |          |
|  | MU-CC480             |          |             | Ø4.8<br>L=4.3   | Cementing cap for Multi-Unit                                        | Titanium        |          |
|  | MU-CP048             |          |             | L=10            | Plastic cylinder for Multi-Unit                                     | Burnout plastic |          |
|  | MU-S0220             |          |             | L=3.8           | Prosthetic screw for Multi-Unit, 25Ncm                              | Titanium        |          |
|  | MU-G0220             |          |             | L=3.8           | Prosthetic gold screw for Multi-Unit, 20Ncm                         | Gold-s          |          |
|  | MU-H4480             |          |             | Ø4.8<br>L=4.3   | Healing cap for Multi-Unit                                          | Titanium        |          |
|  | MU-I0480             |          |             | Ø4.8<br>L=10    | Impression coping for Multi-Unit, open tray,                        | Titanium        |          |

H = Gingival height

C = Crown height

L = Length



|                                                                                     | Cat. No. | Platform | Implant | Dimensions (mm) | Description                                   | Material               | Platform |
|-------------------------------------------------------------------------------------|----------|----------|---------|-----------------|-----------------------------------------------|------------------------|----------|
|    | MU-IT480 |          |         | Ø4.8<br>L=8     | Impression coping, closed tray for Multi-Unit | Titanium               |          |
|    | MU-MG480 |          |         | Ø4.8<br>L=14.5  | Direct gold for Multi-Unit                    | Burnout plastic+gold-a |          |
|    | MU-TO480 |          |         | Ø4.8<br>L=10    | Titanium temporary cylinder for Multi-Unit    | Titanium               |          |
|    | MU-SP102 |          |         | L=10            | Scan Post for Multi-Unit                      | Titanium               |          |
|   | MU-MAN10 |          |         | Ø4.8<br>L=10    | Digital model analog for Multi-Unit           | Titanium               |          |
|  | MT-MURL2 |          |         | L=21            | Straight Multi-Unit long ratchet key          | Stainless steel        |          |
|  | MT-MURS2 |          |         | L=16            | Straight Multi-Unit short ratchet key         | Stainless steel        |          |
|  | MT-MUML2 |          |         | L=31            | Straight Multi-Unit long motor key            | Stainless steel        |          |
|  | MT-MUMS2 |          |         | L=24            | Straight Multi-Unit short motor key           | Stainless steel        |          |

H = Gingival height

C = Crown height

L = Length



Screw Retained Restoration

## Gold Plastic Cylinders



|                                                                                     | Cat. No. | Platform | Implant | Dimensions (mm) | Description                                                     | Material                | Platform                                                                                                                                                                                                                                                          |
|-------------------------------------------------------------------------------------|----------|----------|---------|-----------------|-----------------------------------------------------------------|-------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | CN-GP010 | NP       | V3/C1   | Ø4.3, H=1       | Direct gold plastic cylinder, free rotation, conical connection | Burnout plastic+ Gold-a |    |
|                                                                                     | CS-GP010 | SP       | V3/C1   | Ø4.8, H=1       |                                                                 |                         |                                                                                                                                                                                                                                                                   |
|                                                                                     | CW-GP010 | WP       | C1      | Ø5.5, H=1       |                                                                 |                         |                                                                                                                                                                                                                                                                   |
|  | CN-GPC10 | NP       | C1      | Ø4.3, H=1       | Direct gold plastic cylinder, anti-rotation, conical connection | Burnout plastic+ Gold-a |    |
|                                                                                     | VN-GPC10 | NP       | V3      | Ø4.3, H=1       |                                                                 |                         |                                                                                                                                                                                                                                                                   |
|                                                                                     | CS-GPC10 | SP       | V3/C1   | Ø4.8, H=1       |                                                                 |                         |                                                                                                                                                                                                                                                                   |
|                                                                                     | CW-GPC10 | WP       | C1      | Ø5.5, H=1       |                                                                 |                         |                                                                                                                                                                                                                                                                   |
|  | MN-S0160 | NP       | V3/C1   | L=7.6           | Direct prosthetic screw, 30 Ncm                                 | Titanium                |    |
|                                                                                     | MD-S0220 | SP       | V3/C1   |                 |                                                                 |                         |                                                                                                                                                                                                                                                                   |
|                                                                                     | MD-S0220 | WP       | C1      |                 |                                                                 |                         |                                                                                                                                                                                                                                                                   |
|  | MD-G0220 | SP<br>WP | V3/C1   | L=7.7           | Direct gold prosthetic screw, 20 Ncm                            | Gold-s                  |                                                                                         |

H = Gingival height


C = Crown height

L = Length



Digital Restoration

## EZ-Base System



|                                                                                     | Cat. No. | Platform | Implant | Dimensions (mm)   | Description                                     | Material        | Platform                                                                                                                                                                                                                                                          |
|-------------------------------------------------------------------------------------|----------|----------|---------|-------------------|-------------------------------------------------|-----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | CN-EZB15 | NP       | C1      | Ø4.25, H=1.5, C=5 | EZ-Base abutment, conical connection            | Titanium        |    |
|                                                                                     | VN-EZB15 | NP       | V3      | Ø4.25, H=1.5, C=5 |                                                 |                 |                                                                                                                                                                                                                                                                   |
|                                                                                     | CS-EZB15 | SP       | V3/C1   | Ø4.3, H=1.5, C=5  |                                                 |                 |                                                                                                                                                                                                                                                                   |
|                                                                                     | CW-EZB15 | WP       | C1      | Ø5, H=1.5, C=5    |                                                 |                 |                                                                                                                                                                                                                                                                   |
|  | CN-LEZ15 | NP       | C1      | Ø4.6, H=1.5, C=10 | EZ-Base adjustable abutment, conical connection | Titanium        |    |
|                                                                                     | VN-LEZ15 | NP       | V3      | Ø4.6, H=1.5, C=10 |                                                 |                 |                                                                                                                                                                                                                                                                   |
|                                                                                     | CS-LEZ15 | SP       | V3/C1   | Ø4.8, H=1.5, C=10 |                                                 |                 |                                                                                                                                                                                                                                                                   |
|                                                                                     | CW-LEZ15 | WP       | C1      | Ø5, H=1.5, C=10   |                                                 |                 |                                                                                                                                                                                                                                                                   |
|  | MT-ELM10 |          |         | L=30              | Long motor key for angled screw channel         | Stainless steel |                                                                                                                                                                                                                                                                   |
|  | MT-ESM10 |          |         | L=24              | Short motor key for angled screw channel        | Stainless steel |                                                                                                                                                                                                                                                                   |
|  | MT-ELR10 |          |         | L=22.5            | Long ratchet key for angled screw channel       | Stainless steel |                                                                                                                                                                                                                                                                   |
|  | MT-ESR10 |          |         | L=18.7            | Short ratchet key for angled screw channel      | Stainless steel |                                                                                                                                                                                                                                                                   |
|  | MN-E0160 | NP       |         | L=7.1             | EZ-Base prosthetic screw                        | Titanium        |                                                                                         |
|                                                                                     | MD-E0220 | SP       |         | L=7.6             |                                                 |                 |                                                                                                                                                                                                                                                                   |

H = Gingival height


C = Crown height

L = Length



Digital Restoration

## Ti-Bases



|  | Cat. No. | Platform | Implant | Dimensions (mm)  | Description                                      | Material | Platform |
|--|----------|----------|---------|------------------|--------------------------------------------------|----------|----------|
|  | CN-TB001 | NP       | C1      | Ø4.2, H=0.5, C=4 | Ti-Basis, anti-rotation, conical connection      | Titanium |          |
|  | VN-TB001 | NP       | V3      | Ø4.2, H=0.5, C=4 |                                                  |          |          |
|  | CS-TB001 | SP       | V3/C1   | Ø4.3, H=0.5, C=4 |                                                  |          |          |
|  | CW-TB001 | WP       | C1      | Ø5.7, H=0.5, C=4 |                                                  |          |          |
|  | CN-TB015 | NP       | C1      | Ø4.2, H=1.5, C=4 | Ti-Basis, anti-rotation, conical connection      | Titanium |          |
|  | VN-TB015 | NP       | V3      | Ø4.2, H=1.5, C=4 |                                                  |          |          |
|  | CS-TB015 | SP       | V3/C1   | Ø4.3, H=1.5, C=4 |                                                  |          |          |
|  | CW-TB015 | WP       | C1      | Ø5.7, H=1.5, C=4 |                                                  |          |          |
|  | CN-TB030 | NP       | C1      | Ø4.2, H=3, C=4   | Ti-Basis, anti-rotation, conical connection      | Titanium |          |
|  | VN-TB030 | NP       | V3      | Ø4.2, H=3, C=4   |                                                  |          |          |
|  | CS-TB030 | SP       | V3/C1   | Ø4.3, H=3, C=4   |                                                  |          |          |
|  | CW-TB030 | WP       | C1      | Ø5.7, H=3, C=4   |                                                  |          |          |
|  | CN-IN060 | NP       | C1      | Ø4.2, H=0.5, C=6 | Incisor esthetic Ti abutment, conical connection | Titanium |          |
|  | VN-IN060 | NP       | V3      | Ø4.2, H=0.5, C=6 |                                                  |          |          |
|  | CS-IN060 | SP       | V3/C1   | Ø4.3, H=0.5, C=6 |                                                  |          |          |

H = Gingival height

C = Crown height

L = Length

|                                                                                     | Cat. No. | Platform | Implant | Dimensions (mm)   | Description                                      | Material | Platform                                                                                                                                                                                                                                                          |
|-------------------------------------------------------------------------------------|----------|----------|---------|-------------------|--------------------------------------------------|----------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | CN-IN160 | NP       | C1      | Ø4.2, H=1.5, C=6  | Incisor esthetic Ti abutment, conical connection | Titanium |                                                                                             |
|                                                                                     | VN-IN160 | NP       | V3      | Ø4.2, H=1.5, C=6  |                                                  |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CS-IN160 | SP       | V3/C1   | Ø4.3, H=1.5, C=6  |                                                  |          |                                                                                                                                                                                                                                                                   |
|    | CN-IN360 | NP       | C1      | Ø4.2, H=3, C=6    | Incisor esthetic Ti abutment, conical connection | Titanium |                                                                                             |
|                                                                                     | VN-IN360 | NP       | V3      | Ø4.2, H=3, C=6    |                                                  |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CS-IN360 | SP       | V3/C1   | Ø4.3, H=3, C=6    |                                                  |          |                                                                                                                                                                                                                                                                   |
|    | CN-TB002 | NP       | V3/C1   | Ø4.25, H=0.5, C=4 | Ti-Base, free rotation, conical connection       | Titanium |          |
|                                                                                     | CS-TB002 | SP       | V3/C1   | Ø4.3, H=0.5, C=4  |                                                  |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CW-TB002 | WP       | C1      | Ø5.7, H=0.5, C=4  |                                                  |          |                                                                                                                                                                                                                                                                   |
|    | CN-TBF15 | NP       | V3/C1   | Ø4.25, H=1.5, C=4 | Ti-Base, free rotation, conical connection       | Titanium |          |
|                                                                                     | CS-TBF15 | SP       | V3/C1   | Ø4.3, H=1.5, C=4  |                                                  |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CW-TBF15 | WP       | C1      | Ø5.7, H=1.5, C=4  |                                                  |          |                                                                                                                                                                                                                                                                   |
|    | CN-TBF30 | NP       | V3/C1   | Ø4.25, H=3, C=4   | Ti-Base, free rotation, conical connection       | Titanium |          |
|                                                                                     | CS-TBF30 | SP       | V3/C1   | Ø4.3, H=3, C=4    |                                                  |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CW-TBF30 | WP       | C1      | Ø5.7, H=3, C=4    |                                                  |          |                                                                                                                                                                                                                                                                   |
|  | MN-S0160 | NP       | V3/C1   | L=7.6             | Direct prosthetic screw, 30 Ncm                  | Titanium |    |
|                                                                                     | MD-S0220 | SP       | V3/C1   |                   |                                                  |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | MD-S0220 | WP       | C1      |                   |                                                  |          |                                                                                                                                                                                                                                                                   |


H = Gingival height

C = Crown height


L = Length



## 2 Engagement Types



Free Rotation



Anti-Rotation



Digital Restoration

Scan Post Kit

MK-0057  
Scan post kit



|  | Cat. No. | Platform | Implant | Dimensions (mm) | Description                   | Material | Platform |
|--|----------|----------|---------|-----------------|-------------------------------|----------|----------|
|  | CN-SP102 | NP       | C1      | L=10            | Scan Post, conical connection | Titanium |          |
|  | VN-SP102 | NP       | V3      |                 |                               |          |          |
|  | CS-SP102 | SP       | V3/C1   |                 |                               |          |          |
|  | CW-SP102 | WP       | C1      |                 |                               |          |          |
|  | MU-SP102 |          |         | L=10            | Scan Post for Multi-Unit      | Titanium |          |



Digital Restoration

## Digital Model Analogs



| Cat. No. | Platform | Implant | Dimensions (mm) | Description | Material | Platform |
|----------|----------|---------|-----------------|-------------|----------|----------|
|----------|----------|---------|-----------------|-------------|----------|----------|

|                                                                                     |          |    |       |              |                                             |          |                                                                                                                                                                                                                                                                   |
|-------------------------------------------------------------------------------------|----------|----|-------|--------------|---------------------------------------------|----------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | CN-MAN10 | NP | C1    | Ø3.3, L=10   | Digital model analog,<br>conical connection | Titanium |    |
|                                                                                     | VN-MAN10 | NP | V3    | Ø3.3, L=10   |                                             |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CS-MAN10 | SP | V3/C1 | Ø3.8, L=10   |                                             |          |                                                                                                                                                                                                                                                                   |
|                                                                                     | CW-MAN10 | WP | C1    | Ø4.6, L=10   |                                             |          |                                                                                                                                                                                                                                                                   |
|  | MU-MAN10 |    |       | Ø4.8<br>L=10 | Digital model analog for Multi-Unit         | Titanium |                                                                                                                                                                                                                                                                   |
|  | MM-MAN40 |    |       | Ø4<br>L=10   | CONNECT digital model analog                | Titanium |                                                                                                                                                                                                                                                                   |
|  | MM-MAN57 |    |       | Ø5.7<br>L=12 | CONNECT digital model analog                | Titanium |                                                                                                                                                                                                                                                                   |

H = Gingival height

C = Crown height

L = Length



Digital Restoration

## Titanium Blanks



|  | Cat. No.                                     | Platform             | Implant                 | Dimensions (mm) | Description                                                          | Material        | Platform |
|--|----------------------------------------------|----------------------|-------------------------|-----------------|----------------------------------------------------------------------|-----------------|----------|
|  | CN-IB091<br>VN-IB091<br>CS-IB091<br>CW-IB091 | NP<br>NP<br>SP<br>WP | C1<br>V3<br>V3/C1<br>C1 | Ø9              | Titanium Blank, anti-rotation, conical connection                    | Titanium        |          |
|  | CN-IB121<br>VN-IB121<br>CS-IB121<br>CW-IB121 | NP<br>NP<br>SP<br>WP | C1<br>V3<br>V3/C1<br>C1 | Ø12             | Titanium Blank, anti-rotation, conical connection                    | Titanium        |          |
|  | CN-AGT12<br>VN-AGT12<br>CS-AGT12<br>CW-AGT12 | NP<br>NP<br>SP<br>WP | C1<br>V3<br>V3/C1<br>C1 | Ø12             | Titanium Blank for Amann Girrbach, anti-rotation, conical connection | Titanium        |          |
|  | CN-MAB10<br>VN-MAB10<br>CS-MAB10<br>CW-MAB10 | NP<br>NP<br>SP<br>WP | C1<br>V3<br>V3/C1<br>C1 | Ø5<br>L=13.5    | Milling analog for Titanium Blank, conical connection                | Stainless Steel |          |

H = Gingival height

C = Crown height

L = Length





Cemented Restoration

## Cementable Abutments



|  | Cat. No. | Platform | Implant | Dimensions (mm)  | Description                                            | Material | Platform |
|--|----------|----------|---------|------------------|--------------------------------------------------------|----------|----------|
|  | CN-AN001 | NP       | C1      | Ø4.8, H=1, C=7   | Cementable abutment, conical connection                | Titanium | NP SP WP |
|  | VN-AT001 | NP       | V3      | Ø4.8, H=1, C=7   |                                                        |          |          |
|  | CS-A0010 | SP       | V3/C1   | Ø4.8, H=1, C=7.5 |                                                        |          |          |
|  | VS-AT001 | SP       | V3/C1   | Ø5.8, H=1, C=8   |                                                        |          |          |
|  | CW-AT001 | WP       | C1      | Ø5.5, H=1, C=7.5 |                                                        |          |          |
|  | CN-AT002 | NP       | C1      | Ø4.8, H=2, C=7   | Cementable abutment, conical connection                | Titanium | NP SP WP |
|  | VN-AT002 | NP       | V3      | Ø4.8, H=2, C=7   |                                                        |          |          |
|  | CS-AT002 | SP       | V3/C1   | Ø4.8, H=2, C=7.5 |                                                        |          |          |
|  | VS-AT002 | SP       | V3/C1   | Ø5.8, H=2, C=8   |                                                        |          |          |
|  | CW-AT002 | WP       | C1      | Ø5.5, H=2, C=7.5 |                                                        |          |          |
|  | CN-AN003 | NP       | C1      | Ø4.8, H=3, C=7   | Cementable abutment, conical connection                | Titanium | NP SP WP |
|  | VN-AT003 | NP       | V3      | Ø4.8, H=3, C=7   |                                                        |          |          |
|  | CS-A0030 | SP       | V3/C1   | Ø4.8, H=3, C=7.5 |                                                        |          |          |
|  | VS-AT003 | SP       | V3/C1   | Ø5.8, H=3, C=8   |                                                        |          |          |
|  | CW-AT003 | WP       | C1      | Ø5.5, H=3, C=7.5 |                                                        |          |          |
|  | CN-AN010 | NP       | C1      | Ø4.8, H=1, C=7   | Cementable angulated abutment, 10°, conical connection | Titanium | NP       |
|  | VN-AT101 | NP       | V3      | Ø4.8, H=1, C=7   |                                                        |          |          |
|  | CN-AT102 | NP       | C1      | Ø4.8, H=2, C=7   | Cementable angulated abutment, 10°, conical connection | Titanium | NP       |
|  | VN-AT102 | NP       | V3      | Ø4.8, H=2, C=7   |                                                        |          |          |

H = Gingival height

C = Crown height

L = Length

## Cementable Abutments

|                                                                                     | Cat. No.                         | Platform       | Implant                 | Dimensions (mm)                                          | Description                                            | Material | Platform                                                                                                                                                                                                                                                          |
|-------------------------------------------------------------------------------------|----------------------------------|----------------|-------------------------|----------------------------------------------------------|--------------------------------------------------------|----------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | CN-AN310<br>VN-AT103             | NP<br>NP       | C1<br>V3                | Ø4.8, H=3, C=7<br>Ø4.8, H=3, C=7                         | Cementable angulated abutment, 10°, conical connection | Titanium |                                                                                                                                                                                |
|    | CS-A1510<br>CS-AN151<br>CW-AN151 | SP<br>SP<br>WP | V3/C1<br>V3/C1<br>C1    | Ø4.8, H=1, C=7.5<br>Ø5.8, H=1, C=7.5<br>Ø5.5, H=1, C=7.5 | Cementable angulated abutment, 15°, conical connection | Titanium |                                                                                             |
|    | CS-AT152<br>VS-AT152<br>CW-AT152 | SP<br>SP<br>WP | V3/C1<br>V3/C1<br>C1    | Ø4.8, H=2, C=7.5<br>Ø5.8, H=2, C=7.5<br>Ø5.5, H=2, C=7.5 | Cementable angulated abutment, 15°, conical connection | Titanium |                                                                                             |
|    | CS-A1530<br>VS-AT153<br>CW-AT153 | SP<br>SP<br>WP | V3/C1<br>V3/C1<br>C1    | Ø4.8, H=3, C=7.5<br>Ø5.8, H=3, C=7.5<br>Ø5.5, H=3, C=7.5 | Cementable angulated abutment, 15°, conical connection | Titanium |                                                                                             |
|   | CN-AN020<br>VN-AT201             | NP<br>NP       | C1<br>V3                | Ø4.8, H=1, C=7<br>Ø4.8, H=1, C=7                         | Cementable angulated abutment, 20°, conical connection | Titanium |                                                                                                                                                                                |
|  | CN-AT202<br>VN-AT202             | NP<br>NP       | C1<br>V3                | Ø4.8, H=2, C=7<br>Ø4.8, H=2, C=7                         | Cementable angulated abutment, 20°, conical connection | Titanium |                                                                                                                                                                              |
|  | CN-AN320<br>VN-AT203             | NP<br>NP       | C1<br>V3                | Ø4.8, H=3, C=7<br>Ø4.8, H=3, C=7                         | Cementable angulated abutment, 20°, conical connection | Titanium |                                                                                                                                                                              |
|  | CS-A2510<br>CS-AN251<br>CW-AT251 | SP<br>SP<br>WP | V3/C1<br>V3/C1<br>C1    | Ø4.8, H=1, C=7.5<br>Ø5.8, H=1, C=7.5<br>Ø5.5, H=1, C=7.5 | Cementable angulated abutment, 25°, conical connection | Titanium |                                                                                         |
|  | CS-AT252<br>VS-AT252<br>CW-AT252 | SP<br>SP<br>WP | V3/C1<br>V3/C1<br>C1    | Ø4.8, H=2, C=7.5<br>Ø5.8, H=2, C=7.5<br>Ø5.5, H=2, C=7.5 | Cementable angulated abutment, 25°, conical connection | Titanium |                                                                                         |
|  | CS-A2530<br>VS-AT253<br>CW-AT253 | SP<br>SP<br>WP | V3/C1<br>V3/C1<br>C1    | Ø4.8, H=3, C=7.5<br>Ø5.8, H=3, C=7.5<br>Ø5.5, H=3, C=7.5 | Cementable angulated abutment, 25°, conical connection | Titanium |                                                                                         |
|  | MD-S0200                         | SP<br>WP       | V3/C1                   | L=6.9                                                    | Mini prosthetic screw, 30 Ncm                          | Titanium |                                                                                         |
|  | MN-S0160<br>MD-S0220<br>MD-S0220 | NP<br>SP<br>WP | V3/C1<br>V3/C1<br>V3/C1 | L=7.6                                                    | Direct prosthetic screw, 30 Ncm                        | Titanium |    |

H = Gingival height

C = Crown height

L = Length



Cemented Restoration

## Cementing Posts



|  | Cat. No. | Platform | Implant | Dimensions (mm) | Description                          | Material | Platform |
|--|----------|----------|---------|-----------------|--------------------------------------|----------|----------|
|  | CN-MAC10 | NP       | C1      | Ø4, H=1, C=10   | Cementing post, conical connection   | Titanium |          |
|  | VN-MC101 | NP       | V3      | Ø4, H=1, C=10   |                                      |          |          |
|  | CS-MAC10 | SP       | V3/C1   | Ø4.8, H=1, C=10 |                                      |          |          |
|  | VS-MC101 | SP       | V3/C1   | Ø5.8, H=1, C=10 |                                      |          |          |
|  | CW-SMAC1 | WP       | C1      | Ø5.5, H=1, C=10 |                                      |          |          |
|  | CN-MAC30 | NP       | C1      | Ø4, H=3, C=10   | Cementing post, conical connection   | Titanium |          |
|  | VN-MC103 | NP       | V3      | Ø4, H=3, C=10   |                                      |          |          |
|  | CS-MC103 | SP       | V3/C1   | Ø4.8, H=3, C=10 |                                      |          |          |
|  | VS-MC103 | SP       | V3/C1   | Ø5.8, H=3, C=10 |                                      |          |          |
|  | CW-MC103 | WP       | C1      | Ø5.5, H=3, C=10 |                                      |          |          |
|  | MN-S0160 | NP       | V3/C1   | L=7.6           | Direct prosthetic screw, 30 Ncm      | Titanium |          |
|  | MD-S0220 | SP       |         |                 |                                      |          |          |
|  | MD-S0220 | WP       |         |                 |                                      |          |          |
|  | MD-G0220 | SP<br>WP | V3/C1   | L=7.7           | Direct gold prosthetic screw, 20 Ncm | Gold-s   |          |

H = Gingival height

C = Crown height

L = Length



Cemented Restoration

## CPK - Complete Prosthetic Kit



|  | Cat. No. | Platform | Implant | Dimensions (mm) | Description                                 | Material | Platform |
|--|----------|----------|---------|-----------------|---------------------------------------------|----------|----------|
|  | CK-NPK41 | NP       | C1      | Ø4, H=1, C=4    | Complete Prosthetic Kit, conical connection | Titanium | NP SP    |
|  | VK-NPK41 | NP       | V3      | Ø4, H=1, C=4    |                                             |          |          |
|  | CK-CPK41 | SP       | V3/C1   | Ø4.8, H=1, C=4  |                                             |          |          |
|  | CK-NPK42 | NP       | C1      | Ø4, H=2, C=4    | Complete Prosthetic Kit, conical connection | Titanium | NP SP    |
|  | VK-NPK42 | NP       | V3      | Ø4, H=2, C=4    |                                             |          |          |
|  | CK-CPK42 | SP       | V3/C1   | Ø4.8, H=2, C=4  |                                             |          |          |
|  | CK-NPK43 | NP       | C1      | Ø4, H=3, C=4    | Complete Prosthetic Kit, conical connection | Titanium | NP SP    |
|  | VK-NPK43 | NP       | V3      | Ø4, H=3, C=4    |                                             |          |          |
|  | CK-CPK43 | SP       | V3/C1   | Ø4.8, H=3, C=4  |                                             |          |          |

H = Gingival height

C = Crown height

L = Length



|  | Cat. No.                                                 | Platform                   | Implant                          | Dimensions (mm)                                                                    | Description                                 | Material | Platform |
|--|----------------------------------------------------------|----------------------------|----------------------------------|------------------------------------------------------------------------------------|---------------------------------------------|----------|----------|
|  | CK-NPK44<br>VK-NPK44<br>CK-CPK44                         | NP<br>NP<br>SP             | C1<br>V3<br>V3/C1                | Ø4, H=4, C=4<br>Ø4, H=4, C=4<br>Ø4.8, H=4, C=4                                     | Complete Prosthetic Kit, conical connection | Titanium | NP SP    |
|  | CK-NPK61<br>VK-NPK61<br>CK-CPK61<br>VK-CPK61<br>CK-WPK61 | NP<br>NP<br>SP<br>SP<br>WP | C1<br>V3<br>V3/C1<br>V3/C1<br>C1 | Ø4, H=1, C=6<br>Ø4, H=1, C=6<br>Ø4.8, H=1, C=6<br>Ø5.5, H=1, C=6<br>Ø5.5, H=1, C=6 | Complete Prosthetic Kit, conical connection | Titanium | NP SP WP |
|  | CK-NPK62<br>VK-NPK62<br>CK-CPK62<br>VK-CPK62<br>CK-WPK62 | NP<br>NP<br>SP<br>SP<br>WP | C1<br>V3<br>V3/C1<br>V3/C1<br>C1 | Ø4, H=2, C=6<br>Ø4, H=2, C=6<br>Ø4.8, H=2, C=6<br>Ø5.5, H=2, C=6<br>Ø5.5, H=2, C=6 | Complete Prosthetic Kit, conical connection | Titanium | NP SP WP |
|  | CK-NPK63<br>VK-NPK63<br>CK-CPK63<br>VK-CPK63<br>CK-WPK63 | NP<br>NP<br>SP<br>SP<br>WP | C1<br>V3<br>V3/C1<br>V3/C1<br>C1 | Ø4, H=3, C=6<br>Ø4, H=3, C=6<br>Ø4.8, H=3, C=6<br>Ø5.5, H=3, C=6<br>Ø5.5, H=3, C=6 | Complete Prosthetic Kit, conical connection | Titanium | NP SP WP |
|  | CK-NPK64<br>VK-NPK64<br>CK-CPK64<br>VK-CPK64<br>CK-WPK64 | NP<br>NP<br>SP<br>SP<br>WP | C1<br>V3<br>V3/C1<br>V3/C1<br>C1 | Ø4, H=4, C=6<br>Ø4, H=4, C=6<br>Ø4.8, H=4, C=6<br>Ø5.5, H=4, C=6<br>Ø5.5, H=4, C=6 | Complete Prosthetic Kit, conical connection | Titanium | NP SP WP |
|  | CK-CPK81<br>VK-CPK81                                     | SP<br>SP                   | V3/C1<br>V3/C1                   | Ø4.8, H=1, C=8<br>Ø5.5, H=1, C=8                                                   | Complete Prosthetic Kit, conical connection | Titanium | SP       |
|  | CK-CPK82<br>VK-CPK82                                     | SP<br>SP                   | V3/C1<br>V3/C1                   | Ø4.8, H=2, C=8<br>Ø5.5, H=2, C=8                                                   | Complete Prosthetic Kit, conical connection | Titanium | SP       |
|  | CK-CPK83<br>VK-CPK83                                     | SP<br>SP                   | V3/C1<br>V3/C1                   | Ø4.8, H=3, C=8<br>Ø5.5, H=3, C=8                                                   | Complete Prosthetic Kit, conical connection | Titanium | SP       |
|  | CK-CPK84<br>VK-CPK84                                     | SP<br>SP                   | V3/C1<br>V3/C1                   | Ø4.8, H=4, C=8<br>Ø5.5, H=4, C=8                                                   | Complete Prosthetic Kit, conical connection | Titanium | SP       |
|  | CN-CPK41<br>VN-P4140<br>CS-CPK41                         | NP<br>NP<br>SP             | C1<br>V3<br>V3/C1                | Ø4, H=1, C=4<br>Ø4, H=1, C=4<br>Ø4.8, H=1, C=4                                     | Transgingival abutment, conical connection  | Titanium | NP SP    |
|  | CN-CPK42<br>VN-P4240<br>CS-CPK42                         | NP<br>NP<br>SP             | C1<br>V3<br>V3/C1                | Ø4, H=2, C=4<br>Ø4, H=2, C=4<br>Ø4.8, H=2, C=4                                     | Transgingival abutment, conical connection  | Titanium | NP SP    |
|  | CN-CPK43<br>VN-P4340<br>CS-CPK43                         | NP<br>NP<br>SP             | C1<br>V3<br>V3/C1                | Ø4, H=3, C=4<br>Ø4, H=3, C=4<br>Ø4.8, H=3, C=4                                     | Transgingival abutment, conical connection  | Titanium | NP SP    |

H = Gingival height

C = Crown height

L = Length

# CPK - Complete Prosthetic Kit

|                                                                                     | Cat. No. | Platform | Implant | Dimensions (mm) | Description                                | Material | Platform                                                                                                                                                                                                                                                    |
|-------------------------------------------------------------------------------------|----------|----------|---------|-----------------|--------------------------------------------|----------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | CN-CPK44 | NP       | C1      | Ø4, H=4, C=4    | Transgingival abutment, conical connection | Titanium |                                                                                       |
|                                                                                     | VN-P4440 | NP       | V3      | Ø4, H=4, C=4    |                                            |          |                                                                                                                                                                                                                                                             |
|                                                                                     | CS-CPK44 | SP       | V3/C1   | Ø4.8, H=4, C=4  |                                            |          |                                                                                                                                                                                                                                                             |
|    | CN-CPK61 | NP       | C1      | Ø4, H=1, C=6,   | Transgingival abutment, conical connection | Titanium |    |
|                                                                                     | VN-P6140 | NP       | V3      | Ø4, H=1, C=6,   |                                            |          |                                                                                                                                                                                                                                                             |
|                                                                                     | CS-CPK61 | SP       | V3/C1   | Ø4.8, H=1, C=6  |                                            |          |                                                                                                                                                                                                                                                             |
|                                                                                     | VS-P6155 | SP       | V3/C1   | Ø5.5, H=1, C=6  |                                            |          |                                                                                                                                                                                                                                                             |
|                                                                                     | CW-CPK61 | WP       | C1      | Ø5.5, H=1, C=6  |                                            |          |                                                                                                                                                                                                                                                             |
|    | CN-CPK62 | NP       | C1      | Ø4, H=2, C=6    | Transgingival abutment, conical connection | Titanium |    |
|                                                                                     | VN-P6240 | NP       | V3      | Ø4, H=2, C=6    |                                            |          |                                                                                                                                                                                                                                                             |
|                                                                                     | CS-CPK62 | SP       | V3/C1   | Ø4.8, H=2, C=6  |                                            |          |                                                                                                                                                                                                                                                             |
|                                                                                     | VS-P6255 | SP       | V3/C1   | Ø5.5, H=2, C=6  |                                            |          |                                                                                                                                                                                                                                                             |
|                                                                                     | CW-CPK62 | WP       | C1      | Ø5.5, H=2, C=6  |                                            |          |                                                                                                                                                                                                                                                             |
|    | CN-CPK63 | NP       | C1      | Ø4, H=3, C=6    | Transgingival abutment, conical connection | Titanium |    |
|                                                                                     | VN-P6340 | NP       | V3      | Ø4, H=3, C=6    |                                            |          |                                                                                                                                                                                                                                                             |
|                                                                                     | CS-CPK63 | SP       | V3/C1   | Ø4.8, H=3, C=6  |                                            |          |                                                                                                                                                                                                                                                             |
|                                                                                     | VS-P6355 | SP       | V3/C1   | Ø5.5, H=3, C=6  |                                            |          |                                                                                                                                                                                                                                                             |
|                                                                                     | CW-CPK63 | WP       | C1      | Ø5.5, H=3, C=6  |                                            |          |                                                                                                                                                                                                                                                             |
|    | CN-CPK64 | NP       | C1      | Ø4, H=4, C=6    | Transgingival abutment, conical connection | Titanium |    |
|                                                                                     | VN-P6440 | NP       | V3      | Ø4, H=4, C=6    |                                            |          |                                                                                                                                                                                                                                                             |
|                                                                                     | CS-CPK64 | SP       | V3/C1   | Ø4.8, H=4, C=6  |                                            |          |                                                                                                                                                                                                                                                             |
|                                                                                     | VS-P6455 | SP       | V3/C1   | Ø5.5, H=4, C=6  |                                            |          |                                                                                                                                                                                                                                                             |
|                                                                                     | CW-CPK64 | WP       | C1      | Ø5.5, H=4, C=6  |                                            |          |                                                                                                                                                                                                                                                             |
|  | CS-CPK81 | SP       | V3/C1   | Ø4.8, H=1, C=8  | Transgingival abutment, conical connection | Titanium |                                                                                                                                                                        |
|                                                                                     | VS-P8155 | SP       |         | Ø5.5, H=1, C=8  |                                            |          |                                                                                                                                                                                                                                                             |
|  | CS-CPK82 | SP       | V3/C1   | Ø4.8, H=2, C=8  | Transgingival abutment, conical connection | Titanium |                                                                                                                                                                        |
|                                                                                     | VS-P8255 | SP       |         | Ø5.5, H=2, C=8  |                                            |          |                                                                                                                                                                                                                                                             |
|  | CS-CPK83 | SP       | V3/C1   | Ø4.8, H=3, C=8  | Transgingival abutment, conical connection | Titanium |                                                                                                                                                                        |
|                                                                                     | VS-P8355 | SP       |         | Ø5.5, H=3, C=8  |                                            |          |                                                                                                                                                                                                                                                             |
|  | CS-CPK84 | SP       | V3/C1   | Ø4.8, H=4, C=8  | Transgingival abutment, conical connection | Titanium |                                                                                                                                                                        |
|                                                                                     | VS-P8455 | SP       |         | Ø5.5, H=4, C=8  |                                            |          |                                                                                                                                                                                                                                                             |

H = Gingival height

C = Crown height

L = Length



Overdentures

## LOCKiT Attachment System



|  | Cat. No.                         | Implant              | Platform       | Dimensions (mm)                       | Description                         | Material | Platform |
|--|----------------------------------|----------------------|----------------|---------------------------------------|-------------------------------------|----------|----------|
|  | CN-LOC01<br>CS-LOC01<br>CW-LOC01 | V3/C1<br>V3/C1<br>C1 | NP<br>SP<br>WP | Ø3.75, H=1<br>Ø3.75, H=1<br>Ø4.5, H=1 | LOCKiT abutment, conical connection | Titanium | NP SP WP |
|  | CN-LOC02<br>CS-LOC02             | V3/C1<br>V3/C1       | NP<br>SP       | Ø3.75, H=2<br>Ø3.75, H=2              | LOCKiT abutment, conical connection | Titanium | NP SP    |
|  | CN-LOC03<br>CS-LOC03<br>CW-LOC03 | V3/C1<br>V3/C1<br>C1 | NP<br>SP<br>WP | Ø3.75, H=3<br>Ø3.75, H=3<br>Ø4.5, H=3 | LOCKiT abutment, conical connection | Titanium | NP SP WP |
|  | CN-LOC04<br>CS-LOC04             | V3/C1<br>V3/C1       | NP<br>SP       | Ø3.75, H=4<br>Ø3.75, H=4              | LOCKiT abutment, conical connection | Titanium | NP SP    |
|  | CN-LOC05<br>CS-LOC05<br>CW-LOC05 | V3/C1<br>V3/C1<br>C1 | NP<br>SP<br>WP | Ø3.75, H=5<br>Ø3.75, H=5<br>Ø4.5, H=5 | LOCKiT abutment, conical connection | Titanium | NP SP WP |

H = Gingival height

C = Crown height

L = Length

# LOCKiT Attachment System

|                                                                                     | Cat. No.                      | Implant              | Platform       | Dimensions (mm)                       | Description                                  | Material        | Platform |
|-------------------------------------------------------------------------------------|-------------------------------|----------------------|----------------|---------------------------------------|----------------------------------------------|-----------------|----------|
|    | CK-NLC1<br>CK-SLC1<br>CK-WLC1 | V3/C1<br>V3/C1<br>C1 | NP<br>SP<br>WP | Ø3.75, H=1<br>Ø3.75, H=1<br>Ø4.5, H=1 | LOCKiT abutment kit, conical connection      | Titanium        | NP SP WP |
|    | CK-NLC2<br>CK-SLC2            | V3/C1<br>V3/C1       | NP<br>SP       | Ø3.75, H=2<br>Ø3.75, H=2              | LOCKiT abutment kit, conical connection      | Titanium        | NP SP    |
|    | CK-NLC3<br>CK-SLC3<br>CK-WLC3 | V3/C1<br>V3/C1<br>C1 | NP<br>SP<br>WP | Ø3.75, H=3<br>Ø3.75, H=3<br>Ø4.5, H=3 | LOCKiT abutment kit, conical connection      | Titanium        | NP SP WP |
|    | CK-NLC4<br>CK-SLC4            | V3/C1<br>V3/C1       | NP<br>SP       | Ø3.75, H=4<br>Ø3.75, H=4              | LOCKiT abutment kit, conical connection      | Titanium        | NP SP    |
|   | CK-NLC5<br>CK-SLC5<br>CK-WLC5 | V3/C1<br>V3/C1<br>C1 | NP<br>SP<br>WP | Ø3.75, H=5<br>Ø3.75, H=5<br>Ø4.5, H=5 | LOCKiT abutment kit, conical connection      | Titanium        | NP SP WP |
|  | MT-LOC10                      |                      |                | Ø7.8<br>L=86                          | LOCKiT core tool                             | Stainless steel |          |
|  | MM-LOC1                       |                      |                | Ø5.4<br>H=4.4                         | LOCKiT impression coping                     |                 |          |
|  | MM-LOA40<br>MM-LOA50          | NP/SP<br>WP          |                | Ø4, L=10.2<br>Ø5, L=10.2              | LOCKiT analog                                |                 |          |
|  | MM-LOC02                      |                      |                |                                       | LOCKiT male processing package               |                 |          |
|  | MM-CLE10                      |                      |                | Ø4.75<br>H=1.78                       | LOCKiT male retention cap (clear)            | Nylon           |          |
|  | MM-PIN10                      |                      |                | Ø4.72<br>H=1.78                       | LOCKiT light retention male cap (pink)       | Nylon pink      |          |
|  | MM-BLU10                      |                      |                | Ø4.72<br>H=1.78                       | LOCKiT extra light male retention cap (blue) | Nylon blue      |          |

H = Gingival height

C = Crown height

L = Length



|                                                                                     | Cat. No. | Implant | Platform | Dimensions (mm) | Description                                                | Material           | Platform |
|-------------------------------------------------------------------------------------|----------|---------|----------|-----------------|------------------------------------------------------------|--------------------|----------|
|    | MM-GRE10 |         |          | Ø4.75<br>H=1.78 | LOCKiT extended range male retention cap (green)           | Nylon green        |          |
|    | MM-RED10 |         |          | Ø4.72<br>H=1.78 | LOCKiT extra light extended range male retention cap (red) | Nylon red          |          |
|    | MM-BLA10 |         |          | Ø4.75<br>H=1.9  | LOCKiT male processing cap (black)                         | Polyethylene black |          |
|    | MM-TIT10 |         |          | Ø5.5<br>H=2.25  | LOCKiT metal housing titanium cap                          |                    |          |
|   | MT-RSD19 |         |          | L=19            | LOCKiT short torque ratchet driver                         | Stainless steel    |          |
|  | MT-RSD24 |         |          | L=24            | LOCKiT long torque ratchet driver                          | Stainless steel    |          |
|  | MT-LOC11 |         |          | Ø7.35<br>H=4.1  | LOCKiT abutment holder sleeve                              |                    |          |


H = Gingival height

C = Crown height


L = Length



## 5 Retention Options




Extra Light Retention




Light Retention



Standard Retention



Extra Light Retention



Standard Retention



Overdentures

# OT-Equator Kits



|  | Cat. No.                      | Implant              | Platform       | Dimensions (mm)                   | Description                        | Material | Platform |
|--|-------------------------------|----------------------|----------------|-----------------------------------|------------------------------------|----------|----------|
|  | CK-NOE1<br>CK-SOE1<br>CK-WOE1 | V3/C1<br>V3/C1<br>C1 | NP<br>SP<br>WP | Ø4, H=1<br>Ø4.4, H=1<br>Ø4.4, H=1 | OT-Equator kit, conical connection | Titanium |          |
|  | CK-NOE2<br>CK-SOE2            | V3/C1<br>V3/C1       | NP<br>SP       | Ø4, H=2<br>Ø4.4, H=2<br>Ø4.4, H=2 | OT-Equator kit, conical connection | Titanium |          |
|  | CK-NOE3<br>CK-SOE3<br>CK-WOE3 | V3/C1<br>V3/C1<br>C1 | NP<br>SP<br>WP | Ø4, H=3<br>Ø4.4, H=3<br>Ø4.4, H=3 | OT-Equator kit, conical connection | Titanium |          |
|  | CK-NOE4<br>CK-SOE4            | V3/C1<br>V3/C1       | NP<br>SP       | Ø4, H=4<br>Ø4.4, H=4              | OT-Equator kit, conical connection | Titanium |          |
|  | CK-NOE5<br>CK-SOE5            | V3/C1<br>V3/C1       | NP<br>SP       | Ø4, H=5<br>Ø4.4, H=5              | OT-Equator kit, conical connection | Titanium |          |

H = Gingival height

C = Crown height

L = Length height

# OT-Equator Kits

|                                                                                     | Cat. No. | Implant | Platform | Dimensions (mm) | Description                                  | Material        | Platform |
|-------------------------------------------------------------------------------------|----------|---------|----------|-----------------|----------------------------------------------|-----------------|----------|
|    | ET-IT002 |         |          |                 | OT-Equator handling tool                     |                 |          |
|    | OE-IC002 |         |          |                 | OT-Equator impression coping cap             |                 |          |
|    | OE-RSM60 |         |          |                 | OT-Equator analog                            | Stainless steel |          |
|    | MT-RDL30 |         |          | L=27            | Long driver for 0.05 inch hex.               | Stainless steel |          |
|   | MT-RDS30 |         |          | L=21.3          | Short driver for 0.05 inch hex.              | Stainless steel |          |
|  | OE-RCP01 |         |          |                 | OT-Equator soft retentive cap (pink)         |                 |          |
|  | OE-RCW01 |         |          |                 | OT-Equator standard retentive cap (white)    |                 |          |
|  | OE-RCV01 |         |          |                 | OT-Equator rigid retentive cap (violet)      |                 |          |
|  | OE-RCY01 |         |          |                 | OT-Equator extra soft retentive cap (yellow) |                 |          |
|  | OE-RCB01 |         |          |                 | OT-Equator laboratory retentive cap (black)  |                 |          |
|  | OE-IC001 |         |          |                 | OT-Equator impression transfer cap           |                 |          |

H = Gingival height

C = Crown height

L = Length



Overdentures

## Ball Attachment System



|  | Cat. No.                         | Implant              | Platform       | Dimensions (mm)                   | Description                                | Material                          | Platform |
|--|----------------------------------|----------------------|----------------|-----------------------------------|--------------------------------------------|-----------------------------------|----------|
|  | CN-B1330<br>CS-B1375<br>CW-B1500 | V3/C1<br>V3/C1<br>C1 | NP<br>SP<br>WP | Ø4, H=1<br>Ø4.4, H=1<br>Ø4.4, H=1 | Ball attachment anchor, conical connection | Titanium alloy + Titanium nitride | NP SP WP |
|  | CN-B2330<br>CS-B2375             | V3/C1                | NP<br>SP       | Ø4, H=2<br>Ø4.4, H=2              | Ball attachment anchor, conical connection | Titanium alloy + Titanium nitride | NP SP    |
|  | CN-B3330<br>CS-B3375<br>CW-B3500 | V3/C1<br>V3/C1<br>C1 | NP<br>SP<br>WP | Ø4, H=3<br>Ø4.4, H=3<br>Ø4.4, H=3 | Ball attachment anchor, conical connection | Titanium alloy + Titanium nitride | NP SP WP |
|  | CN-B4330<br>CS-B4375             | V3/C1                | NP<br>SP       | Ø4, H=4<br>Ø4.4, H=4              | Ball attachment anchor, conical connection | Titanium alloy + Titanium nitride | NP SP    |
|  | CS-B5375                         | V3/C1                | SP             | Ø4.4, H=5                         | Ball attachment anchor, conical connection | Titanium alloy + Titanium nitride | SP       |
|  | MB-RS375                         |                      |                | L=15.1                            | Ball attachment analog                     | Stainless steel                   |          |

H = Gingival height

C = Crown height

L = Length


## Ball Attachment System

|                                                                                     | Cat. No. | Implant | Platform | Dimensions (mm) | Description                                            | Material        | Platform |
|-------------------------------------------------------------------------------------|----------|---------|----------|-----------------|--------------------------------------------------------|-----------------|----------|
|    | MB-SF200 |         |          | Ø3.7<br>L=3.5   | Plastic ball cap                                       | Plastic         |          |
|    | MB-GPS10 |         |          | L=2.4           | Gold ball cap                                          | Gold-e          |          |
|    | MB-TTP10 |         |          | Ø3.5<br>L=3.1   | Titanium ball cap                                      | Titanium        |          |
|    | MB-PR010 |         |          | Ø3.8<br>L=2     | Plastic mounting sleeve                                | Plastic         |          |
|    | MB-TBS10 |         |          | Ø2.5            | Titanium ball cap spring, standard                     | Steel           |          |
|  | MB-TBS20 |         |          | Ø2.6            | Titanium ball cap spring, soft                         | Steel           |          |
|  | MB-DMH10 |         |          |                 | Double metal housing                                   |                 |          |
|  | MB-PPC10 |         |          | Ø4.2<br>L=2.6   | Standard plastic cap for metal housing                 | Plastic         |          |
|  | MB-PPR10 |         |          | Ø4.2<br>L=2.6   | Soft plastic cap for metal housing, 800-950g retention | Plastic         |          |
|  | MB-DB235 |         |          | Ø15<br>L=2.6    | Plastic disc for ball attachment                       | Plastic         |          |
|  | MT-RB225 |         |          | L=22            | Ball attachment anchor key                             | Stainless steel |          |
|  | MT-NB225 |         |          | L=22            | Ball attachment anchor key                             | Stainless steel |          |

H = Gingival height


C = Crown height

L = Length



|  | Cat. No. | Implant | Platform | Dimensions (mm)       | Description       | Material        | Platform |
|--|----------|---------|----------|-----------------------|-------------------|-----------------|----------|
|  | MM-TRU10 |         |          | Ø2.8<br>L=50<br>H=2.9 | Titanium rider    | Titanium        |          |
|  | MM-PBU10 |         |          | Ø1.6<br>L=80<br>H=2.3 | Plastic bar       | Plastic         |          |
|  | MM-GB100 |         |          | Ø1.9<br>L=50          | Gold bar          | Gold-e          |          |
|  | MM-GBU10 |         |          | L=50                  | Gold bar          | Gold-e          |          |
|  | MM-PC100 |         |          | Ø1.9<br>L=85          | Plastic round bar | Burnout plastic |          |







## MGUIDE - DIGITAL DENTISTRY

MGUIDE is an innovative guided surgery system, which was designed to provide clinicians with a worry-free solution, and flawlessly incorporates prosthetically driven planning and easy-to-use tools, to deliver the most predictable results.





|                                                                                     | Cat. No. | Dimensions (mm) | Description                            | Material | Sold |
|-------------------------------------------------------------------------------------|----------|-----------------|----------------------------------------|----------|------|
|  | MG-TMP10 |                 | MGUIDE value pack                      |          |      |
|  | MG-TMP20 |                 | MGUIDE value pack for edentulous cases |          |      |



MGUIDE

## MGUIDE Surgical Set



MG-K004   
MGUIDE Surgical Set

|                                                                                     | Cat. No. | Dimensions (mm) | Description                                | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|-----------------|--------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | MG-TP055 | Ø5.5<br>L=26    | MGUIDE tissue punch                        | Stainless steel |   |
|  | MG-BM047 | Ø4.7<br>L=24.5  | MGUIDE bone mill                           | Stainless steel |   |
|  | MG-BM055 | Ø5.5<br>L=24.5  | MGUIDE bone mill                           | Stainless steel |   |
|  | MG-D0624 | Ø2.4<br>L=30.8  | MGUIDE pilot drill for 6mm length implant  | Stainless steel |   |
|  | MG-D0824 | Ø2.4<br>L=32.8  | MGUIDE pilot drill for 8mm length implant  | Stainless steel |   |
|  | MG-D1024 | Ø2.4<br>L=34.8  | MGUIDE pilot drill for 10mm length implant | Stainless steel |   |

L = Length

# MGUIDE Surgical Set

|                                                                                     | Cat. No. | Dimensions (mm)     | Description                                                                  | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|---------------------|------------------------------------------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | MG-D1124 | Ø2.4<br>L=36.3      | MGUIDE pilot drill for 11.5mm length implant                                 | Stainless steel |       |
|    | MG-D1324 | Ø2.4<br>L=37.8      | MGUIDE pilot drill for 13mm length implant                                   | Stainless steel |       |
|    | CG-D0833 | Ø2.3-3.25<br>L=32.8 | MGUIDE drill for Ø3.30 x 8mm length implant, conical connection              | Stainless steel |       |
|    | CG-D1033 | Ø2.3-3.25<br>L=34.8 | MGUIDE drill for Ø3.30 x 10mm length implant, conical connection             | Stainless steel |       |
|   | CG-D1133 | Ø2.3-3.25<br>L=36.3 | MGUIDE drill for Ø3.30 x 11.5mm length implant, conical connection           | Stainless steel |       |
|  | CG-D1333 | Ø2.3-3.25<br>L=37.8 | MGUIDE drill for Ø3.30 x 13mm length implant, conical connection             | Stainless steel |   |
|  | CG-D0839 | Ø2.9-3.85<br>L=32.8 | MGUIDE drill for Ø3.75 and Ø3.90 x 8mm length implant, conical connection    | Stainless steel |   |
|  | CG-D1039 | Ø2.9-3.85<br>L=34.8 | MGUIDE drill for Ø3.75 and Ø3.90 x 10mm length implant, conical connection   | Stainless steel |   |
|  | CG-D1139 | Ø2.9-3.85<br>L=36.3 | MGUIDE drill for Ø3.75 and Ø3.90 x 11.5mm length implant, conical connection | Stainless steel |   |
|  | CG-D1339 | Ø2.9-3.85<br>L=37.8 | MGUIDE drill for Ø3.75 and Ø3.90 x 13mm length implant, conical connection   | Stainless steel |   |
|  | CG-D0843 | Ø3.4-4.25<br>L=32.8 | MGUIDE drill for Ø4.20 and Ø4.30 x 8mm length implant, conical connection    | Stainless steel |   |
|  | CG-D1043 | Ø3.4-4.25<br>L=34.8 | MGUIDE drill for Ø4.20 and Ø4.30 x 10mm length implant, conical connection   | Stainless steel |   |

L = Length

|                                                                                     | Cat. No. | Dimensions (mm)     | Description                                                                  | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|---------------------|------------------------------------------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | CG-D1143 | Ø3.4-4.25<br>L=36.3 | MGUIDE drill for Ø4.20 and Ø4.30 x 11.5mm length implant, conical connection | Stainless steel |       |
|    | CG-D1343 | Ø3.4-4.25<br>L=37.8 | MGUIDE drill for Ø4.20 and Ø4.30 x 13mm length implant, conical connection   | Stainless steel |       |
|    | CG-D0850 | Ø3.9-4.95<br>L=32.8 | MGUIDE drill for Ø5 x 8mm length implant, conical connection                 | Stainless steel |       |
|    | CG-D1050 | Ø3.9-4.95<br>L=34.8 | MGUIDE drill for Ø5 x 10mm length implant, conical connection                | Stainless steel |       |
|   | CG-D1150 | Ø3.9-4.95<br>L=36.3 | MGUIDE drill for Ø5 x 11.5mm length implant, conical connection              | Stainless steel |     |
|  | CG-D1350 | Ø3.9-4.95<br>L=37.8 | MGUIDE drill for Ø5 x 13mm length implant, conical connection                | Stainless steel |   |
|  | MG-DLG55 | Ø5.5<br>L=55        | MGUIDE drill length gauge                                                    | Titanium        |   |
|  | CG-GRN01 | Ø5.5<br>L=30.5      | MGUIDE C1 direct ratchet insertion tool, conical connection, narrow platform | Stainless steel |   |
|  | CG-GMN10 | Ø5.5<br>L=31        | MGUIDE C1 motor insertion tool, conical connection, narrow platform          | Stainless steel |   |
|  | CG-GRN10 | Ø5.5<br>L=26        | MGUIDE C1 ratchet insertion tool, conical connection, narrow platform        | Stainless steel |   |
|  | VG-GRN01 | Ø5.5<br>L=30.5      | MGUIDE V3 direct ratchet insertion tool, conical connection, narrow platform | Stainless steel |   |
|  | VG-GMN10 | Ø5.5<br>L=31        | MGUIDE V3 motor insertion tool, conical connection, narrow platform          | Stainless steel |   |

L = Length



# MGUIDE Surgical Set

|                                                                                     | Cat. No. | Dimensions (mm) | Description                                                                       | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|-----------------|-----------------------------------------------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | VG-GRN10 | Ø5.5<br>L=25.5  | MGUIDE V3 ratchet insertion tool, conical connection, narrow platform             | Stainless steel |       |
|    | CG-GRS01 | Ø5.5<br>L=31    | MGUIDE V3/C1 direct ratchet insertion tool, conical connection, standard platform | Stainless steel |       |
|    | CG-GMS10 | Ø5.5<br>L=32    | MGUIDE V3/C1 motor insertion tool, conical connection, standard platform          | Stainless steel |       |
|    | CG-GRS10 | Ø5.5<br>L=27    | MGUIDE C1 ratchet insertion tool, conical connection, standard platform           | Stainless steel |       |
|    | VG-GRS10 | Ø5.5<br>L=26    | MGUIDE V3 ratchet insertion tool, conical connection, standard platform           | Stainless steel |       |
|  | CG-GRW01 | Ø5.5<br>L=31    | MGUIDE C1 direct ratchet insertion tool, conical connection, wide platform        | Stainless steel |   |
|  | CG-GMW10 | Ø5.5<br>L=32    | MGUIDE C1 motor insertion tool, conical connection, wide platform                 | Stainless steel |   |
|  | CG-GRW10 | Ø5.5<br>L=27    | MGUIDE C1 ratchet insertion tool, conical connection, wide platform               | Stainless steel |   |
|  | CG-TAS55 | Ø5.5<br>L=24    | MGUIDE template anchoring screw, conical connection, standard platform            | Titanium        |   |
|  | MG-FP020 | Ø2<br>L=29      | MGUIDE fixation pin                                                               | Titanium        |   |
|  | MG-DFP20 | Ø2<br>L=35.3    | MGUIDE drill for fixation pin                                                     | Stainless steel |   |
|  | MT-RDL30 | L=27            | Long driver for 0.05 inch hex.                                                    | Stainless steel |   |

|                                                                                   | Cat. No. | Dimensions (mm) | Description             | Material        | Sold                                                                                                                                                                    |
|-----------------------------------------------------------------------------------|----------|-----------------|-------------------------|-----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | MT-RA480 |                 | Ratchet to hex. adapter | Stainless steel |   |
|  | MT-RI030 | L=84            | Monoblock ratchet       |                 |   |

L = Length

# 6 Diameters





MG-K005  
MGUIDE Surgical Set,  
Narrow Sleeve



|                                                                                     | Cat. No. | Dimensions (mm) | Description                                             | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|-----------------|---------------------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | MG-NBM40 | Ø4<br>L=24.5    | MGUIDE Narrow Sleeve bone mill                          | Stainless steel |   |
|  | MG-NTP40 | Ø4<br>L=26      | MGUIDE Narrow Sleeve tissue punch                       | Stainless steel |   |
|  | MG-N0624 | Ø2.4<br>L=30.8  | MGUIDE Narrow Sleeve pilot drill for 6mm length implant | Stainless steel |   |
|  | MG-N0824 | Ø2.4<br>L=32.8  | MGUIDE Narrow Sleeve pilot drill for 8mm length implant | Stainless steel |   |

L = Length



|                                                                                     | Cat. No. | Dimensions (mm)     | Description                                                                                | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|---------------------|--------------------------------------------------------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | MG-N1024 | Ø2.4<br>L=34.8      | MGUIDE Narrow Sleeve pilot drill for 10mm length implant                                   | Stainless steel |       |
|    | MG-N1124 | Ø2.4<br>L=36.3      | MGUIDE Narrow Sleeve pilot drill for 11.5mm length implant                                 | Stainless steel |       |
|    | MG-N1324 | Ø2.4<br>L=37.8      | MGUIDE Narrow Sleeve pilot drill for 13mm length implant                                   | Stainless steel |       |
|    | CG-N0833 | Ø2.3-3.25<br>L=32.8 | MGUIDE Narrow Sleeve drill for Ø3.30 x 8mm length implant, conical connection              | Stainless steel |       |
|   | CG-N1033 | Ø2.3-3.25<br>L=34.8 | MGUIDE Narrow Sleeve drill for Ø3.30 x 10mm length implant, conical connection             | Stainless steel |     |
|  | CG-N1133 | Ø2.3-3.25<br>L=36.3 | MGUIDE Narrow Sleeve drill for Ø3.30 x 11.5mm length implant, conical connection           | Stainless steel |   |
|  | CG-N1333 | Ø2.3-3.25<br>L=37.8 | MGUIDE Narrow Sleeve drill for Ø3.30 x 13mm length implant, conical connection             | Stainless steel |   |
|  | CG-N0839 | Ø2.9-3.85<br>L=32.8 | MGUIDE Narrow Sleeve drill for Ø3.75 and Ø3.90 x 8mm length implant, conical connection    | Stainless steel |   |
|  | CG-N1039 | Ø2.9-3.85<br>L=34.8 | MGUIDE Narrow Sleeve drill for Ø3.75 and Ø3.90 x 10mm length implant, conical connection   | Stainless steel |   |
|  | CG-N1139 | Ø2.9-3.85<br>L=36.3 | MGUIDE Narrow Sleeve drill for Ø3.75 and Ø3.90 x 11.5mm length implant, conical connection | Stainless steel |   |
|  | CG-N1339 | Ø2.9-3.85<br>L=37.8 | MGUIDE Narrow Sleeve drill for Ø3.75 and Ø3.90 x 13mm length implant, conical connection   | Stainless steel |   |
|  | MG-NLG40 | Ø4<br>L=55          | MGUIDE Narrow Sleeve drill length gauge                                                    | Titanium        |   |

L = Length

# MGUIDE Surgical Set, Narrow Sleeve

|                                                                                     | Cat. No. | Dimensions (mm) | Description                                                                             | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|-----------------|-----------------------------------------------------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | CG-NMN10 | Ø4<br>L=31      | MGUIDE Narrow Sleeve C1 motor insertion tool, conical connection, narrow platform       | Stainless steel |       |
|    | CG-NMS10 | Ø4<br>L=32      | MGUIDE Narrow Sleeve V3/C1 motor insertion tool, conical connection, standard platform  | Stainless steel |       |
|    | CG-NRN10 | Ø4<br>L=26      | MGUIDE Narrow Sleeve C1 ratchet insertion tool, conical connection, narrow platform     | Stainless steel |       |
|    | CG-NRN01 | Ø4<br>L=30.5    | MGUIDE Narrow Sleeve C1 direct insertion tool, conical connection, narrow platform      | Stainless steel |       |
|   | CG-NRS10 | Ø4<br>L=27      | MGUIDE Narrow Sleeve C1 ratchet insertion tool, conical connection, standard platform   | Stainless steel |     |
|  | CG-NRS01 | Ø4<br>L=31      | MGUIDE Narrow Sleeve V3/C1 direct insertion tool, conical connection, standard platform | Stainless steel |   |
|  | VG-NMN10 | Ø4<br>L=31      | MGUIDE Narrow Sleeve V3 motor insertion tool, conical connection, narrow platform       | Stainless steel |   |
|  | VG-NRN10 | Ø4<br>L=25.5    | MGUIDE Narrow Sleeve V3 ratchet insertion tool, conical connection, narrow platform     | Stainless steel |   |
|  | VG-NRN01 | Ø4<br>L=30.5    | MGUIDE Narrow Sleeve V3 direct insertion tool, conical connection, narrow platform      | Stainless steel |   |
|  | VG-NRS10 | Ø4<br>L=26      | MGUIDE Narrow Sleeve V3 ratchet insertion tool, conical connection, standard platform   | Stainless steel |   |
|  | CG-NTAS0 | Ø4<br>L=24      | MGUIDE Narrow Sleeve template anchoring screw, conical connection, standard platform    | Titanium        |   |

L = Length



MGUIDE

## MGUIDE 16mm Drills Kit

Please note, products which are part of the MGUIDE drills kit may not be registered or available in every country/region.

MG-KD07  
MGUIDE 16mm Drills Kit



| Cat. No.                                                                                     | Dimensions (mm) | Description                                | Material        | Sold                                                                                                                                                                        |
|----------------------------------------------------------------------------------------------|-----------------|--------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  MG-D0024 | Ø2.2<br>L=16    | MGUIDE marking drill                       | Stainless steel |   |
|  MG-D0624 | Ø2.4<br>L=30.8  | MGUIDE pilot drill for 6mm length implant  | Stainless steel |   |
|  MG-D1624 | Ø2-2.4<br>L=16  | MGUIDE pilot drill for 16mm length implant | Stainless steel |   |


L = Length

MGUIDE 16mm  
Drills Kit

|                                                                                     | Cat. No. | Dimensions (mm)  | Description                                                                   | Material        | Sold                                                                                                                                                                        |
|-------------------------------------------------------------------------------------|----------|------------------|-------------------------------------------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | CG-D1633 | Ø2.3-3.3<br>L=16 | MGUIDE drill for Ø3.30, 16mm length implant, conical connection               | Stainless steel |       |
|    | CG-D1639 | Ø2.9-3.9<br>L=16 | MGUIDE drill for Ø3.90, 16mm length implant, conical connection               | Stainless steel |       |
|    | CG-D1643 | Ø3.4-4.3<br>L=16 | MGUIDE drill for Ø4.30, 16mm length implant, conical connection               | Stainless steel |       |
|    | CG-D1650 | Ø3.9-5<br>L=16   | MGUIDE drill for Ø5, 16mm length implant, conical connection                  | Stainless steel |       |
|   | MG-N0024 | Ø2.2<br>L=16     | MGUIDE Narrow Sleeve marking drill                                            | Stainless steel |     |
|  | MG-N0624 | Ø2.4<br>L=30.8   | MGUIDE Narrow Sleeve pilot drill for 6mm length implant                       | Stainless steel |   |
|  | MG-N1624 | Ø2-2.4<br>L=16   | MGUIDE Narrow Sleeve pilot drill for 16mm length implant                      | Stainless steel |   |
|  | CG-N1633 | Ø2.3-3.3<br>L=16 | MGUIDE Narrow Sleeve drill for Ø3.30, 16mm length implant, conical connection | Stainless steel |   |
|  | CG-N1639 | Ø2.9-3.9<br>L=16 | MGUIDE Narrow Sleeve drill for Ø3.90, 16mm length implant, conical connection | Stainless steel |   |

L = Length





## REGENERATIVE SOLUTIONS

MIS offers a range of regenerative materials for dental use aimed at providing solutions for a wide variety of clinical needs. MIS regenerative products are available in varying volumes and sizes for simple and easy delivery directly to the surgical site. High quality engineering and manufacturing processes ensure safety of all MIS products.





Bone Graft Substitutes

4MATRIX+  
Bone Graft Putty

Please note, 4MATRIX+ products might not be registered or available in every country/region.



Cat. No.

Description



BS-4MXP1

4MATRIX+ bone graft putty, 1cc



BS-4MXP5

4MATRIX+ bone graft putty, 0.5cc





Bone Graft Substitutes

## BONDBONE Bone Graft Cement



Cat. No.

Description



BS-BB005

BONDBONE 0.5cc, 3 drivers



BS-BB010

BONDBONE 1cc, 3 drivers



Bone Graft Substitutes

4BONE BCH  
Bone Graft Particulate



Cat. No.      Description



BS-4BS01      4BONE BCH syringe 0.5cc, granule size: 0.5-1mm



BS-4BS1S      4BONE BCH tube 1cc, granule size: 0.5-1mm



BS-4BS25      4BONE BCH tube 2.5cc, granule size: 1-2mm



Barrier Membrane

4BONE RCM  
Resorbable Collagen  
Membrane



Cat. No. Description

|                                                                                     |          |                             |
|-------------------------------------------------------------------------------------|----------|-----------------------------|
|  | BR-C1525 | 4BONE RCM membrane, 15x25mm |
|  | BR-C2030 | 4BONE RCM membrane, 20x30mm |
|  | BR-C3040 | 4BONE RCM membrane, 30x40mm |






## ADDITIONAL INFORMATION

For your convenience, the following section provides an additional resource about products and processes that may require expanded explanations in the form of visuals, schematics, technical information or documentation.



Additional information


Pilot Drills,  
External Irrigation





Additional information

## Step Drills, External Irrigation





Guidelines for using MIS Restorative Components

| Location                                  | Anterior Maxilla                                                                                                         | Canine, Premolars and Molars                                                                                |                                                                                                                  |                                                                                                                                   |                                                                                                                | Premolars and Molars                                                                                               |
|-------------------------------------------|--------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------|
| Crown/<br>Implant<br>Inclination<br>Ratio | 15° or 25°                                                                                                               | Up to 3.5°                                                                                                  |                                                                                                                  |                                                                                                                                   | 15° or 25°                                                                                                     | Up to 4°                                                                                                           |
| Gingival<br>Height                        | Up to 3mm<br>buccal<br>Up to 4mm<br>palatal                                                                              | Up to 3mm<br>buccal<br>Up to 4mm<br>lingual/palatal                                                         | Adjustable,<br>based on<br>abutment<br>preparation.<br>Heights 1, 3mm                                            | Up to 2mm                                                                                                                         | Up to 3mm<br>buccal<br>Up to 4mm<br>lingual                                                                    | According to<br>gingival height<br>available in<br>heights of<br>1,2,3,4mm                                         |
| Catalog<br>Number                         | CS-A1510<br>CS-A2510                                                                                                     | CS-A0010                                                                                                    | CS-MAC10<br>CS-MC103                                                                                             | CS-GPC10<br>CS-GP010                                                                                                              | CS-AN151<br>CS-AN251                                                                                           | CS-CPK61<br>CS-CPK62<br>CS-CPK63<br>CS-CPK64                                                                       |
| Abutment<br>description                   | Esthetic<br>angulated<br>abutment<br> | Esthetic<br>abutment<br> | Standard post<br>abutment<br> | Screw-retained<br>gold plastic<br>abutment<br> | Angulated<br>abutment<br> | Transgingival<br>abutment<br> |
| Upper view                                |                                       |                          |                               | Bottom view<br>                               |                           |                               |

\*Recommendation only





Additional information

## Alloy Composition

|             | Material           | Titanium Alloy<br>Ti 6Al 4V ELI | Gold-a    | Gold-s   | Gold-e   |
|-------------|--------------------|---------------------------------|-----------|----------|----------|
|             | Application        | C, K                            | L, G, K   | L, G, K  | L, K     |
| Composition | Gold %             | -                               | 60        | 61       | 70       |
|             | Platinum %         | -                               | 19        | 13.5     | 8.50     |
|             | Silver %           | -                               | -         | 16.5     | 13.40    |
|             | Palladium %        | -                               | 20        | -        | -        |
|             | Copper %           | -                               | -         | 9        | 7.50     |
|             | Zinc %             |                                 | -         | Balance  | 0.05     |
|             | Iridium %          |                                 | Balance   | -        | 0.10     |
|             | Titanium (max)     | Balance                         | -         | -        | -        |
|             | Carbon (max)       | 0.08                            | -         | -        | -        |
|             | Iron (max)         | 0.25                            | -         | -        | -        |
|             | Oxygen (max)       | 0.13                            | -         | -        | -        |
|             | Nitrogen (max)     | 0.05                            | -         | -        | -        |
|             | Hydrogen (max)     | 0.013                           | -         | -        | -        |
|             | Aluminum           | 5.5-6.5                         | -         | -        | -        |
|             | Vanadium           | 3.5-4.5                         | -         | -        | -        |
|             | Melting range °C   | 1604-1660                       | 1400-1490 | 950-1050 | 895-1010 |
| Hardness    | Soft annealing HV5 | -                               | 145       | 160      | 170      |
|             | Hardened HV5       | 350-385                         | 240       | 250      | 295      |
|             | Soft °C / min      | -                               | 100/30    | 750/5    | 750/5    |
|             | Hardened °C / min  | -                               | 700/30    | 350/15   | 450/15   |

Legend: C-bonding K-resin L-soldering G-casting onto





Cat. No. Page

B

|          |     |
|----------|-----|
| BR-C1525 | 117 |
| BR-C2030 | 117 |
| BR-C3040 | 117 |
| BS-4BS01 | 116 |
| BS-4BS1S | 116 |
| BS-4BS25 | 116 |
| BS-4MXP1 | 114 |
| BS-4MXP5 | 114 |
| BS-BB005 | 115 |
| BS-BB010 | 115 |

C

|           |    |
|-----------|----|
| C1-08375  | 17 |
| C1-08420  | 19 |
| C1-08500  | 21 |
| C1-10330  | 15 |
| C1-10375  | 17 |
| C1-10420  | 19 |
| C1-10500  | 21 |
| C1-11330  | 15 |
| C1-11375  | 17 |
| C1-11420  | 19 |
| C1-11500  | 21 |
| C1-13330  | 15 |
| C1-13375  | 17 |
| C1-13420  | 19 |
| C1-13500  | 21 |
| C1-16330  | 15 |
| C1-16375  | 17 |
| C1-16420  | 19 |
| C1-16500  | 21 |
| C1B-08375 | 23 |
| C1B-08420 | 23 |
| C1B-08500 | 23 |
| C1B-10330 | 23 |
| C1B-10375 | 23 |
| C1B-10420 | 23 |
| C1B-10500 | 23 |
| C1B-11330 | 23 |
| C1B-11375 | 23 |
| C1B-11420 | 23 |
| C1B-11500 | 23 |

Cat. No. Page

|           |     |
|-----------|-----|
| C1B-13330 | 23  |
| C1B-13375 | 23  |
| C1B-13420 | 23  |
| C1B-13500 | 23  |
| C1B-16330 | 23  |
| C1B-16375 | 23  |
| C1B-16420 | 23  |
| C1B-16500 | 23  |
| CG-D0833  | 102 |
| CG-D0839  | 102 |
| CG-D0843  | 102 |
| CG-D0850  | 103 |
| CG-D1033  | 102 |
| CG-D1039  | 102 |
| CG-D1043  | 102 |
| CG-D1050  | 103 |
| CG-D1133  | 102 |
| CG-D1139  | 102 |
| CG-D1143  | 103 |
| CG-D1150  | 103 |
| CG-D1333  | 102 |
| CG-D1339  | 102 |
| CG-D1343  | 103 |
| CG-D1350  | 103 |
| CG-D1633  | 110 |
| CG-D1639  | 110 |
| CG-D1643  | 110 |
| CG-D1650  | 110 |
| CG-GMN10  | 103 |
| CG-GMS10  | 104 |
| CG-GMW10  | 104 |
| CG-GRN01  | 103 |
| CG-GRN10  | 103 |
| CG-GRS01  | 104 |
| CG-GRS10  | 104 |
| CG-GRW01  | 104 |
| CG-GRW10  | 104 |
| CG-N0833  | 107 |
| CG-N0839  | 107 |
| CG-N1033  | 107 |
| CG-N1039  | 107 |
| CG-N1133  | 107 |
| CG-N1139  | 107 |
| CG-N1333  | 107 |
| CG-N1339  | 107 |
| CG-N1633  | 110 |
| CG-N1639  | 110 |
| CG-NMN10  | 108 |
| CG-NMS10  | 108 |
| CG-NRN01  | 108 |
| CG-NRN10  | 108 |
| CG-NRS01  | 108 |
| CG-NRS10  | 108 |
| CG-NTAS0  | 108 |
| CG-TAS55  | 104 |

Cat. No. Page

|          |    |
|----------|----|
| CK-CPK41 | 86 |
| CK-CPK42 | 86 |
| CK-CPK43 | 86 |
| CK-CPK44 | 87 |
| CK-CPK61 | 87 |
| CK-CPK62 | 87 |
| CK-CPK63 | 87 |
| CK-CPK64 | 87 |
| CK-CPK81 | 87 |
| CK-CPK82 | 87 |
| CK-CPK83 | 87 |
| CK-CPK84 | 87 |
| CK-NLC1  | 90 |
| CK-NLC2  | 90 |
| CK-NLC3  | 90 |
| CK-NLC4  | 90 |
| CK-NLC5  | 90 |
| CK-NOE1  | 92 |
| CK-NOE2  | 92 |
| CK-NOE3  | 92 |
| CK-NOE4  | 92 |
| CK-NOE5  | 92 |
| CK-NPK41 | 86 |
| CK-NPK42 | 86 |
| CK-NPK43 | 86 |
| CK-NPK44 | 87 |
| CK-NPK61 | 87 |
| CK-NPK62 | 87 |
| CK-NPK63 | 87 |
| CK-NPK64 | 87 |
| CK-SLC1  | 90 |
| CK-SLC2  | 90 |
| CK-SLC3  | 90 |
| CK-SLC4  | 90 |
| CK-SLC5  | 90 |
| CK-SOE1  | 92 |
| CK-SOE2  | 92 |
| CK-SOE3  | 92 |
| CK-SOE4  | 92 |
| CK-SOE5  | 92 |
| CK-WLC1  | 90 |
| CK-WLC3  | 90 |
| CK-WLC5  | 90 |
| CK-WOE1  | 92 |
| CK-WOE3  | 92 |
| CK-WPK61 | 87 |
| CK-WPK62 | 87 |
| CK-WPK63 | 87 |
| CK-WPK64 | 87 |
| CM-N1480 | 73 |
| CM-N2480 | 73 |
| CM-N3480 | 73 |
| CM-N5480 | 73 |
| CM-S1480 | 73 |
| CM-S2480 | 73 |

| Cat. No. | Page | Cat. No. | Page | Cat. No. | Page |
|----------|------|----------|------|----------|------|
| CM-S3480 | 73   | CN-LOC03 | 89   | CS-CPK64 | 88   |
| CM-S4480 | 73   | CN-LOC04 | 89   | CS-CPK81 | 88   |
| CM-S5480 | 73   | CN-LOC05 | 89   | CS-CPK82 | 88   |
| CM-SO171 | 74   | CN-MAB10 | 82   | CS-CPK83 | 88   |
| CM-SO172 | 74   | CN-MAC10 | 85   | CS-CPK84 | 88   |
| CM-SO301 | 74   | CN-MAC30 | 85   | CS-EZB15 | 77   |
| CM-SO302 | 74   | CN-MAN10 | 81   | CS-GP010 | 76   |
| CM-W1480 | 73   | CN-PF330 | 65   | CS-GPC10 | 76   |
| CM-W3480 | 73   | CN-RSM10 | 65   | CS-I0375 | 64   |
| CM-WO171 | 74   | CN-SP102 | 80   | CS-IB091 | 82   |
| CM-WO172 | 74   | CN-TB001 | 78   | CS-IB121 | 82   |
| CM-WO301 | 74   | CN-TB002 | 79   | CS-IC481 | 64   |
| CM-WO302 | 74   | CN-TB015 | 78   | CS-IC485 | 65   |
| CN-AGT12 | 82   | CN-TB030 | 78   | CS-IN060 | 78   |
| CN-AN001 | 83   | CN-TBF15 | 79   | CS-IN160 | 79   |
| CN-AN003 | 83   | CN-TBF30 | 79   | CS-IN360 | 79   |
| CN-AN010 | 83   | CN-TC013 | 66   | CS-IO391 | 64   |
| CN-AN020 | 84   | CN-TC023 | 66   | CS-IO395 | 64   |
| CN-AN310 | 84   | CN-TC033 | 67   | CS-IO485 | 64   |
| CN-AN320 | 84   | CN-TCI23 | 67   | CS-IT300 | 65   |
| CN-AT002 | 83   | CN-TCI33 | 67   | CS-LEZ15 | 77   |
| CN-AT102 | 83   | CN-TPA10 | 67   | CS-LOC01 | 89   |
| CN-AT202 | 84   | CN-TPH50 | 67   | CS-LOC02 | 89   |
| CN-B1330 | 94   | CS-A0010 | 83   | CS-LOC03 | 89   |
| CN-B2330 | 94   | CS-A0030 | 83   | CS-LOC04 | 89   |
| CN-B3330 | 94   | CS-A1510 | 84   | CS-LOC05 | 89   |
| CN-B4330 | 94   | CS-A1530 | 84   | CS-MAB10 | 82   |
| CN-C2040 | 68   | CS-A2510 | 84   | CS-MAC10 | 85   |
| CN-C3040 | 68   | CS-A2530 | 84   | CS-MAN10 | 81   |
| CN-CPK41 | 87   | CS-AGT12 | 82   | CS-MC103 | 85   |
| CN-CPK42 | 87   | CS-AN151 | 84   | CS-PF375 | 65   |
| CN-CPK43 | 87   | CS-AN251 | 84   | CS-RSM10 | 65   |
| CN-CPK44 | 88   | CS-AT002 | 83   | CS-SP102 | 80   |
| CN-CPK61 | 88   | CS-AT152 | 84   | CS-TB001 | 78   |
| CN-CPK62 | 88   | CS-AT252 | 84   | CS-TB002 | 79   |
| CN-CPK63 | 88   | CS-B1375 | 94   | CS-TB015 | 78   |
| CN-CPK64 | 88   | CS-B2375 | 94   | CS-TB030 | 78   |
| CN-EZB15 | 77   | CS-B3375 | 94   | CS-TBF15 | 79   |
| CN-GP010 | 76   | CS-B4375 | 94   | CS-TBF30 | 79   |
| CN-GPC10 | 76   | CS-B5375 | 94   | CS-TC013 | 66   |
| CN-I0330 | 64   | CS-C1540 | 68   | CS-TC023 | 66   |
| CN-IB091 | 82   | CS-C1557 | 68   | CS-TC033 | 67   |
| CN-IB121 | 82   | CS-C2040 | 68   | CS-TCI23 | 67   |
| CN-IC401 | 64   | CS-C2057 | 68   | CS-TCI33 | 67   |
| CN-IC405 | 65   | CS-C3040 | 68   | CS-TPA10 | 67   |
| CN-IN060 | 78   | CS-C3057 | 68   | CS-TPH50 | 67   |
| CN-IN160 | 79   | CS-C4040 | 68   | CT-BGP32 | 47   |
| CN-IN360 | 79   | CS-C4057 | 68   | CT-BGP35 | 47   |
| CN-IO331 | 64   | CS-CPK41 | 87   | CT-BGP46 | 47   |
| CN-IO335 | 64   | CS-CPK42 | 87   | CT-BP320 | 47   |
| CN-IO405 | 64   | CS-CPK43 | 87   | CT-BP355 | 47   |
| CN-IT300 | 65   | CS-CPK44 | 88   | CT-BP460 | 47   |
| CN-LEZ15 | 77   | CS-CPK61 | 88   | CT-BTC24 | 42   |
| CN-LOC01 | 89   | CS-CPK62 | 88   | CT-BTC30 | 42   |
| CN-LOC02 | 89   | CS-CPK63 | 88   | CT-BTC35 | 42   |

| Cat. No. | Page       |
|----------|------------|
| CT-BTC40 | 42         |
| CT-NLI10 | 34, 43     |
| CT-NLM30 | 32, 43     |
| CT-NLR10 | 34, 43     |
| CT-NLR30 | 34, 43     |
| CT-NSI10 | 34, 44     |
| CT-NSM30 | 32, 43     |
| CT-NSR10 | 34, 44     |
| CT-NSR30 | 32, 43     |
| CT-P2408 | 29, 41, 46 |
| CT-P2410 | 29, 42, 46 |
| CT-P2411 | 29, 42, 46 |
| CT-P2413 | 29, 42, 46 |
| CT-P2416 | 29, 42, 46 |
| CT-SLI10 | 35, 44     |
| CT-SLM30 | 32, 43     |
| CT-SLR10 | 34, 44     |
| CT-SLR30 | 33, 43     |
| CT-SSI10 | 35, 44     |
| CT-SSM30 | 33, 43     |
| CT-SSR10 | 34, 44     |
| CT-SSR30 | 33, 43     |
| CT-TDC30 | 30, 42, 46 |
| CT-TDC35 | 30, 42, 46 |
| CT-TDC40 | 30, 42, 46 |
| CT-TDC45 | 30         |
| CT-UN001 | 55         |
| CT-US001 | 55         |
| CT-UW001 | 55         |
| CT-WLI10 | 35, 44     |
| CT-WLR10 | 35, 44     |
| CT-WLR10 | 44         |
| CT-WSI10 | 35, 44     |
| CT-WSR10 | 35, 44     |
| CW-AGT12 | 82         |
| CW-AN151 | 84         |
| CW-AT001 | 83         |
| CW-AT002 | 83         |
| CW-AT003 | 83         |
| CW-AT152 | 84         |
| CW-AT153 | 84         |
| CW-AT251 | 84         |
| CW-AT252 | 84         |
| CW-AT253 | 84         |
| CW-B1500 | 94         |
| CW-B3500 | 94         |
| CW-C1540 | 68         |
| CW-C1557 | 68         |
| CW-C2040 | 68         |
| CW-C2057 | 68         |
| CW-C3040 | 68         |
| CW-C3057 | 68         |
| CW-C4040 | 68         |
| CW-C4057 | 68         |
| CW-CPK61 | 88         |

| Cat. No. | Page |
|----------|------|
| CW-CPK62 | 88   |
| CW-CPK63 | 88   |
| CW-CPK64 | 88   |
| CW-EZB15 | 77   |
| CW-GP010 | 76   |
| CW-GPC10 | 76   |
| CW-IO500 | 64   |
| CW-IB091 | 82   |
| CW-IB121 | 82   |
| CW-IC551 | 64   |
| CW-IC555 | 65   |
| CW-IO501 | 64   |
| CW-IO505 | 64   |
| CW-IO555 | 64   |
| CW-IT300 | 65   |
| CW-LEZ15 | 77   |
| CW-LOC01 | 89   |
| CW-LOC03 | 89   |
| CW-LOC05 | 89   |
| CW-MAB10 | 82   |
| CW-MAN10 | 81   |
| CW-MC103 | 85   |
| CW-PF500 | 65   |
| CW-RSM10 | 65   |
| CW-SMAC1 | 85   |
| CW-SP102 | 80   |
| CW-TB001 | 78   |
| CW-TB002 | 79   |
| CW-TB015 | 78   |
| CW-TB030 | 78   |
| CW-TBF15 | 79   |
| CW-TBF30 | 79   |
| CW-TC013 | 66   |
| CW-TC023 | 66   |
| CW-TC033 | 67   |
| CW-TCI23 | 67   |
| CW-TPA10 | 67   |

|          |            |
|----------|------------|
| E        |            |
| ET-IT002 | 40, 56, 93 |

M

|           |                        |
|-----------|------------------------|
| MB-DB235  | 95                     |
| MB-DMH10  | 95                     |
| MB-GPS10  | 95                     |
| MB-PPC10  | 95                     |
| MB-PPR10  | 95                     |
| MB-PR010  | 95                     |
| MB-RS375  | 94                     |
| MB-SF200  | 95                     |
| MB-TBS10  | 95                     |
| MB-TBS20  | 95                     |
| MB-TTP10  | 95                     |
| MD-E0220  | 77                     |
| MD-G0213  | 65                     |
| MD-G0220  | 76, 85                 |
| MD-S0200  | 84, 65, 67, 76, 79, 85 |
| MG2-06160 | 25                     |
| MG2-08160 | 25                     |
| MG2-10160 | 25                     |
| MG2-12160 | 25                     |
| MG-BM047  | 101                    |
| MG-BM055  | 101                    |
| MG-D0024  | 109                    |
| MG-D0624  | 101, 109               |
| MG-D0824  | 101                    |
| MG-D1024  | 101                    |
| MG-D1124  | 102                    |
| MG-D1324  | 102                    |
| MG-D1624  | 109                    |
| MG-DFP20  | 104                    |
| MG-DLG55  | 103                    |
| MG-FP020  | 104                    |
| MG-N0024  | 110                    |
| MG-N0624  | 106, 110               |
| MG-N0824  | 106                    |
| MG-N1024  | 107                    |
| MG-N1124  | 107                    |
| MG-N1324  | 107                    |
| MG-N1624  | 110                    |
| MG-NBM40  | 106                    |
| MG-NLG40  | 107                    |
| MG-NTP40  | 106                    |
| MG-TMP10  | 100                    |
| MG-TMP20  | 100                    |
| MG-TP055  | 101                    |
| MK-0022   | 58                     |
| MK-0023   | 52                     |
| MK-0024   | 53                     |
| MK-BC101  | 48                     |

[illegible]

|          |    |
|----------|----|
| OE-IC001 | 93 |
| OE-IC002 | 93 |

Cat. No.      Page

|          |    |
|----------|----|
| OE-RCB01 | 93 |
| OE-RCP01 | 93 |
| OE-RCV01 | 93 |
| OE-RCW01 | 93 |
| OE-RCY01 | 93 |
| OE-RSM60 | 93 |

V

|           |     |
|-----------|-----|
| V3-08390  | 9   |
| V3-08430  | 11  |
| V3-08500  | 13  |
| V3-10330  | 7   |
| V3-10390  | 9   |
| V3-10430  | 11  |
| V3-10500  | 13  |
| V3-11330  | 7   |
| V3-11390  | 9   |
| V3-11430  | 11  |
| V3-11500  | 13  |
| V3-13330  | 7   |
| V3-13390  | 9   |
| V3-13430  | 11  |
| V3-13500  | 13  |
| V3-16330  | 7   |
| V3-16390  | 9   |
| V3-16430  | 11  |
| V3-16500  | 13  |
| V3B-08390 | 23  |
| V3B-08430 | 23  |
| V3B-08500 | 23  |
| V3B-10330 | 23  |
| V3B-10390 | 23  |
| V3B-10430 | 23  |
| V3B-10500 | 23  |
| V3B-11330 | 23  |
| V3B-11390 | 23  |
| V3B-11430 | 23  |
| V3B-11500 | 23  |
| V3B-13330 | 23  |
| V3B-13390 | 23  |
| V3B-13430 | 23  |
| V3B-13500 | 23  |
| V3B-16330 | 23  |
| V3B-16390 | 23  |
| V3B-16430 | 23  |
| V3B-16500 | 23  |
| VG-GMN10  | 103 |

Cat. No.      Page

|          |     |
|----------|-----|
| VG-GRN01 | 103 |
| VG-GRN10 | 104 |
| VG-GRS10 | 104 |
| VG-NMN10 | 108 |
| VG-NRN01 | 108 |
| VG-NRN10 | 108 |
| VG-NRS10 | 108 |
| VK-CPK61 | 87  |
| VK-CPK62 | 87  |
| VK-CPK63 | 87  |
| VK-CPK64 | 87  |
| VK-CPK81 | 87  |
| VK-CPK82 | 87  |
| VK-CPK83 | 87  |
| VK-CPK84 | 87  |
| VK-NPK41 | 86  |
| VK-NPK42 | 86  |
| VK-NPK43 | 86  |
| VK-NPK44 | 87  |
| VK-NPK61 | 87  |
| VK-NPK62 | 87  |
| VK-NPK63 | 87  |
| VK-NPK64 | 87  |
| VN-AGT12 | 82  |
| VN-AT001 | 83  |
| VN-AT002 | 83  |
| VN-AT003 | 83  |
| VN-AT101 | 83  |
| VN-AT102 | 83  |
| VN-AT103 | 84  |
| VN-AT201 | 84  |
| VN-AT202 | 84  |
| VN-AT203 | 84  |
| VN-EZB15 | 77  |
| VN-GPC10 | 76  |
| VN-IB091 | 82  |
| VN-IB121 | 82  |
| VN-IC400 | 65  |
| VN-IC401 | 64  |
| VN-IC405 | 65  |
| VN-IN060 | 78  |
| VN-IN160 | 79  |
| VN-IN360 | 79  |
| VN-IO331 | 64  |
| VN-IO335 | 64  |
| VN-IO401 | 64  |
| VN-IO405 | 64  |
| VN-LEZ15 | 77  |
| VN-MAB10 | 82  |
| VN-MAN10 | 81  |
| VN-MC101 | 85  |
| VN-MC103 | 85  |
| VN-P4140 | 87  |
| VN-P4240 | 87  |
| VN-P4340 | 87  |

Cat. No.      Page

|          |    |
|----------|----|
| VN-P4440 | 88 |
| VN-P6140 | 88 |
| VN-P6240 | 88 |
| VN-P6340 | 88 |
| VN-P6440 | 88 |
| VN-PF340 | 65 |
| VN-RSM10 | 65 |
| VN-SP102 | 80 |
| VN-TB001 | 78 |
| VN-TB015 | 78 |
| VN-TB030 | 78 |
| VN-TCI23 | 67 |
| VN-TCI33 | 67 |
| VN-TPA10 | 67 |
| VN-TPH50 | 67 |
| VS-AT001 | 83 |
| VS-AT002 | 83 |
| VS-AT003 | 83 |
| VS-AT152 | 84 |
| VS-AT153 | 84 |
| VS-AT252 | 84 |
| VS-AT253 | 84 |
| VS-MC101 | 85 |
| VS-MC103 | 85 |
| VS-P6155 | 88 |
| VS-P6255 | 88 |
| VS-P6355 | 88 |
| VS-P6455 | 88 |
| VS-P8155 | 88 |
| VS-P8255 | 88 |
| VS-P8355 | 88 |
| VS-P8455 | 88 |
| VS-TC013 | 66 |
| VS-TC023 | 66 |
| VS-TC033 | 67 |
| VS-TCI23 | 67 |
| VS-TCI33 | 67 |











IT'S SIMPLE  
WHEN YOU  
HAVE ALL THE  
OPTIONS



**MIS<sup>®</sup> | C1**

**THE CONNECTION FOR PREDICTABLE BIOLOGY. MAKE IT SIMPLE**

The MIS C1 conical connection implant system offers uncompromising accuracy, with high initial and biological stability and a safe, yet simple, procedure. A consistent concave emergence profile of the C1 prosthetic components improves soft tissue esthetic results. Learn more about the C1 implant and MIS at: [www.mis-implants.com](http://www.mis-implants.com)

