

# NATIONAL PROMOTION POLAND

SPRING/SUMMER 2021

Valid from 15<sup>th</sup> of May  
until 31<sup>st</sup> of July 2021

**MEDIF**  
DENTISTRY

# Introducing...

# HuFriedyGroup

The Best In Practice

## Who We Are

We are a combination of companies that have come together under one roof to deliver the products, services, and solutions you need to be the best in practice. With our passion for infection prevention solutions, premium instrumentation, and service and support programs, HuFriedyGroup will help you perfect your craft and continue to deliver the highest level of care to each patient.

## What Can You Expect From HuFriedyGroup?

HuFriedyGroup is here to provide the best dental experience to the industry. Our unique combination of products, services, and support enable dental professionals like you to reduce risk, improve efficiencies, drive compliance, and enhance patient and staff safety.

### World-Class Products

Our unwavering passion and unrivaled expertise in the dental industry is evident in the instruments and infection prevention products you trust each and every day.

### Compliance & Protection

We provide solutions that help you to ensure a safe environment for patients and staff.

### Service & Support

Through our community platforms, customer service and educational programs, we provide you with the information and tools you need to perfect your craft.

## Proud members of HuFriedyGroup:



# DIAGNOSTIC INSTRUMENTS HU-FRIEDY®



## EXD5

Explorer #5 hdl #31

Sugg. Retail zł 135

Promo zł 88



## EXD3CH

Explorer #3 cow horn

Sugg. Retail zł 135

Promo zł 88



## EXD11/126

Explorer Dominion #11/12  
hdl #6

Sugg. Retail zł 142

Promo zł 92



## PCPUNC157

Probe UNC #UNC15 hdl  
#7 Qulix 1-2-3-...15

Sugg. Retail zł 198

Promo zł 129



## PCP11

Probe #11 hdl #30 Qulix  
3-6-8-11

Sugg. Retail zł 168

Promo zł 109



## PQ2N6

Bifurcation Probe  
Nabers #2N hdl #6

Sugg. Retail zł 299

Promo zł 194



## DPU17

Dressing Pliers Universal  
#17

Sugg. Retail zł 97

Promo zł 58



## DP2

Dressing Pliers College  
#2

Sugg. Retail zł 117

Promo zł 76



## DP18

Dressing Pliers #18

Sugg. Retail zł 128

Promo zł 77



## MIR5HD12E

High Definition Mirror #5  
single-sided, European,  
12 pcs./pack

Sugg. Retail zł 467

Promo zł 280



## MIR4HD12E

High Definition Mirror #4  
single-sided, European,  
12 pcs./pack

Sugg. Retail zł 467

Promo zł 280



## MHE6

Mirror Handle hdl #6  
European Thread

Sugg. Retail zł 120

Promo zł 72



# NO SHARPENING... DO YOU REALLY BELIEVE IT?



## SG5/695E2

Curette Gracey #5/6 hdl #9 oral/labial, yellow

**Sugg. Retail**

zł 292

**Promo**

zł 190

Valid for all #9 HANDLE



## EVEREDGE™ 2.0 INSTRUMENTS SCALERS AND CURETTES

SCALERS EVEREDGE™ WERE SHARP..

SCALERS EVEREDGE™ 2.0 ARE EVEN  
MORE SHARPENED.

Using an enhanced finishing process, Hu-Friedy's professional artisans are creating edges that are scientifically measured to be over 60% sharper than the next closest competitor. Less force means easier, more comfortable scaling for the clinician and patient.

\*Products tested: Hu-Friedy EverEdge™ 2.0 H7, American Eagle XP H6-7, PDT H6-7, LM EM/EN H6-7, and G M Hartzell H6-7. Hu-Friedy cooperated with CATRA, an independent specialized blade laboratory, to conduct this testing. For more information, contact David Platt, Hu-Friedy Mfg. Co., LLC, 3232 N. Rockwell St, Chicago, IL 60618.



For more information  
on **EVEREDGE™ 2.0** Instruments,  
visit our site [HuFriedy.eu](https://www.HuFriedy.eu)







## IMP204SDT

Implant Scaler Titanium  
#204SD

Sugg. Retail zł 344

Promo zł 224



## IMPLG1/2T

Implant Scaler Titanium  
#Langer 1/2

Sugg. Retail zł 344

Promo zł 224



## IMPM1/2T

Implant Scaler Titanium  
#Mini-Five 1/2

Sugg. Retail zł 344

Promo zł 224



## IMPM13/14T

Implant Scaler Titanium  
#Mini-Five 13/14

Sugg. Retail zł 344

Promo zł 224



## IMPM11/12T

Implant Scaler Titanium  
#Mini-Five 11/12

Sugg. Retail zł 344

Promo zł 224



## IMEDIMPLM

Titanium Implant Maintenance Kit  
(IMS™ Infinity Series™ Cassette 1/4 DIN, purple + 5 Instruments)

Sugg. Retail zł 2.053

Promo zł 1.335

# RESTORATIVE INSTRUMENTS HU-FRIEDY®



## **EXC17H**

Excavator #17 hdl #41 heavy

Sugg. Retail zł 198

**Promo zł 119**



## **EXC18H**

Excavator #18H hdl #41 heavy

Sugg. Retail zł 198

**Promo zł 119**



## **EXC19H**

Excavator #19 hdl #41 heavy

Sugg. Retail zł 198

**Promo zł 119**



## **CIGFT16**

Plastic Filling Instrument  
Goldstein #1, hdl #6 Flexi-Thin

Sugg. Retail zł 229

**Promo zł 137**



## **CIGFT26**

Plastic Filling Instrument  
Goldstein #2, hdl #6 Flexi-Thin

Sugg. Retail zł 229

**Promo zł 137**



## **CIGFT36**

Plastic Filling Instrument  
Goldstein #3, hdl #6 Flexi-Thin

Sugg. Retail zł 229

**Promo zł 137**



## **CIGFT46**

Plastic Filling Instrument  
Goldstein #4, hdl #6 Flexi-Thin

Sugg. Retail zł 229

**Promo zł 137**



## **CIGFT5**

Plastic Filling Instrument  
Goldstein #5, hdl #41 Flexi-Thin

Sugg. Retail zł 198

**Promo zł 119**



## **CIGFT6**

Plastic Filling Instrument  
Goldstein #6, hdl #41 Flexi-Thin

Sugg. Retail zł 198

**Promo zł 119**



## **CIGFTMIN16**

Plastic Filling Instrument  
Goldstein #1, hdl #6  
mini Flexi-Thin

Sugg. Retail zł 229

**Promo zł 137**



## **CIGFTMIN36**

Plastic Filling Instrument  
Goldstein #3, hdl #6  
mini Flexi-Thin

Sugg. Retail zł 229

**Promo zł 137**



## **CIGFTMIN46**

Plastic Filling Instrument  
Goldstein #4, hdl #6  
mini Flexi-Thin

Sugg. Retail zł 229

**Promo zł 137**


# RESTORATIVE INSTRUMENTS HU-FRIEDY®

## AKRO-FLEX™ COMPOSITE INSTRUMENT

### CREATE BEAUTIFUL RESTORATIONS


Akro-Flex™ Composite Instrument is extremely flexible which allows clinicians to effortlessly manipulate composite materials during aesthetic restorations.

**The inclusion of a ductile material allows Akro-Flex™ Composite Instrument to perform as a solid brush.**



#### HYPER-THIN PROFILE

Ultra thin working ends reach narrow interproximal spaces with ease. The flexible, versatile working ends allow for better visibility as compared to traditional composite instruments.



#### UNPARALLELED FLEXIBILITY

By incorporating Nickel Titanium, a material known for its super elasticity, Akro-Flex™ Composite Instrument acts as a solid brush. The resilient working ends are excellent when creating fine anatomical detail with delicate, artistic strokes.



#### ERGONOMIC HANDLE

The smooth, lightweight handle offers increased control due to the large diameter. It creates an ergonomically friendly option that provides maximum comfort and helps reduce hand fatigue. Reducing hand fatigue can increase the longevity of a clinician's career.

### AKRO-FLEX

Akro-Flex™  
Composite Instrument

**Sugg. Retail**

zł 385

**Promo**

zł 269



WATCH OUR  
**TUTORIAL VIDEO** FOR AKRO-FLEX™  
COMPOSITE INSTRUMENT ON YOUTUBE



Patent pending.



# RESTORATIVE INSTRUMENTS BY WERESTORE.IT<sup>†</sup>

## CALIBRA™ COMPOSITE INSTRUMENT HU-FRIEDY®

### FOR CALIBRATED RESTORATIONS

Calibra™ instrument is an **XTS (AlTiN)<sup>††</sup>** coated composite instrument which helps to provide the proper anatomical curvature to the dentin material layer of a restoration and leaves the ideal space for the final layer of enamel material. The **Calibra™ “IN”** point is used to define the dentinal skeleton, while the **Calibra™ “OUT”** point helps to reduce the chances of excess composite material being placed during the final layer of the restoration. Calibra Composite Instrument acts as a **guide**, helping to prevent over-contouring at the preparation margin while assisting in reducing excess composite thickness build up as well as **modeling** the final layer of composite material.

### TNCALIBRA

XTS PFI WeRestore.it  
Calibra™ Instrument

NEW

Sugg. Retail

Promo

zł 333

zł 233

HOW  
CALIBRA™ IN  
POINT WORKS



HOW  
CALIBRA™ OUT  
POINT WORKS



Calibra™ Composite Instrument consists of two terminals: “IN” which defines the dentinal skeleton. “OUT” which allows to verify the thickness and shape of the last layer without excess. This instrument is a caliper to check the thickness (a practical alternative to the sagittal silicone index) and a modelling tool for the last layer.

<sup>†</sup> WeRestore.it is not owned by or affiliated with Hu-Friedy. <sup>††</sup> Data on file, available upon request.



### KOL CORNER

Instruments made  
in collaboration with  
**Dr. Gaetano Paolone**  
and **Dr. Salvatore**  
**Scolavino**, Opinion  
leaders Hu-Friedy.

CALIBRA™  
OUT POINT

CALIBRA™  
IN POINT






**CALIBRA™**  
COMPOSITE  
INSTRUMENT BY  
WERESTORE.IT  
WATCH OUR  
PRACTICAL VIDEO

# RESTORATIVE INSTRUMENTS BY WERESTORE.IT

## HU-FRIEDY 3SSENTIAL KIT HU-FRIEDY®




**ANTERIOR**

**TNCIPAAX**  
XTS PFI Essential  
WeRestore.it Anterior

Sugg. Retail zł 262
<b>Promo zł 170</b>




**POSTERIOR**

**TNEXBKR1X**  
XTS PFI Essential  
WeRestore.it Posterior

Sugg. Retail zł 262
<b>Promo zł 170</b>

**SPATULA**

**CVKOT1X**  
XTS PFI Essential  
WeRestore.it Spatula

Sugg. Retail zł 262
<b>Promo zł 170</b>


## THE ESSENTIAL RESTORATIVE KIT INCLUDING "PLUG&SCULPT" INSTRUMENT

The Hu-Friedy 3SSENTIAL KIT by Werestore.it is the most straightforward kit of restorative ever made, with just 3 essential instruments at once: Anterior, Posterior and Spatula, from the Hu-Friedy Black Line Collection. The three instruments are synergic and interchangeable and the word "Essential" refers to the easiness of all restorative procedures focused on the essential steps to create high esthetic composite restorations while using this kit.

The revolutionary concept of "Plug&Sculpt" of the Posterior is very simple: one tip is used to plug the composite into the cavity and the other one can sculpt composite in additive and subtractive modelling techniques. The Anterior, instead, is called by the doctors "solid brush" since it models and spreads the composite very easily, just like a brush.

### RESTOREITDIN

Set of 3 instruments and cassette



#### Ce kit comprend:

IMS™ Infinity Series™ Cassette 1/4 DIN, 5 instruments, black	<b>IME14DIN57</b>
XTS PFI Essential WeRestore.it Spatula	<b>CVKOT1X</b>
XTS PFI Essential WeRestore.it Posterior	<b>TNEXBKR1X</b>
XTS PFI Essential WeRestore.it Anterior	<b>TNCIPAAX</b>

DESCRIPTION	CODE	Sugg. Retail	Promo
3SSENTIAL Kit	<b>RESTOREITDIN</b>	zł 1.280	<b>zł 832</b>

# LUXATING ELEVATORS HU-FRIEDY®



## **EL2S**

Luxating Elevator hdl #510,  
straight 2 mm

Sugg. Retail zł 519

**Promo zł 363**



## **EL2C**

Luxating Elevator hdl #510,  
curved 2 mm

Sugg. Retail zł 519

**Promo zł 363**



## **EL3S**

Luxating Elevator hdl #510,  
straight 3 mm

Sugg. Retail zł 519

**Promo zł 363**



## **EL3C**

Luxating Elevator hdl #510,  
curved 3 mm

Sugg. Retail zł 519

**Promo zł 363**



## **EL4S**

Luxating Elevator hdl #510,  
straight 4 mm

Sugg. Retail zł 519

**Promo zł 363**



## **EL4C**

Luxating Elevator hdl #510,  
curved 4 mm

Sugg. Retail zł 519

**Promo zł 363**



# SURGICAL CURETTES HU-FRIEDY®



## CL846

Surgical Curette Lucas #84,  
hdl #6 2,2 mm

Sugg. Retail zł 441

Promo zł 265



## CL856

Surgical Curette Lucas #85,  
hdl #6 2,8 mm

Sugg. Retail zł 441

Promo zł 265



## CL866

Surgical Curette Lucas #86,  
hdl #6 2,8 mm

Sugg. Retail zł 441

Promo zł 265



## CL876

Surgical Curette Lucas #87,  
hdl #6 3,5 mm

Sugg. Retail zł 441

Promo zł 265



## CM2/4S6

Surgical Curette Molt #2/4,  
short hdl #6

Sugg. Retail zł 444

Promo zł 267

**DID YOU KNOW?** All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like HU-FRIEDY® Instrument Lubricant Spray (ILS), will prevent rust, corrosion and stiff joints to ensure smooth operation. **All hinged instruments should be sterilised in the open position.**



# NEEDLE HOLDERS, SCISSORS & TISSUE PLIERS HU-FRIEDY®

## NEEDLE HOLDERS HU-FRIEDY®



**NHMH**  
Needle Holder Mayo-Hegar,  
grooved, 1-0 - 4-0 Needle, 16 cm

Sugg. Retail zł 251

**Promo zł 163**



**NHCW**  
Needle Holder Crile-Wood,  
grooved, 3-0 - 6-0, 15 cm

Sugg. Retail zł 229

**Promo zł 149**



**NH5038**  
Needle Holder Crile-Wood  
#5038, PermaSharp™, 15 cm

Sugg. Retail zł 762

**Promo zł 495**



**NH5024R**  
Needle Holder Castroviejo  
#5024R, round handle,  
PermaSharp™, straight, 18 cm

Sugg. Retail zł 2.213

**Promo zł 1.439**



**NH5024SIM**  
Needle Holder Simion #5024R,  
round handle, PermaSharp™,  
straight, 18 cm

Sugg. Retail zł 2.213

**Promo zł 1.439**

### SUPER CUT™

Cutting blade honed/serrated to „razor” sharpness with opposing blade finished with fine serrations to avoid tissue slippage.

### PERMA SHARP™

Featuring tungsten carbide inserts bonded to surgical-grade stainless steel on both cutting edges for smooth, controlled cutting and maximum durability.

## SCISSORSS HU-FRIEDY®



**SCORT**  
Scissors Cortellini curved 12 cm

Sugg. Retail zł 307

**Promo zł 200**



**S16**  
Scissors Goldman-Fox #16  
curved serrated 12,5 cm

Sugg. Retail zł 255

**Promo zł 166**



**S14**  
Scissors LaGrange #14 2-x  
angled 11,5 cm

Sugg. Retail zł 340

**Promo zł 221**



**S18SC**  
Scissors Iris #18 fine angled  
SuperCut™ 11,5 cm

Sugg. Retail zł 444

**Promo zł 289**

## TISSUE PLIERS HU-FRIEDY®



**TP41**  
Tissue Pliers Adson #41  
anatomical straight 12 cm

Sugg. Retail zł 154

**Promo zł 100**



**TP42**  
Tissue Pliers Adson #42  
surgical 1x2-straight 12 cm

Sugg. Retail zł 154

**Promo zł 100**

# PERIODONTAL KNIVES AND SINUS LIFT INSTRUMENTS HU-FRIEDY®

## PERIODONTAL KNIVES HU-FRIEDY®



**TKN1**  
Tunneling Knife #1 hdl #6

Sugg. Retail zł 538

**Promo zł 323**



**TKN2**  
Tunneling Knife #2 hdl #6

Sugg. Retail zł 538

**Promo zł 323**



**PFIWDS1MK**  
Tunneling knife hdl #6

Sugg. Retail zł 287

**Promo zł 172**



**IMP65786**  
Sinuslift Instrument  
Kramer-Nevins #LG/SM  
hdl #6, 2x angled

Sugg. Retail zł 553

**Promo zł 387**



**IMP6524S6**  
Sinuslift Instrument short hdl #6

Sugg. Retail zł 553

**Promo zł 387**



**IMP6577SC6**  
Sinuslift Instrument  
Kramer-Nevins #LG/SM  
hdl #6, strongly curved

Sugg. Retail zł 553

**Promo zł 387**



# ENDODONTICS HU-FRIEDY®

## RUBBER DAM PUNCH HU-FRIEDY®



**RDPDA**  
Rubber Dam Punch Ainsworth, 17 cm

Sugg. Retail zł 1.150

**Promo zł 690**

## RUBBER DAM FORCEPS HU-FRIEDY®



**RDF4E**  
Rubber Dam Clamp Forceps #4 European

Sugg. Retail zł 684

**Promo zł 410**



**EXDG16**  
Explorer Endo #DG16 hdl #31

Sugg. Retail zł 176

**Promo zł 106**

## RUBBER DAM CLAMPS HU-FRIEDY®



For more information  
on **RUBBER DAM CLAMPS**  
visit [HuFriedy.eu](https://www.HuFriedy.eu)



\*Valid for all the stainless steel  
Rubber Dam Clamps Line

**RDCM00**  
Rubber Dam Clamps #00 UK anterior

Sugg. Retail zł 90

**Promo zł 63\***

# IMS™ INSTRUMENT MANAGEMENT SYSTEM HU-FRIEDY®

## IMS™ INFINITY SERIES™ CASSETTES



Description	CODE	Sugg. Retail	Promo
1 IMS™ Infinity Series™ Cassette DIN, 18 instruments, Surgical, Dimensions: 290 x 184 x 34 mm	<b>IMEDIN18_S</b>	zł 1.064	<b>zł 692</b>
2 IMS™ Infinity Series™ Cassette 1/4 DIN, 5 instruments, Dimensions: 72 x 184 x 34 mm	<b>IME14DIN5_</b>	zł 654	<b>zł 425</b>
3 IMS™ Infinity Series™ Cassette 1/2 DIN, 10 instruments, Dimensions: 152 x 184 x 34 mm	<b>IME12DIN1_</b>	zł 732	<b>zł 476</b>

## SIGNATURE SERIES™ TUBS

### IMS-140\_

Together, IMS™ Signature Series™ Tubs and Cassettes create the most efficient and organised system for managing instruments and consumable products.

#### Include:

- Tub Signature
- Tub Liner
- Tub Divider (package of 8)
- Single Tub Cup w/Cover
- Double Tub Cup w/Cover
- Divided Slide Tray
- Locking Tub Cover




Sugg. Retail	Promo
zł 474	<b>zł 308</b>

### CHOOSE THE COLOR.

Just add the color coding number to the last part of the partcode (replace "\_" with the number)

Red	0	Lavender	6
Grey	1	Black	7
White	2	Blue	8
Orange	3	Green	9
Purple	4	Ocean Blue	11
Yellow	5		



For more information on Hu-Friedy products or to place an order, call at 00800 48 37 43 39 or visit [HuFriedy.eu](https://HuFriedy.eu)  
For more information on Crosstex products, call +39 05 24 52 74 53 or visit [crosstex.com](https://crosstex.com)  
For more information on Omnia products or to place an order, call +39 05 24 52 74 53 or visit [omniaspa.eu](https://omniaspa.eu)



Biuro Obsługi Klienta  
+48 22 338 70 50

[www.medif.com](https://www.medif.com)

As some of the products included in this catalogue might not be available in some countries, we invite you to verify the availability of the product with your local authorised dealer or to contact our customer care at the free phone number 00800 48 37 43 39.  
All prices indicated are intended. VAT included.  
Hu-Friedy Mfg. Co., LLC • European Headquarters • Lyoner Str. 9 • 60528 Frankfurt am Main, Germany • [HuFriedy.eu](https://HuFriedy.eu)  
All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.  
©2021 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-071/0521