


Kits & Collections	A602	Explorers	A629
Additional Collections & Kits	A604	Excavators	A630
Rubber Dam Clamps	A607	Spreaders	A631
Anterior Clamps	A607	Nickel Titanium	A633
Premolar Clamps	A608	Pluggers	A635
General Molar Clamps	A609	Post & Silver Point Removal Forceps	A639
Upper Molar Clamps	A610	Rulers	A641
Lower Molar Clamps	A611	Locking Pliers	A641
Rubber Dam Kits	A613		
Rubber Dam Accessories	A615		
Dental Dams	A617		
Microsurgical	A618		
Velvart Mirrors	A618		
Micro Mirrors	A619		
Mirror Handles	A619		
Curettes	A619		
Explorers	A621		
Pluggers	A622		
Periosteal Elevators	A624		
Excavators	A624		
Condensers/Burnishers	A625		



Endodontic

Our extensive line of endodontic hand instruments and accessories, many developed through our ongoing partnerships with Key Opinion Leaders, function as a seamless extension of your skill. From routine root canals to complex microsurgery, our endodontic dental instruments deliver exceptional response, durability and integrity for a variety of endodontic procedures. The result is outstanding patient care.



Kits & Collections

Our Endodontic portfolio is highlighted by an array of Endodontic Kits and Collections. Specifically curated for a slew of Endodontic procedures with help from our Key Opinion Leaders, these Kits come with all products needed to be the Best in Practice.



Learn more about our Key Opinion Leader [Bertrand Khayat](#).



Dr. Bertrand Khayat Endodontic Surgical Kit	IMKHAYAT1
Large Cassette - 20 Instrument, Blue	IM4208
#5 Scalpel Handle	10-130-05E
Mini Scalpel Blade/Mirror Handle Stainless Steel	SHDPV
Micro Mirror Round 3mm	MM4
Khayat Endo Microsurgical Plugger	PLGBK6
Khayat Endo Plastic Filling Instrument	PFI8ABK6
#85 Lucas DE Surgical Curette, #6 Handle	CL856
Velvart Endo Plugger	PLGVEL56
Perma Sharp™ Tissue Pliers	TP5072
PR-3S Prichard Periosteal, #6 Handle	PPR3S6
Columbia 13/14 Curette, #6 Handle	SC13/146
34/35 Mini Jacquette Scaler, #6 Handle	SJ34/356
89/92 Discoid-Cleoid Carver, #6 Handle	CD89/926
DG16/17 DE Endo Explorer, #6 Handle	EXDG16/176
Cawood-Minnesota Retractor	CRM2
Perm Sharp™ Castro Needle Holder, Curve Round	NH5021RC
#14 LaGrange Scissors, Super Cut	S14SC

Basic Endodontic Kit	ENDOKIT
#6 Cone Socket Mirror Handle	MH6
#4 Front Surface Rhodium Single-Sided Mirror	MIR4
1 Endo Locking Pliers	EPL1
DG16 DE Endo Explorer, Black Line	EXDG16X
31L DE Endo Excavator, Black Line	EXC31LX
D11NT SE Root Canal Spreader	RCS11NT
1/3 DE Root Canal Plugger, Black Line	RCP1/3X
5/7 DE Root Canal Plugger, Black Line	RCP5/7X
9/11 DE Root Canal Plugger, Black Line	RCP9/11X
1 DE Glick Plastic Plugger, Black Line	RCPGL1X



Chivian-Arens Endo Surgery Kit	ENDO-SXKIT
Mini CF II Amalgam Carrier	A1975
#85 Lucas DE Surgical Curette, #6 Handle	CL856
87 Lucas DE Surgical Curette	CL87
#2/4 Molt DE Curette	CM2/4
#1 DE Explorer	EXD1
Mathieu-Kocher Perma Sharp™ Needle Holder	NH5074
#1 Hourigan Periosteal	PH1
#149 Periosteal Elevator	P149
Columbia 13/14 Curette, #6 Handle	SC13/14
Indiana University 17/18 Curette	SIU17/18
18 Curved Iris Scissors	S18
Modified Serrated Minnesota Retractor	CRMT
#23/UNC15 Expro	XP23/UNC15
#5 Scalpel Handle, Round Straight	10-130-05

Dr. Thomas Von Arx Kit	KIT-VON
#5 Scalpel Handle	10-130-05E
Micro Mirror 3 mm x 9 mm	MM8
Mini Scalpel Blade/Mirror Handle Stainless Steel	SHDPV
Apical Explorer, #6 Handle	MEX16
Titanium Micro Tissue Pliers, Straight	TPASTMBH
#84 Lucas DE Surgical Curette, #6 Handle	CL846
#86 Lucas DE Surgical Curette, #6 Handle	CL866
Buser Periosteal, #6 Handle	PPBUSER6
#3/9A Periosteal, #6 Handle	PPR3/9A6
8A Composite Instrument	TNPF18A
5 Novatech Long/Fluted	CSNT5
Crile-Wood Perma Sharp™ Needle Holder	NH5038
Kelly Perma Sharp™ Scissors, 6.25" Curved	S5001
Large Left Micro-Condenser/Burnisher	MCBLL
Medium Left Micro-Condenser/Burnisher	MCBLM
Medium Right Micro-Condenser/Burnisher	MCBRM
Small Left Micro-Condenser/Burnisher	MCBLS
Small Right Micro-Condenser/Burnisher	MCBRS
Large Right Micro-Condenser/Burnisher	MCBRL
5 Abou-Rass Angled Burnisher	MAR-B5
Hartmann Hemostat, 4" Straight	HHS
Kocher Hemostat, 5.5" Curved	HKC
Oral Surgery Signature Series Cassette - 20 Instrument, Gray	IMOS201NC8
Hinged Instrument Clips, 2 Pack	IM1000



Dr. Ralan Wong Kit	KIT-WONG
#4 Front Surface Rhodium Single-Sided Mirror	MIR4
#6 Cone Socket Mirror Handle (QTY: 2)	MH6
#5 Front Surface Rhodium Single-Sided Mirror	MIR5
Qulix SE Probe, 3-6-9-12, #6 Handle	PCP126
DG16 DE Endo Explorer	EXDG16
18L Locking Dressing Plier	DP18L
#9 Molt Periosteal Elevator	P9
#23 Seldin Periosteal Retractor	P23
Gracey 1/2 Curette, #9 Handle	SG1/29E2
Woodson Periosteal	PFIWDS1
31L DE Endo Excavator	EXC31L
33L DE Endo Excavator	EXC33L
#10 Miller DE Surgical Curette	CM10
#11 Miller DE Surgical Curette	CM11
#12 Miller DE Surgical Curette	CM12
#42 Adson 1x2 Tissue Pliers	TP42
#18 Iris Scissors, Curved	S18
Perma Sharp™ Hegar-Baumgar Needle Holder	NH5034
Large Cassette - 16 Instrument, Blue	IM4168

Additional Collections & Kits

Abou-Rass Microsurgical Endodontic Kit	ABOU-RASS MICROSURGICAL ENDODONTIC KIT
Large Cassette - 16 Instrument, Blue	IM4168
Hinged Instrument Clips, 2 Pack	IM1000
#9 Molt Periosteal Elevator	P9
#1 Henahan Retractor, Left	TRH1
3 Abou-Rass Endo Excavator DE	MAR-C3
3 Abou-Rass Apical Surgical Curette SE	MAR-SC3
1 Abou-Rass Endo Explorer DE	MAR-EX1
DG16/17 DE Endo Explorer	EXDG16/17
1 Abou-Rass Straight Apical Condenser	MAR-SAC1
1 Abou-Rass Angulated Apical Condenser	MAR-AAC1
3 Abou-Rass Straight Burnisher	MAR-SB3
3 Abou-Rass Angled Burnisher	MAR-B3
#50L Fox Tissue Pliers w/Lock	TP50L
1 Velvart	MMVEL1
2 Velart	MMVEL2
Metzenbaum Perma Sharp™ Scissors, Curved/Blunt	S5055
Crile-Wood Perma Sharp™ Needle Holder	NH5038
Perma Sharp™ Halsey Needle Holder, Serrated	NH5037
4 Kramer-Nevins DE Curette	SKN4

Hahn Endodontic Microsurgical Kit	HAHN ENDODONTIC MICROSURGICAL KIT
Large Cassette - 16 Instrument, Blue	IM4168
Castro Microsurgical Scissors	SPV
Micro Cooley SinusLine Curved Tissue Plier	TPSLCOCMN
Microsurgical Needle Holder	NHDPV
Mini Scalpel Blade/Mirror Handle	SHDPV
Periosteal PH2MBH	PH26M
Pr. Spahr Papilla Periosteal, #6 Handle	PPSPAHR6
Mini Scalpel Blade 69 FRE Box of 12	MB69
Labanca Periosteal, #6 Handle	PPLABANCA6

Velvart Endodontic Microsurgical Kit	VELVART ENDODONTIC MICROSURGICAL KIT
Large Cassette - 16 Instrument, Blue	IM4168
Micro Mirror Round 3mm	MM4
Micro Mirror Round, 5mm	MM5
Micro Mirror Round, 2mm x 7mm	MM7
Micro Mirror 3 mm x 9 mm	MM8
Microsurgical Needle Holder	NHDPV
Mini Scalpel Blade/Mirror Handle Stainless Steel	SHDPV
Mini Scalpel Blade 67 Box of 12	MB67
Mini Scalpel Blade 69 FRE Box of 12	MB69
Castro Microsurgical Scissors	SPV
Anatomic Microsurgical Tissue Pliers	TPDAPV
Microsurgical Tissue Pliers, Curved	TPDSCPV
Microsurgical Tissue Pliers, Straight	TPDSSPV

Machtou Endodontic Pre-Treatment Kit	MACHTOU ENDODONTIC PRE-TREATMENT KIT
Large Cassette - 16 Instrument, Blue	IM4168
Hinged Instrument Clips, 2 Pack	IM1000
IMS™ Bur Cushion, Short Lid	IMS-1372S
Articulation Paper Pliers, European	APF2E
3S Hollenback, #6 Handle	CVHL3S6
U17 Utility Pick-up Dressing Pliers	DPU17
32L DE Endo Excavator, #6 Handle	EXC32L6
#6L SE Explorer	EXS6XL
1 Kelly Hemostat, Straight	H1
#5 Front Surface Rhodium Single-Sided Mirror	MIR5
#6 Simple Stem European Mirror Handle	MHE6
Qulix SE UNC Probe, 1-15, #6 Handle	PCPUNC156
2 Loesche, #6 Handle	PFIGML26
Rubber Dam Frame, Adult	RDAF6
Rubber Dam Clamp #1	RDCM1
Rubber Dam Clamp #212	RDCM212
Rubber Dam Clamp #26N	RDCM26N
Rubber Dam Clamp #4	RDCM4
Rubber Dam Clamp #8	RDCM8
Rubber Dam Clamp #9	RDCM9
Rubber Dam Clamp #W8A	RDCMW8A
Rubber Dam Forceps #4	RDF4
#17 Iris Scissors, Straight	S17
204SD Posterior Sickle Scaler, #6 Handle	S204SD6
Cook-Waite Anesthetic Aspirating Syringe, 1.8ml	SYRCWE

Machtou Endodontic Treatment Kit	MACHTOU ENDODONTIC TREATMENT KIT
Large Cassette - 16 Instrument, Blue	IM4168
IMS™ Large Parts Box	IMS-1273
A6 Rigid, #6 Handle	CSA6
#1 Dressing Pliers (Perry)	DP1
DG16 DE Endo Explorer, #7 Handle	EXDG167
5 European Front Surface Mouth Mirror, 12 pack	M5C
#6 Simple Stem European Mirror Handle	MHE6
10 SE Root Canal Plugger	RCP10
11 SE Root Canal Plugger	RCP11
8 SE Root Canal Plugger	RCP8
9 SE Root Canal Plugger	RCP9
1/2 Machtou Root Canal Plugger	RCPM1/2
3/4 Machtou Root Canal Plugger	RCPM3/4
D11T Root Canal Spreader	RSCD11T
#17 Iris Scissors, Straight	S17



Endodontic Infinity Kit	ENDODONTIC INFINITY KIT
Large Cassette with Accessory Area - 12 Instrument, Gray	IMN4121 or IMEDIN131
#6 Cone Socket Mirror Handle	MH6
#4 Front Surface Rhodium Single-Sided Mirror	MIR4
DG16 DE Endo Explorer, #6 Handle	EXDG166
1 Endo Locking Pliers	EPL1
31L DE Endo Excavator, #6 Handle	EXC31L6
D11NT SE Root Canal Spreader	RCSD11NT
1/3 DE Root Canal Plugger	RCP1/3
5/7 DE Root Canal Plugger	RCP5/7
9/11 DE Root Canal Plugger	RCP9/11
1 DE Glick Plastic Plugger	RCPGL1
Miller Articulating Paper Forceps	APF2
Cook-Waite Anesthetic Aspirating Syringe, 1.8ml	SYRCW
Rubber Dam Clamp #1	RDCM1
Rubber Dam Clamp #2	RDCM2
Rubber Dam Clamp #9	RDCM9
Rubber Dam Clamp #26N	RDCM26N
IMS™ Endodontic Stand - Holds 24	IMS-1275
IMS™ Small Parts Box	IMS-1271
IMS™ Bur Cushion, Short Lid	IMS-1372S
Signature Series™ Air/Water Clip	IM1005

Endodontic Infinity Kit 2	ENDODONTIC INFINITY KIT 2
Large Cassette with Accessory Area - 12 Instrument, Gray	IMN4121 OR IMEDIN131
#6 Cone Socket Mirror Handle	MH6
#4 Front Surface Rhodium Single-Sided Mirror	MIR4
#5 DE Explorer	EXD5
Qulix SE Probe, 3-6-9-12	PCP12
DG16 DE Endo Explorer	EXDG16
1 Endo Locking Pliers	EPL1
31L DE Endo Excavator	EXC31L
D11T SE Root Canal Spreader	RCSD11T
5/7 DE Root Canal Plugger	RCP5/7
1 DE Glick Plastic Plugger	RCPGL1
IMS™ Bur Cushion, Short Lid	IMS-1372S
IMS™ Small Parts Box	IMS-1271
Rubber Dam Clamp #1	RDCM1
Rubber Dam Clamp #2	RDCM2
Rubber Dam Clamp #9	RDCM9
Rubber Dam Clamp #26N	RDCM26N
Miller Articulating Paper Forceps	APF2

Berutti Pasqualini Microsurgical Endo Kit	MBPKIT
Berutti Pasqualini Explorer/Debriber	MBP-EX1
Berutti Pasqualini Elevator	MBP-E1
Berutti Pasqualini Plugger	MBP-PLG1
Berutti Pasqualini Plugger	MBP-PLG2
Berutti Pasqualini Plugger	MBP-PLG3
Berutti Pasqualini Curette	MBP-C1
Berutti Pasqualini Spatula/Burnisher	MBP-CB1
Retractor Columbia Cheek	CRC
#23 Seldin Periosteal Retractor	P23
Micro Tissue Plier, Straight 18cm	TPSLCOSM
Goldman-Fox, Straight, Super-Cut, Black Line	S20SCX
FSP Perma Sharp™ Needle Holder	NH5052
2 Kelly Hemostat, Curved	H2
Mini Scalpel Blade/Mirror Handle	SHDPV
Mini Scalpel Blade/Mirror Handle Stainless Steel	SHDPV
Mini Scalpel Blade 69 FRE Box of 12	MB69
Micro Mirror Round 3mm	MM4


Rubber Dam Clamps

Our Rubber Dam Clamps are made from proprietary Satin Steel™ for corrosion resistance, superior strength, and a fine balance between flexibility and memory. The matte finish on all 48 styles absorbs and dissipates light instead of reflecting it back, so you can experience a comfortable field of vision that is easier on the eyes — to help you be the Best in Practice.



Easily identify the perfect [Rubber Dam Clamp](#) for your procedure with our interactive [Pick Your Perfect Clamp Tool](#).

Anterior Clamps



Rubber Dam Clamp #6

Slightly offset jaws for labial caries on anterior teeth.

Black Line [RDCM6X](#)

Satin Steel™ Alloy [RDCM6](#)



Rubber Dam Clamp #9

Offset flat jaws and stiff bow for anterior teeth.

Black Line [RDCM9X](#)

Satin Steel™ Alloy [RDCM9](#)



Rubber Dam Clamp #9S

Offset flat serrated jaws for anterior teeth.

Black Line [RDCM9SX](#)

Satin Steel™ Alloy [RDCM9S](#)



Rubber Dam Clamp #212

Flat jaws for labial caries on anterior teeth.

Black Line [RDCM212X](#)

Satin Steel™ Alloy [RDCM212](#)



Rubber Dam Clamps #212SA

Slightly offset staggered jaws for odd shaped anteriors.

Black Line [RDCM212SAX](#)

Satin Steel™ Alloy [RDCM212SA](#)



Premolar Clamps



Rubber Dam Clamp #0

Flat jaws and high bow for longer premolars.

Black Line [RDCM0X](#)

Satin Steel™ Alloy [RDCM0](#)



Rubber Dam Clamp #1

Slightly festooned jaws for upper premolars.

Black Line [RDCM1X](#)

Satin Steel™ Alloy [RDCM1](#)



Rubber Dam Clamp #1A

Deeply festooned jaws for canines and premolars.

Black Line [RDCM1AX](#)

Satin Steel™ Alloy [RDCM1A](#)



Rubber Dam Clamp #2

Flat jaws for premolars.

Black Line [RDCM2X](#)

Satin Steel™ Alloy [RDCM2](#)



Rubber Dam Clamp #2A

Flat jaws for larger premolars.

Black Line [RDCM2AX](#)

Satin Steel™ Alloy [RDCM2A](#)



Rubber Dam Clamp #2AS

Flat serrated jaws for larger premolars.

Black Line [RDCM2ASX](#)

Satin Steel™ Alloy [RDCM2AS](#)



Rubber Dam Clamp #W2

Flat jaws for premolars.

Black Line [RDCMW2X](#)

Satin Steel™ Alloy [RDCMW2](#)



Rubber Dam Clamp #27N

Slightly festooned jaws for premolars.

Black Line [RDCM27NX](#)

Satin Steel™ Alloy [RDCM27N](#)



Rubber Dam Clamp #29

Flat jaws for lower premolars.

Black Line [RDCM29X](#)

Satin Steel™ Alloy [RDCM29](#)



Rubber Dam Clamp #209

Small, flat jaws with stiff bow for lower premolars.

Black Line [RDCM209X](#)

Satin Steel™ Alloy [RDCM209](#)



Rubber Dam Clamp #00

Flat jaws and high bow for small lower canines and first premolars.

Black Line	RDCM00X
Satin Steel™ Alloy	RDCM00

General Molar Clamps



Rubber Dam Clamp #10

Flat jaws designed for small upper left and lower right molars; paired with RDCM11.

Black Line	RDCM10X
Satin Steel™ Alloy	RDCM10



Rubber Dam Clamp #11

Flat jaws designed for small upper right and lower left molars; paired with RDCM10.

Black Line	RDCM11X
Satin Steel™ Alloy	RDCM11



Rubber Dam Clamp #14

Deeply festooned jaws for partially erupted or irregularly shaped molars.

Black Line	RDCM14X
Satin Steel™ Alloy	RDCM14



Rubber Dam Clamp #14A

Deeply festooned jaws for partially erupted or irregularly shaped molars; larger than RDCM14.

Black Line	RDCM14AX
Satin Steel™ Alloy	RDCM14A



Rubber Dam Clamp #W14A

Deeply festooned jaws for partially erupted or irregularly shaped molars.

Black Line	RDCMW14AX
Satin Steel™ Alloy	RDCMW14A



Rubber Dam Clamp #18

Flat jaws and stiff bow for molars.

Black Line	RDCM18X
Satin Steel™ Alloy	RDCM18



Rubber Dam Clamps



Rubber Dam Clamp #27

Flat jaws for large molars.

Black Line [RDCM27X](#)

Satin Steel™ Alloy [RDCM27](#)



Rubber Dam Clamp #56

Slightly festooned jaws for large molars.

Black Line [RDCM56X](#)

Satin Steel™ Alloy [RDCM56](#)



Rubber Dam Clamp #56S

Slightly festooned serrated jaws for large partially erupted or irregularly shaped molars.

Black Line [RDCM56SX](#)

Satin Steel™ Alloy [RDCM56S](#)

Upper Molar Clamps



Rubber Dam Clamp #4

Flat narrow jaws for small upper molars.

Black Line [RDCM4X](#)

Satin Steel™ Alloy [RDCM4](#)



Rubber Dam Clamp #W4

Flat narrow jaws for small upper molars.

Black Line [RDCMW4X](#)

Satin Steel™ Alloy [RDCMW4](#)



Rubber Dam Clamp #5

Slightly festooned jaws for large upper molars.

Black Line [RDCM5X](#)

Satin Steel™ Alloy [RDCM5](#)



Rubber Dam Clamp #8

Deeply festooned jaws for upper molars.

Black Line [RDCM8X](#)

Satin Steel™ Alloy [RDCM8](#)



Rubber Dam Clamp #26N

Flat jaws for upper molars.

Black Line **RDCM26NX**

Satin Steel™ Alloy **RDCM26N**



Rubber Dam Clamp #30

Offset jaws for upper right molars; paired with RDCM31.

Black Line **RDCM30X**

Satin Steel™ Alloy **RDCM30**



Rubber Dam Clamp #31

Offset jaws for upper left molars; paired with RDCM30.

Black Line **RDCM31X**

Satin Steel™ Alloy **RDCM31**



Rubber Dam Clamp #201

Slightly festooned jaws for upper molars.

Black Line **RDCM201X**

Satin Steel™ Alloy **RDCM201**



Rubber Dam Clamp #205

Slightly festooned jaws for large upper molars.

Black Line **RDCM205X**

Satin Steel™ Alloy **RDCM205**

Lower Molar Clamps



Rubber Dam Clamp #24

Flat jaws for upper right and lower left molars; paired with RDCM25.

Black Line **RDCM24X**

Satin Steel™ Alloy **RDCM24**



Rubber Dam Clamp #25

Flat jaws for upper left and lower right molars; paired with RDCM24.

Black Line **RDCM25X**

Satin Steel™ Alloy **RDCM25**



Rubber Dam Clamps



Rubber Dam Clamp #3

Flat narrow jaws for small lower molars.

Black Line **RDCM3X**

Satin Steel™ Alloy **RDCM3**



Rubber Dam Clamp #W3

Flat narrow jaws for small lower molars.

Black Line **RDCMW3X**

Satin Steel™ Alloy **RDCMW3**



Rubber Dam Clamp #7

Flat jaws for lower molars.

Black Line **RDCM7X**

Satin Steel™ Alloy **RDCM7**



Rubber Dam Clamp #7A

Slightly reversed festooned jaws for large lower molars.

Black Line **RDCM7AX**

Satin Steel™ Alloy **RDCM7A**



Rubber Dam Clamp #W7

Flat jaws for lower molars.

Black Line **RDCMW7X**

Satin Steel™ Alloy **RDCMW7**



Rubber Dam Clamp #12A

Serrated offset jaws for lower right molars; paired with RDCM13A.

Black Line **RDCM12AX**

Satin Steel™ Alloy **RDCM12A**



Rubber Dam Clamp #13A

Serrated offset jaws for lower left molars; paired with RDCM12A.

Black Line **RDCM13AX**

Satin Steel™ Alloy **RDCM13A**



Rubber Dam Clamp #26

Slightly reversed festooned jaws for small lower molars.

Black Line **RDCM26X**

Satin Steel™ Alloy **RDCM26**



Rubber Dam Clamp #28

Flat jaws for lower molars.

Black Line **RDCM28X**

Satin Steel™ Alloy **RDCM28**



Learn more about our **Rubber Dam Clamps and Accessories.**

Rubber Dam Kits

Our Rubber Dam Clamps are not only available separately, they are also available in sets and kits. Popular and commonly used Rubber Dam Clamps are bundled together for ease and convenience. Rubber Dam Clamp Kits options are also available that include accessories.





Black Line Rubber Dam Clamp Kit	RDKITX
Rubber Dam Clamp #13A, Black Line	RDCM13AX
Rubber Dam Clamp #12A, Black Line	RDCM12AX
Rubber Dam Clamp #27N, Black Line	RDCM27NX
Rubber Dam Clamp #W8A, Black Line	RDCMW8AX
Rubber Dam Clamp #212, Black Line	RDCM212X
Rubber Dam Clamp #00, Black Line	RDCM00X
Rubber Dam Clamp #9, Black Line	RDCM9X
Rubber Dam Clamp Organizing Board	RDCOB

Satin Steel™ Rubber Dam Clamp Set	RDCSET
Rubber Dam Clamp #0, Satin Steel™ Alloy	RDCM0
Rubber Dam Clamp #1, Satin Steel™ Alloy	RDCM1
Rubber Dam Clamp #14A, Satin Steel™ Alloy	RDCM14A
Rubber Dam Clamp #2, Satin Steel™ Alloy	RDCM2
Rubber Dam Clamp #212, Satin Steel™ Alloy	RDCM212
Rubber Dam Clamp #7, Satin Steel™ Alloy	RDCM7
Rubber Dam Clamp #8, Satin Steel™ Alloy	RDCM8
Rubber Dam Clamp #8A, Satin Steel™ Alloy	RDCM8A
Rubber Dam Clamp Organizing Board	RDCOB




Rubber Dam Clamp Set 7	RDCSET7
Rubber Dam Clamp #00, Satin Steel™ Alloy	RDCM00
Rubber Dam Clamp #2, Satin Steel™ Alloy	RDCM2
Rubber Dam Clamp #7, Satin Steel™ Alloy	RDCM7
Rubber Dam Clamp #8, Satin Steel™ Alloy	RDCM8
Rubber Dam Clamp #8A, Satin Steel™ Alloy	RDCM8A
Rubber Dam Clamp #9, Satin Steel™ Alloy	RDCM9
Rubber Dam Clamp #14A, Satin Steel™ Alloy	RDCM14A



Rubber Dam Kit	RDKIT
Rubber Dam Forceps	RDF
Rubber Dam Punch	RDP
Rubber Dam Frame, Adult	RDAF6
Universal Rubber Dam Template	RDTUNIV
FDI Rubber Dam Template	RDTFDI
Rubber Dam Clamp Organizing Board	RDCOB
Rubber Dam Clamp #0, Satin Steel™ Alloy	RDCM0
Rubber Dam Clamp #1, Satin Steel™ Alloy	RDCM1
Rubber Dam Clamp #2, Satin Steel™ Alloy	RDCM2
Rubber Dam Clamp #7, Satin Steel™ Alloy	RDCM7
Rubber Dam Clamp #8, Satin Steel™ Alloy	RDCM8
Rubber Dam Clamp #8A, Satin Steel™ Alloy	RDCM8A
Rubber Dam Clamp #14A, Satin Steel™ Alloy	RDCM14A
Rubber Dam Clamp #212, Satin Steel™ Alloy	RDCM212



Rubber Dam Accessories




Universal Rubber Dam Template

6" x 6" (15.2cm x 15.2cm)

Adult template shown.

Pedodontic template also included.

RTUNIV



FDI Rubber Dam Template

6" x 6" (15.2cm x 15.2cm)

Adult template shown.

Pedodontic template also included.


RDTFDI



Rubber Dam Punch, Ainsworth

6.75" (17.5cm)


RDPA



Rubber Dam Punch

6.75" (16cm)

RDP




Rubber Dam Clamp Organizing Board

4.62" x 2.4" x .375" (11.7 cm x 6.1 cm x 1 cm)

Store, sterilize and organize with eight pegs that securely hold your clamps. Clamps not included.

RDCOB



Rubber Dam Frame

Child - 5" (13cm)


RDCF5



Rubber Dam Frame

Adult - 6" (15cm)

RDAF6



Rubber Dam Forceps #4

7" (17.5cm)

RDF4



Rubber Dam Forceps

6.5" (17cm)

RDF

Dental Dams

Our Dental Dam are manufactured from medium modulus natural latex and are powder-free – less than 5mg/sheet (ASTM D6124). These Dental Dams are low protein, extractable protein content is below FDA allowable claim limit of 50µg/g per sheet (ASTM D5712).



Learn more about our [Dental Dams](#).



Dental Dam

Non-Latex, Powder-Free/Low Protein, Peppermint
6"x 6" - 15 per box

Medium	19500
--------	--------------



Dental Dam

Latex, Powder-Free/Low Protein, Unflavored
5"x 5" - 52 per box

Thin	19500
------	--------------

Medium	19100
--------	--------------

Heavy	19101
-------	--------------

6" x 6" - 36 per box

Thin	19302
------	--------------

Medium	19300
--------	--------------

Heavy	19301
-------	--------------



Dental Dam

Latex, Powder-Free/Low Protein, Mint
5" x 5" - 52 per box

Thin	19202
------	--------------

Medium	19200
--------	--------------

Heavy	19201
-------	--------------

6" x 6" - 36 per box

Thin	19402
------	--------------

Medium	19400
--------	--------------

Heavy	19401
-------	--------------



Microsurgical

Made from highly polished stainless steel, come in round or rectangular shapes, are single or double-ended, and are offered in a variety of sizes.



Learn more about our line of [Microsurgical Instruments](#).

Velvart Mirrors



1 Small Rectangular

2 x 7 mm

Round, 1 Ring

MM1



2 Medium Rectangular

3x9 mm

Round, 2 Rings

MM2



3 Large Rectangular

4 x 12 mm

Round, 3 Rings

MM3



Small Round

2.5 mm

Round

MMRS



Large Round

4.5 mm

Round

MMR



Micro Mirrors



Micro Mirror Round 3mm
MM4



Micro Mirror Round 5mm
MM5



Micro Mirror 2 mm x 7 mm
MM7



Micro Mirror 3 mm x 9 mm
MM8

Mirror Handles



1 Mirror Handle
MH6692K



7 Mini Blade Scalpel Handle
10-130-70



8 Mini Blade Scalpel Handle
10-130-08



Microsurgical Scalpel Handle, Stainless Steel
SHDPV



Microsurgical Scalpel Handle, Titanium
SHTMBH

Curettes



1 Velvart
3 x 6 mm; 2 x 6 mm
#6

MMVEL1



2 Velvart
3.5 mm; 3.5 mm
#6

MMVEL2



1 Abou-Rass Apical Surgical Curette SE

1.5 mm

#40

MAR-SC1



2 Abou-Rass Apical Surgical Curette SE

2.5 mm

#40

MAR-SC2



3 Abou-Rass Apical Surgical Curette SE

3 mm

#40

MAR-SC3



4 Abou-Rass Apical Surgical Curette SE

4 mm

#40

MAR-SC4




1 Al-Fouzan Curette

1.0 mm

#6

MCAF1




2 Al-Fouzan Curette

1.2 mm

#6

MCAF2




3 Al-Fouzan Curette

1.4 mm

#6

MCAF3



Berutti Pasqualini Curette

0.5 mm

Black Line

MBP-C1

Refer to [page A606](#) for a detailed listing of
Berutti Pasqualini Microsurgical Endo Kit



Explorers



1 Abou-Rass Endo Explorer DE

#33


MAR-EX1



3 Abou-Rass Calibrated Endo Explorer DE

#1

MAR-EX3



2 Abou-Rass Endo Explorer SE

#41


MAR-EX2



4 Abou-Rass Endo Explorer SE

#40

MAR-EX4



5 Abou-Rass Endo Explorer SE

#40

MAR-EX5



6 Abou-Rass Endo Explorer SE

#40

MAR-EX6



1 Al-Fouzan Endo Explorer

#41

EXAF1



2 AI-Fouzan Endo Explorer

#41

EXAF2



1 AI-Fouzan Endo Explorer, Left

#41

MEXAF1L



1 AI-Fouzan Endo Explorer, Right

#41

MEXAF1R



2 AI-Fouzan Endo Explorer, Left

#41

MEXAF2L



2 AI-Fouzan Endo Explorer, Right

#41

MEXAF2R



Berutti Pasqualini Explorer/Debriber

Black Line

MBP-EX1



Micro DE Explorer

#41

MEX1

Pluggers



Berutti Pasqualini Plugger

0.5 mm & 0.7 mm

Black Line

MBP-PLG1



Berutti Pasqualini Plugger

0.5 mm

Black Line

MBP-PLG2



Berutti Pasqualini Plugger

0.7 mm

Black Line

MBP-PLG3



1 Retro Filling DE Plugger

0.80 mm

#41

PLGRF1



2 Retro Filling DE Plugger

1.10 mm

#41

PLGRF2



Velvart Endo Plugger

0.25 mm / 0.33 mm

#41

PLGVEL56



Khayat Endo Plastic Filling Instrument

#6 Satin Steel

PF18ABK6

- Refer to [page A602](#) for a detailed listing of
- [Dr. Bertrand Khayat Endodontic Surgical Kit](#)



Left Apical Micro-Retro Filling Instrument

#41

MRFL



Right Apical Micro-Retro Filling Instrument

#41

MRFR



Periosteal Elevators



1 Al-Fouzan Periosteal

0.80 mm

#6

MPAF1




2 Al-Fouzan Periosteal

1.15 mm

#6

MPAF2



3 Al-Fouzan Periosteal

1.35 mm

#6

MPAF3



Berutti Pasqualini Elevator

1.4 mm

Black Line

MBP-E1

Excavators



1 Abou-Rass Endo Excavator DE

1 mm

#41

MAR-C1



2 Abou-Rass Endo Excavator DE

2 mm

#41

MAR-C2



3 Abou-Rass Endo Excavator DE

3.2 mm

#41

MAR-C3



4 Abou-Rass Endo Excavator DE

4 mm

#41

MAR-C4



Condensers/Burnishers



.5 Abou-Rass Straight Apical Condenser

0.5 mm/0.6 mm

#41

MAR-SAC.5



1 Abou-Rass Straight Apical Condenser

0.8 mm/0.8 mm

#41

MAR-SAC1



1.5 Abou-Rass Straight Apical Condenser

1.5 mm/2.0 mm

#41

MAR-SAC1.5



2 Abou-Rass Straight Apical Condenser

1.5 mm/2.0 mm

#41

MAR-SAC2



3 Abou-Rass Straight Apical Condenser

1.5 mm/3.0 mm

#41

MAR-SAC3



.5 Abou-Rass Angulated Apical Condenser

0.5 mm/0.6 mm

#41

MAR-AAC.5



1 Abou-Rass Angulated Apical Condenser

0.8 mm/0.8 mm

#41

MAR-AAC1



1.5 Abou-Rass Angulated Apical Condenser

1.5 mm/1.5 mm

#41

MAR-AAC1.5



2 Abou-Rass Angulated Apical Condenser

1.5 mm/2.0 mm

#41

MAR-AAC2



3 Abou-Rass Angulated Apical Condenser

1.5 mm/3.0 mm

#41

MAR-AAC3



1 Abou-Rass Straight Burnisher

1.5 mm

#40

MAR-SB1



2 Abou-Rass Straight Burnisher

1.75 mm

#40

MAR-SB2



3 Abou-Rass Straight Burnisher

2 mm

#40

MAR-SB3



XX 4 Abou-Rass Straight Burnisher

2.25 mm

#40

MAR-SB4



5 Abou-Rass Straight Burnisher

2.5 mm

#40

MAR-SB5



1 Abou-Rass Angled Burnisher

0.5 mm

#40

MAR-B1



2 Abou-Rass Angled Burnisher

1 mm

#40

MAR-B2



3 Abou-Rass Angled Burnisher

1.25 mm

#40

MAR-B3



4 Abou-Rass Angled Burnisher

1.75 mm

#40

MAR-B4



Al-Fouzan Burnisher, Left

0.78 mm

#41

MBPFIAFL



Al-Fouzan Burnisher, Right

0.78 mm

#41

MBPFIAFR



5 Abou-Rass Angled Burnisher

2 mm

#40

MAR-B5



Berutti Pasqualini Spatula/Burnisher

0.55 mm

Black Line

MBP-CB1




Large Left Micro-Condenser/Burnisher

0.1 mm / 1.4 mm

#41

MCBLL




Large Right Micro-Condenser/Burnisher

0.1 mm / 1.4 mm

#41

MCBRL



Medium Right Micro-Condenser/Burnisher

0.75 mm / 1.4 mm

#41

MCBRM



Medium Right Micro-Condenser/Burnisher

0.75 mm / 1.4 mm

#41

MCBLM




Small Right Micro-Condenser/Burnisher

0.50 mm / 1.4 mm

#41

MCBRS



Small Left Micro-Condenser/Burnisher

0.50 mm / 1.4 mm

#41

MCBLS



Universal Small Micro-Condenser/Burnisher

0.30 mm / 1.4 mm

#41

MCBUS



Universal Medium Micro-Condenser/Burnisher

0.75 mm / 1.4 mm

#41

MCBUM



Universal Large Micro-Condenser/Burnisher

0.90 mm / 1.4 mm

#41

MCBUL

Explorers

Our endodontic explorers are made from a special blend of alloy which provides the flexibility, strength, tactile sensitivity and resilience clinicians need when exploring and assessing the treatment site. They are available in several double-ended designs and handles.



Learn more about our [Explorers](#).



DG16 DE Endo Explorer

Black Line	EXDG16X
#6	EXDG166
#7	EXDG167
#31	EXDG16
Thick	EXDG16T



DG16/17 DE Endo Explorer

Black Line	EXDG16/17X
#31	EXDG16/17

DG16/23 DE Endo Explorer

Black Line	EXDG16/23X
#31	EXDG16/23



Excavators

Our endodontic excavators feature extended-length shanks and are meticulously handcrafted from the highest grade stainless steel. They help the clinician access the pulp chamber for easy and efficient removal of dentin.



Learn more about our [Excavators](#).



2 Glick DE Endo Excavator

1.5 mm

#41

[EXCGL2](#)



31L DE Endo Excavator

1 mm

Black Line


[EXC31LX](#)

#6

[EXC31L6](#)

#41

[EXC31L](#)



31W DE Oval Spoon Endo Excavator

1.3 mm

Black Line


[EXC31WX](#)

#6

[EXC31W6](#)

#41

[EXC31W](#)



32L DE Endo Excavator

1.5 mm

Black Line


[EXC32LX](#)

#6

[EXC32L6](#)

#41

[EXC32L](#)



33L DE Endo Excavator

2 mm

Black Line

[EXC33LX](#)

#6

[EXC33L6](#)

#41

[EXC33L](#)

Spreaders

We manufacture a wide variety of spreaders, all handcrafted from the finest stainless steel to precise specifications. Our spreaders help the clinician gain access into the most difficult canals for optimal compaction of filling material during lateral condensation.



Learn more about our [Spreaders](#).



GP1 SE Root Canal Spreader

Black Line	RCSGP1X
#6	RCSGP16
#41	RCSGP1



1S Wakai SE Root Canal Spreader

#41	RCSW1S
-----	------------------------



2S Wakai SE Root Canal Spreader

#41	RCSW2S
-----	------------------------



D11S SE Root Canal Spreader

#41	RCSD11S
-----	-------------------------



D11 SE Root Canal Spreader

Black Line	RCSD11X
#6	RCSD116
#41	RCSD11



D11T SE Root Canal Spreader

Black Line	RCSD11TX
#6	RCSD11T6
#41	RCSD11T



D11TS SE Root Canal Spreader

Black Line **RCSD11TSX**

#41 **RCSD11TS**



MA57 SE Root Canal Spreader

Black Line **RCSMA57X**



ISO Size 30 Endo Spreader

24 mm - Blue

#41 **RCS30**



ISO Size 40 Endo Spreader

24 mm - Black

#41 **RCS40**



ISO Size 50 Endo Spreader

28 mm - Yellow

#41 **RCS50**



ISO Size 60 Endo Spreader

28 mm - Blue

#41 **RCS60**

Nickel Titanium

Nickel Titanium, or NiTi, points provide these spreaders with extra flexibility and “memory” while the ergonomic, color coded Satin Steel™ Handles provide maximum comfort, control and easy identification. Choose from nine popular designs.



D11NT SE Root Canal Spreader

.30 mm

#7

RCSD11NT



D11TNT SE Root Canal Spreader

.23 mm

#7

RCSD11TNT



D11TSNT SE Root Canal Spreader

.25 mm

#7

RCSD11TSNT



25SNT SE Root Canal Spreader

.25 mm

#7

RCS25SNT



25NT SE Root Canal Spreader

.40 mm

#7

RCS25NT



40NT SE Root Canal Spreader

.40 mm

#7

RCS40NT



40SNT SE Root Canal Spreader

.20 mm

#7

RCS40SNT



MA57NT SE Root Canal Spreader

.35 mm

#7

RCSMA57NT



3NT SE Root Canal Spreader

.35 mm

#7

RCS3NT



3 SE Root Canal Spreader

.35 mm

#41

RCS3




GP2 SE Root Canal Spreader

.35 mm

#41

RCSGP2




0 SE Posterior Spreader

.25 mm

#41

RCSOP




00 SE Posterior Spreader

.55 mm

#41

RCS00P




MA57 SE Root Canal Spreader

.20 mm

#41

RCSMA57



GP3 SE Root Canal Spreader

.25 mm

#41

RCSGP3

Pluggers

Our variety of pluggers give the clinician very different working ends for enhanced maneuvering and versatility.



Learn more about our [Pluggers](#).



ISO Sized 30 Plugger

.30 mm - Blue

#41

[RCP30](#)



ISO Sized 40 Plugger

.40 mm - Black

#41

[RCP40](#)




ISO Sized 50 Plugger

.50 mm - Yellow

#41

[RCP50](#)



ISO Sized 60 Plugger

.60 mm - Blue

#41

[RCP60](#)



1/3 DE Root Canal Plugger

.40 mm / .45 mm

Black Line

[RCP1/3X](#)

#41

[RCP1/3](#)



5/7 DE Root Canal Plugger


.50 mm / .75 mm

Black Line

[RCP5/7X](#)

#41

[RCP5/7](#)



9/11 DE Root Canal Plugger

1.00 mm / 1.15 mm

Black Line







[RCP9/11X](#)


#41

[RCP9/11](#)



Pluggers





2 Luks SE Root Canal Plugger

0.50 mm

#41

RCPL2




3 Luks SE Root Canal Plugger

0.75 mm

#41

RCPL3




4 Luks SE Root Canal Plugger

1.00 mm

#41

RCPL4



8 SE Root Canal Plugger

0.40 mm

#41

RCP8A




8-1/2 SE Root Canal Plugger

0.45 mm

#41

RCP8-1/2A




9 SE Root Canal Plugger

0.55 mm

#41

RCP9A




9-1/2 SE Root Canal Plugger

0.65 mm

#41

RCP9-1/2A




10 SE Root Canal Plugger

0.75 mm

#41

RCP10A




10-1/2 SE Root Canal Plugger

0.90 mm

#41

RCP10-1/2A




11 SE Root Canal Plugger

1.05 mm

#41

RCP11A




11-1/2 SE Root Canal Plugger

1.15 mm

#41

RCP11-1/2A



12 SE Root Canal Plugger

1.30 mm

#41

RCP12A




8 SE Root Canal Plugger

0.40 mm

#41

RCP8



8-1/2 SE Root Canal Plugger

0.45 mm

#41

RCP8-1/2




9 SE Root Canal Plugger

0.55 mm

#41

RCP9




9-1/2 SE Root Canal Plugger

0.65 mm

#41

RCP9-1/2



10 SE Root Canal Plugger

0.75 mm

#41

RCP10



10-1/2 SE Root Canal Plugger

0.90 mm

#41

RCP10-1/2




11 SE Root Canal Plugger

1.05 mm

#41

RCP11




11-1/2 SE Root Canal Plugger

1.15 mm

#41

RCP11-1/2



12 SE Root Canal Plugger

1.30 mm

#41

RCP12

Post & Silver Point Removal Forceps

Made from high quality stainless steel, our selection of Endodontic Forceps are designed for optimal balance, leverage, and control.



Learn more about our [Post & Silver Point Removal Forceps](#).



Steiglitz Straight

5-1/2" (14 cm)

RFSTR



Steiglitz (Archer) 45°

5-1/4" (13-1/2 cm)

RF45



Steiglitz 90°

4-3/4" (12 cm)

RF90



Peet 45°

4-3/4" (12 cm)

RFP



Rulers



6" Stainless Steel Ruler

6" (15.24 cm)

CLR6

Locking Pliers

These Pliers include a locking feature for easier control of gutta percha points. Locking Pliers are made of high quality stainless steel and are handcrafted to precise specifications.



1 Endo Locking Pliers

6-1/2" (16-1/2 cm)

EPL1