

Cassettes & Accessories**B102**

Infinity Series™ Cassettes

B102

Signature Series™ Cassettes

B110

Resin Cassettes

B116

Accessories & Organization Products

B117

Cassette Replacement Parts

B122

Procedural Kits**B132**

ADA CDR Coding Collection Series

B145


Procedural Tubs**B148**

Organizational Tubs

B148

Tub Accessories

B150



Instrument Management

Our cassette-based Instrument Management System (IMS) allows you to move instruments from cleaning to chairside efficiently and safely through reduced handling of contaminated instruments. By removing multiple time-consuming and potentially dangerous steps in the process, IMS™ Cassettes can save your facility five to ten minutes per procedure, reduce the chance of sharps injuries, and protect your bottom line through instrument integrity and increased productivity - all while promoting compliance to infection prevention standards of care.



Cassettes & Accessories

IMS™ Cassettes are manufactured from high quality, electroplished stainless steel shells to protect against corrosion while providing a long-lasting, professional appearance. IMS™ Cassettes feature medical-grade silicone rails, robust and ergonomic latches with positive locking mechanism, and are available in multiple sizes and colors that can be configured to best fit the needs of your facility.



Learn more about our [Cassette](#) portfolio.

Infinity Series™ Cassettes

Infinity Series™ Cassettes feature a revolutionary, silicone rail design and patented hole pattern to provide more access to instruments during reprocessing and increase compatibility with today's cleaning equipment. Smooth round corners promote draining and drying, while the beveled edge allows for a more ergonomic grasp. The patented latch design provides ergonomic, single-handed opening and closing and easy identification if the cassette is locked/unlocked. Available in 11 standard colors.

Exam Cassette



Exam Cassette - 5 Instrument

Holds 5 instruments

3" x 8" x 1.25" (76 mm x 203 mm x 32 mm)

	Red	IMN6050
	Gray	IMN6051
	White	IMN6052
	Orange	IMN6053
	Purple	IMN6054
	Yellow	IMN6055
	Lavender	IMN6056
	Black	IMN6057
	Blue	IMN6058
	Green	IMN6059
	Ocean Blue	IMN60511



Small Cassettes



Small Cassette - 8 Instrument		
<i>Holds 8 instruments with an accessory area.</i>		
5.5" x 8" x 1.25" (140 mm x 203 mm x 32 mm)		
 Red		IMN5080
 Gray		IMN5081
 White		IMN5082
 Orange		IMN5083
 Purple		IMN5084
 Yellow		IMN5085
 Lavender		IMN5086
 Black		IMN5087
 Blue		IMN5088
 Green		IMN5089
 Ocean Blue		IMN50811



Small Cassette - 10 Instrument		
<i>Holds 10 instruments. No accessory area.</i>		
5.5" x 8" x 1.25" (140 mm x 203 mm x 32 mm)		
 Red		IMN5100
 Gray		IMN5101
 White		IMN5102
 Orange		IMN5103
 Purple		IMN5104
 Yellow		IMN5105
 Lavender		IMN5106
 Black		IMN5107
 Blue		IMN5108
 Green		IMN5109
 Ocean Blue		IMN51011



Learn more about our [Instrument Management System](#) and request a consultation today.



12 Instrument Signa-Stat Cassette



Signa-Stat Cassette - 12 Instrument		
Holds 12 instruments with an accessory area. Compatible with SCiCan® STATIM® Autoclave or Midmark® M3 Autoclave. 6.5" X 10.5" X 1.25" (165 mm x 267 mm x 32 mm)		
	Red	IMN8120
	Gray	IMN8121
	White	IMN8122
	Orange	IMN8123
	Purple	IMN8124
	Yellow	IMN8125
	Lavender	IMN8126
	Black	IMN8127
	Blue	IMN8128
	Green	IMN8129
	Ocean Blue	IMN81211

Large Cassettes



Large Cassette With Accessory Area		
Holds 12 or 16 instruments with large accessory area. 8" X 11" X 1.25" (203 mm x 280 mm x 32 mm)		
	12 Instrument	16 Instrument
	Red IMN4120	IMN4160
	Gray IMN4121	IMN4161
	White IMN4122	IMN4162
	Orange IMN4123	IMN4163
	Purple IMN4124	IMN4164
	Yellow IMN4125	IMN4165
	Lavender IMN4126	IMN4166
	Black IMN4127	IMN4167
	Blue IMN4128	IMN4168
	Green IMN4129	IMN4169
	Ocean Blue IMN41211	IMN41611

IMN416x, 16 Instrument Cassette shown. Also available in 12 instrument cassette with larger accessory area (IMN412x).


Large Cassette - 20 Instrument
Holds 20 instruments. No accessory area.

8" X 11" X 1.25" (203 mm x 280 mm x 32 mm)

● Red	IMN4200
● Gray	IMN4201
○ White	IMN4202
● Orange	IMN4203
● Purple	IMN4204
● Yellow	IMN4205
● Lavender	IMN4206
● Black	IMN4207
● Blue	IMN4208
● Green	IMN4209
● Ocean Blue	IMN42011

Oral Surgery Cassettes

Oral Surgery Cassette

Accommodates 12 hand instruments, 4 elevators, 3 hinged instruments in the lid and a large accessory area. Hinged instruments should be in open position. Fits in most 10" chamber or larger sterilizers, most instrument washers, and larger ultrasonic baths.

8" X 14.5" X 1.25" (203 mm x 368 mm x 32 mm)

● Red	IMNOS0
● Gray	IMNOS1
○ White	IMNOS2
● Orange	IMNOS3
● Purple	IMNOS4
● Yellow	IMNOS5
● Lavender	IMNOS6
● Black	IMNOS7
● Blue	IMNOS8
● Green	IMNOS9
● Ocean Blue	IMNOS11



Double-Decker™ Cassettes




Double-Decker™ Cassette - 8 Instrument
Two-Tier cassettes, holds 4 instruments on each side for a total of 8 instruments.
3" X 8" X 1.5" (76 mm x 203 mm x 38 mm)

● Red	IMN9080
● Gray	IMN9081
○ White	IMN9082
● Orange	IMN9083
● Purple	IMN9084
● Yellow	IMN9085
● Lavender	IMN9086
● Black	IMN9087
● Blue	IMN9088
● Green	IMN9089
● Ocean Blue	IMN90811



Double-Decker™ Cassette - 14 Instrument
Two-Tier cassettes, holds 7 instruments on each side for a total of 14 instruments.
4.5" X 8" X 1.5" (114 mm x 203 mm x 38 mm)

● Red	IMN9140
● Gray	IMN9141
○ White	IMN9142
● Orange	IMN9143
● Purple	IMN9144
● Yellow	IMN9145
● Lavender	IMN9146
● Black	IMN9147
● Blue	IMN9148
● Green	IMN9149
● Ocean Blue	IMN91411


Double-Decker™ Casette - 20 Instrument

Two-Tier cassettes, holds 10 instruments on each side for a total of 20 instruments. Features an accessory area with a built-in needle cap holder and slot for a winged syringe. 6" X 9.7" X 1.5" (152 mm x 246 mm x 38 mm)

● Red	IMN9200
● Gray	IMN9201
○ White	IMN9202
● Orange	IMN9203
● Purple	IMN9204
● Yellow	IMN9205
● Lavender	IMN9206
● Black	IMN9207
● Blue	IMN9208
● Green	IMN9209
● Ocean Blue	IMN92011

Ultrasonic Insert/Handpiece Cassettes

The Small Ultrasonic Insert/Handpiece Cassette holds 2 ultrasonic inserts or 1 handpiece while the Large Ultrasonic Insert/Handpiece Cassette holds up to 3 high-speed handpieces or 5 ultrasonic inserts. Both cassettes feature silicone finger mat for protection of delicate items.


Ultrasonic Insert/Handpiece Cassette

	Small 1.5" X 8" X 1.25" (38 mm x 203 mm x 32 mm)	Large 3" X 8" X 1.25" (76 mm x 203 mm x 32 mm)
● Red	IMN30HU0	IMN60H0
● Gray	IMN30HU1	IMN60H1
○ White	IMN30HU2	IMN60H2
● Orange	IMN30HU3	IMN60H3
● Purple	IMN30HU4	IMN60H4
● Yellow	IMN30HU5	IMN60H5
● Lavender	IMN30HU6	IMN60H6
● Black	IMN30HU7	IMN60H7
● Blue	IMN30HU8	IMN60H8
● Green	IMN30HU9	IMN60H9
● Ocean Blue	IMN30HU11	IMN60H11



Small Double-Decker™ Ortho Cassette



Double-Decker™ Small Orthodontic Cassette

Holds up to 3 pliers, a Mathieu needle holder, and 2 hand instruments. Hinged instruments should be in open position. 4.5" X 8" X 1.5" (114 mm x 203 mm x 38 mm)

● Red	IMN9140-O
● Gray	IMN9141-O
○ White	IMN9142-O
● Orange	IMN9143-O
● Purple	IMN9144-O
● Yellow	IMN9145-O
● Lavender	IMN9146-O
● Black	IMN9147-O
● Blue	IMN9148-O
● Green	IMN9149-O
● Ocean Blue	IMN91411-O

Thin Ortho Cassettes



Large Thin Orthodontic Cassette

Holds up to 7 pliers and 4 hand instruments. Hinged instruments should be in open position. 8" X 11" X .75" (203 mm x 280 mm x 19 mm)

● Purple	IMNORTHL4
● Blue	IMNORTHL8
● Green	IMNORTHL9



Small Thin Orthodontic Cassette

Holds up to 4 pliers and 3 hand instruments. Hinged instruments should be in open position. 7" X 8" X .75" (178 mm x 203 mm x 19 mm)

● Purple	IMNORTHS4
● Blue	IMNORTHS8
● Green	IMNORTHS9



Radiography Cassette



Radiography Cassette

Radiography Cassette. Holds heat sterilizable x-ray accessories.

5.5" X 8" X 1.25" (140 mm x 203 mm x 32 mm)

IMN7000



Signature Series™ Cassettes

Signature Series™ Cassettes feature high quality, electroplated stainless steel, long lasting, medical grade silicone rails, and robust latches and are available in 8 standard colors.

Exam Cassettes



Exam Cassette - 5 Instrument

Holds 5 instruments.
3" X 8" X 1.25" (76 mm x 203 mm x 32 mm)

● Red	IM6050
● Gray	IM6051
○ White	IM6052
● Orange	IM6053
● Yellow	IM6055
● Lavender	IM6056
● Blue	IM6058
● Green	IM6059

Small Cassettes

8 instrument Small Cassettes include an accessory area while the 10 instrument Small Cassette does not include an accessory area.



8 instrument cassette shown. Also available in 10 instrument cassette with no accessory area.

Small Cassette

	8 Instrument 5.5" X 8" X 1.25" (140 mm x 203 mm x 32 mm)	10 Instrument 5.5" X 8" X 1.25" (140 mm x 203 mm x 32 mm)
● Red	IM5080	IM5100
● Gray	IM5081	IM5101
○ White	IM5082	IM5102
● Orange	IM5083	IM5103
● Yellow	IM5085	IM5105
● Lavender	IM5086	IM5106
● Blue	IM5088	IM5108
● Green	IM5089	IM5109



12 Instrument Signa-Stat Cassette



Signa-Stat Cassette - 12 Instrument

Holds 12 instruments with an accessory area. Compatible with SCiCan® STATIM® Autoclave or Midmark® M3 Autoclave.
6.5" X 10.5" X 1.25" (165 mm x 267 mm x 32 mm)

● Red	IM8120
● Gray	IM8121
○ White	IM8122
● Orange	IM8123
● Yellow	IM8125
● Lavender	IM8126
● Blue	IM8128
● Green	IM8129

Large Cassettes

Large Cassettes that hold 12 and 16 instruments include a large accessory area.

The Large Cassette that holds 20 instruments does not include an accessory area.



Large Cassette


8" X 11" X 1.25" (203 mm x 280 mm x 32 mm)

	12 Instrument	16 Instrument	20 Instrument
● Red	IM4120	IM4160	IM4200
● Gray	IM4121	IM4161	IM4201
○ White	IM4122	IM4162	IM4202
● Orange	IM4123	IM4163	IM4203
● Yellow	IM4125	IM4165	IM4205
● Lavender	IM4126	IM4166	IM4206
● Blue	IM4128	IM4168	IM4208
● Green	IM4129	IM4169	IM4209

16 instrument cassette shown above. Also available in 12 and 20 instrument cassettes



8 Instrument Double-Decker™ Cassette



Double-Decker™ Cassette - 8 Instrument

Two-Tier cassettes, holds 4 instruments on each side for a total of 8 instruments.

3" X 8" X 1.5" (76 mm x 203 mm x 38 mm)

Red	IM9080
Gray	IM9081
White	IM9082
Orange	IM9083
Yellow	IM9085
Lavender	IM9086
Blue	IM9088
Green	IM9089

14 Instrument Double-Decker™ Cassette



Double-Decker™ Cassette - 14 Instrument

Two-Tier cassettes, holds 7 instruments on each side for a total of 14 instruments.

4.5" X 8" X 1.5" (114 mm x 203 mm x 38 mm)

Red	IM9140
Gray	IM9141
White	IM9142
Orange	IM9143
Yellow	IM9145
Lavender	IM9146
Blue	IM9148
Green	IM9149

20 Instrument Double-Decker™ Cassette



Double-Decker™ Cassette - 20 Instrument

Two-Tier cassettes, holds 10 instruments on each side for a total of 20 instruments. Features an accessory area with a built-in needle cap holder and slot for a winged syringe.

6" X 9.7" X 1.5" (152 mm x 246 mm x 38 mm)

Red	IM9200
Gray	IM9201
White	IM9202
Orange	IM9203
Yellow	IM9205
Lavender	IM9206
Blue	IM9208
Green	IM9209



Oral Surgery Cassette



Oral Surgery Cassette

Accommodates 12 hand instruments, 4 elevators, 3 hinged instruments in the lid and a large accessory area. Hinged instruments should be in open position. Fits in most 10" chamber or larger sterilizers, most instrument washers, and larger ultrasonic baths. 8" X 14.5" X 1.25" (203 mm x 368 mm x 32 mm)

● Red	IMOS0
● Gray	IMOS1
○ White	IMOS2
● Orange	IMOS3
● Yellow	IMOS5
● Lavender	IMOS6
● Blue	IMOS8
● Green	IMOS9

Ultrasonic Insert/Handpiece Cassettes

Small Ultrasonic Insert/Handpiece Cassette holds 2 ultrasonic inserts or 1 handpiece while the Large Ultrasonic Insert/Handpiece Cassette holds up to 3 high-speed handpieces or 5 ultrasonic inserts. Both options feature silicone finger mats for protection of delicate items.



Ultrasonic Insert/Handpiece Cassette

	Small 1.5" X 8" X 1.25" (38 mm x 203 mm x 32 mm)	Large 3" X 8" X 1.25" (76 mm x 203 mm x 32 mm)
● Red	IM30HU0	IM60H0
● Gray	IM30HU1	IM60H1
○ White	IM30HU2	IM60H2
● Orange	IM30HU3	IM60H3
● Yellow	IM30HU5	IM60H5
● Lavender	IM30HU6	IM60H6
● Blue	IM30HU8	IM60H8
● Green	IM30HU9	IM60H9



Thin Ortho Cassettes



Large Thin Orthodontic Cassette

Holds up to 7 pliers and 4 hand instruments. Hinged instruments should be in open position.

8" X 11" X .75" (203 mm x 280 mm x 19 mm)

● Blue

IMORTHLG8



Small Thin Orthodontic Cassette

Holds up to 4 pliers and 3 hand instruments. Hinged instruments should be in open position.

7" X 8" X .75" (178 mm x 203 mm x 19 mm)

● Blue

IMORTHSM8

Small Double-Decker™ Ortho Cassette



Double-Decker™ Small Orthodontic Cassette

Holds up to 3 pliers, a Mathieu needle holder, and 2 hand instruments. Hinged instruments should be in open position.

4.5" X 8" X 1.5" (114 mm x 203 mm x 38 mm)

● Red

IM9140-OR

● Gray

IM9141-OR

○ White

IM9142-OR

● Orange

IM9143-OR

● Yellow

IM9145-OR

● Lavender

IM9146-OR

● Blue

IM9148-OR

● Green

IM9149-OR



Large Double-Decker™ Ortho Cassette



Double-Decker™ Large Orthodontic Cassette

Holds up to 6 pliers, 1 Mathieu needle holder, and 3 hand instruments. No accessory area. Hinged instruments should be in open position. 6" X 8" X 1.5" (152 mm x 203 mm x 38 mm)

● Red	IM5200-OR
● Gray	IM5201-OR
○ White	IM5202-OR
● Orange	IM5203-OR
● Yellow	IM5205-OR
● Lavender	IM5206-OR
● Blue	IM5208-OR
● Green	IM5209-OR

Piezo Handpiece Cassettes


Small Piezo Handpiece Cassette holds 1 Piezo Handpiece while the Large Piezo Handpiece Cassette holds 2 Piezo Handpieces.



Piezo Handpiece Cassette

	Single 1.5" X 8" X 1.25" (38 mm x 203 mm x 32 mm)	Double 3" X 8" X 1.25" (76 mm x 203 mm x 32 mm)
● Blue	IM30IQ8	IM60IQ8

Piezo Wrench Cassette



Piezo Wrench Cassette

Holds 3 Piezo Wrenches. 5.63" X 1.95" X 1.77" (143 mm x 50 mm x 45 mm)

● Blue	IMI38
--------	--------------



Resin Cassettes

IMS™ Resin Cassettes are constructed of high-grade medical resin. Soft silicone instrument rails gently hold instruments in place. Resin cassettes cool down quickly after sterilization and are lightweight.



Small Resin Cassette

Holds 7, 12 or 16 instruments.

Ideal for exam and prophylactic procedures.

4.5" X 8" X 1.5" (114 mm x 203 mm x 38 mm)

7 Instrument	IMS-3107
12 Instrument	IMS-3112
16 Instrument	IMS-3116



Large Resin Cassette, No Accessory Area

Holds 18 or 14 instruments. No accessory area.

8" X 11" X 1.5" (203 mm x 279 mm x 38 mm)

18 Instrument	IMS-1112
24 Instrument	IMS-1124



Large Resin Cassette with Accessory Area

Holds 12 or 16 instruments with accessory area.

8" X 11" X 1.5" (203 mm x 279 mm x 38 mm)

12 Instrument	IMS-1112
16 Instrument	IMS-1116



Accessories & Organization Products

We offer many cassette accessories and organizational products available to support smaller items in a procedure as well as color-coding and tracking needs.



IMS™ Parts Box

Designed to keep small items such as mirror heads, rubber dam clamps and matrix band retainers organized during the sterilization process.

Large
3.5" X 1.12" X 1" (89 mm x 29 mm x 25 mm)

IMS-1273

Small
1.56" X 1.25" X 1" (40 mm x 32 mm x 25 mm)

IMS-1271

Barcode Plates

Metal plate can have barcode label affixed or can be lasermarked with text. Barcodes cannot be lasermarked onto the metal.



Signature Series™ Barcode Plate

.88" X 1.88" X .047" (22 mm x 48 mm x 1.2 mm)

IMS-1378



Infinity Series™ Barcode Plate

.88" X 1.75" X .047" (22 mm x 44 mm x 1.2 mm)

IMS-1379



IMS™ Endodontic Stand

Hinged, double-lid stand holds endodontic files and reamers. Lid is etched with a ruler to measure stop positions easily from either direction.

2.6" X 2.25" X 1.2" (67 mm x 57 mm x 30 mm)

Holds 24	IMS-1275
Holds 48	IMS-1275D



Air/Water Clip


Hold air/water syringe tips securely in place throughout cleaning and sterilization. Clips snap directly onto cassettes. Ask your rep for options to permanently install a/w clips to your cassettes.

For Infinity Series™ Cassette	IM1005N
For Signature Series™ Cassette	IM1005
For Infinity Series™ and Signature Series™ Double-Decker™ Cassette	IM1006



IMS™ Resin Cassette Clip

Large	IMS-1278
Small	IMS-1279



Hinged Instrument Clips

Secure scissors, hemostats, and other hinged instruments in the open position for optimal cleaning and sterilization with hinged instrument clips. Ask your rep for options to permanently install hinged instrument clips to your cassettes.

2 Pack	IM1000
6 Pack	IM1002

**Needle Cap Holders**

Needle cap holds attach to the cassette to help maintain the cap in upright position for one-handed needle recapping. Ask your rep for options to permanently install needle cap holders to your cassettes.

2 Pack **IM1001**

6 Pack **IM1003**

Combination Pack **IM1004**

Includes: 6 snap-in needle cap holders
and 6 snap-in hinged instrument clips

**Rubber Dam Clamp Organizing Board**

Rubber Dam Clamp Organizing Board. Holds up to 8 clamps during reprocessing. Fits into cassette accessory of most IMS™ Cassettes.

4.6" X 2.4" X .375" (117 mm x 61 mm x 10 mm) **RDCOB**

Bur Cushions

IMS™ Bur Cushions feature a hollow, rectangular, silicone base with a stainless steel lid designed to hold burs securely in place during cleaning and sterilization. Hollow design allows drainage and proper cleaning. Silicone base allows for different size burs. Please check bur manufacturer instructions for use to ensure ability to be reprocessed.

**IMS™ 1/2 Bur Cushion, Short Lid**

1.75" X 1.3" X .9" (44 mm x 33 mm x 22 mm) **IMS-1372SH**

IMS™ 1/2 Bur Cushion, Tall Lid

1.75" X 2.44" X .9" (44 mm x 62 mm x 22 mm) **IMS-1372TH**

IMS™ Bur Cushion, Short Lid

3.3" X 1.3" X .9" (84 mm x 33 mm x 22 mm) **IMS-1372S**

IMS™ Bur Cushion, Tall Lid

3.3" X 2.44" X .9" (84 mm x 62 mm x 22 mm) **IMS-1372T**



IMS™ Color Code Instrument Rings - Large

Silicone rings used for identifying and organizing instruments. Helpful for maintaining proper instrument order in cassettes or identifying by procedure or clinician. Large rings compatible with #6, #7, #8, and #9 handles.
50 per jar

Gray	IMS-1281L
White	IMS-1282L
Orange	IMS-1283L
Brown	IMS-1284L
Yellow	IMS-1285L
Red	IMS-1286L
Green	IMS-1287L
Blue	IMS-1288L
Black	IMS-1289L
Pink	IMS-12810L
Lavender	IMS-12811L
Ocean Blue	IMS-12812L
Purple	IMS-12813L
Assorted Colors (90 per jar)	IMS-1280L

IMS™ Color Code Instrument Rings - Small

Silicone rings used for identifying and organizing instruments. Helpful for maintaining proper instrument order in cassettes or identifying by procedure or clinician. Small rings compatible with #2 and #4 handles.
50 per jar

Gray	IMS-1281
White	IMS-1282
Orange	IMS-1283
Brown	IMS-1284
Yellow	IMS-1285
Red	IMS-1286
Green	IMS-1287
Blue	IMS-1288
Black	IMS-1289
Pink	IMS-12810
Lavender	IMS-12811
Ocean Blue	IMS-12812
Purple	IMS-12813
Assorted Colors (90 per jar)	IMS-1280



IMS™ Color Tabs

*Silicone color tabs fit on cassette lid for color-code identification by procedure or clinician.
5 per package*

 Red	IMS-1290N
 Gray	IMS-1291N
 White	IMS-1292N
 Orange	IMS-1293N
 Purple	IMS-1294N
 Yellow	IMS-1295N
 Lavender	IMS-1296N
 Black	IMS-1297N
 Blue	IMS-1298N
 Green	IMS-1299N
 Pink	IMS-12910N
 Ocean Blue	IMS-12911N
 Brown	IMS-12913N

Assorted Colors (13 per package)	IMS-12912N
----------------------------------	----------------------------



Cassette I.D. Labels and Marker

Specially designed polyester pressure-sensitive labels for cassette identification. Use the permanent indelible ink marker to write procedural specific information on I.D. labels or on Monitor Tape.

IMS™ Cassette I.D. Labels (50 per package)	IMS-1234
--	--------------------------

IMS™ I.D. Marker	IMS-1235
------------------	--------------------------

- For instrument **Grips**, see [page A209](#)
- under Instrument Maintenance.



Cassette Replacement Parts

There are times where instrument rails or other small cassette components need replaced after a long time of use. Replacement items are available to keep your cassettes looking like new.

Infinity Series™ Replacement Rail Packs

Replacement packs provide both the compression and instrument rails for your Infinity Series™ Cassettes. It is best to replace them all at the same time for consistency. To identify which pack you need, look at the lid of your Infinity Series™ Cassette. If you have two embossed stripes across the lid/infinity pattern, you need to order the replacement pack for the embossed lid.



Infinity Series™ Small Handpiece Cassette - Rail Replacement Set

● Red	PK-30HU0RL
● Gray	PK-30HU1RL
○ White	PK-30HU2RL
● Orange	PK-30HU3RL
● Purple	PK-30HU4RL
● Yellow	PK-30HU5RL
● Lavender	PK-30HU6RL
● Black	PK-30HU7RL
● Blue	PK-30HU8RL
● Green	PK-30HU9RL
● Ocean Blue	PK-30HU11R



Plain Lid



Embossed lid

Infinity Series™ Large 12 Instrument Cassette - Rail Replacement Set

	Plain Lid	Embossed lid
● Red	PK-4120RL	PKE-4120R
● Gray	PK-4121RL	PKE-4121R
○ White	PK-4122RL	PKE-4122R
● Orange	PK-4123RL	PKE-4123R
● Purple	PK-4124RL	PKE-4124R
● Yellow	PK-4125RL	PKE-4125R
● Lavender	PK-4126RL	PKE-4126R
● Black	PK-4127RL	PKE-4127R
● Blue	PK-4128RL	PKE-4128R
● Green	PK-4129RL	PKE-4129R
● Ocean Blue	PK-41211RL	PKE-41211R




Plain Lid

Embossed lid

Infinity Series™ Large 16 Instrument Cassette - Rail Replacement Set

	Plain Lid	Embossed lid
Red	PK-4160RL	PKE-4160R
Gray	PK-4161RL	PKE-4161R
White	PK-4162RL	PKE-4162R
Orange	PK-4163RL	PKE-4163R
Purple	PK-4164RL	PKE-4164R
Yellow	PK-4165RL	PKE-4165R
Lavender	PK-4166RL	PKE-4166R
Black	PK-4167RL	PKE-4167R
Blue	PK-4168RL	PKE-4168R
Green	PK-4169RL	PKE-4169R
Ocean Blue	PK-41611RL	PKE-41611R



Plain Lid

Embossed lid

Infinity Series™ Large 20 Instrument Cassette - Rail Replacement Set

	Plain Lid	Embossed lid
Red	PK-4200RL	PKE-4200R
Gray	PK-4201RL	PKE-4201R
White	PK-4202RL	PKE-4202R
Orange	PK-4203RL	PKE-4203R
Purple	PK-4204RL	PKE-4204R
Yellow	PK-4205RL	PKE-4205R
Lavender	PK-4206RL	PKE-4206R
Black	PK-4207RL	PKE-4207R
Blue	PK-4208RL	PKE-4208R
Green	PK-4209RL	PKE-4209R
Ocean Blue	PK-42011RL	PKE-42011R

**Infinity Series™ Large Thin Ortho Cassette - Rail Replacement Set**

Purple	PK-ORTL4L
Blue	PK-ORTL8L
Green	PK-ORTL9L



Infinity Series™ Small Thin Ortho Cassette - Rail Replacement Set

● Purple	PK-ORTS4L
● Blue	PK-ORTS8L
● Green	PK-ORTS9L



Infinity Series™ Small 8 Instrument Cassette - Rail Replacement Set

● Red	PK2-5080RL
● Gray	PK2-5081RL
○ White	PK2-5082RL
● Orange	PK2-5083RL
● Purple	PK2-5084RL
● Yellow	PK2-5085RL
● Lavender	PK2-5086RL
● Black	PK2-5087RL
● Blue	PK2-5088RL
● Green	PK2-5089RL
● Ocean Blue	PK2-50811R



Infinity Series™ Small 10 Instrument Cassette - Rail Replacement Set

● Red	PK2-5100RL
● Gray	PK2-5101RL
○ White	PK2-5102RL
● Orange	PK2-5103RL
● Purple	PK2-5104RL
● Yellow	PK2-5105RL
● Lavender	PK2-5106RL
● Black	PK2-5107RL
● Blue	PK2-5108RL
● Green	PK2-5109RL
● Ocean Blue	PK2-51011R

**Infinity Series™ Exam Cassette - 5 Instrument - Rail Replacement Set**

● Red	PK2-6050RL
● Gray	PK2-6051RL
○ White	PK2-6052RL
● Orange	PK2-6053RL
● Purple	PK2-6054RL
● Yellow	PK2-6055RL
● Lavender	PK2-6056RL
● Black	PK2-6057RL
● Blue	PK2-6058RL
● Green	PK2-6059RL
● Ocean Blue	PK2-60511R

**Infinity Series™ Large Handpiece Cassette - Rail Replacement Set**

● Red	PK-60H0RL
● Gray	PK-60H1RL
○ White	PK-60H2RL
● Orange	PK-60H3RL
● Purple	PK-60H4RL
● Yellow	PK-60H5RL
● Lavender	PK-60H6RL
● Black	PK-60H7RL
● Blue	PK-60H8RL
● Green	PK-60H9RL
● Ocean Blue	PK-60H11RL



Plain Lid



Embossed lid

Infinity Series™ Signa Stat Cassette - 12 Instrument Rail Replacement Set

	Plain Lid	Embossed lid
● Red	PK-8120RL	PKE-8120R
● Gray	PK-8121RL	PKE-8121R
○ White	PK-8122RL	PKE-8122R
● Orange	PK-8123RL	PKE-8123R
● Purple	PK-8124RL	PKE-8124R
● Yellow	PK-8125RL	PKE-8125R
● Lavender	PK-8126RL	PKE-8126R
● Black	PK-8127RL	PKE-8127R
● Blue	PK-8128RL	PKE-8128R
● Green	PK-8129RL	PKE-8129R
● Ocean Blue	PK-81211RL	PKE-81211R



Infinity Series™ 8 Instrument Double-Decker™ Cassette - Rail Replacement Set

● Red	PK-9080RL
● Gray	PK-9081RL
○ White	PK-9082RL
● Orange	PK-9083RL
● Purple	PK-9084RL
● Yellow	PK-9085RL
● Lavender	PK-9086RL
● Black	PK-9087RL
● Blue	PK-9088RL
● Green	PK-9089RL
● Ocean Blue	PK-90811RL



Infinity Series™ 14 Instrument Double-Decker™ Cassette - Rail Replacement Set

● Red	PK-9140RL
● Gray	PK-9141RL
○ White	PK-9142RL
● Orange	PK-9143RL
● Purple	PK-9144RL
● Yellow	PK-9145RL
● Lavender	PK-9146RL
● Black	PK-9147RL
● Blue	PK-9148RL
● Green	PK-9149RL
● Ocean Blue	PK-91411RL




Infinity Series™ Double-Decker™ Small Orthodontic Cassette - Rail Replacement Set

● Red	PK-9140OR
● Gray	PK-9141OR
○ White	PK-9142OR
● Orange	PK-9143OR
● Purple	PK-9144OR
● Yellow	PK-9145OR
● Lavender	PK-9146OR
● Black	PK-9147OR
● Blue	PK-9148OR
● Green	PK-9149OR
● Ocean Blue	PK-91411OR



Infinity Series™ 20 Instrument Double-Decker™ Cassette - Rail Replacement Set
Embossed lid

● Red	PK-9200RL
● Gray	PK-9201RL
○ White	PK-9202RL
● Orange	PK-9203RL
● Purple	PK-9204RL
● Yellow	PK-9205RL
● Lavender	PK-9206RL
● Black	PK-9207RL
● Blue	PK-9208RL
● Green	PK-9209RL
● Ocean Blue	PK-92011RL



Plain Lid



Embossed lid

Infinity Series™ Oral Surgery Cassette - Rail Replacement Set

	Plain Lid	Embossed lid
● Red	PK-OS0RAIL	PKE-OS0R
● Gray	PK-OS1RAIL	PKE-OS1R
○ White	PK-OS2RAIL	PKE-OS2R
● Orange	PK-OS3RAIL	PKE-OS3R
● Purple	PK-OS4RAIL	PKE-OS4R
● Yellow	PK-OS5RAIL	PKE-OS5R
● Lavender	PK-OS6RAIL	PKE-OS6R
● Black	PK-OS7RAIL	PKE-OS7R
● Blue	PK-OS8RAIL	PKE-OS8R
● Green	PK-OS9RAIL	PKE-OS9R
● Ocean Blue	PK-OS11RAI	PKE-OS11R



Signature Series™ Replacement Rails & Packs

Replacement packs provide both the compression and instrument rails for your Signature Series™ Cassettes. It is best to replace them all at the same time for consistency.



Signature Series™ Small Handpiece Cassette - Rail Replacement Set

● Red	IM30HU0RL
● Gray	IM30HU1RL
○ White	IM30HU2RL
● Orange	IM30HU3RL
● Yellow	IM30HU5RL
● Lavender	IM30HU6RL
● Blue	IM30HU8RL
● Green	IM30HU9RL



Signature Series™ Small 8 Instrument Cassette - Rail Replacement Set

● Red	IM5080RAIL
● Gray	IM5081RAIL
○ White	IM5082RAIL
● Orange	IM5083RAIL
● Yellow	IM5085RAIL
● Lavender	IM5086RAIL
● Blue	IM5088RAIL
● Green	IM5089RAIL



Signature Series™ Small 10 Instrument Cassette - Rail Replacement Set

● Red	IM5100RAIL
● Gray	IM5101RAIL
○ White	IM5102RAIL
● Orange	IM5103RAIL
● Yellow	IM5105RAIL
● Lavender	IM5106RAIL
● Blue	IM5108RAIL
● Green	IM5109RAIL

**Signature Series™ Large 16 Instrument Cassette - Rail Replacement Set**

● Red	IM4160RAIL
● Gray	IM4161RAIL
○ White	IM4162RAIL
● Orange	IM4163RAIL
● Yellow	IM4165RAIL
● Lavender	IM4166RAIL
● Blue	IM4168RAIL
● Green	IM4169RAIL

**Signature Series™ Large 20 Instrument Cassette - Rail Replacement Set**

● Red	IM4200RAIL
● Gray	IM4201RAIL
○ White	IM4202RAIL
● Orange	IM4203RAIL
● Yellow	IM4205RAIL
● Lavender	IM4206RAIL
● Blue	IM4208RAIL
● Green	IM4209RAIL

**Signature Series™ Large Double-Decker™ - Ortho Cassette - Rail Replacement Set**

● Red	IM5200-ORR
● Gray	IM5201-ORR
○ White	IM5202-ORR
● Orange	IM5203-ORR
● Yellow	IM5205-ORR
● Lavender	IM5206-ORR
● Blue	IM5208-ORR
● Green	IM5209-ORR

**Signature Series™ Small Double-Decker™ - Ortho Cassette - Rail Replacement Set**

● Red	IM9140ORL
● Gray	IM9141ORL
○ White	IM9142ORL
● Orange	IM9143ORL
● Yellow	IM9145ORL
● Lavender	IM9146ORL
● Blue	IM9148ORL
● Green	IM9149ORL



Signature Series™ Large Handpiece Cassette - Rail Replacement Set

Red	IM60H0RAIL
Gray	IM60H1RAIL
White	IM60H2RAIL
Orange	IM60H3RAIL
Yellow	IM60H5RAIL
Lavender	IM60H6RAIL
Blue	IM60H8RAIL
Green	IM60H9RAIL



Signature Series™ Exam Cassette - 5 Instrument - Rail Replacement Set

Red	IM6050RAIL
Gray	IM6051RAIL
White	IM6052RAIL
Orange	IM6053RAIL
Yellow	IM6055RAIL
Lavender	IM6056RAIL
Blue	IM6058RAIL
Green	IM6059RAIL



Signature Series™ 8 Instrument Double-Decker™ Cassette - Rail Replacement Set

Red	IM9080RAIL
Gray	IM9081RAIL
White	IM9082RAIL
Orange	IM9083RAIL
Yellow	IM9085RAIL
Lavender	IM9086RAIL
Blue	IM9088RAIL
Green	IM9089RAIL



Signature Series™ 14 Instrument Double-Decker™ Cassette - Rail Replacement Set

Red	IM9140RAIL
Gray	IM9141RAIL
White	IM9142RAIL
Orange	IM9143RAIL
Yellow	IM9145RAIL
Lavender	IM9146RAIL
Blue	IM9148RAIL
Green	IM9149RAIL

**Signature Series™ 20 Instrument Double-Decker™ Cassette - Rail Replacement Set**

● Red	IM9200RAIL
● Gray	IM9201RAIL
○ White	IM9202RAIL
● Orange	IM9203RAIL
● Yellow	IM9205RAIL
● Lavender	IM9206RAIL
● Blue	IM9208RAIL
● Green	IM9209RAIL

**Signature Series™ Signa Stat Cassette - 12 Instrument - Rail Replacement Set**

● Red	IM8120RAIL
● Gray	IM8121RAIL
○ White	IM8122RAIL
● Orange	IM8123RAIL
● Yellow	IM8125RAIL
● Lavender	IM8126RAIL
● Blue	IM8128RAIL
● Green	IM8129RAIL

**Signature Series™ Oral Surgery Cassette - Rail Replacement Set**

● Red	IMOS0RAIL
● Gray	IMOS1RAIL
○ White	IMOS2RAIL
● Orange	IMOS3RAIL
● Yellow	IMOS5RAIL
● Lavender	IMOS6RAIL
● Blue	IMOS8RAIL
● Green	IMOS9RAIL

**Signature Series™ Thin Ortho Cassette - Rail Replacement Set**

● Blue (Small)	IMORTS8RL
● Blue (Large)	IMORTL8RL



Procedural Kits

The IMS™ Collection Series Kits feature a compilation of the most popular procedural set-ups. One part code for a full procedure consisting of a cassette and instruments is provided at a reduced price vs a la carte ordering. An ideal solution for recent dental graduates.



Learn more about our [Procedural Kits](#).



Exam Standard Set-Up

3" x 8" x 1.25" (76 mm x 203 mm x 32 mm)

	Infinity Series™ Set-Up IMNEXAM	Signature Series™ Set-Up IMEXAM
Exam Cassette - 5 Instrument, Yellow	IMN6055	IM6055
#5 Front Surface Rhodium Single- Sided Mirror (qty: 3)	MIR5/3	MIR5/3
#6 Cone Socket Mirror Handle	MH6	MH6
#5 DE Explorer, #6 Handle	EXD56	EXD56
Qulix SE Probe, 3-6-9-12, #6 Handle	PCP126	PCP126
Dressing Pliers, Utility Pick-Up	DPU17	DPU17



**Total Restorative Standard Set-Up
(Composite, Crown & Bridge, Amalgam)**
8" x 11" x 1.25" (203 mm x 280 mm x 32 mm)

	Infinity Series™ Set-Up IMNRESTOR	Signature Series™ Set-Up IMRESTOR
Large 16 Instrument Cassette, Yellow	IMN4165	IM4165
#6 Cone Socket Mirror Handle	MH6	MH6
#5 Front Surface Rhodium Single-Sided Mirror (qty: 3)	MIR5/3	MIR5/3
#23/CP-12 Expro, #6 Handle	XP23/126	XP23/126
18 Excavator, #6 Handle	EXC186	EXC186
113 Gingival Cord Packer, Serrated, #6 Handle	GCP1136	GCP1136
Calcium Hydroxide Placer, #6 Handle	PICH6	PICH6
0/1 Serrated Marquette Condenser, #6 Handle	PLG0/16	PLG0/16
1/2 Serrated Plugger, #6 Handle	PLG1/26	PLG1/26
Large/ Regular	AC5202	AC5202
3/6 Carver, #6 Handle	CD3/66	CD3/66
8 Wiland, #6 Handle	CVWI86	CVWI86
1/2 Hollenback, #6 Handle	CVHL1/26	CVHL1/26
Interproximal Carver, XTS™ Handle	TNCVIPC	TNCVIPC
Goldstein Flexi-Thin #2, XTS™ Handle	TNCIGFT2	TNCIGFT2
Goldstein Flexi-Thin, Extra-Flex #3, XTS™ Handle	TNCIGFT3	TNCIGFT3
27/29 Burnisher, #6 Handle	BB27/296	BB27/296
24 Flexible Spatula, #6 Handle	CS246	CS246
#2 Dressing Pliers (College)	DP2	DP2
Miller Articulating Paper Forceps	APF2	APF2
18 Iris Scissors, Curved	S18	S18
Cook-Waite Anesthetic Aspirating Syringe	SYRCW	SYRCW
Amalgam Well	WA	WA
IMS™ Bur Cushion, Short Lid	IMS-1272S	IMS-1372S
Hinged Instrument Clips, 2 Pack	IM1000	IM1000
Air/Water Clip	IM1005N	IM1005



Prophy Standard Set-Up

5.5" x 8" x 1.25" (140 mm x 203 mm x 32 mm)

	Infinity Series™ Set-Up IMNPROPHY	Signature Series™ Set-Up IMPROPHY
10 Instrument Cassette, Lavender	IMN5106	IM5106
#6 Cone Socket Mirror Handle	MH6	MH6
#5 Front Surface Rhodium Single- Sided Mirror (QTY: 3)	MIR5/3	MIR5/3
#11/12 ODU DE Explorer, #6 Handle	EXD11/126	EXD11/126
Qulix SE Probe, 3-6-9-12, #6 Handle	PCP126	PCP126
H6/7 Sickle Scaler, #9 Handle	SH6/79E2	SH6/79E2
Gracey 11/12 Standard Curette, #9 Handle	SG11/129E2	SG11/129E2
Gracey 13/14 Standard Curette, #9 Handle	SG13/149E2	SG13/149E2
Gracey 13/14 Standard Curette, #9 Handle	SC13/149E2	SC13/149E2
#3 Ceramic Stone	SS3C	SS3C
Air/Water Clip	IM1005N	IM1005

**Perio Standard Set-Up**

4.5" x 8" x 1.5" (114 mm x 203 mm x 38 mm)

	Infinity Series™ Set-Up	Signature Series™ Set-Up
	IMNPERIO	IMPERIO
Double-Decker™ Cassette - 14 Instrument, Blue	IMN9148	IM9148
#6 Cone Socket Mirror Handle	MH6	MH6
#5 Front Surface Rhodium Single- Sided Mirror (qty: 3)	MIR5/3	MIR5/3
#11/12 ODU After Five DE Explorer, #6 Handle	EXD11/12A6	EXD11/12A6
Qulix SE Probe, 3-6-9-12, #6 Handle	PCP126	PCP126
H6/7 Sickle Scaler, #9 Handle	SH6/79E2	SH6/79E2
204S Posterior Sickle Scaler, #9 Handle	S204S9E2	S204S9E2
Gracey 1/2 Mini Five™ Curette, #9 Handle	SAS1/29E2	SAS1/29E2
Gracey 7/8 After Five™ Curette, #9 Handle	SRPG7/89E2	SRPG7/89E2
Gracey 11/12 After Five™ Curette, #9 Handle	SRPG11/129E2	SRPG11/129E2
Gracey 13/14 After Five™ Curette, #9 Handle	SRPG13/149E2	SRPG13/149E2
Columbia 4R/4L Curette, #9 Handle	SC4R/4L9E2	SC4R/4L9E2
#3 Ceramic Stone	SS3C	SS3C
Air/Water Clip	IM1006	IM1006



Composite/Crown Prep Standard Set-Up

8" x 11" x 1.25" (203 mm x 280 mm x 32 mm)

	Infinity Series™ Set-Up	Signature Series™ Set-Up
	IMNCMPOSIT	IMCOMPOSIT
16 Instrument Cassette, Blue	IMN4168	IM4168
#6 Cone Socket Mirror Handle	MH6	MH6
#5 Front Surface Rhodium Single- Sided Mirror (qty: 3)	MIR5/3	MIR5/3
#5 DE Explorer, #6 Handle	EXD56	EXD56
#18 DE Excavator, #6 Handle	EXC186	EXC186
#113 Gingival Cord Packer, Serrated, #6 Handle	GCP1136	GCP1136
Interproximal Carver, XTS™ Handle	TNCVIPC	TNCVIPC
Goldstein Flexi-Thin #2, XTS™ Handle	TNCIGFT2	TNCIGFT2
Goldstein Flexi-Thin, Extra-Flex #3, XTS™ Handle	TNCIGFT3	TNCIGFT3
24 Flexible Spatula, #6 Handle	CS246	CS246
#2 Dressing Pliers (College)	DP2	DP2
Miller Articulating Paper Forceps	APF2	APF2
Cook-Waite Anesthetic Aspirating Syringe	SYRCW	SYRCW
18 Iris Scissors, Curved	S18	S18
IMS™ Bur Cushion, Short Lid	IMS-1372S	IMS-1372S
Hinged Instrument Clips, 2 Pack	IM1000	IM1000
Air/Water Clip	IM1005N	IM1005

**Endo Standard Set-Up**

8" x 11" x 1.25" (203 mm x 280 mm x 32 mm)

	Infinity Series™ Set-Up	Signature Series™ Set-Up
	IMNENDO	IMENDO
Large Cassette - 12 Instrument, Gray	IMN4121	IM4121
#6 Cone Socket Mirror Handle	MH6	MH6
#4 Front Surface Rhodium Single-Sided Mirror, Cone Socket (qty: 3)	MIR4/3	MIR4/3
DG16 DE Endo Explorer, #6 Handle	EXDG166	EXDG166
1 Endo Locking Pliers	EPL1	EPL1
31L DE Endo Excavator, #6 Handle	EXC31L6	EXC31L6
D11NT SE Root Canal Spreader	RCSD11NT	RCSD11NT
1/3 DE Root Canal Plugger	RCP1/3	RCP1/3
5/7 DE Root Canal Plugger	RCP5/7	RCP5/7
9/11 DE Root Canal Plugger	RCP9/11	RCP9/11
1 DE Glick Plastic Plugger	RCPGL1	RCPGL1
Miller Articulating Paper Forceps	APF2	APF2
Cook-Waite Anesthetic Aspirating Syringe	SYRCW	SYRCW
Rubber Dam Forceps #4	RDF	RDF
Rubber Dam Clamp #1	RDCM1	RDCM1
Rubber Dam Clamp #2	RCDM2	RCDM2
Rubber Dam Clamp #9	RDCM9	RDCM9
Rubber Dam Clamp #26N	RDCM26N	RDCM26N
IMS Endodontic Stand - Holds 24	IMS-1275	IMS-1275
IMS Small Parts Box	IMS-1271	IMS-1271
IMS Bur Cushion, Short Lid	IMS-1372S	IMS-1372S
Air/Water Clip	IM1005N	IM1005



Suture Removal Standard Set-Up

5.5" x 8" x 1.25" (140 mm x 203 mm x 32 mm)

	Infinity Series™ Set-Up	Signature Series™ Set-Up
	IMNSUTURE	IMSUTURE
Small Cassette - 10 Instrument, Orange	IMN5103	IM5103
#6 Cone Socket Mirror Handle	MH6	MH6
#5 Front Surface Rhodium Single- Sided Mirror (qty: 3)	MIR5/3	MIR5/3
#5 DE Explorer, #6 Handle	EXD56	EXD56
#13S Suture Scissors	S13S	S13S
2 College Dressing Pliers	DP2	DP2
Air/Water Clip	IM1005N	IM1005
Hinged Instrument Clips, 2 Pack	IM1000	IM1000



Cement (Crown Seat) Standard Set-Up

5.5" x 8" x 1.25" (140 mm x 203 mm x 32 mm)

	Infinity Series™ Set-Up	Signature Series™ Set-Up
	IMNCEMENT	IMCEMENT
Small Cassette - 10 Instrument, Red	IMN5100	IM5100
#6 Cone Socket Mirror Handle	MH6	MH6
#5 Front Surface Rhodium Single- Sided Mirror (qty: 3)	MIR5/3	MIR5/3
#5 DE Explorer, #6 Handle	EXD56	EXD56
8 Wiland, #6 Handle	CVWI86	CVWI86
24 Flexible Spatula, #6 Handle	CS246	CS246
Miller Articulating Paper Forceps	APF2	APF2
Cook-Waite Anesthetic Aspirating Syringe	SYRCW	SYRCW
2 College Dressing Pliers	DP2	DP2
IMS™ 1/2 Bur Cushion, Short Lid	IMS-1372SH	IMS-1372SH
Air/Water Clip	IM1005N	IM1005N

**Crown & Bridge Standard Set-Up**

8" x 11" x 1.25" (203 mm x 280 mm x 32 mm)

	Infinity Series™ Set-Up	Signature Series™ Set-Up
	IMNCROWN	IMCROWN
Large Cassette With Accessory Area - 12 Instrument, Green	IMN4169	IM4169
#6 Cone Socket Mirror Handle	MH6	MH6
#5 Front Surface Rhodium Single- Sided Mirror (qty: 3)	MIR5/3	MIR5/3
#5 DE Explorer, #6 Handle	EXD56	EXD56
Qulix SE Probe, 3-6-9-12, #6 Handle	PCP126	PCP126
8 Wiland, #6 Handle	CVWI86	CVWI86
18 Excavator, #6 Handle	EXC186	EXC186
113 Gingival Cord Packer, Serrated, #6 Handle	GCP1136	GCP1136
Goldstein Flexi-Thin #2, XTS™ Handle	TNCIGFT2	TNCIGFT2
24 Flexible Spatula, #6 Handle	CS246	CS246
2 College Dressing Pliers	DP2	DP2
Miller Articulating Paper Forceps	APF2	APF2
Cook-Waite Anesthetic Aspirating Syringe	SYRCW	SYRCW
18 Iris Scissors, Curved	S18	S18
IMS Bur Cushion, Short Lid	IMS-1372S	IMS-1372S
Hinged Instrument Clips, 2 Pack	IM1000	IM1000
Air/Water Clip	IM1005N	IM1005



Amalgam Standard Set-Up

8" x 11" x 1.25" (203 mm x 280 mm x 32 mm)

	Infinity Series™ Set-Up	Signature Series™ Set-Up
	IMNAMALGAM	IMAMALGAM
Large Cassette With Accessory Area - 12 Instrument, White	IMN4162	IM4162
#6 Cone Socket Mirror Handle	MH6	MH6
#5 Front Surface Rhodium Single- Sided Mirror (qty: 3)	MIR5/3	MIR5/3
#5 DE Explorer, #6 Handle	EXD56	EXD56
18 Excavator, #6 Handle	EXC186	EXC186
Calcium Hydroxide Placer, #6 Handle	PICH6	PICH6
Large/ Regular Amalgam Carrier	AC5202	AC5202
0/1 Serrated Marquette Condenser, #6 Handle	PLG0/16	PLG0/16
1/2 Serrated Plugger, #6 Handle	PLG1/26	PLG1/26
3/6 Carver, #6 Handle	CD3/66	CD3/66
89/92 Carver, #6 Handle	CD89/926	CD89/926
1/2 Hollenback, #6 Handle	CVHL1/26	CVHL1/26
21B Burnisher, #6 Handle	BB21B6	BB21B6
27/29 Burnisher, #6 Handle	BB27/296	BB27/296
2 College Dressing Pliers	DP2	DP2
Miller Articulating Paper Forceps	APF2	APF2
Cook-Waite Anesthetic Aspirating Syringe	SYRCW	SYRCW
Amalgam Well	WA	WA
IMS Bur Cushion, Short Lid	IMS-1372S	IMS-1372S
Air/Water Clip	IM1005N	IM1005

**Sinus Lift Standard Set-Up**

8" x 11" x 1.25" (203 mm x 280 mm x 32 mm)

	Infinity Series™ Set-Up	Signature Series™ Set-Up
	IMNSINUS	IMSINUS
Large Cassette - 12 Instrument, Gray	IMN4121	IM4121
Hinged Instrument Clips, 2 Pack	IM1000	IM1000
Buser Periosteal, #6 Handle	PPBUSER6	PPBUSER6
Palti Bone Graft Scoop	IMP6522	IMP6522
Bone Tamper	BONT	BONT
Bone Rasp	RASP	RASP
Sinus Retractor Misch 10	SRM10	SRM10
1 Sinus Curette	SINC1	SINC1
2 Sinus Curette	SINC2	SINC2
Bone Grafting Syringe	SYBG	SYBG
#1 Mead Immunity Mallet	MAL1	MAL1
Crile-Wood Perma Sharp® Needle Holder	NH5038	NH5038
2/4 Molt DE Curette	CM2/4	CM2/4
#UNC12 Colorvue™ Probe Kit (1-2-3-4-5- 6-7-8-9-10-11-12)	PCVNCKIT6	PCVNCKIT6
23 Seldin Periosteal Retractor	P23	P23
Bone Packer DE	BPACK	BPACK



Ortho Banding/Bonding Standard Set-Up 8" x 11" x 1.25" (203 mm x 280 mm x 32 mm)	IMORTHOBB
Double-Decker™ Large Orthodontic Cassette, Orange	IM5203-OR
#6 Cone Socket Mirror Handle	MH6
#5 Front Surface Rhodium Single-Sided Mirror (qty: 3)	MIR5/3
Ligature Director/Band Pusher	678-912
1/2 Hollenback, #6 Handle	CVHL1/26
#18 Dressing Pliers	DP18
Universal Cut & Hold Distal End Cutter, Long Handle	678-101L
Micro Mini Pin & Ligature Cutter	678-107
Weingart Pliers, Long Handle	678-201L
Direct Bond Bracket Holder	678-212
Mathieu Pliers	678-330

Available in Signature Series only



Ortho Archwire Adjustment Standard Set-Up 6" x 8" x 1.5" (152 mm x 203 mm x 38 mm)	IMORTHOAA
Double-Decker™ Large Orthodontic Cassette, Green	IM5209-OR
#6 Cone Socket Mirror Handle	MH6
#5 Front Surface Rhodium Single-Sided Mirror (qty: 3)	MIR5/3
Ligature Director/Scaler	678-905
Universal Cut & Hold Distal End Cutter, Long Handle	678-101L
Hard Wire Cutter, 15°	678-103
Micro Mini Pin & Ligature Cutter	678-107
Weingart Pliers, Long Handle	678-201L
Direct Bond Bracket Holder	678-212
Bird Beak Pliers	678-304
Three Jaw Pliers	678-302
Mathieu Pliers	678-330

Available in Signature Series only



Basic Extraction Standard Set-Up 8" x 11" x 1.25" (203 mm x 280 mm x 32 mm)	IMEXTRACTB
Special Large Signature Series 16 Instrument Cassette, 2 Elevator Rails, 4 Posted Hinged Instrument Clips, Orange	IM4163E2C4
Posted Hinged Instrument Clips (2 packs of 2)	I1000
Cone Socket Mirror Handle, Satin Steel Handle	MH6
4 Front Surface Cone Socket Mirror, Single	MIR4
5 Round Scalpel Handle	10-130-05
9 Mold Periosteal Elevator, Satin Steel Handle	P96
13/14 Heidbrink DE Root Pick, Satin Steel Handle	EHB13/146
86 Lucas Surgical Curette, Satin Steel Handle	CL866
1X Miller-Colburn Bone File	BF1X6
2 Dressing Pliers (college)	DP2
Retractor U of Minnesota	CRM
Crile-Wood Needle Holder	NHCW
34 Seldin Elevator	E34S
301 Apexo Elevator	E301
42 Adson 1/X Tissue Pliers	TP42
Anesthetic Syringe Cook-Waite	SYRCW
1 Kelly Scissor, Curved	S1
Child Mouth Prop (1 pack of 2)	MP52

Available in Signature Series only



Full Extraction Standard Set-Up

8" x 14.5" x 1.25" (203 mm x 368 mm x 32 mm)

IMEXTRACTF

Special Signature Series Oral Surgery
Cassette, 2 Elevator Rails, 4 Posted Hinged
Instrument Clips, 6 Low Profile Rails,
Lavender

IMOS6E2C4L

Posted Hinged Instrument Clips
(2 packs of 2)

I1000

Cone Socket Mirror Handle,
Satin Steel Handle

MH6

4 Front Surface Cone Socket Mirror, Single

MIR4

5 Round Scalpel Handle

10-30-05

9 Mold Periosteal Elevator,
Satin Steel Handle

P96

13/14 Heidbrink DE Root Pick,
Satin Steel Handle

EHB13/146

86 Lucas Surgical Curette,
Satin Steel Handle

CL866

2/4 Molt DE Curette, Satin Steel Handle

CM2/46

1X Miller-Colburn Bone File

BF1X6

2 Dressing Pliers (college)

DP2

Retractor U of Minnesota

CRM

Crile-Wood Perma Sharp Needle Holder

NHCW

34S Seldin Elevator

E34S

E301 Apexo Elevator

E301

42 Adson 1/X Tissue Pliers

TP42

Anesethetic Syringe Cook-Waite

SYRCW

1 Kelly Scissor, Curved

S1

Child Mouth Prop (1 pack of 2)

MP52

15 Cryer Forceps

F150

151 Cryer Forceps

F151

Available in Signature Series only



ADA CDR Coding Collection Series

We have created a collection of procedural instrument kits that follow American Dental Association (ADA) — Current Dental Terminology (CDT) Coding. Start pairing patients and insurance codes more effectively to move from standardized care to personalized care. Each kit, uniquely color-coded, offers the highest quality, technology, and ergonomic instrumentation within a small, space-efficient cassette for ensured equipment compatibility and safety.



Evaluation Kit 3" x 8" x 1.25" (76 mm x 203 mm x 32 mm)	IM-D0120/80
Exam Cassette - 5 Instrument, Blue	IM6058
#6 Cone Socket Mirror Handle	MH6
#5 HD Double-Sided Mirror	MIR5DSHD
#5 DE Explorer, #6 Handle	EXD56
Qulix SE UNC Probe, 1-12, #6 Handle	PCPUNC126
#11/12 ODU After Five DE Explorer, #6 Handle	EXD11/12A6
#2 Color Coded Nabers DE Probe 3-6-9-12, #6 Handle	PQ2N6



Prophylaxis Kit 3" x 8" x 1.25" (76 mm x 203 mm x 32 mm)	IM-D1110
Exam Cassette - 5 Instrument, Green	IM6059
Plastic Test Stick	PTS-IND
Nevi™ 2 Posterior Sickle Scaler, #9 Handle	SCNEVI29E2
Nevi™ 4 Posterior Sickle Scaler, #9 Handle	SCNEVI49E2
Barnhart 5/6 Curette, #9 Handle	SBH5/69E2
Gracey 15/16 Curette, #9 Handle	SG15/169E2
Gracey 13/14 Standard Curette, #9 Handle	SG13/149E2



Scaling And Root Planing Kit 3" x 8" x 1.5" (76 mm x 203 mm x 38 mm)	IM-D4341/42
Double-Decker™ Cassette - 8 Instrument, Red	IM9080
Plastic Test Stick	PTS-IND
Langer 3/4 After Five™ Curette, #9 Handle	SL3/4AF9E2
Nevi™ 2 Posterior Sickle Scaler, #9 Handle	SCNEVI29E2
Nevi™ 4 Posterior Sickle Scaler, #9 Handle	SCNEVI49E2
Gracey 1/2 Mini Five™ Curette, #9 Handle	SAS1/29E2
Gracey 7/8 After Five™ Curette, #9 Handle	SRPG7/89E2
Gracey 15/16 After Five™ Curette, #9 Handle	SRPG15/169E2
Gracey 13/14 After Five™ Curette, #9 Handle	SRPG13/149E2
Gracey 17/18 Curette, #9 Handle	SG17/189E2



Pedo Propy Kit

3" x 8" x 1.25" (76 mm x 203 mm x 32 mm)

IM-D1120

Exam Cassette - 5 Instrument, Lavender	IM6056
#6 Cone Socket Mirror Handle	MH6
#5 HD Double-Sided Mirror	MIR5DSHD
#5 DE Explorer, #6 Handle	EXD56
Gracey 13/14 Standard Curette, #9 Handle	SC13/149E2
Nevi™ 3 Posterior Sickle Scaler, #9 Handle	SCNEVI39E2
Nevi™ 1 Anterior Sickle Scaler, #9 Handle	SCNEVI19E2
Plastic Test Stick	PTS-IND



Presence Of Gingival Inflammation Kit

3" x 8" x 1.25" (76 mm x 203 mm x 32 mm)

IM-D4346

Exam Cassette - 5 Instrument, Yellow	IM6055
Plastic Test Stick	PTS-IND
Nevi™ 3 Posterior Sickle Scaler, #9 Handle	SCNEVI39E2
Nevi™ 4 Posterior Sickle Scaler, #9 Handle	SCNEVI49E2
Gracey 1/2 Mini Five™ Curette, #9 Handle	SAS1/29E2
Gracey 15/16 After Five™ Curette, #9 Handle	SRPG15/169E2
Gracey 13/14 After Five™ Curette, #9 Handle	SRPG13/149E2



Debridement For Comprehensive Evaluation Kit

3" x 8" x 1.25" (76 mm x 203 mm x 32 mm)

IM-D4355

Exam Cassette - 5 Instrument, Orange	IM6053
Plastic Test Stick	PTS-IND
Nevi™ 4 Posterior Sickle Scaler, #9 Handle	SCNEVI4E92
McCall 17/18 Curette, #9 Handle	SM17/189E2
Gracey 1/2 Mini Five™ Curette, #9 Handle	SAS1/2R9E2
Gracey 15/16 After Five™ Curette, #9 Handle	SRP15/16R9E2
Gracey 13/14 After Five™ Curette, #9 Handle	SRP13/14R9E2



Periodontal Maintenance Kit	IM-D4910
3" x 8" x 1.5" (76 mm x 203 mm x 38 mm)	
Double-Decker™ Cassette - 8 Instrument, Gray	IM9081
Plastic Test Stick	PTS-IND
Langer 3/4 Curette, #9 Handle	SL3/49E2
Nevi™ 2 Posterior Sickle Scaler, #9 Handle	SCNEVI29E2
Nevi™ 4 Posterior Sickle Scaler, #9 Handle	SCNEVI49E2
Pattison Gracey Lite 1/2 Curette, #9 Handle, Gray	SG1/29APE2
Pattison Gracey Lite 7/8 Curette, #9 Handle, Green	SG7/89APE2
Pattison Gracey Lite 11/12 Curette, #9 Handle, Orange	SG11/129APE2
Pattison Gracey Lite 13/14 Curette, #9 Handle, Blue	SG13/149APE2
Gracey 15/16 Curette, #9 Handle	SG15/169E2



Implant Maintenance Kit	IM-D6080
3" x 8" x 1.25" (76 mm x 203 mm x 32 mm)	
Exam Cassette - 5 Instrument, White	IM6052
Colorvue™ Probe Kit 12 University of North Carolina	PCVPROMONC
Titanium Implant Scaler 204SD	IMP204SDT
Titanium Implant Scaler Langer 1/2	IMPLG1/2T
Titanium Implant Scaler Mini Five™ 1/2	IMPM1/2T
Implacare™ II H6/7 Scaler	IMPLH6/76
Plastic Test Stick	PTS-IND



Procedural Tubs

Together, procedure tubs and cassettes are an efficient and organized system for managing instruments and disposables by procedure. The time savings achieved by using IMS™ Cassettes and procedural tubs allows staff to spend more time on value added activities in the facility. Tub colors are designed to match the IMS™ Cassette rail colors and are treated with Sanafor® P0-5 antimicrobial* to protect the product by helping prevent the growth of stain and odor causing bacteria.

Organizational Tubs

Organizational tubs are available as complete tubs with accessories or can be ordered separately. Tubs available in 11 colors to match IMS™ Cassette rail colors.



Signature Series™ Complete Tub	
Each Tub includes: Tub Liner, Dividers, Single Tub Cup With Cover, Double Tub Cup With Cover, Divided Slide Tray, Locking Tub Cover 12-3/4" X 11-3/16" X 4" (323.85 mm x 284.162 mm x 101.6 mm)	
● Red	IMS-1400
● Gray	IMS-1401
○ White	IMS-1402
● Orange	IMS-1403
● Purple	IMS-1404
● Yellow	IMS-1405
● Lavender	IMS-1406
● Black	IMS-1407
● Blue	IMS-1408
● Green	IMS-1409
● Ocean Blue	IMS-14011

*These antimicrobial properties are built in to protect the product. This product does not protect users against bacteria, viruses, germs or other disease organisms. Clear components and tub liner do not contain antimicrobial product protection.



Signature Series™ Tub <i>Tub available in 11 color options. Lid ordered separately.</i> 12-1/4" x 10-7/8" x 2-3/4" (311.15 mm x 276.225 mm x 69.85 mm)	
 Red	IMS-1410
 Gray	IMS-1411
 White	IMS-1412
 Orange	IMS-1413
 Purple	IMS-1414
 Yellow	IMS-1415
 Lavender	IMS-1416
 Black	IMS-1417
 Blue	IMS-1418
 Green	IMS-1419
 Ocean Blue	IMS-14111

Quick Tip: Red Signature Series™ Tubs can also be used as a solution to transport your contaminated instruments and small cassettes.



Locking Tub Cover	IMS-1420
<i>Locking tub cover protects items from airborne materials during use, transportation, and storage.</i> 12-3/4" x 11-3/16" x 2" (323.85 mm x 284.162 mm x 50.8 mm)	

Tub Racks



Adjustable 3 Tub Rack	IMS-1433
<i>Fits most standard storage cabinets. Holds a maximum of 3 tubs with lids or 6 cassettes. Non-slip feet. Sturdy metal construction. No screws, tools nuts or bolts needed.</i> Min 7-1/2" Max 15-1/2" x 10-1/2" x 13" (Min 190.5 mm Max 393.7 mm x 266.7 mm x 330.2 mm)	



Adjustable 4 Tub Rack	IMS-1434
<i>Holds a maximum of 4 tubs with lids or 8 cassettes. Sturdy metal construction. Non-slip feet.</i> Min 7-1/2" Max 15-1/2" x 10-1/2" x 18-3/4" (Min 190.5 mm Max 393.7 mm x 266.7 mm x 476.25 mm)	



Tub Accessories

Several accessories are available to support organization of disposables needed for each procedure.



Tub Insert

Features 10 compartments for easy organization of tub materials.
Sits on rim of tub allowing for storage underneath.

9-3/4" x 10-3/4" x 11/32" (247.65 mm x 273.05 mm x 0.0129 mm)

IMS-1423



Divided Slide Tray

Stores cotton rolls, cotton squares, clamps, syringe tips and more.

9-3/4" x 5-3/8" x 11/16" (247.65 mm x 136.525 mm x 17.463 mm)

IMS-1421




Flat Slide Tray

Can be used to store discs, strips, tofflemire, articulating paper, etc.

9-3/4" x 5-3/16" x 1-3/8" (247.65 mm x 131.762 mm x 34.925 mm)

IMS-1422



Dividers

Includes 8 dividers.

12-3/4" x 11-3/16" X 1-3/8" (323.85 mm x 284.162 mm x 34.925 mm)

IMS-1427



Tub Liner

Non-slip liner prevents small items from slipping under divider sections.

9-5/8" x 8-7/8" X 1/8" (244.475 mm x 225.425 mm x 3.175 mm)

IMS-1428

**Capsule Clip-On Holder**

Allows quick access to composite material during procedure. Conveniently clips on any side of tub. Holds up to 8 capsules, 4 syringe tubes, and capsule/gun dispenser.

8" x 2-1/2" x 3/4" (203.2 mm x 63.5 mm x 19.05 mm)

IMS-1429

**Syringe Clip-On Holder**

Allows quick access to composite material during procedure. Conveniently clips on any side of tub. Holds up to 15 syringe tubes and 3 applicator tips.

8" x 2-5/8" x 3/4" (203.2 mm x 66.675 mm x 19.05 mm)

IMS-1430

**10-Unit Syringe Add On**

An add-on for the Syringe 20 unit stand. Provides additional space for syringe storage. Holds up to 10 syringe tubes.

8-1/2" x 2-3/8" x 3/4" (215.9 mm x 60.325 mm x 19.05 mm)

IMS-1432

**20-Unit Syringe Stand**

Stands up for easy access to syringes. Lays flat to fit inside a tub or drawer. Holds up to 20 syringe tubes.

9" x 4-3/4" x 2-1/2" (228.6 mm x 120.65 mm x 63.5 mm)

IMS-1431

Tub Cups & Covers

Used to store brushes, clamps, cotton pellets, syringe tips, prophy paste, prophy angle tips, pumice, implant tools, etc. Long tub Cup Cover provides five storage compartments. Dividers are removable for customization. Perfect for storing capsules, discs, abrasive strips, wedges, matrix bands, etc.

**Single Tub Cup With Cover**

2-1/4" x 1-5/8" x 1-1/8"
(57.15 mm x 41.275 mm
x 28.575 mm)

IMS-1424

**Double Tub Cup With Cover**

3-3/8" x 2-1/4" x 1-1/8"
(85.725 mm x 57.15 mm
x 28.575 mm)

IMS-1425

**Long Tub Cup With Cover**

8-1/3" x 2-1/4" x 1-1/8"
(211.582 mm x 57.15 mm
x 28.575 mm)

IMS-1426