


Handle Options	A102	Specialty Scalers	A129
Procedural Kits	A104	Scaler and Curette Sets	A130
Aerosol Mitigation	A104	McCall Scalers and Curettes	A130
Gracey Curettes	A105	Goldman-Fox Scalers and Curettes	A131
Standard Gracey Curettes	A105	Cattoni Scalers	A132
Gracey After Five™ Curettes	A108		
Gracey Mini Five™ Curettes	A110	Files and Hoes	A133
Gracey Micro Mini Five™ Curettes	A112	Hoe Scalers	A133
Pattison Gracey Lite Curettes	A113	File Scalers	A134
Vision Curvettes	A114	Chisels	A135
BioGent Curettes	A115	DiamondTec™ File Scalers	A135
Universal Curettes	A116	Empowerment Kit	A136
Langer Curettes	A116	Empowerment Scalers and Curettes	A136
Columbia Curettes	A118		
Barnhart Curettes	A119		
Furcation Curettes	A120		
Standard Universal Curettes	A121		
Sickle Scalers	A123		
Nevi™ Scalers	A123		
Anterior Sickle Scalers	A124		
Posterior Sickle Scalers	A126		
Jacquette Sickle Scalers	A128		



Hand Scaling

Our hand scalers are made from the finest grade stainless steel to ensure superior performance and longevity. With the largest selection of tip designs in the industry, we also offer a variety of innovative and ergonomic handle designs including the new line of Harmony™ Ergonomic Scalers and Curettes designed with TrueFit™ Technology. And, whether your choice is an EverEdge™ 2.0 Blade with superior edge retention characteristics or the classic Immunity Steel™ Blade, you can be assured that your instruments are meticulously crafted to allow you to be the Best in Practice.



Handle Options

We offer a wide variety of handle options for our hand scaling instruments giving clinicians a choice between styles, materials, and textures. Take a look below at our EverEdge™ 2.0 Scaler and Curette Handles as well as our Immunity Steel™ Scalers and Curettes.



Learn more about
[EverEdge™ 2.0 Technology.](#)

EverEdge™ 2.0 Handles

Harmony™ Ergonomic Handle



#9 Handle



#C8 ResinEight™ Colors Handle



Immunity Steel™ Handles

#6 Satin Steel™ Handle



#7 Satin Steel™ Colours Handle



#8 ResinEight™ Handle



#2 Octagonal Handle





HuFriedyGroup Introduces **HARMONY**TM

SCIENTIFICALLY PROVEN ERGONOMIC SCALERS AND CURETTES

Knowing nearly 70% of RDHs experience hand, wrist, and/or arm pain,¹ HuFriedyGroup proactively developed an innovative and truly ergonomic hand scaling solution. Building off the success of EverEdgeTM 2.0, the NEW HarmonyTM Ergonomic Scalers and Curettes designed with TrueFitTM Technology are the latest addition to the best-in-class hand instrument portfolio.



Learn How to Scale in Perfect Harmony at
HuFriedyGroup.com/Harmony

1. Compared to other leading scaler designs. Data on file. Available upon request.



Procedural Kits

Procedural kits are a compilation of clinicians' most popular procedural set-ups designed for ease of ordering.

Aerosol Mitigation

Many studies have shown that the same clinical outcomes are achieved with hand scaling as with ultrasonics, allowing you to provide enhanced patient care while mitigating aerosol risks. If you are prioritizing hand scaling take advantage of our Aerosol Mitigation Kit, available with or without the HVESolo™ Disposable Evacuation Tips.



Perio Recovery Kit	RECPERIO-KIT
H6/7 Sickle Scaler, #9 Handle	SH6/79E2
204S Posterior Sickle Scaler, #9 Handle	S204S9E2
Barnhart 5/6 Curette, #9 Handle	SBH5/69E2
Gracey 11/12 Standard Curette, #9 Handle	SG11/129E2
Gracey 13/14 Standard Curette, #9 Handle	SG13/149E2

Perio Recovery Kit with HVESolo™ Evacuation Tips	RECPSOLO-KIT
H6/7 Sickle Scaler, #9 Handle	SH6/79E2
204S Posterior Sickle Scaler, #9 Handle	S204S9E2
Barnhart 5/6 Curette, #9 Handle	SBH5/69E2
Gracey 11/12 Standard Curette, #9 Handle	SG11/129E2
Gracey 13/14 Standard Curette, #9 Handle	SG13/149E2
HVESolo™ Disposable Evacuation Tips (2.5" length)	1808S



Gracey Curettes

Gracey curettes are area-specific curettes with blades having a single cutting edge, offset at 70° to the shank, to provide a perfect working angulation to the tooth surface.



Learn more about our
[Gracey Curettes](#).

Standard Gracey Curettes

Standard Gracey Curettes allow for deep scaling, root planing and periodontal debridement. Gracey Rigid curettes have thicker, stronger and less flexible shanks than Standard Graceys to allow for removal of moderate to heavy calculus. Gracey Extra Rigid Curettes are ideal for tenacious calculus removal. Standard, Rigid and Extra Rigid Gracey Curettes all have identical blade dimensions.



Gracey 1/2 Curette	Standard	Rigid
Harmony	SG1/2XE2	SG1/2RXE2
#9	SG1/29E2	SG1/2R9E2
#C8, Gray	SG1/2C8E2	SG1/2RC8E2
#6	SG1/26	SG1/2R6
#7	SG1/27	SG1/2R7
#8	SG1/28	SG1/2R8
#2	SG1/2	SGR1/2R



Gracey 3/4 Curette	Standard	Rigid
Harmony	SG3/4XE2	SG3/4RXE2
#9	SG3/49E2	SG3/4R9E2
#C8, Gray	SG3/4C8E2	SG3/4RC8E2
#6	SG3/46	SG3/4R6
#7	SG3/47	SG3/4R7
#8, Gray	SG3/48	SG3/4R8
#2	SG3/4	SGR3/4R



Gracey 5/6 Curette	Standard	Rigid	Extra Rigid
Harmony	SG5/6XE2	SG5/6RXE2	N/A
#9	SG5/69E2	SG5/6R9E2	N/A
#C8, Yellow	SG5/6C8E2	SG5/6RC8E2	N/A
#6	SG5/66	SG5/6R6	SG5/6ER6
#7	SG5/67	SG5/6R7	SG5/6ER7
#8, Gray	SG5/68	SG5/6R8	SG5/6ER8
#2	SG5/6	SGR5/6R	N/A



Gracey Currettes • Standard Gracey Currettes



Gracey 7/8 Curette	Standard	Rigid	Extra Rigid
Harmony	SG7/8XE2	SG7/8RXE2	N/A
#9	SG7/89E2	SG7/8R9E2	N/A
#C8, Green	SG7/8C8E2	SG7/8RC8E2	N/A
#6	SG7/86	SG7/8R6	SG7/8ER6
#7	SG7/87	SG7/8R7	SG7/8ER7
#8, Gray	SG7/88	SG7/8R8	SG7/8ER8
#2	SG7/8	SGR7/8R	N/A



Gracey 9/10 Curette	Standard	Rigid
Harmony	SG9/10XE2	SG9/10RXE2
#9	SG9/109E2	SG9/10R9E2
#C8, Green	SG9/10C8E2	N/A
#6	SG9/106	SG9/10R6
#7	SG9/107	SG9/10R7
#8, Gray	SG9/108	SG9/10R8
#2	SG9/10	SGR9/10R



Gracey 11/12 Curette	Standard	Rigid	Extra Rigid
Harmony	SG11/12XE2	SG11/12RXE2	N/A
#9	SG11/129E2	SG11/12R9E2	N/A
#C8, Orange	SG11/12C8E2	SG11/12RC8E2	N/A
#6	SG11/126	SG11/12R6	SG11/12ER6
#7	SG11/127	SG11/12R7	SG11/12ER7
#8, Gray	SG11/128	SG11/12R8	SG11/12ER8
#2	SG11/12	SGR11/12R	N/A



Gracey 13/14 Curette	Standard	Rigid	Extra Rigid
Harmony	SG13/14XE2	SG13/14RXE2	N/A
#9	SG13/149E2	SG13/14R9E2	N/A
#C8, Blue	SG13/14C8E2	SG13/14RC8E2	N/A
#6	SG13/146	SG13/14R6	SG13/14ER6
#7	SG13/147	SG13/14R7	SG13/14ER7
#8, Gray	SG13/148	SG13/14R8	SG13/14ER8
#2	SG13/14	SGR13/14R	N/A



Gracey 15/16 Curette	Standard	Rigid	Extra Rigid
Harmony	SG15/16XE2	SG15/16RXE2	N/A
#9	SG15/169E2	SG15/16R9E2	N/A
#C8, Orange	SG15/16C8E2	SG15/16RC8E2	N/A
#6	SG15/166	SG15/16R6	SG15/16ER6
#7	SG15/167	SG15/16R7	N/A
#8, Gray	SG15/168	SG15/16R8	SG15/16ER8
#2	SG15/16	SGR15/16R	N/A



Gracey 17/18 Curette	Standard	Rigid
Harmony	SG17/18XE2	SG17/18RXE2
#9	SG17/189E2	SG17/18R9E2
#C8, Blue	SG17/18C8E2	SG17/18RC8E2
#6	SG17/186	SG17/18R6
#7	SG17/187	SG17/18R7
#8, Gray	SG17/188	SG17/18R8
#2	SG17/18	SGR17/18R



Gracey 11/14 Curette	Standard	Rigid
Harmony	SG11/14XE2	N/A
#9	SG11/149E2	N/A
#C8, Red	SG11/14C8E2	N/A
#6	SG11/146	N/A
#7	SG11/147	N/A
#8, Gray	SG11/148	N/A
#2	SG11/14	SG11/14R



Gracey 12/13 Curette	Standard	Rigid
Harmony	SG12/13XE2	N/A
#9	SG12/139E2	N/A
#C8, Red	SG12/13C8E2	N/A
#6	SG12/136	N/A
#7	SG12/137	N/A
#8, Gray	SG12/138	N/A
#2	SG12/13	SG12/13R



Gracey After Five™ Currettes

Gracey After Five™ Currettes feature a terminal shank that is 3mm longer and a blade that is 10% thinner than the Standard Gracey Curette. The thinner blade allows for less tissue distention when accessing deep pockets.



Gracey 1/2 After Five™ Curette	Standard	Rigid
Harmony	SRPG1/2XE2	SRP1/2RXE2
#9	SRPG1/29E2	SRP1/2R9E2
#C8, Gray	SRPG1/2C8E2	SRP1/2RC8E2
#6	SRPG1/26	SRP1/2R6
#7	SRPG1/27	SRP1/2R7
#8, Gray	SRPG1/28	SRP1/2R8
#2	SRPG1/2	SRPG1/2R



Gracey 3/4 After Five™ Curette	Standard	Rigid
Harmony	SRPG3/4XE2	SRP3/4RXE2
#9	SRPG3/49E2	SRP3/4R9E2
#C8, Gray	SRPG3/4C8E2	SRP3/4RC8E2
#6	SRPG3/46	SRP3/4R6
#7	SRPG3/47	SRP3/4R7
#8, Gray	SRPG3/48	SRP3/4R8
#2	SRPG3/4	SRPG3/4R



Gracey 5/6 After Five™ Curette	Standard	Rigid	Extra Rigid
Harmony	SRPG5/6XE2	SRP5/6RXE2	N/A
#9	SRPG5/69E2	SRP5/6R9E2	N/A
#C8, Yellow	SRPG5/6C8E2	SRP5/6RC8E2	N/A
#6	SRPG5/66	SRP5/6R6	SRP5/6E6
#7	SRPG5/67	SRP5/6R7	SRP5/6E7
#8, Gray	SRPG5/68	SRP5/6R8	SRP5/6E8
#2	SRPG5/6	SRPG5/6R	N/A



Gracey 7/8 After Five™ Curette	Standard	Rigid	Extra Rigid
Harmony	SRPG7/8XE2	SRP7/8RXE2	N/A
#9	SRPG7/89E2	SRP7/8R9E2	N/A
#C8, Green	SRPG7/8C8E2	SRP7/8RC8E2	N/A
#6	SRPG7/86	SRP7/8R6	SRP7/8E6
#7	SRPG7/87	SRP7/8R7	SRP7/8E7
#8, Gray	SRPG7/88	SRP7/8R8	SRP7/8E8
#2	SRPG7/8	SRPG7/8R	N/A



Gracey 11/12 After Five™ Curette	Standard	Rigid	Extra Rigid
Harmony	SRPG11/12XE2	SRP11/12RXE2	N/A
#9	SRPG11/129E2	SRP11/12R9E2	N/A
#C8, Orange	SRPG11/12C8E2	SRP11/12RC8E2	N/A
#6	SRPG11/126	SRP11/12R6	SRP11/12E6
#7	SRPG11/127	SRP11/12R7	SRP11/12E7
#8, Gray	SRPG11/128	SRP11/12R8	SRP11/12E8
#2	SRPG11/12	SRPG11/12R	N/A



Gracey 13/14 After Five™ Curette	Standard	Rigid	Extra Rigid
Harmony	SRPG13/14XE2	SRP13/14RXE2	N/A
#9	SRPG13/149E2	SRP13/14R9E2	N/A
#C8, Blue	SRPG13/14C8E2	SRP13/14RC8E2	N/A
#6	SRPG13/146	SRP13/14R6	SRP13/14E6
#7	SRPG13/147	SRP13/14R7	SRP13/14E7
#8, Gray	SRPG13/148	SRP13/14R8	SRP13/14E8
#2	SRPG13/14	SRPG13/14R	N/A



Gracey 15/16 After Five™ Curette	Standard	Rigid
Harmony	SRPG15/16XE2	SRP15/16RXE2
#9	SRPG15/169E2	SRP15/16R9E2
#C8, Orange	SRPG15/16C8E2	SRP15/16RC8E2
#6	SRPG15/166	SRP15/16R6
#7	SRPG15/167	SRP15/16R7
#8, Gray	SRPG15/168	SRP15/16R8
#2	SRPG15/16	SRPG15/16R



Gracey 11/14 After Five™ Curette	
Harmony	SRPG11/14XE2
#9	SRPG11/149E2
#C8, Red	SRPG11/14C8E2
#6	SRPG11/146
#7	SRPG11/147
#2	SRPG11/14



Gracey 12/13 After Five™ Curette	
Harmony	SRPG12/13XE2
#9	SRPG12/139E2
#C8, Red	SRPG12/13C8E2
#6	SRPG12/136
#7	SRPG12/137
#8, Gray	SRPG12/138
#2	SRPG12/13



Gracey Mini Five™ Currettes

Gracey Mini Five™ Currettes feature a terminal shank that is 3mm longer and a blade that is 50% shorter and 10% thinner than a Standard Gracey Curette. Gracey Mini Five™ Currettes are ideal for deep, narrow pockets.

Hand Scaling



Gracey 1/2 Mini Five™ Curette	Standard	Rigid
Harmony	SAS1/2XE2	SAS1/2RXE2
#9	SAS1/29E2	SAS1/2R9E2
#C8, Gray	SAS1/2C8E2	SAS1/2RC8E2
#6	SAS1/26	SAS1/2R6
#7	SAS1/27	SAS1/2R7
#8, Gray	SAS1/28	SAS1/2R8
#2	SAS1/2	SAS1/2R



Gracey 3/4 Mini Five™ Curette	Standard	Rigid
Harmony	SAS3/4XE2	SAS3/4RXE2
#9	SAS3/49E2	SAS3/4R9E2
#C8, Gray	SAS3/4C8E2	SAS3/4RC8E2
#6	SAS3/46	SAS3/4R6
#7	SAS3/47	SAS3/4R7
#8, Gray	SAS3/48	SAS3/4R8
#2	SAS3/4	SAS3/4R



Gracey 5/6 Mini Five™ Curette	Standard	Rigid
Harmony	SAS5/6XE2	SAS5/6RXE2
#9	SAS5/69E2	SAS5/6R9E2
#C8, Yellow	SAS5/6C8E2	SAS5/6RC8E2
#6	SAS5/66	SAS5/6R6
#7	SAS5/67	SAS5/6R7
#8, Gray	SAS5/68	SAS5/6R8
#2	SAS5/6	SAS5/6R



Gracey 7/8 Mini Five™ Curette	Standard	Rigid
Harmony	SAS7/8XE2	SAS7/8RXE2
#9	SAS7/89E2	SAS7/8R9E2
#C8, Green	SAS7/8C8E2	SAS7/8RC8E2
#6	SAS7/86	SAS7/8R6
#7	SAS7/87	SAS7/8R7
#8, Gray	SAS7/88	SAS7/8R8
#2	SAS7/8	SAS7/8R



Gracey 11/12 Mini Five™ Curette	Standard	Rigid
Harmony	SAS11/12XE2	SAS11/12RXE2
#9	SAS11/129E2	SAS11/12R9E2
#C8, Orange	SAS11/12C8E2	SAS11/12RC8E2
#6	SAS11/126	SAS11/12R6
#7	SAS11/127	SAS11/12R7
#8, Gray	SAS11/128	SAS11/12R8
#2	SAS11/12	SAS11/12R



Gracey 13/14 Mini Five™ Curette	Standard	Rigid
Harmony	SAS13/14XE2	SAS13/14RXE2
#9	SAS13/149E2	SAS13/14R9E2
#C8, Blue	SAS13/14C8E2	SAS13/14RC8E2
#6	SAS13/146	SAS13/14R6
#7	SAS13/147	SAS13/14R7
#8, Gray	SAS13/148	SAS13/14R8
#2	SAS13/14	SAS13/14R



Gracey 15/16 Mini Five™ Curette	Standard	Rigid
Harmony	SAS15/16XE2	SAS15/16RXE2
#9	SAS15/169E2	SAS15/16R9E2
#C8, Orange	SAS15/16C8E2	SAS15/16RC8E2
#6	SAS15/166	SAS15/16R6
#7	SAS15/167	SAS15/16R7
#8, Gray	SAS15/168	SAS15/16R8
#2	SAS15/16	SAS15/16R



Gracey 11/14 Mini Five™ Curette	
Harmony	SAS11/14XE2
#9	SAS11/149E2
#C8, Red	SAS11/14C8E2
#6	SAS11/146
#7	SAS11/147
#8, Gray	SAS11/148
#2	SAS11/14



Gracey 12/13 Mini Five™ Curette	
Harmony	SAS12/13XE2
#9	SAS12/139E2
#C8, Red	SAS12/13C8E2
#6	SAS12/136
#7	SAS12/137
#8, Gray	SAS12/138
#2	SAS12/13

Gracey Micro Mini Five™ Cures

Gracey Micro Mini Five™ Cures feature an elongated shank to access deep, narrow periodontal pockets and a blade that is 50% shorter and 28% thinner than the Standard Gracey Curette.



Gracey 1/2 Micro Mini Five™ Curette	
Harmony	SMS1/2XE2
#9	SMS1/29E2



Gracey 7/8 Micro Mini Five™ Curette	
Harmony	SMS7/8XE2
#9	SMS7/89E2



Gracey 11/12 Micro Mini Five™ Curette	
Harmony	SMS11/12XE2
#9	SMS11/129E2



Gracey 13/14 Micro Mini Five™ Curette	
Harmony	SMS13/14XE2
#9	SMS13/149E2



Pattison Gracey Lite Currettes

The Anna Pattison Currettes are modified from the original Gracey Currettes. These new instruments have shorter, thinner, smaller blades honed at a 60° angle with a modified rigid shank specially designed to allow easier access and insertion under tight tissue.



Learn more about our key opinion leader [Anna Pattison](#).



Pattison Gracey Lite 1/2 Curette

#9, Gray

[SG1/29APE2](#)



Pattison Gracey Lite 7/8 Curette

#9, Green

[SG7/89APE2](#)



Pattison Gracey Lite 11/12 Curette

#9, Orange

[SG11/129APE2](#)



Pattison Gracey Lite 13/14 Curette

#9, Blue

[SG13/149APE2](#)



Pattison Gracey Lite 11/14 Curette

#9, Orange

[SG11/149APE2](#)



Pattison Gracey Lite 12/13 Curette

#9, Orange

[SG12/139APE2](#)



Vision Curvettes

Vision Curvettes are modified Gracey Cures that feature 5mm and 10mm shank markings to help visually assess pocket depth; and a blade identification mark (+) on the handle for easy determination of the cutting edge. The short, curved blades are ideal for root surfaces and furcations.

Hand Scaling



Sub-0 Vision Curvette
#6 **SGCSUB-06**



1/2 Vision Curvette
#6 **SGC1/26**



11/12 Vision Curvette
#6 **SGC11/126**



13/14 Vision Curvette
#6 **SGC13/146**



BioGent Currettes

Designed by Drs. Antonella Labriola and Pierpaolo Cortellini BioGent Currettes feature reduced blade width for non-traumatic and efficient access to the periodontium, modified angulation for better angulation and color coding for efficient time management.



Labriola-Cortellini Gracey 1/2 BioGent Curette

#9

SG1/29LCE2



Labriola-Cortellini Gracey 7/8 BioGent Curette

#9

SG7/89LCE2



Labriola-Cortellini Gracey 11/12 BioGent Curette

#9

SG11/129LCE2



Labriola-Cortellini Gracey 13/14 BioGent Curette

#9

SG13/149LCE2



Universal Currettes

Universal Currettes have two parallel cutting edges that meet at a rounded toe. The back surface is rounded giving a semi-circular cross section. Due to this rounded back surface, currettes are ideal for the removal of subgingival debris.



Learn more about our [Universal Currettes](#).

Langer Currettes

The Langer Curette design combines the benefits of a universal curette blade with two cutting edges and the shank designs of area specific currettes. The Langer Currettes were originally designed to be used in a set of three (1/2, 3/4, 5/6), to adapt to most tooth surfaces, and to be suitable for scaling as well as periodontal surgery.



Langer 1/2 Curette	Standard	Rigid
Harmony	SL1/2XE2	N/A
#9	SL1/29E2	N/A
#C8, Lavender	SL1/2C8E2	N/A
#6	SL1/26	SL1/2R6
#7	SL1/27	SL1/2R7
#8, Gray	SL1/28	N/A
#2	SL1/22	N/A

Langer 1/2 After Five™ Curette	
A 3 mm longer terminal shank and a thinner blade for improved access	
#9	SL1/2AF9E2
#6	SL1/2AF6

Langer 1/2 Mini Five™ Curette	
Longer terminal shank, standard diameter, and blade half the length of the After Five™ Currettes	
#9	SL1/2MF9E2
#6	SL1/2MF6



Langer 3/4 Curette	Standard	Rigid
Harmony	SL3/4XE2	N/A
#9	SL3/49E2	N/A
#C8, Lavender	SL3/4C8E2	N/A
#6	SL3/46	SL3/4R6
#7	SL3/47	SL3/4R7
#8, Gray	SL3/48	N/A
#2	SL3/42	N/A

Langer 3/4 Mini Five™ Curette

Longer terminal shank, standard diameter, and blade half the length of the After Five™ Curettes

#9	SL3/4MF9E2
#6	SL3/4MF6



Langer 5/6 Curette	Standard	Rigid
Harmony	SL5/6XE2	N/A
#9	SL5/69E2	N/A
#C8, Lavender	SL5/6C8E2	N/A
#6	SL5/66	SL5/6R6
#7	SL5/67	N/A
#8, Gray	SL5/68	N/A

Langer 5/6 Mini Five™ Curette

Longer terminal shank, standard diameter, and blade half the length of the After Five™ Curettes

#9	SL5/6MF9E2
#6	SL5/6MF6

Langer 5/6 After Five™ Curette

#6	SL5/6AF6
----	--------------------------



Langer 17/18 Curette	Standard	Rigid
Harmony	SL17/18XE2	N/A
#9	SL17/189E2	N/A
#C8, Lavender	SL17/18C8E2	N/A
#6	SL17/186	SL17/18R6
#7	SL17/187	N/A
#8, Gray	SL17/188	N/A



Columbia Currettes

Columbia Currettes feature two parallel cutting edges that meet at a rounded toe. Anterior, posterior and universal designs are available.



Columbia 2R/2L Curette

Harmony	SC2R/2LXE2
#9	SC2R/2L9E2
#C8, Lavender	SC2R/2LC8E2
#6	SC2R/2L6
#7	SC2R/2L7
#8, Gray	SC2R/2L8
#2	SC2R/2L



Columbia 4R/4L Curette

Harmony	SC4R/4LXE2
#9	SC4R/4L9E2
#C8, Lavender	SC4R/4LC8E2
#6	SC4R/4L6
#7	SC4R/4L7
#8, Gray	SC4R/4L8
#2	SC4R/4L



Columbia 13/14 Curette	Standard	Rigid
Harmony	SC13/14XE2	N/A
#9	SC13/149E2	N/A
#C8, Lavender	SC13/14C8E2	N/A
#6	SC13/146	SC13/14R6
#7	SC13/147	SC13/14R7
#8, Gray	SC13/148	SC13/14R8
#2	SC13/14	SC13/14R



Columbia 13/14 Curved Tip Curette

#2	SC13/14C
----	-----------------



Barnhart Curettes

Barnhart Universal Curettes are designed with two cutting edges and a rounded toe for use on all posterior surfaces. Barnhart Curettes have a less pronounced blade curvature than other universal curettes combined with a long shank. Available in standard, rigid and USC versions.

Hand Scaling



Barnhart 1/2 Curette	Standard	Rigid
Harmony	SBH1/2XE2	SBH1/2RXE2
#9	SBH1/29E2	SBH1/2R9E2
#C8, Lavender	SBH1/2C8E2	SBH1/2RC8E2
#6	SBH1/26	SBH1/2R6
#7	SBH1/27	SBH1/2R7
#8, Gray	SBH1/28	SBH1/2R8
#2	SBH1/2	SBH1/2R



Barnhart 5/6 Curette	Standard	Rigid
Harmony	SBH5/6XE2	SBH5/6RXE2
#9	SBH5/69E2	SBH5/6R9E2
#C8, Lavender	SBH5/6C8E2	SBH5/6RC8E2
#6	SBH5/66	SBH5/6R6
#7	SBH5/67	SBH5/6R7
#8, Gray	SBH5/68	SBH5/6R8
#2	SBH5/6	SBH5/6R



Barnhart USC 1/2 Curette	
Harmony	SBH1/2USCXE2
#9	SBH1/2USC9E2



Barnhart USC 5/6 Curette	
Harmony	SBH5/6USCXE2
#9	SBH5/6USC9E2



Furcation Currettes


Furcation Currettes are designed for scaling root concavities and furcations. The DeMarco Currettes are excellent for stain removal on anterior lingual surfaces.



#1 DeMarco Furcation Curette

#6


SDM16



#2 DeMarco Furcation Curette

#6

SDM26




Quetin Buccal-Lingual Furcation Curette

.09 mm

#6

SQBL16




Quetin Buccal-Lingual Furcation Curette

1.3 mm

#6

SQBL26




Quetin Mesial-Distal Furcation Curette

.09 mm

#6

SQMD16




Quetin Mesial-Distal Furcation Curette

1.3 mm

#6

SQMD26



Quetin Labial-Lingual Scaler

#6

SQLL36




Standard Universal Curettes



Younger-Good 7/8 Curette

Harmony	SYG7/8XE2
#9	SYG7/89E2
#C8, Lavender	SYG7/8C8E2
#6	SYG7/86
#7	SYG7/87
#8, Gray	SYG7/88
#2	SYG7/8



Indiana University 13/14 Curette

#6	SIU13/146
#2	SIU13/14



Indiana University 17/18 Curette

#6	SIU17/186
#7	SIU17/187
#8, Gray	SIU17/188
#2	SIU17/18



Loma Linda 10/11 Curette	Standard	Rigid	Slim
#6	SLL10/116	SLL10/11R6	N/A
#7	SLL10/117	N/A	N/A
#8, Gray	SLL10/118	N/A	SLL10/11S8
#2	SLL10/112	N/A	N/A



Bunting 1/2 Curette

#6	SBN1/26
#2	SBN1/2



Bunting 5/6 Curette

#2	SBN5/6
----	------------------------



Ratcliff 3/4 Curette	
Harmony	SUCR3/4XE2
#9	SUCR3/49E2
#6	SUCR3/46
#8, Gray	SUCR3/48



#411/412 Curette	
#2	S411/412



#413/414 Curette	
#2	S413/414



ME-1 (10-11) Curette	
#2	SME1



Mallery 1/2 Curette	
#2	SMAL1/2

Rule 3/4 Curette	
#2	SR3/4

Goldman 31/31 Curette	
#2	SHG31/32



Sickle Scalers

Sickle Scalers have two cutting edges that meet at a pointed tip. The back surface is pointed giving a triangular cross-section. The blade can be straight or curved. Sickle Scalers are designed primarily for the removal of supragingival debris.

Nevi™ Scalers

Nevi™ Scalers feature unique angles and distinct characteristics to allow for improved access and reach into interproximal areas.




Nevi™ 1 Anterior Sickle Scaler

Harmony	SCNEVI1XE2
#9	SCNEVI19E2
#C8, Lavender	SCNEVI1C8E2
#6	SCNEVI1
#8, Gray	SCNEVI18
#7	SCNEVI17



Nevi™ 1/H5 Anterior Sickle/Discoild Scaler

Harmony	SCNEVI1XE2
#9	SCNEVI19E2
#C8, Lavender	SCNEVI1C8E2
#6	SCNEVI1
#8, Gray	SCNEVI18
#7	SCNEVI17



Nevi™ 2 Posterior Sickle Scaler

Harmony	SCNEVI2XE2
#9	SCNEVI29E2
#C8, Red	SCNEVI2C8E2
#6	SCNEVI2
#7	SCNEVI27
#8, Gray	SCNEVI28



Nevi™ 3 Posterior Sickle Scaler

Harmony	SCNEVI3XE2
#9	SCNEVI39E2
#C8, Red	SCNEVI3C8E2



Nevi™ 4 Posterior Sickle Scaler

Harmony	SCNEVI4XE2
#9	SCNEVI49E2
#C8, Red	SCNEVI4C8E2
#6	SCNEVI46



Nevi™ 4 Slim Posterior Sickle Scaler

Harmony	SCNEVI4LXE2
#9	SCNEVI4L9E2
#C8, Red	SCNEVI4LC8E2



Anterior Sickle Scalers

Anterior Sickle Scalers have a straight or slightly angled shank. The blades can be curved, straight, or combined.



Learn more about our
[Anterior Sickle Scalers](#).

Hand Scaling



H6/7 Sickle Scaler	
Harmony	SH6/7XE2
#9	SH6/79E2
#C8, Lavender	SH6/7C8E2
#6	SH6/76
#7	SH6/77
#8, Gray	SH6/78
#2	SH6/7



H5/33 Sickle Scaler	
Harmony	SH5/33XE2
#9	SH5/339E2
#C8, Lavender	SH5/33C8E2
#6	SH5/336
#7	SH5/337
#8, Gray	SH5/338
#2	SH5/33



U15/30 Towner-Jacquette Sickle Scaler	
Harmony	SU15/30XE2
#9	SU15/309E2
#C8, Purple	SU15/30C8E2
#6	SU15/306
#7	SU15/307
#8, Gray	SU15/308
#2	SU15/30



U15/33 Towner-Jacquette Sickle Scaler	
Harmony	SU15/33XE2
#9	SU15/339E2
#C8, Purple	SU15/33C8E2
#6	SU15/336
#7	SU15/337
#8, Gray	SU15/338
#2	SU15/33



H5/U15 Hygienist-Towner Sickle Scaler	
#6	SH5/U156
#2	SH5/U15



U15 Towner Sickle, Single-End	
#6	SU156
#7	SU157
#8, Gray	SU158
#24	SU151



Goldman 6/7 Sickle Scaler

#6	SHG6/76
#7	SHG6/77
#8, Gray	SHG6/78
#2	SHG6/7



Nebraska 128 Scaler

Harmony	SN128DEXE2
#9	SN128DE9E2
#C8, Lavender	SN128DEC8E2
#6	SN1286
#7	SN1287
#8, Gray	SN1288
#24	SN128



Morse 0/00 Scaler

Harmony	SM0/00XE2
#9	SM0/009E2
#C8, Lavender	SM0/00C8E2
#6	SM0/006
#7	SM0/007
#8, Gray	SM0/008
#2	SM0/00



USC 128 Scaler

#6	SUSC1286
#24	SUSC128



#1 Springfield Scaler, Single-End

#6	SSP16
#24	SSP1



Younger-Good 15 Scaler, Single-End

#6	SYG156
#1	SYG15



#1 Springfield Scaler, Single-End

#24	SWH2
-----	----------------------

Posterior Sickle Scalers

Posterior Sickle Scalers are designed with contra-angled shanks for areas difficult to access, such as posterior interproximal surfaces. The blades can be curved or straight.



Learn more about our [Posterior Sickle Scalers](#).



204 Posterior Sickle Scaler

#6	S2046
#2	S204



204S Posterior Sickle Scaler

Harmony	S204SXE2
#9	S204S9E2
#C8, Red	S204SC8E2
#6	S204S6
#7	S204S7
#8, Gray	S204S8
#2	S204S




204SD Posterior Sickle Scaler

Harmony	S204SDXE2
#9	S204SD9E2
#C8, Red	S204SDC8E2
#6	S204SD6
#7	S204SD7
#8, Gray	S204SD8
#2	S204SD2



204IFUW Sickle Scaler	
Harmony	SIUFW204XE2
#9	SIUFW2049E2
#C8, Red	SIUFW204C8E2
#6	SIUFW2046
#7	SIUFW2047
#8, Gray	SIUFW2048
#2	SIUFW204




Springfield 2/3 Scaler	
#6	SSP2/36
#2	SSP2/3



#135 Scaler	
Harmony	SN135XE2
#9	SN1359E2
#C8, Red	SN135C8E2
#6	SN1356
#7	SN1357
#8, Gray	SN1358



Taylor T2/T3 Scaler	
Harmony	SN135XE2
#9	SN1359E2
#7	SN1357
#8, Gray	SN1358



Remington 3/4 Scaler	
#6	SRM3/46
#2	SRM3/4



#129 Scaler	
#6	S1296
#2	S129




Jacquette Sickle Scalers

Jacquette Sickle Scalers have two cutting edges on a straight blade that end at a sharp point. The back surface of the blade is pointed giving a triangular cross-section.



30/33 Standard/Mini Jacquette Scaler

Harmony	SJ30/33XE2
#9	SJ30/339E2
#C8, Lavender	SJ30/33C8E2
#6	SJ30/336
#7	SJ30/337
#8, Gray	SJ30/338
#2	SJ30/33



31/32 Standard Jacquette Scaler

Harmony	SJ31/32XE2
#9	SJ31/329E2
#C8, Red	SJ31/32C8E2
#6	SJ31/326
#7	SJ31/327
#8, Gray	SJ31/328
#2	SJ31/32



34/35 Mini Jacquette Scaler

Harmony	SJ34/35XE2
#9	SJ34/359E2
#C8, Red	SJ34/35C8E2
#6	SJ34/356
#7	SJ34/357
#8, Gray	SJ34/358
#2	SJ34/35



H3/H4 Jacquette Scaler

#6	SH3/46
#2	SH3/4



14/15 Posterior Jacquette Scaler

#6	SJ14/156
#2	SJ14/15



Specialty Scalers



CI #1 Scaler

#6	SCI16
#24	SCI1



CI 2/3 DE Scaler

#6	SCI2/36
#7	SCI2/37
#2	SCI2/3



Crane-Kaplan #6 Scaler

#6	SCK66
#2	SCK6



Ball #1 Scaler, Single-End

#24	SB1
-----	------------



Ball #4 Scaler, Single-End

#24	SB4
-----	------------



Ball #2/3 Scaler

#6	SB2/36
#2	SB2/3



Ball #5/6 Scaler

#2	SB5/6
----	--------------



Scaler and Curette Sets

We offer a wide variety of Scaler and Curette Sets including McCall, Goldman-Fox, and Cattoni. All scalers and curettes included in these sets are available in multiple handle options.

McCall Scalers and Curettes

McCall Scalers and Curettes are posterior instruments featuring dual cutting edges and available with either a rounded or pointed toe.



McCall 11/12 Scaler	
#6	SM11/126
#7	SM11/127
#8, Gray	SM11/128
#2	SM11/12



McCall 11/12A Scaler	
#6	SM11/12A6
#2	SM11/12A



McCall 11/12 Thin Scaler	
#6	S11/12S6



McCall 13/14 Curette	
Harmony	SM13/14SXE2
#9	SM13/149E2
#C8, Lavender	SM13/14C8E2
#6	SM13/146
#7	SM13/147
#8, Gray	SM13/148
#2	SM13/14



McCall 13/14 Pointed Curette	
Harmony	SM13/14SXE2
#9	SM13/14S9E2
#C8, Lavender	SM13/14SC8E2
#6	SM13/14S6
#8	SM13/14S8
#7	SM13/14S7
#2	SM13/14S


McCall 13/14 Pointed Curette, Extended Blade

#6 **SMT13/146**
McCall 17/18 Pointed Curette

#9 **SM17/18PT9E2**

#6 **SM17/18PT6**

#2 **SM17/18PT**

McCall 17/18 Curette

Harmony **SM17/18XE2**

#9 **SM17/189E2**

#C8, Lavender **SM17/18C8E2**

#6 **SM17/186**

#7 **SM17/187**

#8, Gray **SM17/188**

#2 **SM17/18**

Goldman-Fox Scalers and Curettes

Goldman-Fox Scalers have robust, elongated shanks for the removal of moderate to heavy calculus deposits.


Goldman-Fox #1 Universal Scaler

#6 **SGF16**
Goldman-Fox #2 Anterior Scaler

#6 **SGF26**

Goldman-Fox #3 Posterior Scaler

#6 **SGF36**
Goldman-Fox #4 Anterior Scaler

#6 **SGF46**

Goldman-Fox #21 Scaler/Spoon

#6 **SGF216**

Goldman-Fox/Jacquette Mini Scaler

#6 **SGF21S/33**



Cattoni Scalers

Cattoni Scalers are sets designed for optimal removal of moderate to heavy supra and subgingival calculus.



Cattoni 103/106 Scaler

#6 [S103/1066](#)



Cattoni 104/105 Scaler

#6 [S104/1056](#)



Cattoni 107/108 Scaler

#6 [S107/1086](#)

#7 [S107/1087](#)

#8 [S107/1088](#)



Cattoni 107/108 Slim Scaler

#2 [S107/108S](#)



Files and Hoes

Files and Hoes are specialty instruments used to dislodge or crush and remove heavy calculus deposits.

Hoe Scalers

Hoe Scalers are used to dislodge heavy supragingival calculus using a pull stroke.



Orban 6/7 Mesial-Distal Hoe Scaler

#6	SO6/76
#2	SO6/7



Orban 8/9 Buccal-Lingual Hoe Scaler

#6	SO8/96
#2	SO8/9



Centennial Anterior Sickle/Hoe Scaler

#9	SDH1009E2
#C8, Lavender	SDH100C8E2



File Scalers

File Scalers have multiple cutting edges and are used to crush and remove heavy hard deposits using a pull stroke.

Hand Scaling



Orban 10/11 Buccal-Lingual File	
#6	FO10/116
#2	FO10/11



Orban 12/13 Mesial-Distal File	
#6	FO12/136
#2	FO12/13



Orban 10/11 Buccal-Lingual File	
#6	FH3/76
#2	FH3/7



Hirschfeld 5/11 Mesial-Distal File	
#6	FH5/116
#2	FH5/11



Hirschfeld 9/10 Anterior File	
#6	FH9/106
#2	FH9/10



Chisels

Chisels are designed to dislodge bridges of calculus, using a push stroke, on anterior mandibular teeth.



Cushing #61 Chisel, Single-End

#24

SCU61



Cushing #62 Chisel, Single-End

#24

SCU62



Darby-Perry 5K/6K Chisel

#2

SDP5/6K



Darby-Perry 9K/11K Chisel

#2

SDP9/11K

DiamondTec™ File Scalers

DiamondTec™ File Scalers offer 360° of diamond coated coverage. The absence of blade edges allows use of all sides of the working end. Ideal for furcations, root depressions, mesial/distal line angles and developmental grooves.



DiamondTec™ Mesial Distal File Scaler

#2

SDCM/D7



DiamondTec™ Nabors File Scaler

#2

SDCN7



Empowerment Kit

The perfect set of scalers and curettes to enhance the effectiveness of your power scaling treatment protocols.



Learn more about the [Empowerment Kit](#).

Empowerment Kit	POWCOMPKIT
Small Cassette - 8 Instrument, Purple	IMN5084
Gracey 11/12 Thin Curette, #9 Handle, Orange	SG11/12PC9E2
Nevi™ 4 Thin Scaler, #9 Handle, Red	SCNEVI4LPC9E2
Gracey 1/2 Thin Curette, #9 Handle, Grey	SG1/2PC9E2
Younger-Good 7/8 Thin Scaler, Black	SYG7/8PC9E2
Gracey 13/14 Thin Curette, #9 Handle, Blue	SG13/14PC9E2
Nevi™ 1 Thin Scaler, #9 Handle, Purple	SCNEVI1LPC9E2
Gracey 7/8 Thin Curette, #9 Handle, Green	SG7/8PC9E2
#11/12 ODU DE Explorer	EXD11/126



Empowerment Scalers and Curettes

Empowerment Scalers and Curettes feature blades that are an average of 10% thinner for easier access to tight interproximal and subgingival spaces. EverEdge™ 2.0 Empowerment Scalers and Curettes are color-coded for easy identification.



Nevi™ 1 Thin Scaler
#9, Purple [SCNEVI1LPC9E2](#)



Nevi™ 4 Thin Scaler
#9, Red [SCNEVI4LPC9E2](#)



Younger-Good 7/8 Thin Scaler

#9, Black

[SYG7/8PC9E2](#)



Gracey 1/2 Thin Curette

#9, Gray

[SG1/2PC9E2](#)



Gracey 1/2 After Five Thin Curette

#9, Gray

[SRPG1/2PC9E2](#)



Gracey 7/8 Thin Curette

#9, Green

[SG7/8PC9E2](#)



Gracey 11/12 Thin Curette

#9, Orange

[SG11/12PC9E2](#)



Gracey 13/14 Thin Curette

#9, Blue

[SG13/14PC9E2](#)



Gracey 15/16 After Five Thin Curette

#9, Orange

[SRPG15/16PC9E2](#)



Langer 3/4 Thin Curette

#9, Black

[SL3/4PC9E2](#)